## en MENTOR™ BECKER™ EXPANDERS/BREAST IMPLANTS

August 2025,

LAB100947443v2

|    |  |    |   |
|----|--|----|---|
| ar | MENTOR™ BECKER™ موسعات/غرسات الثدي من شركة                 | it | Espansori/protesi mammarie MENTOR™ BECKER™      |
| bg | Експандери/гърдни импланти MENTOR™ BECKER™                 | lt | „MENTOR™ BECKER™“ plėtikliai / krūtų implantai  |
| cz | Expandéry / prsní implantáty MENTOR™ BECKER™               | lv | MENTOR™ BECKER™ paplašinātāji/krūšu implanti    |
| da | MENTOR™ BECKER™ vævsudvidere/brystimplantater              | nl | MENTOR™ BECKER™-expanders/-borstimplantaten     |
| de | MENTOR™ BECKER™ – Expander/Brustimplantate                 | no | MENTOR™ BECKER™-vevsutvidere/brystimplantater   |
| el | Διατατήρες/εμφυτεύματα μαστού MENTOR™ BECKER™              | pl | Ekspandery/implanty piersi MENTOR™ BECKER™      |
| es | Dispositivos expansores/implantes mamarios MENTOR™ BECKER™ | pt | Expansores/implantes mamários MENTOR™ BECKER™   |
| et | MENTOR™ BECKER™ ekspanderid/rinnaimplantaadid              | ro | Dilatatoare/implanturi mamare MENTOR™ BECKER™   |
| fi | MENTOR™ BECKER™ -venyttimet/-rintaimplantit                | sk | Expandéry/prsné implantáty MENTOR™ BECKER™      |
| fr | Expanseurs/prothèses mammaires MENTOR™ BECKER™             | sl | Razširjevalniki/prsni vsadki MENTOR™ BECKER™    |
| hr | Ekspanderi / implantati za dojke MENTOR™ BECKER™           | sr | MENTOR™ ekspander/implantati dojke BECKER™      |
| hu | MENTOR™ BECKER™ tágitók/mellimplantátumok                  | sv | MENTOR™ BECKER™ expanderare/bröstimplantat      |
| is | MENTOR™ BECKER™ stækkari/brjóstaígræði                     | tr | MENTOR™ BECKER™ Genişleticiler/Meme İmplantları |

---

Only to be **implanted** in an appropriate medical environment by appropriately **trained medical doctors** who are **qualified** or **accredited** in accordance with **national law**

---

ENGLISH

# MENTOR™ BECKER™ EXPANDERS/BREAST IMPLANTS

## PATIENT INFORMATION LEAFLET

August 2025, LAB100947443v2

Breast implantation is an elective procedure and the patient must be well counseled on the risk-benefit relationship.

The surgeon should provide each prospective patient with the following information. This information may also be found in the patient information website

<http://ic.jnjmedicaldevices.com>.

### 1. Device Description

The BECKER™ Expander/Breast Implants is constructed of silicone. The device combines the advantages of tissue expanders with the feel of a gel breast implant. The devices are indicated for temporary overexpansion.

MENTOR™ offers two types of options: SILTEX™ Round BECKER™ 25 and SILTEX™ CONTOUR PROFILE™ BECKER™.

The BECKER™ 25 uses 25 % gel volume. The BECKER™ 25 also uses Cohesive I gel filling.

The BECKER™ 35 uses 35 % gel volume. The BECKER™ 35 also uses Cohesive II gel filling.

BECKER™ Expander/Breast Implant Family (postoperative adjustment of volume)

### 2. Device Materials

The BECKER™ Expander/Breast Implants have been tested for safe use. The Expander/Breast implant that the patient comes in contact with, are 100 % medical implant grade silicone. The TRUE-LOCK™ Connector that the patient comes into contact with is made of polysulfone.

| Device Materials                       | Implant Component                         | Patient Contacting | Component weight in grams of largest CP device <sup>1</sup> |
|--|---|--------------------|---|
| Silicone Elastomer                     | Outer Silicone Shell (with texture layer) | Yes                | ~40grams <sup>2</sup>                                       |
| Silicone Elastomer                     | Patch assembly                            | Yes                | ~4.2 grams <sup>2</sup>                                     |
| Silicone Elastomer                     | Outer valve                               | Yes                | ~ 3 grams <sup>2</sup>                                      |
| Silicone Elastomer with Titanium Oxide | Standard and Micro Dome                   | Yes                | Standard Dome (6 grams)<br>Micro Dome (<2 grams)            |
| Polysulfone                            | True Lock Connector Option                | Yes                | < 0.3 grams   |

<sup>1</sup> This represents the largest BECKER™ CP Device (Cat. No. 324-1605) which presents the heaviest components. It is approximately equal to largest Round device (354-8000D)

<sup>2</sup> Actual values will have a range, as component thicknesses can vary (per specification)

### 3. Information for Safe Use

#### 3.1. General Safety Instructions for patient

Breast implants are for women 18 and older who want breast augmentation or reconstruction.

Breast implants surgery is not for women who have any of the following conditions:

- Pregnancy or nursing mothers

- Autoimmune disease
- Conditions that would increase surgical risk\*
- A history of sensitivity to foreign materials
- Repeated attempts and failures at breast augmentation or reconstruction
- An unwillingness to undergo any further surgery for revision
- Unrealistic expectations,
  - inappropriate attitude or motivation
  - lack of understanding of the risks involved with the surgery and implants.
- Untreated breast cancer or premalignant breast disease\*\*
- A recent history of breast abscess
- Painful breast inflammation or breast tumor
- Uncontrolled diabetes or other disease that will impact healing
- Unsuitable tissue or compromised vascularity
- Any condition that would lead to significant postoperative complications

\* Discuss with your surgeon any possible conditions.

\*\* Speak to your breast surgeon and oncologist about your cancer treatment plan.

### 3.2. Postoperative Warnings, Precautions, or Measures

#### Breast Self-Examination Measures

You should learn to perform breast self-examinations monthly. You should also be shown how to distinguish the implant from your breast tissue.

#### Warnings

- Silicone gel can leak or “bleed” through the undamaged shell into the capsule and adjacent breast tissue. The long-term effects of such “bleed” are unknown. You should be aware of this potentiality.
- The use of microwave diathermy is NOT recommended. It may cause negative side effects to your surgical area.
- You should be aware that any abnormal stress or trauma to the breast could result in damage to the implant.
- MENTOR™ strongly recommends AGAINST the treatment of capsule firmness by forceful external stress (such as closed capsulotomy). It has been reported to cause patient injury and/or implant damage.
- Breast implants may interfere with breastfeeding.
- Patients with connective tissue disease may be more likely to experience some complications while healing.
- Some patients with breast implants have reported neurological and/or rheumatological symptoms.
- As with any surgery, smoking may interfere with the healing process after breast implant surgery

#### Precautions

Advise your healthcare provider about your existing implants if

- You have an autoimmune disease
- You have a weakened immune system
- You are currently taking drugs that weaken the body’s natural resistance to disease

Advise your Healthcare Provider if you are planning any of the following:

- Surgery in the breast area
- Chemotherapy
- Radiation\*
- Mammography\*\*
- MRI\*\*\*
- Breastfeeding\*\*\*\*
- Reoperation\*\*\*\*\*

\*Breast implants can possibly be damaged by radiation treatment. Radiation therapy may increase the chance of capsular contracture and/or change how your breasts look.

\*\* Breast implants can interfere with the early detection of breast cancer/calcium deposits

\*\*\* MENTOR™ BECKER™ Expander/Breast Implant with injection domes and stainless-steel connection system is MR Conditional.

\*\*\*\* Breast implants may interfere with breastfeeding

\*\*\*\*\*Breast implants are not lifetime devices. Additional surgery to your breast and/or implant will be likely in the future.

**Risk and Complications**

| Risks                     | Measures for Reducing Risks  |
|---------------------------|--|
| Infection*                | <ul style="list-style-type: none"> <li>• MENTOR™ device and packaging are sterile when shipped.</li> <li>• Precautioning that infection prior to implant should be treated.</li> <li>• Contraindication for implanting where there is an active infection</li> <li>• Using aseptic technique when filling the implant with saline.</li> </ul>  |
| Hematoma formation        | <ul style="list-style-type: none"> <li>• Warning to the surgeon that careful stoppage of blood flow is important to prevent hematoma.</li> <li>• Warning to the surgeon that recommends that if excessive bleeding persist, the bleeding should be controlled. Then proceed with the surgery.</li> </ul>   |
| Delayed wound healing     | <ul style="list-style-type: none"> <li>• Not implanting where there is a history of compromised wound healing.</li> <li>• Not implanting when there is uncontrolled diabetes or other disease which impacts healing.</li> <li>• Warning that smoking may interfere with the healing process.</li> <li>• * See measures for reducing the risk of infection.</li> </ul>  |
| Capsular Contracture (CC) | <ul style="list-style-type: none"> <li>• Warning that radiation therapy may cause Capsular Contracture.</li> <li>• Warning that Capsular Contracture may occur around any implant placed in contact with soft tissues.</li> <li>• Local complications or infections after surgery may increase the likelihood of Capsular Contracture.</li> </ul>  |
| Chest wall deformity      | <p><b>Implant Selection</b> guidelines:</p> <ul style="list-style-type: none"> <li>• The implant should not be too small or too large in comparison to the patient’s chest-wall dimensions.</li> <li>• Available tissue must provide adequate coverage of the implant.</li> <li>• Submuscular placement of the implant may be preferable in patients with:                             <ul style="list-style-type: none"> <li>➢ Thin</li> <li>➢ Poor quality tissue.</li> </ul> </li> <li>• A well-defined dry pocket of adequate size and symmetry must be created to allow the implant to be placed flat on a smooth surface.</li> </ul> |
| Enlarged lymph nodes      | <ul style="list-style-type: none"> <li>• Precautioning that contamination may occur when adding or removing fluid from the implant.</li> <li>• Using aseptic techniques.</li> <li>• Rinsing clean surgical instruments and gloves to remove impurities</li> <li>• * See measures for reducing the risk of infection.</li> </ul>  |
| Rupture                   | <ul style="list-style-type: none"> <li>• Implant meets ISO safety standards.</li> <li>• Warning to the surgeon that extreme care should be taken not to damage the prosthesis with surgical instruments.</li> <li>• Warning the patient that any abnormal stress or trauma to the breast could result in damage to the implant.</li> </ul>   |

| Risks  | Measures for Reducing Risks  |
|--|--|
| Changes in Nipple and Breast Sensation/Breast Pain | <p>Any surgery that undermines the skin of the breast can result in hypersensitivity or loss of sensation in the nipple-areolar region. These changes can vary in degree and may be temporary or permanent. Changes in nipple/breast sensation may, on occasion, affect sexual response or comfort while nursing. These changes are believed to be a result of nerve damage or stretching of the nerves from the surgery. There is no specific treatment for this condition.</p> <p>Surgeon considers the size of the implant and the firmer nature and higher profile of the SILTEX™ shell when choosing optimum incision size and surgical approach.</p> <p><b>Implant Selection</b> guidelines:</p> <ul style="list-style-type: none"> <li>• The implant should not be too small or too large in comparison to the patient’s chest-wall dimensions.</li> <li>• Available tissue must provide adequate coverage of the implant.</li> <li>• Submuscular placement of the implant may be preferable in patients with <ul style="list-style-type: none"> <li>➢ Thin</li> <li>➢ Poor quality tissue.</li> </ul> </li> <li>• A well-defined dry pocket of adequate size and symmetry must be created to allow the implant to be placed flat on a smooth surface.</li> <li>• * See measures for reducing the risk of infection.</li> </ul> |
| Possible reaction to silicone                      | <ul style="list-style-type: none"> <li>• Implant meets safety standards.</li> <li>• Questioning patients on: <ul style="list-style-type: none"> <li>➢ medical history involving allergic reactions to implant materials or</li> <li>➢ filling agents.</li> </ul> </li> </ul>   |
| Gel Bleed  | <ul style="list-style-type: none"> <li>• Implant meets safety standards.</li> <li>• Warning the surgeon in the IFU of the risk of rupture: that <i>extreme care should be taken not to damage the prosthesis with surgical instruments</i></li> <li>• Warning the patient that any abnormal stress or trauma to the breast could result in damage to the implant.</li> </ul>   |
| Granulomas   | <ul style="list-style-type: none"> <li>• Implant meets safety standards.</li> <li>• Warning the surgeon in the IFU of the risk of rupture: that extreme care should be taken not to damage the prosthesis with surgical instruments.</li> <li>• Warning the patient that any abnormal stress or trauma to the breast could result in damage to the implant.</li> </ul>   |
| Extrusion of the Implant                           | <ul style="list-style-type: none"> <li>• Warning not to use microwave diathermy in patients with breast implants.</li> <li>• Warning that the incidence of extrusion of the prosthesis had been shown to increase when placed in injured areas that are: <ul style="list-style-type: none"> <li>➢ scarred</li> <li>➢ heavily irradiated</li> <li>➢ burned tissue</li> <li>➢ crushed bone areas or</li> <li>➢ severely reduced by surgery.</li> </ul> </li> </ul>   |
| Wrinkling  | <p>Precautioning that surgeon use current techniques for</p> <ul style="list-style-type: none"> <li>• Measuring the patients</li> <li>• sizing the implant and</li> <li>• conducting the surgery.</li> </ul> <p>The surgeon does not underfill or underfill the device.</p>  |
| Thrombosed veins                                   | <p>Warning the surgeon that excessive inflation of the device may result in tissue necrosis/thrombosis</p> <p>Warning that thrombosis may result from the placement of foreign objects in the body.</p>  |

| Risks                                   | Measures for Reducing Risks  |
|---|--|
| Hypertrophic scarring                   | IFU advising the surgeon that: <ul style="list-style-type: none"> <li>• Hypertrophic scars are generally confined to the original site.</li> <li>• Hypertrophic scars respond well to scar revision treatment.                             <ul style="list-style-type: none"> <li>➤ This may include steroid injections to break down the collagen.</li> <li>➤ This may also include surgery to revise the position, direction or line of the scar.</li> </ul> </li> </ul> |
| Anaplastic Large Cell Lymphoma BIA-ALCL | The IFU describes that you should seek a plastic surgeon to evaluate for BIA-ALCL if you experience: <ul style="list-style-type: none"> <li>• a seroma in the breast more than 12 month after implantation</li> <li>• a new breast mass, or</li> <li>• disease of the lymph nodes.</li> </ul>  |

If you are experiencing any breast pain or discomfort, contact your breast surgeon. Discuss any questions, concerns, or potential side effects with your physician.

**Breast Implant-Associated Anaplastic Large Cell Lymphoma (BIA-ALCL)**

Breast implants have been associated with the development of a cancer of the immune system (BIA-ALCL).

This cancer occurs more commonly in patients with textured breast implants than smooth implants, although rates are not well defined.

If you have breast implants, you have a very small, but increased risk of developing breast implant associated anaplastic large cell lymphoma, or BIA-ALCL. BIA-ALCL is not breast cancer.

It is a rare type of non-Hodgkin's lymphoma (cancer of the immune system). In most cases, BIA-ALCL is found in the scar tissue and fluid near the implant.

However, in some cases, it can spread throughout the body. In the cases that have spread beyond the scar tissue and fluid near the implant, rare cases of death have been reported.

You should seek evaluation by a plastic surgeon to evaluate for BIA-ALCL if you experience delayed onset seroma in the breast (delayed means developing more than 12 months after implantation), a new breast mass or lymphadenopathy (disease of the lymph nodes).

If you have breast implants, you should monitor them and follow your routine medical care. It is not necessary to remove your breast implants if you have no symptoms and you have not been diagnosed with BIA-ALCL.

**4. Users and Consumer Exclusions**

The intended users are surgeons who are familiar with the currently available techniques for performing breast tissue expansion and breast implantation surgery. **Persons that do not meet this criterion are excluded from the use of the device.**

Consumers (patients) who may receive this device are females who will undergo breast augmentation (non-medical purpose) or breast reconstruction (medical purpose), or who require tissue expansion and replacement of missing breast tissue.

The use of this prosthesis is contraindicated in patients who have any of the following conditions:

- An active infection anywhere in the body.
- Recent history of breast abscess.
- Diffuse painful cystic mastitis or breast tumor.
- A history of compromised wound healing.
- A compromised immune system.
- Persistent or recurrent breast cancer.
- Pregnancy or nursing mothers.
- Lupus (e.g., SLE and DLE).
- Scleroderma (e.g., progressive systemic sclerosis).
- Uncontrolled diabetes or other disease which impacts healing.
- Unsuitable tissue due to radiation damage on the chest wall, tight thoracic skin grafts, or radical resection of the pectoralis major muscle.
- Compromised vascularity.
- A history of sensitivity to foreign materials or repeated attempts and failures at breast augmentation or reconstruction.
- Any anatomic or physiologic abnormality which could lead to significant postoperative complications

- An unwillingness to undergo any further surgery for revision.
- Psychological instability such as inappropriate attitude or motivation, or a lack of understanding of the risks involved with the surgical procedure and prosthesis.

**Persons that are contraindicated are excluded from the use of the device.**

### 5. Expected Device Lifetime and Follow-up

Breast implantation may not be a one-time surgery. After the surgery, doctor visits and additional surgery may be needed. Implant removal with or without replacement may be needed as well. Breast implants may not last a lifetime. A rupture can occur at any time following implantation. Studies have indicated that the gel breast implant can have an expected average lifetime of over 10 years.

You should have periodic imaging (e.g., MRI, ultrasound) of your implants to screen for implant rupture. This does not replace additional imaging that may be required, such as screening mammography for breast cancer. Even if you have no symptoms, you should have your first ultrasound or MRI at 5-6 years after your initial implant surgery. Then every 2-3 years thereafter.

If you have symptoms at any time or uncertain ultrasound results for breast implant rupture, an MRI is recommended.

Additional doctor visits would be needed to monitor if rupture of the implant has occurred. Speak to your surgeon about regular basis monitoring screening.

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it directly to your healthcare provider or MENTOR™ (Mentor-PE@its.jnj.com or 1-866-250-5115)

### 6. Implant ID Card

You should have received an International Implant Card from your surgeon. This card will have the catalog and serial number of your implants, along with other information. Carry the card with you and show it to doctors or other healthcare providers when you visit them. It will help them treat you appropriately and protect your breast implants during any medical treatment you need in the future. Your doctor should keep a copy of the Implant ID Card with your medical records.

### 7. Patient Information Website

You can find any updates about the device on this website. Your implant card will contain information about your specific device. This can also be found in your medical record

### 8. Volume and size of the device


**Information on the volume and size of the breast implant that you will receive is provided on the implant card following the surgery.**

### 9. Harmonized Standards and Common Specifications Applied

- COMMISSION IMPLEMENTING REGULATION (EU) 2022/2346 of 1 December 2022: laying down common specifications for the groups of products without an intended medical purpose listed in Annex XVI to Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices
- EN ISO 13485 Medical device – Quality management system – Requirements for regulatory purposes
- EN ISO 14971 Medical devices – Application of risk management to medical devices
- EN ISO 14630 Non-active surgical implants – General requirements
- EN 62366 Medical Devices- Application of usability engineering to medical devices
- EN ISO 10993-1 Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management system
- EN 556-1 Sterilization of medical devices – Requirements for medical devices to be designated “STERILE” – Part 1: Requirements for terminally sterilized medical devices
- EN ISO 11737-1 Sterilization of medical devices - Microbiological methods - Part 1: Determination of a population of microorganisms on products
- EN ISO 11737-2 Sterilization of medical devices - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
- EN ISO 11607-1 Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems

- EN ISO 11607-2 Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes
- EN ISO 14644-1 Cleanrooms and associated controlled environments - Part 1: Classification of air cleanliness
- EN ISO 14644-2 Cleanrooms and associated controlled environments - Part 2: Specifications for testing and monitoring to prove continued compliance with ISO 14644-1 (ISO 14644-2:2000)
- EN ISO 14644-5 Cleanrooms and associated controlled environments - Part 5: Operations (ISO 14644-5:2004)

## 10. Symbols Glossary

|   |  |   |                          |   |       |
|---|--|---|--------------------------|---|-------|
|  | Patient name or patient ID                                       |  | Catalogue number         |  | LEFT  |
|  | Date of Implantation   |  | Manufacturer             |  | RIGHT |
|  | Name and Address of implanting healthcare institution / provider |  | Unique Device Identifier |   |       |
|  | Information website for patients                                 |  | Serial number            |   |       |
|  | Device name  |  | MR Safe                  |   |       |

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© MENTOR™ Worldwide, LLC 2025
MENTOR™ Medical Systems B.V.
Johnson & Johnson Medical Ltd