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## Referenzen

1. Anaplastic Large Cell Lymphoma (ALCL) In Women with Breast Implants: Preliminary FDA Findings and Analyses: Center for Devices and Radiological Health, U.S. Food and Drug Administration; January 2011 [cited 2017 26 February]. Available from: <https://www.fda.gov/MedicalDevices/ProductsandMedicalProcedures/ImplantsandProsthetics/BreastImplants/ucm239996.htm>.
2. Anaplastic large cell lymphoma (ALCL). US Food and Drug Administration; 20 January 2016 [cited 2017 26 February]. Available from: <https://www.fda.gov/MedicalDevices/ProductsandMedicalProcedures/ImplantsandProsthetics/BreastImplants/ucm239995.htm>.
3. Lymphome Anaplasique à Grandes Cellules associé aux implants mammaires (LAGC-AIM): Point sur les investigations en cours- Point d'Information.: L'Agence nationale de sécurité du médicament et des produits de santé (ANSM). 06 July 2016 [cited 2016 06 July]. Available from: <http://ansm.sante.fr/S-informer/Points-d-information-Points-d-information/Lymphome-Anaplasique-a-Grandes-Cellules-associe-aux-implants-mammaires-LAGC-AIM-Point-sur-les-investigations-en-cours-Point-d-Information>.
4. Breast Implants: Expert Advisory Panel Advice on Association with Anaplastic Large Cell Lymphoma.: Australian Government Department of Health Therapeutic Goods Administration.; 20 December 2016 [cited 2017 26 February]. Available from: <https://www.tga.gov.au/alert/breast-implants>.
5. Breast Implant-Associated Anaplastic Large Cell Lymphoma (BIA-ALCL): US Food and Drug Administration; 22 March 2017 [cited 2017 22 March]. Available from: <https://www.fda.gov/MedicalDevices/ProductsandMedicalProcedures/ImplantsandProsthetics/BreastImplants/ucm239995.htm>.
6. Clemens MW, Medeiros LJ, Butler CE, Hunt KK, Fanale MA, Horwitz S, et al. Complete Surgical Excision Is Essential for the Management of Patients With Breast Implant-Associated Anaplastic Large-Cell Lymphoma. *J Clin Oncol.* 2016;34(2):160-8. Epub 2015/12/03. doi: 10.1200/jco.2015.63.3412. PubMed PMID: 26628470; PubMed Central PMCID: PMC4872006 online at <http://www.jco.org>. Author contributions are found at the end of this article.
7. Clemens MW, Horwitz SM. NCCN Consensus Guidelines for the Diagnosis and Management of Breast Implant-Associated Anaplastic Large Cell Lymphoma. *Aesthetic surgery journal.* 2017. Epub 2017/02/12. doi: 10.1093/asj/sjw259. PubMed PMID: 28184418.
8. ASPS/ASAPS Update Breast Implant-Associated Anaplastic Large Cell Lymphoma (BIA-ALCL) Quick Facts and FAQs: The American Society For Aesthetic Plastic Surgery & American Society of Plastic Surgeons. Available from: <http://www.surgery.org/downloads/blasts/BIA-ALCL/>.
9. Breast implants - report cases of anaplastic large cell lymphoma (ALCL): Medicines and Healthcare products Regulatory Agency; 10 July 2014 [cited 2017 26 February]. Available from: <https://www.gov.uk/drug-device-alerts/medical-device-alert-breast-implants-report-cases-of-anaplastic-large-cell-lymphoma-alcl>.
10. Information on BIA-ALCL: American Society of Plastic Surgeons; 2017 [cited 2017 26 February]. Available from: <https://www.plasticsurgery.org/patient-safety/information-on-bia-alcl>.
11. Joint ASPS & ASAPS Statement on Breast Implant-Associated ALCL 2016 [cited 2017 26 February]. Available from: <http://www.surgery.org/downloads/private/joint-asps-asaps-statement-on-breast-implant-associated-alcl.pdf>.
12. Frequently Asked Questions (FAQ): A Guide to Breast Implant-Associated Anaplastic Large Cell Lymphoma: American Society of Plastic Surgeons and The Plastic Surgery Foundation; 2016 [cited 2017 22 March]. Available from: <http://www.theapsf.org/Documents/Clinical/PROFILE/profile-faq.pdf>.
13. International Society of Aesthetic Plastic Surgery (ISAPS) 2017. Available from: <http://www.isaps.org>.
14. MemoryShape Post-Approval Cohort Study (formerly Contour Profile Gel Core Study) Final Clinical Study Report. Mentor Worldwide, LLC; 02 June 2015.
15. MemoryGel Core Gel Clinical Study Final Report. Mentor Worldwide, LLC; April 2013.
16. Mentor MemoryShape Post-Approval Continued Access Study (formerly Contour Profile Gel Continued Access Study), Final Report. October 2014.
17. Mentor MemoryGel Breast Implant Large Post Approval Study Re-Op Phase Annual Report. 17 June 2016.
18. Adjunct Study Final Report for Mentor's MemoryGel Silicone Gel-filled Breast Implants. 02 November 2012.
19. Mentor MemoryShape CPG Styles Study: A Study of the Safety of the Contour Profile Gel Breast Implants in Subjects who are Undergoing Primary Breast Augmentation, Primary Breast Reconstruction, or Revision, Final Clinical Study Report. October 2015.

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