

Burden and Epidemiology of Surgical Smoke Evacuation

Value Analysis Brief


EPIDEMIOLOGY AND CHARACTERIZATION OF SURGICAL SMOKE



Surgical smoke is formed when energy-generating devices (electrosurgery units, lasers, powered instruments) raise the intracellular temperature of tissue to at least 100°C (212°F), causing tissue vaporization.¹ Electrocautery smoke may pose potential health risks for the nearly one million surgical staff around the world.²

- A systematic review of potential health implications of surgical smoke for operating room staff reported that the size of particles found in all types of smoke across procedures ranged from 0.05 µm to larger than 25 µm.³ One study of laser corneal surgery noted a mean particle size of 0.22 µm (n = 98).⁴
- The size of an inhaled particle is the most important aspect in determining where it will be deposited in the respiratory tract.^{4,5}
- Particles of 5.0 µm or greater can be deposited on the walls of the nose and trachea. Smaller particles, less than 2.0 µm, may be deposited in the respiratory bronchioles and alveoli.^{4,5}

SIZE AND DEPOSITION OF SURGICAL SMOKE PARTICLES^a



SURGICAL SMOKE TRAVELS FAST AND BECOMES HIGHLY CONCENTRATED IN THE OR



Particles can travel
up to 18 m/s.^{7b}



The concentration of the particles
can rise from 60,000 particles/
cubic foot to over 1 million
particles/ cubic foot within
5 minutes of
electrocautery initiation.⁸



As particles get caught
in air currents, they can
become distributed
in the OR.⁸

Note: OR = operating room

^a Image adapted from Okoshi et al., 2014.⁶ ^b Velocity for smoke particles generated by laser surgery.

BIOLOGICAL AND CHEMICAL CONSTITUENTS OF SURGICAL SMOKE






- Surgical smoke can contain toxic gases, dead and live cellular material - including blood fragments - and, viruses.^c

MORE THAN 150 DIFFERENT CHEMICAL CONSTITUENTS WHICH MAY HAVE IMPACTS ON VARIOUS BODY SYSTEMS HAVE BEEN IDENTIFIED IN SURGICAL SMOKE.⁹

Respiratory	Eyes	CNS (e.g. headache, nausea)	Mutagenic
Acetaldehyde ^{5, 10, 11}	Acetaldehyde ^{5, 10, 11}	Acetylene ^{12, 13}	Acetaldehyde ^{5, 10, 11, 23, 24}
Acrolein ⁵	Acrolein ⁵	Acrylonitrile ^{11-13, 20, 21}	Acrolein ^{5, 23, 24}
Acetonitrile ^{5, 12, 13}	Acrylonitrile ^{11-13, 20, 21}	Benzene ^{5, 11-13, 15, 18, 20}	Acrylonitrile ^{11-13, 20, 21, 23, 24}
Cyclohexanone ¹⁴	Decane ^{14, 15}	Carbon monoxide ^{13, 21, 22}	Benzene ^{5, 11-13, 15, 18, 20, 23, 24}
Decane ^{14, 15}	Formaldehyde ^{5, 13}	Decane ^{14, 15}	Cyclohexanone ^{14, 23, 24}
Formaldehyde ^{5, 13}	Furfural ^{10, 12, 16}	Furfural ^{10, 12, 16}	Formaldehyde ^{5, 13, 23, 24}
Furfural ^{10, 12, 16}	Toluene ^{5, 11-13, 15, 16, 18, 19}	Propylene ¹⁷	Furfural ^{10, 12, 16, 23, 24}
PAH ⁵		Toluene ^{5, 11-13, 15, 16, 18, 19}	PAH ^{5, 23, 24}
Phenol ^{13, 17}			Styrene ^{5, 12, 13, 15, 18, 23, 24}
Pyridine ^{13, 17}			
Styrene ^{5, 12, 13, 15, 18}			
Toluene ^{5, 11-13, 15, 16, 18, 19}			
Xylene ^{5, 13, 15}			

Chemicals are organized according to where/how they have the greatest impacts

BACTERIA AND VIRUSES PREVIOUSLY IDENTIFIED IN SURGICAL SMOKE^{17, 25-27}

 0.010-0.300 µm Hepatitis B virus (HBV)	 0.045 µm Human papillomavirus (HPV)	 0.180 µm Human immunodeficiency virus (HIV)	 0.500 µm Mycobacterium tuberculosis	 0.300-1.500 µm Staphylococcus, Corynebacteriu, Neisseria
---	--	--	--	---

- Bacteria and viruses have been shown to survive and can be regrown from samples collected from surgical smoke.⁶
- Some viruses, such as HPV and HIV, are smaller than the mean particle size found in surgical smoke.^{25, 27}
- The small size of viruses and bacteria can easily pass through the commonly used surgical masks.²⁸

POTENTIAL IMPLICATIONS OF SURGICAL SMOKE

- Exposure to surgical smoke can cause both acute and chronic health effects ranging from eye, nose and throat irritation to emphysema, asthma or chronic bronchitis.^d
- A 2006 study reported on potential risks to staff of surgical smoke. These risks are consistent with reports from healthcare professionals and researchers during the past two decades.²⁰
- As early as 1988, scientists established a causal link between inhaling unfiltered surgical smoke and pulmonary changes, including alveolar congestion and emphysema.³⁰
- Surgical smoke decreases visibility of the laparoscopic surgical field, possibly resulting in procedure delays.³¹



Physical


- Particles that range in size from 0.5 to 5.0 μm are considered to be "lung damaging dust."¹²
- Smoke plumes may cause both acute (e.g., sore eyes, dermatitis) or chronic (e.g., asthma) health effects.^{5, 6, 29, 32-37}
- Common constituents of surgical smoke can also cause neurotoxic symptoms such as drowsiness, headaches, tremor, dizziness, and coma.^{5, 6, 13, 32, 36, 37}
- Smoke plumes can increase risk of pulmonary conditions. A study of surgical residents reported that several developed foreign body sensation (58%) and pharyngeal burning (22%) potentially associated with exposure to electrocautery smoke.^{6, 38}
- Other studies note there is a risk of emphysema, asthma, and chronic bronchitis with exposure to surgical smoke.^{6, 29}



Infection

- Blended current electrosurgery smoke can contain viable bacteria.³⁹
- Viral DNA has been discovered in surgical smoke or plume⁴⁰ and may lead to disease transmission.⁶
- When asked if they were concerned about transmission of infectious disease via surgical smoke, 76% of surveyed dermatology residents responded yes.⁴¹

TIMELINE OF SURGICAL SMOKE RESEARCH




Note: AORN = Association of perioperative Registered Nurses; CDC = Centers for Disease Control and Prevention; CO = Colorado; HPV = Human papillomavirus; LEV = Local Exhaust Ventilation; RI = Rhode Island; WHO = World Health Organization

^d (C2359) NIOSH Study Finds Healthcare Workers' Exposure to Surgical Smoke Still Common. Nov 3, 2015.^{29 †} Indicated by diamonds on the map below.

WIDESPREAD RECOMMENDATIONS BY GLOBAL AND LOCAL ORGANIZATIONS

- Surgical smoke evacuation is strongly recommended by regulatory bodies and industry associations alike, such as OSHA, MHRA, NIOSH, ECRI and AORN.^{5, 37, 42, 44, 46, 49, 50}

EXAMPLE RECOMMENDATIONS:



- MHRA recommended that robust smoke evacuation systems must be used in all theatres that use electrocautery and laser devices.³⁷



- NIOSH recommends evacuation and filtration of smoke produced by surgical procedures and specifies that a smoke evacuator or suction must be within two inches of the surgical site.⁴²




COMMON HISTORICALLY USED METHODS AND THEIR LIMITATIONS



Masks

- Standard surgical masks are ineffective at filtering many substances from surgical smoke.⁵ Most are designed to filter particles > 5 µm.²⁸
- Respirator masks (e.g. N95) are bulky, impede function, and cause discomfort;⁵¹ also need to be fitted.
- Masks worn too loosely or for too long are less effective.⁶



Wall Suction

- Pulls less than 5 cubic feet per meter; only effective in procedures that produce a small amount of smoke.⁵
- Must be used with an inline filter; or else smoke can begin to occlude the smoke particles suction line.⁵
- May be ineffective at removing smoke directly where it is generated.⁵²
- Noisy and disrupts communication between staff.⁵³

Note: ACORN = Australian College of Operating Room Nurses; AfPP = Association of Perioperative Practice; ASLMS = American Society for Laser Medicine and Surgery; AORN = Association of perioperative Registered Nurses; ANSI = American National Standards Institute; BOHS = British Occupational Hygiene Society; CCOHS = Canada's National Centre for Occupational Health and Safety information; CSA = Canadian Standards Association; ECRI = Emergency Care Research Institute; EU-OSHA = European Agency for Safety and Health at Work; IFPN = International Federation of Perioperative Nurses; JCAHO = Joint Commission on Accreditation of Healthcare Organizations; JMS = Japanese Medical Society; MHRA = Medicines and Healthcare products Regulatory Agency; NIOSH = National Institute for Occupation Safety and Health; ORNAC = Operating Room Nurses Association of Canada; OSHA = Occupational Safety and Health Administration; SEORNA = Swedish Operating Nurse Association; SESLHD = South Eastern Sydney Local Health District; Tramontini SLR = Tramontini Systematic Literature Review

HAND-HELD DEVICE FEATURES THAT DIFFER FROM HISTORICALLY USED METHODS

- Studies have shown that the further a smoke evacuation device is from the site of plume generation, the amount of smoke evacuated will decrease significantly, thus allowing residual plume to escape into the air.⁵⁴
 - Using a smoke evacuation pencil design allows plume to be evacuated at the tissue impact site through a vortex motion which promotes greater plume capture.⁵⁴
 - Given their smaller designs, hand-held smoke evacuation devices may not interfere with the operative field.²⁸
 - The MegaVac Plus™ smoke evacuator has a flow rate adjustable to at least 90 liters per minute (lpm) for the HIGH (OPEN) setting and a flow rate from 4 ± 1 lpm to at least 18 lpm for the LOW (LAP) setting,⁵⁵ enabling rapid smoke evacuation in procedures with different requirements.
 - Can sometimes be integrated right into the surgical tool, such as an electrosurgery or electrocautery pencil.⁵³
 - A study demonstrated that a hand-held device was able to capture 99% of surgical smoke when placed one inch from the source.⁵³
-

REFERENCES

1. Spruce L (2018) Back to Basics: Protection From Surgical Smoke: 1.2 www.aornjournal.org/content/cme. AORN J 108 (1): 24-32.
2. Memon AG, Naeem Z, Zaman A, Zahid F (2016) Occupational health related concerns among surgeons. *Int J Health Sci* 10 (2): 279-291.
3. Mowbray N, Ansell J, Warren N, Wall P, Torkington J (2013) Is surgical smoke harmful to theater staff? A systematic review. *Surg Endosc* 27 (9): 3100-3107.
4. Taravella MJ, Viega J, Luiszer F, Drexler J, Blackburn P et al. (2001) Respirable particles in the excimer laser plume. *J Cataract Refract Surg* 27 (4): 604-607.
5. Ulmer BC (2008) The hazards of surgical smoke. *AORN J* 87 (4): 721-734; quiz 735-738.
6. Okoshi K, Kobayashi K, Kinoshita K, Tomizawa Y, Hasegawa S et al. (2015) Health risks associated with exposure to surgical smoke for surgeons and operation room personnel. *Surg Today* 45 (8): 957-965.
7. Nicola J, Nicola E, Vieira R, Braile D, Tanabe M et al. (2002) Speed of particles ejected from animal skin by CO₂ laser pulses, measured by laser Doppler velocimetry. *Phys Med Biol* 47 (5): 847.
8. Brandon HJ, Young LV (1997) Characterization and removal of electrosurgical smoke. *Surg Serv Manag* 3 (3): 14-16.
9. Pierce JS, Lacey SE, Lippert JF, Lopez R, Franke JE (2011) Laser-generated air contaminants from medical laser applications: a state-of-the-science review of exposure characterization, health effects, and control. *J Occup Environ Hyg* 8 (7): 447-466.
10. Takahashi H, Yamasaki M, Hirota M, Miyazaki Y, Moon JH et al. (2013) Automatic smoke evacuation in laparoscopic surgery: a simplified method for objective evaluation. *Surg Endosc* 27 (8): 2980-2987.
11. Zhao C, Kim MK, Kim HJ, Lee SK, Chung YJ et al. (2013) Comparative safety analysis of surgical smoke from transurethral resection of the bladder tumors and transurethral resection of the prostate. *Urology* 82 (3): 744 e9-14.
12. Fan JK, Chan FS, Chu KM (2009) Surgical smoke. *Asian J Surg* 32 (4): 253-257.
13. Barrett WL, Garber SM (2003) Surgical smoke: a review of the literature. Is this just a lot of hot air? *Surg Endosc* 17 (6): 979-987.
14. Al Sahaf OS, Vega-Carrascal I, Cunningham FO, McGrath JP, Bloomfield FJ (2007) Chemical composition of smoke produced by high-frequency electrosurgery. *Ir J Med Sci* 176 (3): 229-232.
15. Choi SH, Kwon TG, Chung SK, Kim TH (2014) Surgical smoke may be a biohazard to surgeons performing laparoscopic surgery. *Surg Endosc* 28 (8): 2374-2380.
16. Hensman C, Baty D, Willis RG, Cuschieri A (1998) Chemical composition of smoke produced by high-frequency electrosurgery in a closed gaseous environment. An in vitro study. *Surg Endosc* 12 (8): 1017-1019.
17. Eickmann U, Falcy M, Fokuhl I, Rügger M, Bloch M et al. (2011) Surgical smoke: Risks and preventive measures. International Section of the ISSA on prevention of occupational risks in health services Hamburg, Germany 5-40.
18. Sagar PM, Meagher A, Sobczak S, Wolff BG (1996) Chemical composition and potential hazards of electrocautery smoke. *Br J Surg* 83 (12): 1792.
19. Lin YW, Fan SZ, Chang KH, Huang CS, Tang CS (2010) A novel inspection protocol to detect volatile compounds in breast surgery electrocautery smoke. *J Formos Med Assoc* 109 (7): 511-516.
20. Alp E, Bijl D, Bleichrodt RP, Hansson B, Voss A (2006) Surgical smoke and infection control. *J Hosp Infect* 62 (1): 1-5.
21. Gianella M, Hahnloser D, Rey JM, Sigrist MW (2014) Quantitative chemical analysis of surgical smoke generated during laparoscopic surgery with a vessel-sealing device. *Surg Innov* 21 (2): 170-179.
22. Beebe DS, Swica H, Carlson N, Palahniuk RJ, Goodale RL (1993) High levels of carbon monoxide are produced by electro-cautery of tissue during laparoscopicolecystectomy. *Anesth Analg* 77 (2): 338-341.
23. International Agency for Research on Cancer (IARC) (Web Page) Agents classified by the IARC monographs, vol. 1109. Updated July 30 2018. Available online at: <https://monographs.iarc.fr/agents-classified-by-the-iarc/>. Accessed: August 27 2018.
24. Health and Safety Executive (Web Page) EH40/2005 Workplace Exposure Limits. Updated Available online at: <https://books.hse.gov.uk/bookstore.asp?ACTION=BOOK&PRODUCTID=9780717664467>. Accessed: August 27 2018.
25. Sandelin A (2017) Operating Room Nurses' work environment due to Surgical Smoke. *EORNA Journal*.
26. Capizzi PJ, Clay RP, Battey MJ (1998) Microbiologic activity in laser resurfacing plume and debris. *Lasers Surg Med* 23 (3): 172-174.

REFERENCES

27. de Boorder T, Verdaasdonk R, Klaessens J (2007) The visualization of surgical smoke produced by energy delivery devices: significance and effectiveness of evacuation systems. In: *Thermal Treatment of Tissue: Energy Delivery and Assessment IV*, p 64400R: International Society for Optics and Photonics.
28. Bree K, Barnhill S, Rundell W (2017) The dangers of electrosurgical smoke to operating room personnel: a review. *Workplace Health Saf* 65 (11): 517-526.
29. Centers for Disease Control and Prevention (CDC) (Web Page) NIOSH Study finds Healthcare Workers' Exposure to Surgical Smoke Still Common. Updated November 3, 2015. Available online at: <https://www.cdc.gov/niosh/updates/upd-11-03-15.html>. Accessed: November 5, 2018.
30. Baggish MS, Baltoyannis P, Sze E (1988) Protection of the rat lung from the harmful effects of laser smoke. *Lasers Surg Med* 8 (3): 248-253.
31. Wu JS, Luttmann DR, Meininger TA, Soper NJ (1997) Production and systemic absorption of toxic byproducts of tissue combustion during laparoscopic surgery. *Surg Endosc* 11 (11): 1075-1079.
32. Abbate C, Giorgianni C, Munao F, Brecciaroli R (1993) Neurotoxicity induced by exposure to toluene. An electrophysiologic study. *Int Arch Occup Environ Health* 64 (6): 389-392.
33. Cometto-Muniz JE, Cain WS (1995) Relative sensitivity of the ocular trigeminal, nasal trigeminal and olfactory systems to airborne chemicals. *Chem Senses* 20 (2): 191-198.
34. Ernstgard L, Gullstrand E, Lof A, Johanson G (2002) Are women more sensitive than men to 2-propanol and m-xylene vapours? *Occup Environ Med* 59 (11): 759-767.
35. Ahaghotu E, Babu RJ, Chatterjee A, Singh M (2005) Effect of methyl substitution of benzene on the percutaneous absorption and skin irritation in hairless rats. *Toxicol Lett* 159 (3): 261-271.
36. Tunsaringkarn T, Siri Wong W, Rungsiyothin A, Nopparatbundit S (2012) Occupational exposure of gasoline station workers to BTEX compounds in Bangkok, Thailand. *Int J Occup Environ Med* 3 (3): 117-125.
37. Medicines and Healthcare Product Regulatory Agency (2008). Guidance on the safe use of lasers, intense light source systems and LEDs in medical, surgical, dental, and aesthetic practices. Available online at: http://webarchive.nationalarchives.gov.uk/20101123113954/http://www.mhra.gov.uk/home/idcplg?IdcService=GET_FILE&dDocName=CON014843&RevisionSelectionMethod=LatestReleased. Accessed: August 9, 2018.
38. Navarro-Meza MC, Gonzalez-Baltazar R, Aldrete-Rodriguez MG, Carmona-Navarro DE, Lopez-Cardona MG (2013) [Respiratory symptoms caused by the use of electrocautery in physicians being trained in surgery in a Mexican hospital]. *Rev Peru Med Exp Salud Publica* 30 (1): 41-44.
39. Schultz L (2015) Can efficient smoke evacuation limit aerosolization of bacteria? *AORN J* 102 (1): 7-14.
40. Christie D, Jefferson P, Ball DR (2005) Diathermy smoke and human health. *Anaesthesia* 60 (6): 632.
41. Chapman LW, Korta DZ, Lee PK, Linden KG (2017) Awareness of surgical smoke risks and assessment of safety practices during electrosurgery among US dermatology residents. *JAMA Dermatol* 153 (5): 467-468.
42. The National Institute for Occupation Safety and Health (1996). Control of smoke from laser/electric surgical procedures. Available online at: <https://www.cdc.gov/niosh/docs/hazardcontrol/pdfs/hc11.pdf>. Accessed: August 9, 2018.
43. Qian Y, Willeke K, Grinshpun SA, Donnelly J, Coffey CC (1998) Performance of N95 respirators: Filtration efficiency for airborne microbial and inert particles. *Am Ind Hyg Assoc J* 59 (2): 128-132.
44. Sehulster L, Chinn RY, Cdc, Hicpac (2003) Guidelines for environmental infection control in health-care facilities. Recommendations of CDC and the Healthcare Infection Control Practices Advisory Committee (HICPAC). *MMWR Recomm Rep* 52 (RR-10): 1-42.
45. Lee SA, Grinshpun SA, Reponen T (2008) Respiratory performance offered by N95 respirators and surgical masks: human subject evaluation with NaCl aerosol representing bacterial and viral particle size range. *Ann Occup Hyg* 52 (3): 177-185.
46. York K, Autry M (2018) Surgical Smoke: Putting the Pieces Together to Become Smoke-Free: 16 www.aornjournal.org/content/cme. *AORN J* 107 (6): 692-703.
47. Association of periOperative Registered Nurses (Web Page) Colorado Second State to Enact Surgical Smoke Evacuation Law. Updated March 28, 2019. Available online at: <https://www.aorn.org/about-aorn/aorn-newsroom/health-policy-news/2019-health-policy-news/co-enacts-surgical-smoke-law>. Accessed: May 15, 2019.
48. Association of periOperative Registered Nurses (Web Page) RI Introduces Smoke Evacuation Bill. Updated March 8, 2017. Available online at: <https://www.aorn.org/government-affairs/my-state/rhode-island/ri-news/ri-introduces-smoke-evacuation-bill>. Accessed: May 15, 2019.
49. Georgesen C, Lipner SR (2018) Surgical Smoke: Risks Assessment and Mitigation Strategies. *J Am Acad Dermatol* 79 (4): 746-755.
50. Occupational Safety and Health Administration (OSHA) (Web Page) Safety and Health Topics: Laser//Electrosurgery Plume. Updated Available online at: <https://www.osha.gov/SLTC/laserelectrosurgerplume/index.html>. Accessed: August 9, 2018.
51. Johnson AT (2016) Respirator masks protect health but impact performance: a review. *J Biol Eng* 10 (1): 4.
52. Hill DS, O'Neill JK, Powell RJ, Oliver DW (2012) Surgical smoke - a health hazard in the operating theatre: a study to quantify exposure and a survey of the use of smoke extractor systems in UK plastic surgery units. *J Plast Reconstr Aesthet Surg* 65 (7): 911-916.
53. Schultz L (2014) An analysis of surgical smoke plume components, capture, and evacuation. *AORN J* 99 (2): 289-298.
54. Ball K (2004) Controlling surgical smoke: A team approach. Information Booklet
55. MEGADYNE (2015) MEGAVAC PLUS™ Smoke Evacuator Operating and Installation Manual. Accessed: November 22, 2019.



ETHICON
PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

Shaping
the future
of surgery

©2019 Ethicon, Inc. All rights reserved. 115417-190529
For complete product details, see Instructions for Use.