

Infection Risk Management

Best Practice BSI Point Prevalence

Data Collection

Hospital Name
City, State

Date

Prepared by:

Objectives for BSI Point Prevalence Program

- To maximize BIOPATCH® Protective Disk with CHG investment and clinical effectiveness through third-party professional vascular access maintenance auditing
- To deliver robust data on vascular access devices broken down by type, dressing, unit and policy compliance
- To provide clinical staff with key takeaways and recommendations, providing insight into strengths and areas of improvement
- To partner with you in your efforts to reduce bloodstream infections

The Joint Commission (TJC)
2018 NPSG.07.04.01

Conduct periodic risk assessments for central line–associated bloodstream infections, monitor compliance with evidence-based practices, and evaluate the effectiveness of prevention efforts. The risk assessments are conducted in time frames defined by the hospital, and this infection surveillance activity is hospital-wide, not targeted¹.

The Changing Healthcare Landscape

This pay-for-performance breakdown shows how your dollars are at stake.¹



DHHS HAI Action Plan ²

2020 Proposal (From 2015 Baseline)

50% ↓ MRSA
(NHSN)


50% ↓
CLABSI

In October 2016, the U.S. Department of Health and Human Services (HHS) announced new targets for the national acute care hospital metrics for the National Action Plan to Prevent Health Care-Associated Infections: Road Map to Elimination (HAI Action Plan). The targets use data from calendar year 2015 as a baseline — and they replace the previous targets that expired in December 2013. ²

1. "Hospitals-Inpatient: Hospital Value Based Purchasing (HVBP), Scoring." https://www.cms.gov/Outreach-and-Education/Medicare-Learning-Network-MLN/MLNProducts/downloads/Hospital_VBPurchasing_Fact_Sheet_ICN907664.pdf. Accessed March 12, 2018.

2. National Action Plan to Prevent Health Care-Associated Infections: Road Map to Elimination. <https://health.gov/hcq/prevent-hai-action-plan.asp>. Accessed March 12, 2018

The Origin of Microorganisms Causing CRBSI¹



Contaminated Catheter Hub
12%

Contaminated Infusate
<1%

Skin Organisms
60%

Undetermined
28%

Dressings Dated



Dressings Dated & Within Policy




Dry Dressings



Intact Dressings



Most BSIs are caused by bacteria from the patients' own skin



HAI – Major Site of Infections	Estimated No.
Pneumonia	157,500
Gastrointestinal illness	123,100
Urinary tract infections	93,300
Primary bloodstream infections	71,900
Surgical site infections from any inpatient surgery	157,500
Other type of infections	118,500
Estimated total number of infections in hospitals	721,800

Reducing Surgical Site Infections: A Review. US National Library of Medicine.

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC2812878/> Accessed 3/4/18

Catheter-related bloodstream infections. US National Library of Medicine.


<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC4093967/> Accessed 3/4/18

HAIs at a Glance. Centers for Disease Control and Prevention. <https://www.cdc.gov/hai/surveillance/index.html>


Accessed 3/4/18

Why Biopatch? Why CHG?


When the Skin is the Source, Every Line Matters



Pre-Prep
Bacteria colonies exist not only on the surface, but below the surface as well, particularly within the hair follicles and sebaceous glands.¹




Post-Prep (immediately following antiseptic application)
Prepping the skin reduces colony counts of bacteria from the surface only — it never completely disinfects the skin.¹




Post-Prep (within 1-2 days following antiseptic application)
Resident bacteria begin to re-colonize the skin surface.¹


Sustained release of CHG



BIOPATCH[®] extends the post-prep environment for up to 7 days²



Pre-Prep



Day 1 – Day 7

1. Hendley JO, Ashe KM. Effect of topical antimicrobial treatment on aerobic bacteria in the stratum corneum of human skin. *Antimicrobial Agents and Chemotherapy*. April 1991;35(4):627-631
 2. Bhende MS, Rothenburger S. In vitro antimicrobial effectiveness of 5 catheter insertion-site dressings. *The Journal of the Association for Vascular Access*. 2007; 12(4):227-231

Guidelines and Standards: CHG Dressings

Publication	Highlights
CDC ¹	<ul style="list-style-type: none"> For patients aged 18 years and older: Chlorhexidine-impregnated dressings with an FDA-cleared label that specifies a clinical indication for reducing catheter-related bloodstream infection (CRBSI) or catheter-associated bloodstream infection (CABSI) are recommended to protect the insertion site of short-term, non-tunneled central venous catheters. (1A)
SHEA ²	<ul style="list-style-type: none"> Use chlorhexidine-containing dressings for CVCs in patients over 2 months of age
INS ³	<ul style="list-style-type: none"> Use chlorhexidine-impregnated dressings over CVADs to reduce infection risk Consider the use of chlorhexidine-impregnated dressings with peripheral arterial catheters as an infection reduction intervention Continuous (implanted) port access has infection rates that are similar to other long-term CVADs
AACN ⁴	<ul style="list-style-type: none"> Apply a chlorhexidine-impregnated sponge to the (CVC, peripheral arterial, or epidural) site Place a chlorhexidine sponge dressing around (Ventricular Assist Device) driveline Apply a chlorhexidine gluconate sponge to the (hemodialysis) VAC insertion site at each dressing change
APIC ⁵	<ul style="list-style-type: none"> The risk for infection is present during the entire dwell time of the catheter. The use of a post-insertion care bundle was associated with a significant reduction
ONS ⁶	<ul style="list-style-type: none"> Use a CHG-impregnated sponge dressing for all (CVC) catheters, including specialty catheters in patients older than 2 months of age Use a CHG-impregnated sponge dressing for any long-term infusion (>4-6 hrs) or if port remains accessed


1. 2017 Updated Recommendations on the Use of Chlorhexidine-Impregnated Dressings for Prevention of Intravascular Catheter-Related Infections Published Nov 1, 2017. <https://www.cdc.gov/infectioncontrol/guidelines/bsi/c-i-dressings/index.html>. 2. Marschall, J., et al. Strategies to Prevent Central Line–Associated Bloodstream Infections in Acute Care Hospitals: 2014 Update. ICHE. 3. Infusion Therapy Standards of Practice, Journal of Infusion Nursing. 2016, V39 (1S). 4. Weigand, DL, ed. *Procedure Manual for High Acuity, Progressive, and Critical Care. 7th Ed*, Philadelphia, PA; Saunders; 2017. 5. Guide to Preventing Central Line Associated Bloodstream Infections. APIC. 2017. 6. Access Device Standards of Practice for Oncology Nursing. ONS 2017

Does Your Insertion Practice Optimize Care & Maintenance?

- Follow Manufacturers' instructions for insertion and use to ensure line is secured for Max CHG delivery




Sample photos



Line Securement does NOT allow for maximum CHG delivery

- Area not accessible for thorough cleaning during dressing changes
- “Hubbed” securement does not allow for 360 degree application of BIOPATCH disk*



Line Secured to allow for maximum CHG Delivery

- “Zero” cm indicator
- Insertion site accessible for thorough cleaning during dressing changes
- Allows for 360 degree application of BIOPATCH disk*

*BIOPATCH placement diagrammed for demonstration purposes only. BIOPATCH should be placed with blue, printed side facing upward

Biopatch Compliance



Gauze



Line Secured for Max CHG Delivery



Sutured Lines



Tubing Labeled



OPTIONAL

Appropriate Protection of Tubing




Disinfection Caps




OPTIONAL

Entry Points of Exogenous Contamination of Vascular Devices




Guidelines & Recommendations: Peripheral IVs (PIV)




Organization	Recommendation
<p>CDC O'Grady, N.P., et al. Guidelines for the Prevention of Intravascular Catheter-Related Infections. AJIC 2011</p>	<p>“As in adults, the use of peripheral venous catheters in pediatric patients might be complicated by phlebitis, infusion extravasation, and catheter infection”</p> <p>“There is no need to replace peripheral catheters more frequently than every 72-96 hours to reduce risk of infection and phlebitis in adults”</p>
<p>SHEA Marshall, J., et al. Strategies to Prevent Central Line-Associated Bloodstream Infections in Acute Care Hospitals: 2014 Update. ICHE</p>	<p>“Peripheral arterial catheters and peripheral venous catheters are not included in most surveillance systems, although they are associated with risk of bloodstream infection independent of CVCs. Future surveillance systems may need to include bloodstream infections associated with these types of catheters”</p>
<p>ONS Access Device Standards of Practice for Oncology Nursing, ONS 2017</p>	<p>“Emerging data suggests that the rate of catheter-related bloodstream infections from peripheral catheters may be higher than once thought”</p>
<p>INS Infusion Therapy Standards of Practice, Journal of Infusion Nursing. 2016, V39 (1S)</p>	<p>“Consider monitoring bloodstream infection rates for peripheral catheters, or vascular catheter associated infections (peripheral) regularly”</p>

The Risk of Bloodstream Infection in Adults With Different Intravascular Devices: A Systematic Review of 200 Published Prospective Studies



“Since almost all the national effort and progress to date to reduce the risk of IVD-related infection have focused on short-term non-cuffed CVCs used in intensive care units, infection control programs must now strive to consistently apply essential control measures and preventive technologies with all types of IVDs.”

Evidence Brief: Clinical and Economic Burden of Peripheral Intravenous Catheter-Associated Complications in a U.S. Hospital Discharge Database¹



Patients with PIV-associated complications have longer LOS, higher costs, and greater risk of death than patients without.

Retrospective Database Analysis, over a 2-year period, from July 2013 to June 2015. Hospital discharge data from the Premier database (more than 700 US hospitals and 6 million patient discharges per year). A large patient population allows for a powered statistical analysis (n=588,375)

1. Lim S, Adams E, Hyde R, Broder M, Chang E, Reddy SR, Tarbox M, Bentley T, Ovington L. Clinical and Economic Burden of Peripheral Intravenous Catheter-Associated Complications in a U.S. Hospital Discharge Database. Poster presented at: 31st Annual Scientific Meeting of the Association for Vascular Access; 2017 Sep 16-19; Phoenix, AZ.

Emerging Evidence

Insertion of an IV catheter is an invasive procedure that introduces multiple risks and potential morbidities, and even mortality, and should be given the respect that it deserves.¹

Clinical Studies Evaluating PIV-BSI		
Author	Year	Publication
Guembe	2017	The Journal of Hospital Infection
Mermel	2017	Clinical Infectious Diseases
Sato	2017	BMC Infectious Diseases
DeVries	2016	Journal of the Association for Vascular Access
Kovacs	2016	American Journal of Infection Control
Heinrich	2013	GMS Hygiene and Infection Control
Mestre	2013	American Journal of Infection Control
Rickard	2013	The Medical Journal of Australia
Dychter	2012	Journal of Infusion Nursing
Hadaway	2012	Journal of Infusion Nursing
Trinh	2011	Infection Control and Hospital Epidemiology
Pien	2010	American Journal of Medicine
Easterlow	2009	Journal of Clinical Nursing
Zing	2009	International Journal of Antimicrobial Agents
Pujol	2007	Journal of Hospital Infection
Maki	2006	Mayo Clin Proc.

1. Helm R. Accepted but Unacceptable: Peripheral IV Catheter Failure Rate. Journal of Infusion Nursing. 2015; 38:3;180-201

Peripheral Data-Site Selection



Dressings Dated



Dressings Dated Within Policy



Dressings Dry



Dressings Intact



Securement Dressing or Device




Protected Peripheral VADs- Biopatch Compliance



OPTIONAL

Multiple Vascular Access Devices



Multiple Vascular Access Devices - Patient Data




OPTIONAL

Knowledge Based Questions

1. Insert question 1 here
2. Insert question 2 here
3. Insert question 3 here
4. Insert question 4 here
5. Insert question 5 here



OPTIONAL



Key Takeaways



Building One Standard of Care for All Vascular Access Devices



Flushing
Protocols

Alcohol
Impregnated
Caps

Stabilization

Scrub the
Hub

Sterile
Transparent
Dressing

Insertion
Bundles

CHG
Skin Prep

Maintenance
Bundles

Protective
Disk with
CHG

BIOPATCH Protective Disk with CHG is a hydrophilic wound dressing that is used to absorb exudate and to cover a wound caused by the use of vascular and non-vascular percutaneous medical devices such as:

- Central Venous Lines
- Arterial Catheters
- Dialysis Catheters
- Drains
- Chest Tubes
- IV Catheters



- Midline Catheters
- Epidural Catheters
- Externally-placed Orthopedic Pins
- Peripherally-inserted Coronary Catheters

It is also intended to reduce local infections, catheter related blood stream infections (CRBSI) and skin colonization of microorganisms commonly related to CRBSI in patients with central venous or arterial catheters



Next Steps

- 100% commitment to providing additional education to help you meet and exceed your institutional goals
- We will partner with you to help support quality patient care, nursing excellence and innovations in professional nursing practice.
- Next Steps: Repeat Point Prevalence _____(tbd)_____



Questions

Guidelines and Standards Appendix

The following slides contain Guidelines and Standards that are commonly referenced for insertion, care and maintenance of vascular access devices.

APPENDIX A: Dressings

Publication	Recommendation
<p data-bbox="222 562 297 601">INS</p> <p data-bbox="65 611 452 701"><i>Infusion Therapy Standards of Practice, Journal of Infusion Nursing. 2016, V39 (1S)</i></p>	<p data-bbox="490 334 1736 425">Label the dressing with the date performed or date to be changed based on organizational policies and procedures</p> <p data-bbox="490 451 1736 542">Replace the catheter site dressing if the dressing becomes damp, loosened, or visibly soiled</p> <p data-bbox="490 554 1804 739">A sterile dressing is applied and maintained on all peripheral, non-tunneled, peripherally inserted central catheters, accessed implanted VADs, and tunneled cuffed catheters, at least until the insertion site is well healed</p> <p data-bbox="490 768 1765 859">Replace dressing on short-term CVC sites at least every 2 days for gauze dressings</p>
<p data-bbox="212 1068 307 1106">CDC</p> <p data-bbox="73 1116 446 1206"><i>O'Grady, N.P., et al. Guidelines for the Prevention of Intravascular Catheter-Related Infections. AJIC 2011</i></p>	<p data-bbox="490 976 1846 1068">Note that a gauze dressing underneath a TSM dressing is considered a gauze dressing and changed at least every 2 days</p> <p data-bbox="490 1093 1736 1185">Replace the catheter site dressing if the dressing becomes damp, loosened, or visibly soiled</p> <p data-bbox="490 1210 1765 1302">Replace dressing on short-term CVC sites at least every 2 days for gauze dressings.</p>

APPENDIX B: CHG Dressings

Publication	Highlights
CDC ¹	<ul style="list-style-type: none"> For patients aged 18 years and older: Chlorhexidine-impregnated dressings with an FDA-cleared label that specifies a clinical indication for reducing catheter-related bloodstream infection (CRBSI) or catheter-associated bloodstream infection (CABSIs) are recommended to protect the insertion site of short-term, non-tunneled central venous catheters. (1A)*
SHEA ²	<ul style="list-style-type: none"> Use chlorhexidine-containing dressings for CVCs in patients over 2 months of age*
INS ³	<ul style="list-style-type: none"> Use chlorhexidine-impregnated dressings over CVADs to reduce infection risk* Consider the use of chlorhexidine-impregnated dressings with peripheral arterial catheters as an infection reduction intervention* Continuous (implanted) port access has infection rates that are similar to other long-term CVADs
AACN ⁴	<ul style="list-style-type: none"> Apply a chlorhexidine-impregnated sponge to the (CVC, peripheral arterial, or epidural) site* Place a chlorhexidine sponge dressing around (Ventricular Assist Device) driveline* Apply a chlorhexidine gluconate sponge to the (hemodialysis) VAC insertion site at each dressing change*
APIC ⁵	<ul style="list-style-type: none"> The risk for infection is present during the entire dwell time of the catheter. The use of a post-insertion care bundle was associated with a significant reduction
ONS ⁶	<ul style="list-style-type: none"> Use a CHG-impregnated sponge dressing for all (CVC) catheters, including specialty catheters in patients older than 2 months of age* Use a CHG-impregnated sponge dressing for any long-term infusion (>4-6 hrs) or if port remains accessed*

1. 2017 Updated Recommendations on the Use of Chlorhexidine-Impregnated Dressings for Prevention of Intravascular Catheter-Related Infections Published Nov 1, 2017. (<https://www.cdc.gov/infectioncontrol/guidelines/bsi/c-i-dressings/index.html>). 2. Marshall, J., et al. Strategies to Prevent Central Line-Associated Bloodstream Infections in Acute Care Hospitals: 2014 Update. ICHE. 3. Infusion Therapy Standards of Practice, Journal of Infusion Nursing. 2016, V39 (1S). 4. Procedure Manual for High Acuity, Progressive, and Critical Care. 7th Ed, AACN 2017. 5. Guide to Preventing Central Line Associated Bloodstream Infections. APIC. 2017. 6. Access Device Standards of Practice for Oncology Nursing, ONS 2017

APPENDIX C: Catheter Securement

Publication	Recommendation
<p>INS <i>Infusion Therapy Standards of Practice, Journal of Infusion Nursing. 2016, V39 (1S)</i></p>	<p>The use of a catheter stabilization device should be considered the preferred alternative to tape or sutures when feasible.²</p>
<p>CDC <i>O'Grady, N.P., et al. Guidelines for the Prevention of Intravascular Catheter-Related Infections. AJIC 2011</i></p>	<p>Use a suture-less securement device to reduce the risk of infection for intravascular catheters.¹</p>

APPENDIX D: Tubing

Publication	Recommendation
<p data-bbox="233 548 311 586">INS</p> <p data-bbox="81 596 467 689"><i>Infusion Therapy Standards of Practice, Journal of Infusion Nursing. 2016, V39 (1S)</i></p>	<p data-bbox="506 522 1765 611">Aseptically attach a new, sterile, compatible covering device to the male luer end of the administration set after each intermittent use.</p> <p data-bbox="591 622 1839 715">Do not attach the exposed male luer end of the administration set to a port on the same set (“looping”).</p>
<p data-bbox="227 772 320 811">CDC</p> <p data-bbox="86 821 461 908"><i>O'Grady, N.P., et al. Guidelines for the Prevention of Intravascular Catheter-Related Infections. AJIC 2011</i></p>	<p data-bbox="506 746 1823 936">In patients not receiving blood, blood products, or fat emulsions, replace administration sets that are continuously used, including secondary sets and add-on devices, no more frequently than 96-hour intervals, but at least every 7 days.</p>

APPENDIX E: Hub Disinfection

Publication	Recommendation
<p data-bbox="208 601 285 639">INS</p> <p data-bbox="67 648 426 739"><i>Infusion Therapy Standards of Practice, Journal of Infusion Nursing. 2016, V39 (1S)</i></p>	<p data-bbox="475 505 1748 596">Perform a vigorous mechanical scrub for manual disinfection of the needleless connector prior to each VAD access and allow it to dry</p> <p data-bbox="475 639 1769 833">Use of passive disinfection caps containing disinfecting agents (e.g., isopropyl alcohol) has been shown to reduce intraluminal microbial contamination and reduce the rates of central line-associated bloodstream infection (CLABSI)</p>
<p data-bbox="59 862 434 901">Safdar N, Maki DG.</p> <p data-bbox="59 911 440 1115"><i>The pathogenesis of catheter-related bloodstream infection with nuncuffed short-term central venous catheters. Int Care Med. 2004; 30:62-67.</i></p>	<p data-bbox="475 968 1570 1006">12% of CRBSI originate from contaminated catheter hubs</p>

APPENDIX F: Central Venous Catheters

Publication	Recommendation
<p>CDC O'Grady, N.P., et al. Guidelines for the Prevention of Intravascular Catheter-Related Infections. AJIC 2011</p>	<p>Prepare clean skin site with > 0.5% chlorhexidine with alcohol</p> <p>Use a chlorhexidine/silver sulfadiazine or minocycline/rifampin - impregnated CVC* (in patients whose catheter is expected to remain in place >5 days if, after successful implementation of a comprehensive strategy to reduce rates of CLABSI, the CLABSI rate is not decreasing)</p> <p>Use ultrasound guidance to place CVCs</p> <p>Use a CVC with the minimum number of necessary ports or lumens</p> <p>Do not administer routine systemic antimicrobial prophylaxis prevent CRBSI</p> <p>Use a subclavian site, Avoid jugular or femoral site</p> <p>Replace the catheter as soon as possible (within 48 hours) When adherence to aseptic technique cannot be ensured and catheters inserted during a medical emergency</p>
<p>CDC Updated Recommendations on the Use of Chlorhexidine-Impregnated Dressings for Prevention of Intravascular Catheter-Related Infections Published Nov 1, 2017</p>	<p>For patients aged 18 years and older: Chlorhexidine-impregnated dressings with an FDA-cleared label that specifies a clinical indication for reducing catheter-related bloodstream infection (CRBSI) or catheter-associated bloodstream infection (CABSI) are recommended to protect the insertion site of short-term, non-tunneled central venous catheters. (1A)</p>
<p>SHEA Marschall, J., et al. Strategies to Prevent Central Line-Associated Bloodstream Infections in Acute Care Hospitals: 2014 Update. ICHE</p>	<p>Factors associated with increased risk of CLABSI includes heavy microbial colonization at the insertion site.</p> <p>Insertion: hand hygiene, max. sterile barrier precautions, avoid fem. vein insertion – select adults, alcoholic CHG solution containing more than 0.5% CHG</p> <p>Daily bathing with CHG for adult ICU pop.</p> <p>Antiseptic- or antimicrobial- impregnated CVCs in adult patients</p> <p>Chlorhexidine-containing dressings for CVCs in patients over 2 mo of age Antiseptic-containing hub/connector cap/port protector to cover connectors</p>

APPENDIX G: Central Venous Catheters *continued*

Publication	Recommendation
<p style="text-align: center;">INS</p> <p>Infusion Therapy Standards of Practice, Journal of Infusion Nursing. 2016, V39 (1S)</p>	<p>Ensure that the intended VAD site is visibly clean prior to application of an antiseptic solution</p> <p>Remove excess hair at the insertion site if needed to facilitate application of VAD dressings (do not shave)</p> <p>Use a standardized insertion check list to be completed by someone other than the inserter</p> <p>Avoid the use of tape or sutures as they are not effective alternatives to an ESD (Engineered Stabilization Device)...Sutures are associated with needle-stick injury, in addition to supporting the growth of biofilm and increasing the risk of catheter-related blood stream infection.</p> <p>Use chlorhexidine-impregnated dressings</p> <p>Recognize the risk of contamination with each manipulation of the infusion system (access with syringes, tubing, etc...)</p> <p>Continuous (implanted) port access has infection rates that are similar to other long-term CVADs</p>
<p style="text-align: center;">AACN</p> <p>Procedure Manual for High Acuity, Progressive, and Critical Care. 7th Ed, AACN 2017</p>	<p>Cleanse the skin, catheter, and stabilizing device with 2% chlorhexidine-based preparation</p> <p>Apply a chlorhexidine-impregnated sponge to the site</p> <p>Assess the catheter site daily and as needed by palpation through an intact dressing</p>
<p style="text-align: center;">APIC</p> <p>Guide to Preventing Central Line Associated Bloodstream Infections. APIC. 2017</p>	<p>Use of an evidence-based bundle of interventions can improve patient outcomes only if the interventions are consistently completed</p> <p>Establishment of a team dedicated to all aspects of intravenous therapy</p> <p>Use of insertion check list</p> <p>The risk for infection is present during the entire dwell time of the catheter</p> <p>The use of a post-insertion care bundle was associated with a significant reduction</p>
<p style="text-align: center;">ONS</p> <p>Access Device Standards of Practice for Oncology Nursing, ONS 2017</p>	<p>Consider use of connectors with ... visible fluid path to assess efficacy of flush technique and a solid, flat surface that easily is disinfected</p> <p>Use a CHG-impregnated sponge dressing for all catheters, including specialty catheters in patients older than 2months of age</p> <p style="padding-left: 20px;">- for any long term infusion (>4-6 hrs) or if port remains accessed</p>

APPENDIX H: Implanted Port

Publication	Recommendation
<p>CDC O'Grady, N.P., et al. <i>Guidelines for the Prevention of Intravascular Catheter-Related Infections</i>. AJIC 2011</p>	<p>Use a chlorhexidine-impregnated sponge dressing for temporary short-term catheters in patients older than 2 months of age if the CLABSI rate is not decreasing despite adherence to basic prevention measures, including education and training, appropriate use of chlorhexidine for skin antisepsis, and MSB</p> <p>No recommendation can be made regarding the frequency for replacing needles to access implantable ports. Unresolved issue</p>
<p>SHEA Marschall, J., et al. <i>Strategies to Prevent Central Line–Associated Bloodstream Infections in Acute Care Hospitals: 2014 Update</i>. ICHE</p>	<p>These recommendations (sic) are not stratified on the basis of catheter type (e.g., tunneled, implanted, cuffed, noncuffed catheter, and dialysis catheter)</p> <p>Use chlorhexidine-containing dressings for CVCs in patients over 2 months of age</p>
<p>ONS <i>Access Device Standards of Practice for Oncology Nursing</i>, ONS 2017</p>	<p>Following chlorhexidine (CHG) skin preparation, use a CHG –impregnated sponge dressing for any long-term infusion exceeding 4-6 hours or if the port remains accessed for intermittent infusion for greater than 4-6 hours</p>
<p>INS <i>Infusion Therapy Standards of Practice</i>, Journal of Infusion Nursing. 2016, V39 (1S)</p>	<p>When used intermittently, ports have a lower incidence of catheter- related bloodstream infection (CR-BSI); however, continuous port access has infection rates that are similar to other long-term CVADs.</p> <p>A sterile dressing is maintained over the access site if the implanted vascular access port remains accessed.</p> <p>Use chlorhexidine-impregnated dressings over CVADs to reduce infection risk when the extraluminal route is the primary source of infection.</p>

APPENDIX I: Hemodialysis Catheters

Publication	Recommendation
<p>CDC Core Interventions for Dialysis Bloodstream Infection Prevention (last update 6/2016)</p>	<p>Apply antibiotic ointment or povidone-iodine ointment to catheter exit sites during dressing change. See information on selecting an antimicrobial ointment for hemodialysis catheter exit sites (selecting an antimicrobial ointment). Use of chlorhexidine-impregnated sponge dressing might be an alternative.</p>
<p>SHEA Marshall, J., et al. <i>Strategies to Prevent Central Line–Associated Bloodstream Infections in Acute Care Hospitals</i>: 2014 Update. ICHE</p>	<p>These recommendations (sic) are not stratified on the basis of catheter type (e.g., tunneled, implanted, cuffed, non-cuffed catheter, and dialysis catheter)</p> <p>Infection prevention and control efforts should include other vulnerable populations, such as patients receiving hemodialysis through catheters, intraoperative patients, and oncology patients.</p> <p>Use chlorhexidine-containing dressings for CVCs in patients over 2 months of age</p>
<p>AACN <i>Procedure Manual for High Acuity, Progressive, and Critical Care</i>. 7th Ed, AACN 2017.</p>	<p>Dressing supplies (... triple-antibiotic ointment or chlorhexidine-impregnated sponge)</p> <p>Procedure Step: Apply a chlorhexidine gluconate sponge to the VAC insertion site at each dressing change</p> <p>Rationale: Lowers the incidence of dialysis catheter–related infections</p>

APPENDIX J: Epidural Catheters

Publication	Recommendation
<p>AACN <i>Procedure Manual for High Acuity, Progressive, and Critical Care. 7th Ed, AACN 2017.</i></p>	<p>Use of povidone-iodine versus chlorhexidine is controversial. The antiseptic solution should be allowed to dry before the initial incision is made. Studies suggest chlorhexidine is neurotoxic.</p> <p>The choice of povidone-iodine or chlorhexidine as an antiseptic agent is an unresolved issue. Both should be allowed to dry completely.</p> <p>Institutional standards may vary with regard to antiseptic choice—follow your institutional standard. Although the package insert for chlorhexidine warns against use before any neuraxial procedures, a large retrospective study and several anesthesiology societies recommend chlorhexidine as an antiseptic. Do not use an alcohol-based preparation. Use aqueous chlorhexidine or povidone-iodine. Preparations with alcohol are neurotoxic to the epidural space. The choice of povidone-iodine or chlorhexidine as an antiseptic agent is controversial. Both should be allowed to dry completely. Studies suggest chlorhexidine is neurotoxic.</p>
<p>ONS <i>Access Device Standards of Practice for Oncology Nursing, ONS 2017</i></p>	<p>“Maintain maximum sterile barrier precautions with mask and sterile gloves with any access or maintenance procedure”</p> <p>“A chlorhexidine impregnated dressing or sponge has been found to significantly reduce the rate of epidural infections”</p>
<p>INS <i>Infusion Therapy Standards of Practice, Journal of Infusion Nursing. 2016, V39 (1S)</i></p>	<p>“Maintain strict aseptic technique while wearing a mask and sterile gloves during any intraspinal access or maintenance procedure.”</p> <p>“Consider the use of chlorhexidine impregnated dressings for patients with an epidural access device. A significant reduction in epidural skin colonization and catheter tip colonization has been demonstrated with their use”</p>

APPENDIX K: Peripheral Arterial Catheters

Publication	Recommendation
<p>CDC O'Grady, N.P., et al. <i>Guidelines for the Prevention of Intravascular Catheter-Related Infections</i>. AJIC 2011</p>	<p>“A minimum of a cap, mask, sterile gloves and a small sterile fenestrated drape should be used during peripheral arterial catheter insertion”</p> <p>“During axillary or femoral artery catheter insertion, maximal sterile barriers precautions should be used”</p> <p>“Do not routinely replace arterial catheters to prevent catheter related infections”</p>
<p>SHEA Marschall, J., et al. <i>Strategies to Prevent Central Line–Associated Bloodstream Infections in Acute Care Hospitals: 2014 Update</i>. ICHE</p>	<p>“Besides central venous catheters (CVCs), peripheral arterial catheters also carry a risk of infection.”</p>
<p>AACN <i>Procedure Manual for High Acuity, Progressive, and Critical Care</i>. 7th Ed, AACN 2017.</p>	<p>“Arterial catheter sites are a source of bloodstream infections, with the femoral site being more heavily associated with colonization compared with other sites. The infective potential of the arterial catheter is equivalent to the short term central venous device regarding colonization and bloodstream infections, and should be assessed together for signs and symptoms of infection”</p> <p>“Apply a chlorhexidine impregnated sponge to the site”</p>
<p>INS <i>Infusion Therapy Standards of Practice</i>, Journal of Infusion Nursing. 2016, V39 (1S)</p>	<p>“Consider the use of chlorhexidine-impregnated dressings with peripheral arterial catheters as an infection reduction intervention”</p>

APPENDIX L: Peripheral IVs (PIV)

Publication	Recommendation
<p>CDC O'Grady, N.P., et al. <i>Guidelines for the Prevention of Intravascular Catheter-Related Infections</i>. AJIC 2011</p>	<p>“As in adults, the use of peripheral venous catheters in pediatric patients might be complicated by phlebitis, infusion extravasation, and catheter infection”</p>
<p>SHEA Marshall, J., et al. <i>Strategies to Prevent Central Line-Associated Bloodstream Infections in Acute Care Hospitals</i>: 2014 Update. ICHE</p>	<p>“Peripheral arterial catheters and peripheral venous catheters are not included in most surveillance systems, although they are associated with risk of bloodstream infection independent of CVCs. Future surveillance systems may need to include bloodstream infections associated with these types of catheters”</p>
<p>ONS <i>Access Device Standards of Practice for Oncology Nursing</i>, ONS 2017</p>	<p>“Emerging data suggests that the rate of catheter-related bloodstream infections from peripheral catheters may be higher than once thought”</p>
<p>INS <i>Infusion Therapy Standards of Practice</i>, Journal of Infusion Nursing. 2016, V39 (1S)</p>	<p>“Consider monitoring bloodstream infection rates for peripheral catheters, or vascular catheter associated infections (peripheral) regularly”</p>

APPENDIX M: Peripheral IV - REMOVAL

Publication	Recommendation
<p data-bbox="164 578 241 614">INS</p> <p data-bbox="54 625 351 715"><i>Infusion Therapy Standards of Practice, Journal of Infusion Nursing. 2016, V39 (1S)</i></p>	<p data-bbox="415 321 1806 514"><i>Remove short peripheral and midline catheters in pediatric and adult patients when clinically indicated based on findings from site assessment and or clinical signs and symptoms of systemic complications (e.g., Bloodstream infection).</i></p> <p data-bbox="415 571 1777 649">Signs and symptoms of complications with or without infusion through the catheter include, but are not limited to the presence of:</p> <ol data-bbox="415 671 1835 1021" style="list-style-type: none"><li data-bbox="415 671 1516 706">1. <i>Any level of pain and or tenderness with or without palpation</i><li data-bbox="415 714 1188 749">2. <i>Changes in color: erythema or blanching</i><li data-bbox="415 756 1188 792">3. <i>Changes in skin temperature: hot or cold</i><li data-bbox="415 799 637 835">4. <i>Edema</i><li data-bbox="415 842 685 878">5. <i>Induration</i><li data-bbox="415 885 1497 921">6. <i>Leakage of fluid or purulent drainage from the puncture site</i><li data-bbox="415 928 1835 1021">7. <i>Other types of dysfunction (e.g., resistance when flushing, absence of the blood return)</i>
<p data-bbox="154 1049 251 1085">CDC</p> <p data-bbox="28 1099 376 1242"><i>O'Grady, N.P., et al. Guidelines for the Prevention of Intravascular Catheter-Related Infections. AJIC 2011</i></p>	<p data-bbox="415 1085 1835 1163">“There is no need to replace peripheral catheters more frequently than every 72-96 hours to reduce risk of infection and phlebitis in adults”</p>

APPENDIX N: Peripheral IV Site Selection

Publication	Recommendation
<p data-bbox="164 768 241 806">INS</p> <p data-bbox="54 816 351 906"><i>Infusion Therapy Standards of Practice, Journal of Infusion Nursing. 2016, V39 (1S)</i></p>	<p data-bbox="415 321 1058 364">Short Peripheral Catheters (adult)</p> <p data-bbox="415 385 1734 478">Use the venous site most likely to last the full length of the prescribed therapy</p> <p data-bbox="415 521 714 564">Consider using:</p> <ul data-bbox="502 585 1831 842" style="list-style-type: none"><li data-bbox="502 585 1831 678">• forearm to increase dwell time, decrease pain during dwell time, promote self-care, and prevent accidental removal and occlusions.<li data-bbox="502 699 1831 842">• veins found on the dorsal and ventral surfaces of the upper extremities, including the metacarpal, cephalic, basilic, and median veins <p data-bbox="415 871 637 913">Avoid using</p> <ul data-bbox="502 935 1850 1163" style="list-style-type: none"><li data-bbox="502 935 1850 1028">• the ventral surface of the wrist due to pain on insertion and possible nerve damage<li data-bbox="502 1049 1445 1092">• areas of flexion and areas of pain on palpation<li data-bbox="502 1113 1671 1163">• compromised areas and sites distal to these compromised <p data-bbox="415 1206 1816 1292">Make no more than 2 attempts at short peripheral intravenous access per clinician, and limit total attempts to no more than 4</p>

APPENDIX O: Peripheral IV Dressings

Publication	Recommendation
<p data-bbox="195 668 272 704">INS</p> <p data-bbox="86 715 382 808"><i>Infusion Therapy Standards of Practice, Journal of Infusion Nursing, 2016, V39 (1S)</i></p>	<p data-bbox="440 382 1843 486">Label the dressing with the date performed or date to be changed based on organizational policies and procedures.</p> <p data-bbox="440 525 1792 686">Perform dressing changes on short peripheral catheters if the dressing becomes damp, loosened, and/or visibly soiled <u>and at least</u> every 5 to 7 days.</p> <p data-bbox="440 715 1682 811">For peripheral catheters, consider two options for catheter stabilization</p> <ul data-bbox="508 843 1850 1053" style="list-style-type: none">an integrated stabilization feature on the catheter hub combined with a bordered polyurethane securement dressing ora standard round hub peripheral catheter in combination with an adhesive ESD (Engineered Stabilization Device).

APPENDIX P: Midline Catheters

Publication	Recommendation
<p>CDC O'Grady, N.P., et al. <i>Guidelines for the Prevention of Intravascular Catheter-Related Infections</i>. AJIC 2011</p>	<p>Use a midline catheter or peripherally inserted central catheter (PICC), instead of a short peripheral catheter, when the duration of IV therapy will likely exceed six days</p> <p>Sterile gloves should be worn for the insertion of arterial, central, and midline catheters</p>
<p>AACN <i>Procedure Manual for High Acuity, Progressive, and Critical Care</i>. 7th Ed, AACN 2017.</p>	<p>Procedure step: "... apply a chlorhexidine impregnated gel dressing or sponge to the site and then cover it with a sterile transparent semipermeable membrane dressing"</p> <p>Rationale: "Decreases catheter-related infections"</p>
<p>INS <i>Infusion Therapy Standards of Practice</i>, Journal of Infusion Nursing. 2016, V39 (1S)</p>	<p>Consider the use of maximal sterile barrier precautions with midline catheter insertion.</p> <p>Remove</p>

APPENDIX Q: Peripheral Arterial Catheters

Publication	Recommendation
<p>CDC O'Grady, N.P., et al. Guidelines for the Prevention of Intravascular Catheter-Related Infections. AJIC 2011</p>	<p>“A minimum of a cap, mask, sterile gloves and a small sterile fenestrated drape should be used during peripheral arterial catheter insertion”</p> <p>“During axillary or femoral artery catheter insertion, maximal sterile barriers precautions should be used”</p> <p>“Do not routinely replace arterial catheters to prevent catheter related infections”</p>
<p>SHEA Marschall, J., et al. Strategies to Prevent Central Line–Associated Bloodstream Infections in Acute Care Hospitals: 2014 Update. ICHE</p>	<p>“Besides central venous catheters (CVCs), peripheral arterial catheters also carry a risk of infection.”</p>
<p>AACN Procedure Manual for High Acuity, Progressive, and Critical Care. 7th Ed, AACN 2017.</p>	<p>“Arterial catheter sites are a source of bloodstream infections, with the femoral site being more heavily associated with colonization compared with other sites. The infective potential of the arterial catheter is equivalent to the short term central venous device regarding colonization and bloodstream infections, and should be assessed together for signs and symptoms of infection”</p> <p>“Apply a chlorhexidine impregnated sponge to the site”</p>
<p>INS Infusion Therapy Standards of Practice, Journal of Infusion Nursing. 2016, V39 (1S)</p>	<p>“Consider the use of chlorhexidine-impregnated dressings with peripheral arterial catheters as an infection reduction intervention”</p>

APPENDIX R: Line Necessity- Central and Peripheral Vascular Access Devices

Publication	Recommendation
<p data-bbox="216 596 291 634">INS</p> <p data-bbox="73 644 436 733"><i>Infusion Therapy Standards of Practice, Journal of Infusion Nursing. 2016, V39 (1S)</i></p>	<p data-bbox="481 554 1837 644">The clinical need for each peripheral and non-tunneled central vascular access device (CVAD) is assessed on a daily basis.</p> <p data-bbox="481 686 1785 776"><i>Remove the short peripheral catheter if it is no longer included in the plan of care or has not been used for 24 hours or more (V)</i></p>
<p data-bbox="206 808 301 845">CDC</p> <p data-bbox="73 856 436 1011">Bloodstream Infection Event. Centers for Disease Control and Prevention website. www.cdc.gov. Accessed 5/16/17</p>	<p data-bbox="481 815 1812 1001">Patients with multiple vascular access devices who develop a Bloodstream Infection may meet the NHSN Central Line Associated Bloodstream Infection (CLABSI) definition if a central line is in place regardless of the causative vascular device.¹</p>