

Enseal™

ENSEAL™ X1 Tissue Sealers Expect more.^{1-3*}



ENSEAL™ X1 CURVED JAW

ENSEAL™ X1 LARGE JAW

*ENSEAL™ X1 Curved Jaw has a longer jaw, longer cut length and wider jaw aperture compared to LigaSure Maryland (LF1937) (p < 0.001). In benchtop testing on porcine arteries, vessels sealed with ENSEAL X1 Curved Jaw had a 22% higher average burst pressure than vessels sealed with LigaSure™ Maryland (LF1937), (1055mmHg vs. 862mmHg, p < 0.001).

References: 1. Ethicon, Project Floyd: Claims Metrology Report, June 2018, PRC079564B (145171-200630)
2. Ethicon, Floyd Relaunch Claims Metrology, June 2020, PRC095763A (145171-200630)
3. Ethicon, Floyd Relaunch Claims Ex-Vivo Sealing, June 2020, PRC094697A (145171-200630)

The ENSEAL™ X1 Tissue Sealers
offer more than LigaSure™



More secure^{1*}

More efficient^{2,3}**

The ENSEAL™ X1 Tissue Sealers are advanced bipolar devices designed for use in open or laparoscopic surgical procedures.[#] We have and will further redesign our Ethicon Energy Advanced Bipolar Portfolio, to provide secure sealing with more intuitive and simplified steps-for-use.^{4,5}

*Preclinical test of distal tip bleeding (ENSEAL® vs. Impact-LF4318) in thick porcine mesentery base (p=0.001). **ENSEAL X1 Curved Jaw Tissue Sealer can capture, seal and transect a longer length of tissue per single activation due to a 16% (or 3.4mm) longer jaw (p < 0.001) and a 19% (or 3.5mm) longer cut length (p < 0.001) compared to LigaSure™ Maryland (LF1937).^{ures} #ENSEAL™ X1 Large Jaw is intended for use in open surgical procedures

REFERENCES: **1.** Ethicon, PSP005173A, Final Report- 2017 Brick Claims: Comparison of Distal Tip Sealing Performance Between ENSEAL™ X1 Large Jaw (NSLX120L) with Algorithm and Ligasure™ Impact (LF4318), Sept 2016, Data on File (154862-201002) **2.** Ethicon, Project Floyd: Claims Metrology Report, June 2018, PRC079564B (145163-200630) **3.** Ethicon, Floyd Relaunch Claims Metrology, June 2020, PRC095763A (145163-200630) **4.** Ethicon, D0C023555A, A Floyd Claims Memo Industrial Design, July 2018, Data on File (155321-201009) **5.** Ethicon, PRC095915A, Curved Design Validation 2 Completion Report, June 2020, Data on File (155321-201009)

ENSEAL™ X1 Curved Jaw

More efficient than LigaSure™ Maryland^{1,2*}

- **Capture more tissue** per bite with a longer jaw and wider jaw aperture^{1,2**}
- 32% stronger grasping^{3#}
- **360° shaft rotation** to enable access to targeted tissue^{4,5}

Curved, tapered tip designed for fine dissection^{3,6}

ENSEAL™ X1 Curved Jaw can capture more tissue per bite with a 16% longer jaw and 9% wider jaw aperture compared to LigaSure™ Maryland^{1,2**}

9%
wider jaw
aperture^{2†}

16%
longer jaw^{1†}

19%
longer cut
length^{1†}

ENSEAL™ X1 Curved Jaw

LigaSure™ Maryland

*ENSEAL™ X1 Curved Jaw Tissue Sealer can capture, seal and transect a longer length of tissue per single activation due to a 16% (or 3.4mm) longer jaw ($p < 0.001$) and a 19% (or 3.5mm) longer cut length ($p < 0.001$) compared to LigaSure™ Maryland (LF1937). **Based on metrology data, ENSEAL™ X1 Curved Jaw Tissue Sealer has a 16% (or 3.4mm) longer jaw than LigaSure™ Maryland (LF1937) ($p < 0.001$) and ENSEAL X1 Curved Jaw Tissue Sealer has a 9% (or 1.15mm) wider jaw aperture than LigaSure™ Maryland (LF1937) ($p < 0.001$). #Grasping force measured as the maximum amount of force required to pull porcine jejunum from the distal tip of device jaws. Comparison of ENSEAL™ X1 Curved Jaw to LigaSure™ Maryland (LF1937) ($p < 0.001$) †Metrology report comparing the jaw aperture of ENSEAL™ X1 Curved Jaw to LigaSure™ Maryland (LF1937) ($p < 0.001$)

REFERENCES: **1.** Ethicon, Project Floyd: Claims Metrology Report, June 2018, PRC079564B (145163-200630, 145041-200629, 152465-200909, 152466-200909) **2.** Ethicon, Floyd Relaunch Claims Metrology, June 2020, PRC095763A (145163-200630, 145041-200629, 145034-200629) **3.** Ethicon, Floyd Relaunch Claims Grasping Force, June 2020, PRC096063A (145160-200630) **4.** As Per Instructions For Use (152441-200909) **5.** Ethicon, DOC023555A, A Floyd Claims Memo - Industrial Design, July 2018, Data on File (152441-200909, 152450-200909) **6.** Ethicon, DOC024721A, Floyd Marketing Claims Memo, Jan 2019, Data on File (152450-200909)

ENSEAL™ X1 Large Jaw

More secure than LigaSure Impact™^{1*}




- Had **less bleeding** at the distal tip vs LigaSure Impact™ in thick tissue^{2,3*}
- Had **41% less lateral thermal spread** than LigaSure™ Impact^{4**}
- The convenient controls of the ENSEAL™ X1 Large Jaw were designed for less hand movement⁵
- Has improved ergonomics which allows a 360° rotation compared to LigaSure Impact™⁵

ENSEAL™ X1 Large Jaw has a larger distal electrode surface area than LigaSure™ Impact.⁶ ENSEAL™ X1 Large Jaw had significantly less bleeding at the distal tip vs LigaSure Impact™ in thick tissue.^{2,**}

ENSEAL™ X1 Large Jaw



LigaSure Impact™



ENSEAL™ X1 Large Jaw



*Preclinical test of distal tip bleeding (ENSEAL® vs. Impact-LF4318) in thick porcine mesentery base (p=0.001). **Preclinical testing on porcine carotids (ENSEAL™ vs. Impact-LF4318) that measured mean max lateral thermal damage via histology (p=0.005)

REFERENCES: 1. Ethicon, PSPO05173A, Final Report- 2017 Brick Claims: Comparison of Distal Tip Sealing Performance Between ENSEAL™ X1 Large Jaw (NSLX120L) with Algorithm and Ligasure™ Impact (LF4318), Sept 2016, Data on File (154862-201002) 2. Ethicon, PSB004548 Final Report Brick Distal Tip Sealing Algo D, Sept 2016, Data on File (116497-190612) 3. Ethicon, PSPO05819, Final Report- 2017 Brick Claims: COMPARISON OF DISTAL TIP SEALING PERFORMANCE BETWEEN ENSEAL™ X1 LARGE JAW (NSLX120L) AND LIGASURE IMPACT (LF4318), Sept 2017, Data on File. (116497-190612) 4. Ethicon, PSB004570 Final Report Brick Thermal Damage vs Impact1, April 2018, Data on File (114378-190516) 5. Ethicon, DOC022221, X1 Large Jaw - better design and ergonomics vs Impact, July 2016, Data on File (130163-200102, 116496-190612) 6. Ethicon, DOC022209A Brick - Superior Distal Tip Sealing Rationale, Aug 2016, Data on File (124114-190923)

Expect more with ENSEAL™ X1 devices


Intelligent energy delivery

Adaptive Tissue Technology, powered by the Ethicon GEN11 Generator, uses an advanced algorithm for intelligent and efficient energy delivery. In ENSEAL™ X1 devices, it continuously:

- **Senses** changes in tissue and device conditions
- **Responds** with the optimal amount of energy

ENSEAL™ X1 Curved Jaw Tissue Sealer produces minimal lateral thermal spread (2.5mm)^{1*}

The Intelligence of Adaptive Tissue Technology



Advancing beyond secure sealing

- ✓ ENSEAL™ X1 Tissue Sealers can seal vessels up to and including 7mm and lymphatics²
- ✓ Vessels sealed with ENSEAL™ X1 Curved Jaw had an average burst pressure of more than 8x normal systolic blood pressure^{3**}
- ✓ Silicone coating minimizes tissue sticking^{4#}

*Mean thermal spread measured via histology on porcine carotid arteries. **In benchtop testing on porcine arteries, average burst pressure was 1055 mmHg.

#Preclinical testing that compared average sticking force (lb) of ENSEAL X1 Large Jaw end-effector coated with and without non-stick silicone (p<0.001).

REFERENCES: 1. Ethicon, PSPO07612A, Evaluation of ENSEAL™ X1 Tissue Sealer, Curved Jaw, in an Acute Porcine Model, Feb 2020, Data on File (145177-200630) 2. As Per Instructions for Use (130176-200102) 3. Ethicon, Floyd Relaunch Claims Ex-Vivo Sealing, June 2020, PRC094697A (145156-200630) 4. Ethicon, PRC074959A, Brick EB2 Silicone versus No Silicone Testing, Sept 2016, Data on File (123391-190913)

ENSEAL™ X1 devices feature ergonomic engineering



- **Intuitive design**¹ Separate seal and cut functionality²
- **Conveniently placed control buttons** are designed for less hand movement³
- ENSEAL™ X1 Tissue Sealers have **360° shaft rotation** to enable access to targeted tissue^{4,5}

DESCRIPTION	PRODUCT CODE	SHAFT LENGTH (cm)	SHAFT DIAMETER (mm)	QUANTITY/SALES UNIT
ENSEAL™ X1 Curved Jaw	NSLX125C	25	5	3
ENSEAL™ X1 Curved Jaw	NSLX137C	37	5	3
ENSEAL™ X1 Curved Jaw	NSLX145C	45	5	3
ENSEAL™ X1 Large Jaw	NSLX120L	20	13	6

- ENSEAL™ X1 Curved Jaw and ENSEAL™ X1 Large Jaw are supplied sterile for single-patient use
- Both devices are compatible with the Ethicon GEN11 Generator (software version 2016-1 or later versions)

How to order

Please contact your local Ethicon Sales Representative.

Customer support

Visit www.jnjmedicaldevices.com for more information about the ENSEAL™ X1 Tissue Sealers.

References: **1.** Ethicon, DOCO23555A, A Floyd Claims Memo - Industrial Design, July 2018, Data on File (155022-201006) **2.** (113402-190503) **3.** Ethicon, DOCO23555 - Floyd Claims Industrial Design Memo, July 2018, Data on File (116498-190612) **4.** As Per Instructions For Use (152441-200909) **5.** Ethicon, DOCO23555A, A Floyd Claims Memo - Industrial Design, July 2018, Data on File (152441-200909)

Please refer always to the Instructions for Use / Package Insert that come with the device for the most current and complete instructions.

The third-party trademarks used herein are the property of their respective owners.

Ethicon Endo-Surgery (Europe) GmbH
Hummelsbütteler Steindamm 71
22851 Norderstedt, Germany

© Ethicon Endo-surgery (Europe) GmbH 2020, 103873-201014 EMEA

ETHICON
PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

Shaping
the future
of surgery

Enseal™