

Harmonic®

Experience standard of **precision and performance** for your procedures



HARMONIC ACE®+ 7 Shears

- ✓ Unites **precision** with **powerful sealing** ability
- ✓ Provides **more multifunctionality**, which may reduce the number of surgical devices needed during a procedure
- ✓ Adaptive Tissue Technology **delivers energy intelligently** and responds to changing tissue conditions

VS.

Olympus Thunderbeat™

- ✗ **Does not offer ultrasonic only mode**
 - Limits multifunctionality
 - Restricted ability to back-score, use cavitation for tissue plane dissection, create otomies and spot coagulate
- ✗ **Does not offer secondary hemostasis and spot coagulation**
 - Advanced Bipolar energy can only be applied with the jaws closed

ETHICON
PART OF THE *Johnson & Johnson* MEDICAL DEVICE COMPANIES

Shaping
the future
of surgery


HARMONIC ACE®+ 7 Shears vs. Olympus Thunderbeat™




Heat and thermal management

The overall temperature profile of the distal shaft of **Thunderbeat** is **hotter**¹

- HARMONIC ACE+ 7 provides superior heat management¹—the distal part of the Thunderbeat shaft is 39% hotter.²
- HARMONIC ACE+ 7 has 27% less lateral thermal damage in MIN mode compared to Thunderbeat.³



Distal shaft temperatures




Seal quality and burst pressure

The overall quality⁴ and **durability of seals produced with Thunderbeat** is **lower**

- HARMONIC ACE®+7 Shears delivers greater first pass vessel sealing than Thunderbeat⁴
- HARMONIC ACE+ 7 median burst pressure is 35% higher than Thunderbeat when sealing vessels 5-7mm in the Advanced Hemostasis mode.⁵

5-7mm vessel burst pressure⁵




Reliability


62.5% Thunderbeat instruments visibly sparked prior to an average of 112 cuts during testing⁶

- HARMONIC ACE+ 7 pad life is as much as 28x greater than that of Thunderbeat.⁷

Pad life testing



When the Thunderbeat pad breaks, you may generate sparking in the jaws



For more information, contact your local Ethicon Product Specialist

¹ In a benchtop study on porcine jejunum, Thunderbeat™ (Seal and Cut mode) exhibited 38.6% higher mean shaft temperature (Celsius) (vs. Ethicon ACE devices with Adaptive Tissue Technology at Max Power Level 5 (p<0.001). (C1935)

² In a benchtop study on porcine jejunum, Thunderbeat™ (Seal and Cut mode) exhibited 38.6% higher mean shaft temperature (Celsius) (vs. Ethicon ACE devices with Adaptive Tissue Technology at Max Power Level 5 (p<0.001). (C1934)

³ In a preclinical study on ≤5mm porcine carotids in MIN mode vs. TB-0535FC (Seal and Cut Mode). Mean damage compared and determined via histology (p<0.001). (C1896)

⁴ In benchtop test on 5-7mm porcine carotids comparing seal leaks at transection, HARMONIC ACE®+7 in Advanced Hemostasis mode (2/152 failures) versus Thunderbeat™ (5/78 failures) (p=0.046). Data on file (PRCO64872B). (C1638)

⁵ In a benchtop study with 5-7mm porcine carotids that compared median burst pressure, HARMONIC ACE®+ 7 vs. Thunderbeat (Seal and Cut mode) (p=0.05). (C1871)

⁶ In a benchtop study on porcine jejunum that compared HARMONIC ACE®+ (MAX mode only) vs. THUNDERBEAT (TB-0535FC) (n=8 for each group). Devices used per IFU until failure of pad. Cuts to failure for THUNDERBEAT: 57°, 70°, 90°, 140°, 150, 205°, 257, 270 (* denotes sparking at failure). No HARMONIC ACE®+ devices sparked. Data on file, Ethicon (PRCO64789 & SCNO33130C). (C1583)

⁷ In a benchtop study on porcine jejunum that compared HARMONIC ACE®+ (MAX mode) vs. THUNDERBEAT (TB-0535FC; Seal and Cut mode) (n=8 for each group). Devices used per IFU until failure of pad (means: 1080 (+/-649) vs 155 (+/-82.3) (p=0.0005), respectively). Pad life for ACE+ was 112 cuts, vs 4 cuts for Thunderbeat via Weibull reliability analysis (99% reliability, 95% confidence: 112/4=28). The probability that ACE+ will last longer than Thunderbeat is 99.92% (p=0.0008), via reliability and life comparison in Reliasoft. Data on file, Ethicon (PRCO64789 & SCNO33130C). (C1678)