

# Summary

## Catheter Ablation or Medical Therapy to Delay Progression of Atrial Fibrillation: The Randomized Controlled Atrial Fibrillation Progression Trial (ATTEST)

Kuck K-H, Lebedev DS, Mikhaylov EN, Romanov A, Gellér L, Kalējs O, Neumann T, Davtyan K, On YK, Popov S, Bongiorno MG, Schlüter M, Willems S, Ouyang F. *Europace*. Published online December 17, 2020.  
[doi:10.1093/europace/euaa298](https://doi.org/10.1093/europace/euaa298).

## STUDY QUESTION

Is RF catheter ablation more effective than AAD therapy alone in delaying the progression to persistent AF in patients with symptomatic paroxysmal AF?

## METHODOLOGY

### DESIGN:

PROSPECTIVE, MULTICENTER,  
RANDOMIZED, CONTROLLED,  
2-ARM OPEN-LABEL STUDY



**29**  
CENTERS

**255**  
PATIENTS  
≥60 years old with  
paroxysmal AF

### INTERVENTIONS

#### Standard-of-Care AF Management



**PVI with irrigated RF  
catheters**  
(n=128)

**VS**



**AAD alone**  
(n=127)

## OUTCOMES

### PRIMARY OUTCOME



**PROGRESSION TO PERSISTENT AF/AT**

### ADDITIONAL OUTCOMES



**FACTORS ASSOCIATED WITH  
DISEASE PROGRESSION**

## RESULTS

**Progression to persistent AF/AT was significantly  
lower in the RF ablation group at 3 years**



**2.4%**  
95% CI: 0.6–9.4%

**VS**



**17.5%**  
95% CI: 10.7–27.8%

**P=0009**

**Factors associated with persistent AF/AT  
disease progression**

**10x**

**Less likely** to develop persistent AF/AT  
for patients treated with **RF ablation**  
Hazard ratio: 0.107 (95% CI: 0.024–0.47; *P*=0.0031)

**4x**

**More likely** to develop persistent AF/AT  
for **patients ≥65 years old**  
Hazard ratio: 3.87 (95% CI: 0.88–17.00; *P*=0.0727)

## CONCLUSION

**RF ablation is superior** to guideline-directed AAD therapy alone  
in **delaying progression** from paroxysmal to **persistent AF**.

## OBJECTIVE

To evaluate if RF ablation is more effective than AAD therapy alone in delaying the progression to persistent AF over the course of 3 years in patients  $\geq 60$  years old with paroxysmal AF.

## METHODS

### Experimental Design

<b>STUDY DESIGN</b>	<ul style="list-style-type: none"> <li>Prospective, multicenter, randomized, controlled, 2-arm open-label study</li> </ul>
<b>STUDY PERIOD AND LOCATION</b>	<ul style="list-style-type: none"> <li>February 2012 to May 2018</li> <li>29 centers in 13 countries</li> </ul>
<b>INCLUSION CRITERIA</b>	<ul style="list-style-type: none"> <li>255 patients <math>\geq 60</math> years old with paroxysmal AF <math>\geq 2</math> years</li> <li><math>\geq 2</math> AF episodes over the 6 months preceding enrollment</li> <li>Failed 1-2 AADs</li> <li>HATCH score 1-4</li> </ul>
<b>PATIENT FOLLOW-UP</b>	<ul style="list-style-type: none"> <li>Standard-of-care monitoring at 3 months, 6 months, 1 year, 2 years, and 3 years</li> </ul>
<b>PRIMARY ENDPOINT</b>	<ul style="list-style-type: none"> <li>First documented occurrence of persistent AF/AT after a 90-day blanking period assessed at 3 years<sup>a</sup></li> </ul>
<b>OTHER OUTCOMES</b>	<ul style="list-style-type: none"> <li>Factors associated with disease progression</li> <li>Time to recurrent AF/AT</li> </ul>

### Interventions Per Standard-of-Care Practice



- Medications were managed for all patients in both groups per current treatment guidelines.
- Crossover from the AAD to the RF group was allowed for treatment-compliant patients who experienced severe side effects or AF symptoms when no other drug options were available.
  - To minimize crossover, noncompliant patients or those taking concurrent incompatible medications underwent medication management and reevaluation before the crossover was authorized.
- Due to budget constraints following slow enrollment, the trial was terminated prematurely by the sponsor in February 2018 after the second planned interim analysis; termination was independent of study outcomes.

AAD, antiarrhythmic drug; AF, atrial fibrillation; AT, atrial tachycardia; ECG, electrocardiogram; FU, follow-up; HATCH score, hypertension=1, age >75=1, transient ischemic attack or stroke=2, chronic obstructive pulmonary disease=1, Heart failure=2; PVI, pulmonary vein isolation; RF, radiofrequency; TTM, transtelephonic monitoring.


<sup>a</sup> Persistent AF was defined as episodes lasting for >7 consecutive days or requiring termination by cardioversion after 48 hours.

<sup>b</sup> RF ablation was performed with noncontact force-sensing and contact force-sensing catheters (NAVISTAR<sup>®</sup> THERMOCOOL<sup>®</sup> Catheter, THERMOCOOL<sup>®</sup> SF NAV Catheter, Bi-Directional Navigation Catheter, THERMOCOOL SMARTTOUCH<sup>®</sup> Bi-Directional Navigation Catheter, THERMOCOOL SMARTTOUCH<sup>®</sup> Uni-Directional Navigation Catheter, NAVISTAR RMT THERMOCOOL<sup>®</sup> Catheter) in conjunction with a 3-dimensional electroanatomic mapping system (CARTO<sup>®</sup> 3, CARTO<sup>®</sup> XP, or CARTO<sup>®</sup> RMT).

<sup>c</sup> Transtelephonic monitoring transmission per specified schedule or whenever patients experienced symptoms. Once AF/AT was identified, daily TTM was initiated for 7 consecutive days.

## RESULTS

### Patient Disposition and Demographics



Demographics and baseline characteristics were generally similar between treatment groups.

PARAMETER	ITT POPULATION	
	RF ABLATION (n=128)	AAD (n=127)
Age, mean (SD), y	67.8 (4.8)	67.6 (4.6)
Male, No. (%)	54 (42.2)	53 (41.7)
Months since first experience of AF, median (range)	51.2 (19–625)	49.8 (25–366)
AF/AT episodes during prior 6 months, median (range)	6.5 (2–180)	6.0 (0–180)
Lone AF	38 (29.7)	39 (30.7)
HATCH score, mean (SD)	1.5 (0.9)	1.7 (0.9)
Left atrial diameter, mean (SD), mm	42.1 (6.1)	43.4 (5.6)
Left ventricular ejection fraction, mean (SD), %	61.8 (5.8)	62.3 (5.2)
AAD Class I/III at baseline	61 (47.7)	69 (54.3)
<b>Medical history, No. (%)</b>		
Hypertension	120 (93.8)	123 (96.9)
Hyperlipidemia/dyslipidemia	67 (52.3)	67 (52.8)
Left ventricular hypertrophy	26 (20.3)	23 (18.1)
Congestive heart failure	24 (18.8)	27 (21.3)
Atrial flutter	15 (11.7)	10 (7.9)
Diabetes	13 (10.2)	14 (11.0)
Transient ischemic attack/stroke	12 (9.4)	8 (6.3)
Cardiomyopathy <sup>a</sup>	6 (4.7)	2 (1.6)
Renal insufficiency	3 (2.3)	4 (3.1)

AAD, antiarrhythmic drug; AF, atrial fibrillation; AT, atrial tachycardia; HATCH score, hypertension=1, age >75=1, transient ischemic attack or stroke=2, chronic obstructive pulmonary disease=1, Heart failure=2; RF, radiofrequency; ITT, intention to treat; PP, per protocol.

<sup>a</sup>Includes ischemic, nonischemic dilated, and hypertrophic obstructive cardiomyopathy.

## RESULTS


### Intervention and Treatment

- Complete PVI was achieved in all patients who underwent RF ablation.
  - Noncontact force-sensing THERMOCOOL® Catheters were the most commonly used catheters (58.8%), followed by contact force-sensing THERMOCOOL SMARTTOUCH® Catheters (31.4%).
- Excluding beta-blockers and calcium channel blockers, 53 of 127 patients in the AAD group (41.7%) and 30 of 128 patients in the RF ablation group (23.4%) were initiated with new Class I/III AADs during the study.

### Primary Efficacy


**The superiority of RF treatment over AAD alone in delaying persistent AF/AT progression was observed at 1-year follow-up and strengthened by the 3-year follow-up.**

Kaplan-Meier Estimate of Persistent AF/AT Progression (ITT Population)



- The primary efficacy analysis on the ITT population was confirmed for the PP population, As-Treated population, and when using the 2017 Heart Rhythm Society expert consensus definition of persistent AF.<sup>1</sup>

Kaplan-Meier Estimate of Persistent AF/AT Progression at 3 Years




AAD, antiarrhythmic drug; AF, atrial fibrillation; AT, atrial tachycardia; HRS, Heart Rhythm Society; ITT, intention to treat; PP, per protocol; PVI, pulmonary vein isolation; RF, radiofrequency.

## RESULTS

### Factors Associated with AF Progression

Treatment modality and age were associated with AF progression in the PP population, suggesting that early RF ablation can delay disease progression.


#### Cox Model With Multiple Baseline Covariates and Treatment as the Time-Dependent Covariate



### Time to Recurrent AF/AT

The incidence of recurrent AF/AT was consistently lower with RF ablation than with AAD treatment from 6 months through the the end of the study.

#### Kaplan-Meier Estimate of Time to Recurrent AF/AT (ITT Population)



AAD, antiarrhythmic drug; AF, atrial fibrillation; AT, atrial tachycardia; ITT, intention to treat; PP, per protocol; RF, radiofrequency.

## CONCLUSION



ATTEST demonstrated that **RF catheter ablation**—as part of standard-of-care AF management including AADs—is **superior to guideline-directed AAD therapy alone** in **delaying the progression to persistent AF** in patients with paroxysmal AF.

- Patients treated with **RF ablation** were **significantly less likely** to **develop persistent AF or persistent AT** than patients treated with AADs.
- **Patients  $\geq 65$  years** were **more likely to progress to persistent AF/AT** than patients  $< 65$ .
- These results suggest that **early RF ablation** may be an **effective treatment strategy** for **delaying AF progression**.

## REFERENCE

1. Calkins H, Hindricks G, Cappato R, et al. 2017 HRS/EHRA/ECAS/APHRS/SOLAECE expert consensus statement on catheter and surgical ablation of atrial fibrillation. *Heart Rhythm* 2017;14:e275–e444.

AAD, antiarrhythmic drug; AF, atrial fibrillation; AT, atrial tachycardia; RF, radiofrequency.

THERMOCOOL® Navigation Catheters are indicated for the treatment of drug refractory recurrent symptomatic paroxysmal atrial fibrillation, when used with CARTO® 3 Systems (excluding NAVISTAR™ RMT THERMOCOOL™ Catheter).

In the US, Biosense Webster, Inc. THERMOCOOL® Non-Navigational Catheters are indicated for the treatment of Type I Atrial Flutter in patients 18 years of age or older.

Important information: Prior to use, refer to the instructions for use supplied with the device for indications, contraindications, side effects, warnings and precautions. Caution: US law restricts this device to sale by or on the order of a physician.

### **Biosense Webster, Inc.**

31 Technology Drive  
Irvine, CA 92618 USA  
Tel: +1-909-839-8500  
Tel: +1-800-729-9010

[www.biosensewebster.com](http://www.biosensewebster.com)

