LAB100947444v2

en **GEL BREAST IMPLANTS**

MENTOR™ SILTEX™ ROUND GEL BREAST IMPLANTS

MENTOR™ SMOOTH ROUND GEL BREAST IMPLANTS

MENTOR™ CPG™ GEL BREAST IMPLANTS

August 2025,

LAB100947444v2

ar	MENTOR™ غرسات الثدي السيليكونية CPG™ من شركة	it	Protesi mammarie riempite con gel MENTOR™ CPG™
bg	CPG™ гел гръдни импланти MENTOR™	lt	„MENTOR™“ CPG™ geliniai krūtų implantai
cz	Gelové prsní implantáty MENTOR™ CPG™	lv	MENTOR™ CPG™ apaļjie gēla krūšu implanti
da	MENTOR™ CPG™ gelfyldte brystimplantater	nl	MENTOR™ CPG™ gelborstimplantaten
de	MENTOR™ CPG™ Gel-Brustimplantate	no	MENTOR™ CPG™ gelbrystimplantater
el	MENTOR™ CPG™ Gel Breast Implants	pl	Żelowe implanty piersi CPG™ MENTOR™
es	Implantes mamarios de gel CPG™ MENTOR™	pt	Implantes mamários de gel MENTOR™ CPG™
et	MENTOR™ CPG™ geeliga rinnaimplantaadid	ro	Implanturi mamare cu gel MENTOR™ CPG™
fi	MENTOR™ CPG™ geelitäytteiset rintaimplantit	sk	Gélové prsné implantáty MENTOR™ CPG™
fr	Implants mammaires en gel MENTOR™ CPG™	sl	Gelasti prsni vsadki MENTOR™ CPG™
hr	Implantati za dojke MENTOR™ CPG™ gelom	sr	MENTOR™ CPG™ gel implantati dojke
hu	MENTOR™ CPG™ gél mellimplantátumok	sv	MENTOR™ CPG™ gelfyllda bröstimplantat
is	MENTOR™ CPG™ hlaupbrjóstaígræði	tr	MENTOR™ CPG™ Yuvarlak Jel Meme İmplantları

Only to be **implanted** in an appropriate medical environment by appropriately **trained medical doctors** who are **qualified** or **accredited** in accordance with **national law**

ENGLISH

GEL BREAST IMPLANTS

MENTOR™ SILTEX™ ROUND GEL BREAST IMPLANTS

MENTOR™ SMOOTH ROUND GEL BREAST IMPLANTS

MENTOR™ CPG™ GEL BREAST IMPLANTS

PATIENT INFORMATION LEAFLET

August 2025, LAB100947444v2

1. Device Description

MENTOR™ Gel Breast Implants is constructed of silicone. MENTOR™ offers Gel Breast Implants with two types of shell surfaces: SILTEX™ and smooth surfaced. The SILTEX™ shell is textured. All MENTOR™ Gel Breast Implants contain silicone gel which is cohesive. MENTOR™ offers Gel Breast Implants with a variety of cohesivity levels of the filling material. The devices are available in a round shape with different projections and in several contour shapes with different heights and projections. The volume indicated on the product label is the fill volume of the gel.

2. Device Materials

MENTOR™ Gel-Filled Breast Implants have been tested for safe use. Materials to which patients can come in contact with, are 100% medical grade silicone.

The shells of the MENTOR™ Gel-Filled Breast Implants are constructed from medical-grade silicone elastomer. The shells are filled with 100% medical-grade, cohesive silicone gel.

3. Information for Safe Use

3.1 General Safety Instructions for patient

Breast implants are for women 18 and older who want breast augmentation or reconstruction.

Breast implants surgery is not for women who have any of the following conditions:

- Pregnancy or nursing mothers.
- Autoimmune disease
- Conditions that would increase surgical risk*
- A history of sensitivity to foreign materials
- Repeated attempts and failures at breast augmentation or reconstruction.
- An unwillingness to have further surgery
- Unrealistic expectations
 - inappropriate attitude or motivation
 - lack of understanding of the risks involved with the surgery and implants.
- Untreated breast cancer or premalignant breast disease *

* Discuss with your surgeon any possible conditions

**Speak to your breast surgeon and oncologist about your cancer treatment plan

3.2 Postoperative Warnings, Precautions, or Measures

Breast Self-Examination Measures

You should learn to perform breast self-examinations monthly. You should also be shown how to distinguish the implant from their breast tissue

Warnings

- Silicone gel can leak or “bleed” through the undamaged shell into the capsule and adjacent breast tissue. The long-term effects of such “bleed” are unknown. You should be aware of this potentiality.
- The use of microwave diathermy is NOT recommended. It may cause negative side effects to your surgical area.
- You should be aware that any abnormal stress or trauma to the breast could result in damage to the implant.
- MENTOR™ strongly recommends AGAINST the treatment of capsule firmness by forceful external stress (such as closed capsulotomy). It has been reported to cause patient injury and/or implant damage.
- As with any surgery, smoking may interfere with the healing process after breast implant surgery

Precautions

Advise your Healthcare Provider about your existing implants if

- you have an autoimmune disease
- have a weakened immune system
- currently taking drugs that weaken the body's natural resistance to disease.

Advise your Healthcare Provider if you are planning any of the following:

- Surgery in the breast area
- Chemotherapy
- Radiation*
- MRI
- Breastfeeding**
- Reoperation***.
- Mammography****

* Radiation treatment may possibly cause damage to the breast implant. Radiation therapy may increase the chance of capsular contracture and/or change how your breasts look

** Breast implants may interfere with breastfeeding

*** Breast implants are not lifetime devices. You will likely undergo implant removal, with or without replacement, in the future. Changes to your breast are irreversible after explant.

**** Breast implants can interfere with the detection of cancer/calcium deposits

All surgical procedures have a small risk of complications such as bleeding, potential drug interactions and the risks associated with anaesthesia. Please discuss these with your surgeon.

Risks and Complications

Risks	Measures for Reducing Risks
Infection*	<ul style="list-style-type: none">• MENTOR™ device and packaging are sterile when shipped.• Precautioning that infection prior to implant should be treated.• Contraindication for implanting where there is an infection in the body.• Using aseptic technique when filling the implant with saline.
Hematoma formation	<ul style="list-style-type: none">• Warning to the surgeon that careful stoppage of blood flow is important to prevent hematoma.• Warning to the surgeon that recommends that if excessive bleeding persist, the bleeding should be controlled. Then proceed with the surgery.
Delayed wound healing	<ul style="list-style-type: none">• Not implanting where there is a history of compromised wound healing.• Warning that smoking may interfere with the healing process.• * See measures for reducing the risk of infection.
Capsular Contracture (CC)	<ul style="list-style-type: none">• Warning that radiation therapy may cause Capsular Contracture.• Warning that Capsular Contracture may occur around any implant placed in contact with soft tissues.

Risks	Measures for Reducing Risks
Chest wall deformity	<p>Implant Selection guidelines:</p> <ul style="list-style-type: none"> • The implant should not be too small or too large in comparison to the patient’s chest-wall dimensions. • Available tissue must provide adequate coverage of the implant. • Submuscular placement of the implant may be preferable in patients with: <ul style="list-style-type: none"> ➤ Thin ➤ Poor quality tissue. • A well-defined dry pocket of adequate size and symmetry must be created to allow the implant to be placed flat on a smooth surface.
Enlarged lymph nodes	<ul style="list-style-type: none"> • Using aseptic techniques. • Rinsing clean surgical instruments and gloves to remove impurities. • * See measures for reducing the risk of infection.
Rupture	<ul style="list-style-type: none"> • Implant meets ISO safety standards. • Warning to the surgeon that extreme care should be taken not to damage the prosthesis with surgical instruments. • Warning the patient that any abnormal stress or trauma to the breast could result in damage to the implant.
Changes in Nipple and Breast Sensation/Breast Pain	<p>Any surgery that undermines the skin of the breast can result in hypersensitivity or loss of sensation in the nipple-areolar region. These changes can vary in degree and may be temporary or permanent. Changes in nipple/breast sensation may, on occasion, affect sexual response or comfort while nursing. These changes are believed to be a result of nerve damage or stretching of the nerves from the surgery. There is no specific treatment for this condition.</p> <p>Surgeon considers the size of the implant and the firmer nature and higher profile of the SILTEX™ shell when choosing optimum incision size and surgical approach.</p> <p>Implant Selection guidelines:</p> <ul style="list-style-type: none"> • The implant should not be too small or too large in comparison to the patient’s chest-wall dimensions. • Available tissue must provide adequate coverage of the implant. • Submuscular placement of the implant may be preferable in patients with <ul style="list-style-type: none"> ➤ Thin ➤ Poor quality tissue. • A well-defined dry pocket of adequate size and symmetry must be created to allow the implant to be placed flat on a smooth surface. • * See measures for reducing the risk of infection.
Possible reaction to silicone	<ul style="list-style-type: none"> • Implant meets safety standards. • Questioning patients on: <ul style="list-style-type: none"> ➤ medical history involving allergic reactions to implant materials or ➤ filling agents.
Gel Bleed	<ul style="list-style-type: none"> • Implant meets safety standards. • Warning the surgeon in the IFU of the risk of rupture: that <i>extreme care should be taken not to damage the prosthesis with surgical instruments</i> • Warning the patient that any abnormal stress or trauma to the breast could result in damage to the implant.
Granulomas	<ul style="list-style-type: none"> • Implant meets safety standards. • Warning the surgeon in the IFU of the risk of rupture: that <i>extreme care should be taken not to damage the prosthesis with surgical instruments.</i> • Warning the patient that any abnormal stress or trauma to the breast could result in damage to the implant.

Risks	Measures for Reducing Risks
Extrusion of the Implant	<ul style="list-style-type: none"> • Warning not to use microwave diathermy in patients with breast implants. • Warning that the incidence of extrusion of the prosthesis had been shown to increase when placed in injured areas that are: <ul style="list-style-type: none"> ➤ scarred ➤ heavily irradiated ➤ burned tissue ➤ crushed bone areas or ➤ severely reduced by surgery.
Wrinkling	Precautioning that surgeon use current techniques for <ul style="list-style-type: none"> • Measuring the patients • sizing the implant and • conducting the surgery. The surgeon does not underfill or overfill the device.
Thrombosed veins	Warning that thrombosis may result from the placement of foreign objects in the body.
Hypertrophic scarring	IFU advising the surgeon that: <ul style="list-style-type: none"> • Hypertrophic scars are generally confined to the original site. • Hypertrophic scars respond well to scar revision treatment. <ul style="list-style-type: none"> ➤ This may include steroid injections to break down the collagen. ➤ This may also include surgery to revise the position, direction or line of the scar.
Anaplastic Large Cell Lymphoma BIA-ALCL	The IFU describes that you should seek a plastic surgeon to evaluate for BIA-ALCL if you experience: <ul style="list-style-type: none"> • a seroma in the breast more than 12 month after implantation • a new breast mass, or • disease of the lymph nodes.

If you are experiencing any breast pain or discomfort, contact your breast surgeon. Discuss any questions, concerns, or potential side effects with your physician.

Breast Implant-Associated Anaplastic Large Cell Lymphoma (BIA-ALCL)

Breast implants have been associated with the development of a cancer of the immune system (BIA-ALCL). This cancer occurs more commonly in patients with textured breast implants than smooth implants, although rates are not well defined. If you have breast implants, you have a very small, but increased risk of developing breast implant associated anaplastic large cell lymphoma, or BIA-ALCL. BIA-ALCL is not breast cancer. It is a rare type of non-Hodgkin's lymphoma (cancer of the immune system). In most cases, BIA-ALCL is found in the scar tissue and fluid near the implant. However, in some cases, it can spread throughout the body. In the cases that have spread beyond the scar tissue and fluid near the implant, rare cases of death have been reported.

You should seek evaluation by a plastic surgeon to evaluate for BIA-ALCL if you experience delayed onset seroma in the breast (delayed means developing more than 12 months after implantation), a new breast mass or lymphadenopathy (disease of the lymph nodes).

If you have breast implants, you should monitor them and follow your routine medical care. It is not necessary to remove your breast implants if you have no symptoms and you have not been diagnosed with BIA-ALCL.

4. Users and Consumer Exclusions

The intended users are surgeons who are familiar with the currently available techniques for performing breast tissue expansion and breast implantation surgery. Persons that do not meet this criterion are excluded from the use of the device.

Consumers (patients) who may receive this device are females who will undergo breast augmentation (non-medical purpose) or breast reconstruction (medical purpose), or who require tissue expansion and replacement of missing breast tissue.

The use of this prosthesis is contraindicated in patients who have any of the following conditions:

- Pregnancy or nursing mothers.
- Lupus (e.g., SLE and DLE).
- Scleroderma (e.g., progressive systemic sclerosis).
- Currently has a condition that could compromise or complicate wound healing (except reconstruction patients).
- Infection or abscess anywhere in the body.
- Demonstrates tissue characteristics, which are clinically incompatible with implant (e.g., tissue damage resulting from radiation, inadequate tissue, or compromised vascularity).
- Possesses any condition, or is under treatment for any condition which, in the opinion of the consulting physician(s), may constitute an unwarranted surgical risk.
- Anatomic or physiologic abnormality, which could lead to significant postoperative complications.
- A history of sensitivity to foreign materials or repeated attempts and failures at breast augmentation or reconstruction.
- An unwillingness to undergo any further surgery for revision.
- Unrealistic expectations, such as inappropriate attitude or motivation, or a lack of understanding of the risks involved with the surgical procedure and implants.
- Premalignant breast disease without a subcutaneous mastectomy.
- Untreated or inappropriately treated breast malignancy, without mastectomy.

Persons that are contraindicated are excluded from the use of the device

5. **Expected Device Lifetime and Follow-up**

Breast implantation may not be a one-time surgery. After the surgery, doctor visits and additional surgery may be needed. Implant removal with or without replacement may be needed as well. Breast implants may not last a lifetime. A rupture can occur at any time following implantation. Studies have indicated that the gel breast implant can have an expected average lifetime of over 10 years.

You should have periodic imaging (e.g., MRI, ultrasound) of your implants to screen for implant rupture. This does not replace additional imaging that may be required, such as screening mammography for breast cancer. Even if you have no symptoms, you should have your first ultrasound or MRI at 5-6 years after your initial implant surgery. Then every 2-3 years thereafter. If you have symptoms at any time or uncertain ultrasound results for breast implant rupture, an MRI is recommended.

Additional doctor visits would be needed to monitor if rupture of the implant has occurred. Speak to your surgeon about regular basis monitoring screening.

If, during the use of this device or as a result of its use, a serious incident has occurred, please report it directly to your healthcare provider or MENTOR™ (Mentor-PE@its.jnj.com or 1-866-250-5115)

6. **Implant ID Card**

You should have received an International Implant Card from your surgeon. This card will have the catalog and serial number of your implants, along with other information. Carry the card with you and show it to doctors or other healthcare providers when you visit them. It will help them treat you appropriately and protect your breast implants during any medical treatment you need in the future. Your doctor should keep a copy of the Implant ID Card with your medical records.

7. **Patient Information Website**

You can find any updates about the device on this website. Your implant card will contain information about your specific device. This can also be found in your medical record

8. **Volume and Sizes**

Information on the volume and size of the breast implant that you will receive is provided on the implant card following the surgery.

9. Harmonized Standards and Common Specifications Applied

- COMMISSION IMPLEMENTING REGULATION (EU) 2022/2346 of 1 December 2022: laying down common specifications for the groups of products without an intended medical purpose listed in Annex XVI to Regulation (EU) 2017/745 of the European Parliament and of the Council on medical devices
- EN ISO 13485 Medical device – Quality management system – Requirements for regulatory purposes
- EN ISO 14971 Medical devices – Application of risk management to medical devices
- EN ISO 14630 Non-active surgical implants – General requirements
- EN 62366 Medical Devices- Application of usability engineering to medical devices
- EN ISO 10993-1 Biological evaluation of medical devices – Part 1: Evaluation and testing within a risk management system
- EN 556-1 Sterilization of medical devices – Requirements for medical devices to be designated “STERILE” – Part 1: Requirements for terminally sterilized medical devices
- EN ISO 11737-1 Sterilization of medical devices - Microbiological methods - Part 1: Determination of a population of microorganisms on products
- EN ISO 11737-2 Sterilization of medical devices - Microbiological methods - Part 2: Tests of sterility performed in the definition, validation and maintenance of a sterilization process
- EN ISO 11607-1 Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems
- EN ISO 11607-2 Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes
- EN ISO 14644-1 Cleanrooms and associated controlled environments - Part 1: Classification of air cleanliness
- EN ISO 14644-2 Cleanrooms and associated controlled environments - Part 2: Specifications for testing and monitoring to prove continued compliance with ISO 14644-1 (ISO 14644-2:2000)
- EN ISO 14644-5 Cleanrooms and associated controlled environments - Part 5: Operations (ISO 14644-5:2004)

10. Symbols Glossary

	Patient name or patient ID		Catalogue number		LEFT
	Date of Implantation		Manufacturer		RIGHT
	Name and Address of implanting healthcare institution / provider		Unique Device Identifier		
	Information website for patients		Serial number		
	Device name		MR Safe		

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© MENTOR Worldwide, LLC 2025
MENTOR Medical Systems B.V.
Johnson & Johnson Medical Ltd