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is

en

ATTUNE™ Medial Stabilized Fixed Bearing Knee

1. Implant name (MD)

ATTUNE™ Medial Stabilized Fixed Bearing Knee

2. Implant name (MD) and number (REF)

The name and number of your knee replacement implant are listed on your patient implant card.

This leaflet has information about your implant. It does not contain all the information and if you have any questions talk to your medical team (doctor).

All medical implants have risks and benefits. Follow your doctor's advice even if it is different to this leaflet.

Please read this leaflet carefully and keep it in a safe place.

If a doctor (or medical team member) asks about your implant, please show them your implant card.

3. Implant Description

Your knee implant is intended to replace your natural knee joint. Like your natural knee joint, the implant works with the tissue and bone around your knee.

- Total knee replacement: this is where the ends of the shin bone and thigh bone are replaced. Part of your kneecap may also be replaced.
- Revision knee replacement: this is where part(s) of your knee replacement is changed during an additional operation.

Your doctor will choose the implants to meet your medical needs.

4. Implant Materials

To identify the material information for your implant match the **MATERIAL** symbol on your implant card to the material below.

Additional trace substances may be present due to the manufacturing process. Please contact the address below if additional information is required. Refer to the implant card for the material (matl) used in your device.

Material identification on implant card label (Material:)	Material Description:	Material nominally contains the following substances
P2	Ultra-High Molecular Weight Polyethylene	Polyethylene (99%) Less than 1% of the antioxidant Pentaerythritol tetrakis [(3-(3,5-di tertiary butyl-4-hydroxyphenyl) propionate)]

Name	Material identification on implant label (Material:)	Material Description:	Percentage (w/w) depending on the size
ATTUNE Medial Stabilized (MS) Tibial Insert	P2	Ultra-High Molecular Weight Polyethylene	100%

5. Information for Safe Use

As part of your operation you should have received a set of instructions regarding exercises, therapies and any limitations on your activities.

It is very important that you follow your doctor's instructions about how to recover and restart activities. Make sure you attend all appointments.

Healing takes time and your doctor will provide information on what to expect.

Not following your doctor's advice may result in complications and the need for additional operations.

Images of your implant are often needed after your operation. X-rays and/or CT scans of your implant will not damage your implant or increase your risk of complications. If you are asked to go for an MRI scan please show your implant card to the doctor and MRI technician. The MRI scan may interact with your implant and may cause complications if your technician is unaware of your implant. Your doctor can locate additional information for this MR conditional device in the Instructions for Use (IFU) provided with the device. This will allow them to manage your MRI scan safely.

6. Possible Side Effects/Risks:

Information about the side effects of your operation are available from your doctor. All operations have risks.

The risk of a serious complication is low. There is a risk that you may require additional operations or treatments for a variety of reasons.

Please talk to your doctor if you have concerns about changes in or near your knee joint or leg.

Possible risks may include:

- Pain, dislocation of the joint, stiffness or uneven leg length.
- Implant loosening change in position or wear.
- Bone loss, fracture or perforation.
- Infection or delayed wound healing.
- Allergic response or adverse biological response.
- Implant breaks, bends, cracks or separates.
- Permanent or short-term nerve damage or tissue damage.
- Heart attack, blood clots, or stroke.

These risks may require additional operations or treatments or may affect the other joints or parts of your body. Another operation may not be possible in all patients. This list does not include all risks.

Your doctor can further explain the risks of your operation.

7. Expected Implant Lifetime and Follow Up

Your implant is designed to remain in your body.

Having a knee joint replaced is a major operation. Most people have a good result, but results vary and you may have complications. Factors such as your anatomy (your size and shape), medical condition, lifestyle (for example, weight and activity level) and surgery affect the outcome.

Implant lifetime is the time from implantation of your knee replacement to when part (or all) of your knee replacement is removed from your body.

Some countries or regions have knee replacement registries which monitor the performance of knee implants. This registry data was analysed and shows that for a person's first total knee operation generally 95% (95 out of 100) last more than 10 years. This means that at 10 years about 5% (5 out of 100) of patients may have had additional operations to remove part(s). For revision knee replacement generally 80% last more than 10 years. Your implant meets this general implant lifetime.

Your implant lifetime may be longer or shorter than the average. It is not possible to predict if you will have complications that may require additional operations. Your surgeon has access to data published for your implant and will be able to provide more information based on your specific needs. Make sure you attend all your medical appointments.

The lifetime of this implant is studied and the results can be found on this website: ec.europa.eu/tools/eudamed. Read Part B of the Summary of Safety and Clinical Performance report. Please note the EUDAMED website will only be live after the EUDAMED database is live.


If this leaflet is updated information will be available on this website: ic.jnjmedicaldevices.com

Symbols

Patient name

Date of implantation

Name and Address of the Implanting Health Institution/Health Care Provider

Name and Address of the manufacturer

Information website for patients

Lot Number/Batch Code – identifies a group of products made at the same timepoint

Unique Device Identifier – identifies your specific implant information

Catalogue number

Implant (medical device) name

Material Composition of Implant

P2

Ultra-High Molecular Weight Polyethylene

DePuy Ireland UC
Loughbeg, Ringaskiddy
Co. Cork
Ireland
Tel: +353 214914000
Fax: +353 214914199

Issued 2024/04
© DePuy Synthes 2024. All rights reserved.

This document is valid only on the date printed. If unsure of the print date, please re-print to ensure use of the latest revision of the PIL (available at www.ic.jnjmedicaldevices.com). The onus resides with the user to ensure that the most up-to-date PIL is used.

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