


FEMORAL NECK SYSTEM (FNS)

Value Analysis Brief




FEMORAL NECK SYSTEM (FNS)

Clinical & Economic Burden

Epidemiology of Femoral Neck Fractures

A majority of hip fractures occur in the femoral neck.¹ In older people, femoral neck fractures (FNF) are the most common traumatic injuries which can lead to severe disability.² Patients with FNF have been on the rise in the last three decades, and the incidence rate of FNF will be continuously increasing in the coming 30 years.³ North America has the highest incidence of FNF in the world at 201 (per 100,000) per year in men and 511 (per 100,000) in women.⁴



Economic Burden

An estimated 1.66 million FNFs have been reported worldwide, and this figure is expected to approach 6.26 million by 2050.⁵

In the United States, the economic burden of hip fractures is amongst the top 20 most expensive diagnoses, with approximately \$5.96 billion spent annually.⁵

The increasingly aging population of the world will face the challenge of coping with the growing number of FNFs and associated economic burden to the healthcare system.⁴

Unmet Need

High Rate of Reoperation
Reoperation rate up to **34%**⁸

The complex block contains two circular icons. The left icon shows a person lying in a hospital bed with a medical monitor. The right icon shows two people: one standing and one sitting in a wheelchair. Between the icons is the text 'High Rate of Reoperation' and 'Reoperation rate up to 34%' with a superscripted 8.

Significant Loss of Function

Decreased mobility & independence in function, health, quality of life and higher rates of institutionalization.⁷

20-60% required assistance for various self-care tasks 1 and 2 years after fracture.⁷

Increased Cost of Treatment
High risk of **reoperation & mortality**⁸ in the elderly which may lead to overall higher treatment cost.⁸

The complex block features a circular icon of a hand holding a stack of money with a dollar sign. To the right of the icon is the text 'Increased Cost of Treatment' and 'High risk of reoperation & mortality' with a superscripted 8, followed by 'in the elderly which may lead to overall higher treatment cost.' with a superscripted 8.

Need For an Improved Solution

Surgical treatment options for FNFs consist of internal fixation, hemiarthroplasty (HA) and total hip arthroplasty (THA).

Internal fixation, including multiple cannulated screws (MCS) or sliding hip screws (SHS), is often a method chosen for patients with non-displaced fractures.⁹ Unfortunately, each of these options have their drawbacks.

Overall rates of reoperation



Mechanical failure rates



Mean operating time




Mean length of stay



■ MCSs □ SHS


Common Clinical Complications[‡]



① Varus Collapse with MCS

MCS have been shown to lack the mechanical stability of sliding hip screws, as they do not provide a fixed angle with additional fixation into the femoral shaft.^{10,11}


This lack of stability is often associated with higher rates of reoperation, which can be as high as 13% due to mechanical failure.^{10,11}



② High Infection Rate with SHS

While SHS offer greater stability when compared to MCS, they require a more invasive approach for implant insertion due to the size of the implant and surgical technique.

This may ultimately result in a larger drop in hemoglobin levels and longer hospital stays.^{9,12,13}



③ Pain Due to Lateral Implant

Lateral protrusion can either occur when the implant moves laterally while the femoral neck is shortening during fracture healing, or when the side plate protrudes from the side of the hip.

In either case, it often results in lateral thigh pain.

Rates of lateral protrusion have been shown to be as high as 5.3% and 3.6% for MCSs and SHSs respectively.^{12,14}

Current evidence suggests²⁻⁷ there is a **need for a solution** that combines the **angular stability of SHSs** with the **minimal invasiveness of MCSs** while **reducing lateral thigh pain and procedural complexity.**

[‡] Percentages are taken from the results of the cited literature. Other publications may report different results.

Femoral Neck System Solution



① Stability^{15,16}

The FNS was designed to provide higher mechanical stability than MCS.

These FNS design features are intended to reduce varus collapse and rotational failures, potentially reducing reoperations due to mechanical instability to a similar level as SHSs.

BENEFITS:

These FNS design features are intended to reduce varus collapse and rotational failures, potentially reducing reoperations due to mechanical instability to a similar level as sliding hip screws.



② Minimally Invasive¹⁷

With its compact design, the FNS was created to minimize implant footprint on the bone. Furthermore, the FNS was designed to reduce the length of incision necessary for implant insertion when compared to a SHS system.

FNS may help reduce blood loss and length of stay, potentially reducing reoperations due to invasiveness to a similar level as MCS.


BENEFITS:

FNS may help reduce blood loss and length of stay, potentially reducing reoperations due to invasiveness to a similar level as MCSs.



71% reduction in footprint compared to SHS¹⁷

60% reduction in incision size compared to SHS¹⁷




③ Reduced Protrusion¹⁷

The bolt design allows the FNS to freely glide within the barrel of the base plate. This allows for 20 mm of controlled collapse of the head fragment, with no lateral protrusion for the first 15 mm.

BENEFITS:

This FNS design feature is intended to reduce incidences of lateral thigh pain.



*Pre-clinical test data is not necessarily indicative of clinical performance

Clinical Value

Lower Complications & Improved Efficacy

He et al.¹⁸ compared the clinical efficacy of the FNS and cannulated screws (CS) in the treatment of FNF in young adults. The FNS had lower rates of overall complications compared to CS (25.0% [9/36] vs. 6.1% [2/33], p=0.032).

Moreover, compared with CS, the FNS reduced the number of intraoperative fluoroscopies (10.58 ± 1.89 vs 18.33 ± 3.82, p<0.000), radiation exposure to medical staff and patients, and short-term complications including femoral neck shortening (3% vs 8.3%) and bone nonunion (0% vs 5.6%).¹⁸

Clinical Efficacy

Cintean et al. compared the outcome of elderly patients with multiple comorbidities with non-displaced FNFs treated with a hemiarthroplasty or the FNS. For the treatment of non-displaced[#] FNFs, FNS use was associated with faster recovery with shorter hospitalization time (9.8 ± 3.8 vs 15.1 ± 5.1 days, P < 0.001) compared to patients treated with hemiarthroplasty.¹⁹

In addition, the study by Zhou et al. compared the early clinical effects of the FNS and three cannulated screws for the treatment of patients with unstable FNFs. Results of the study showed that compared to cannulated screws, FNS

FNS Compared to MCS^{18,20,21}



Lower rates of overall complications



Lower rates of non-union



Lower rates of femoral neck shortening



Fewer intraoperative fluoroscopies



Lower rates short-term complications

demonstrated lower postoperative Visual Analog Scale (VAS) scores (3.13 ± 1.07 vs 3.77 ± 1.04, P = 0.018), and earlier time to walking without crutches (5.23 ± 1.33 vs 6.03 ± 1.45, P < 0.001).^{20§}

Further, Stoffel et al. prospectively assessed the rate of early implant-related complications in FNF osteosynthesis using the FNS. Patients treated with FNS regained their preinjury function and quality of life within one year.^{21*}

FNS demonstrated better clinical efficacy such as lower postoperative Visual Analog Scale (VAS) scores^{20§}



* Harris hip score (HHS), Pre-injury =100 (90.2–100.0), 12 months= 94.4 (70.9–97.0); EQ-5D-5L index score, Pre-injury = 1.0 (0.83–1.00), 12 months= 0.9 (0.79–1.00)

nondisplaced (Garden Type 1 and Type 2) versus displaced (Garden Type 3 and Type 4).²²

§ The Visual Analog Scale (VAS) is a validated, subjective measure of acute and chronic pain related to surgical intervention. A higher score indicates greater pain intensity.

Economic Value



Reduced Operative Time

FNS is indicated for patients with FNF. It provides angular and rotational stability.^{15,17,23} Due to the minimized size of the implant, FNS facilitates insertion through reduced incision.¹⁷

Vazquez et al. compared the operation time for three cannulated screws, dynamic hip screws (DHS), and FNS for non-displaced (Garden I and II) FNF osteosynthesis. Mean operation time was significantly shorter with FNS (43.3 ± 10.1 minutes) compared to three cannulated screws (66.9 ± 25.4 minutes) or DHS system (70.7 ± 20.0 minutes, P < 0.001).²⁴

In a study comparing the FNS to a group of other implants which included a compression hip screw (CHS), Hansson Twin Hook (HTH), three 6.5 mm cannulated compression screws (CCS), & Hansson Pins (HP), the average surgical time was significantly shorter for FNS. (42 ± 13 min vs 53 ± 21 min; P = 0.032).²⁵

Significantly Reduced Operative Time With FNS²⁴



SAVINGS OF
\$35k
 PER PATIENT FOR EACH
REVISION SURGERY
AVOIDED²⁸

Rate of Reoperation Reduced

Compared to other implants*, patients treated with FNS have a significantly lower incidence of reoperation (0% vs 22%, P = 0.023).²⁵

Revision procedures resulting from the failed fixation of the hip can be detrimental to the patient, increasing the risk of mortality, decreasing the ability for patients to return to their original state, and resulting in a two to threefold rise in average cost of treatment.^{26,27}



The **Femoral Neck System** was designed with the aim of **reducing the incidence of costly reoperations and complications**.²⁹ This may provide the opportunity for **cost savings** for the health care system.

*Other implants (compression hip screw: CHS, Hansson Twin Hook: HTH, three 6.5mm cannulated compression screw: CCS, and Hansson pins: HP)

References

1. Femoral Neck Fracture: Types, Symptoms, Treatment, and Recovery. Healthline. Published December 11, 2018. Accessed January 13, 2023. <https://www.healthline.com/health/femoral-neck-fracture>
2. Bäcker HC, Wu CH, Maniglio M, Wittekindt S, Hardt S, Perka C. Epidemiology of proximal femoral fractures. *J Clin Orthop Trauma*. 2021;12(1):161-165. doi:10.1016/j.jcot.2020.07.001
3. Peng P, Xiao F, He X, et al. Global Research Status and Trends of Femoral Neck Fracture Over the Past 27 Years: A Historical Review and Bibliometric Analysis. *Front Surg*. 2022;9. Accessed January 13, 2023. <https://www.frontiersin.org/articles/10.3389/fsurg.2022.875040>
4. Filipov O. EPIDEMIOLOGY AND SOCIAL BURDEN OF THE FEMORAL NECK FRACTURES. *J IMAB – Annu Proceeding Sci Pap*. 2014;20(4):516-518. doi:10.5272/jimab.2014204.516
5. Al-Mohrej OA, Alshaalan FN, Aldakhil SS, Rahman WA. One-Year Mortality Rates Following Fracture of the Femoral Neck Treated With Hip Arthroplasty in an Aging Saudi Population: A Trauma Center Experience. *Geriatr Orthop Surg Rehabil*. 2020;11:2151459320922473. doi:10.1177/2151459320922473
6. Adeyemi A, Delhougne G. Incidence and Economic Burden of Intertrochanteric Fracture: A Medicare Claims Database Analysis. *JB JS Open Access*. 2019;4(1):e0045. doi:10.2106/JBJS.OA.18.00045
7. Dyer SM, Crotty M, Fairhall N, et al. A critical review of the long-term disability outcomes following hip fracture. *BMC Geriatr*. 2016;16(1):158. doi:10.1186/s12877-016-0332-0
8. Overmann AL, Richards JT, O'Hara NN, D'Alleyr and JC, Slobogean GP. Outcomes of elderly patients with nondisplaced or minimally displaced femoral neck fractures treated with internal fixation: A systematic review and meta-analysis. *Injury*. 2019;50(12):2158-2166. doi:10.1016/j.injury.2019.09.039
9. Ma JX, Kuang MJ, Xing F, et al. Sliding hip screw versus cannulated cancellous screws for fixation of femoral neck fracture in adults: A systematic review. *Int J Surg Lond Engl*. 2018;52:89-97. doi:10.1016/j.ijssu.2018.01.050
10. Zhang LL, Zhang Y, Ma X, Liu Y. Multiple cannulated screws vs. dynamic hip screws for femoral neck fractures: A meta-analysis. *Orthopade*. 2017;46(11):954-962. doi:10.1007/s00132-017-3473-8
11. Stiasny J, Dragan S, Kulej M, Martynkiewicz J, Płochowski J, Dragan SŁ. Comparison analysis of the operative treatment results of the femoral neck fractures using side-plate and compression screw and cannulated AO screws. *Ortop Traumatol Rehabil*. 2008;10(4):350-361.
12. Murphy DK, Randell T, Brennan KL, Probe RA, Brennan ML. Treatment and displacement affect the reoperation rate for femoral neck fracture. *Clin Orthop*. 2013;471(8):2691-2702. doi:10.1007/s11999-013-3020-9
13. Rockwood and Green's Fractures in Adults. Accessed January 16, 2023. <https://www.wolterskluwer.com/en/solutions/ovid/rockwood-and-greens-fractures-in-adults-839>
14. Enocson A, Lapidus LJ. The vertical hip fracture – a treatment challenge. A cohort study with an up to 9 year follow-up of 137 consecutive hips treated with sliding hip screw and antirotation screw. *BMC Musculoskeletal Disord*. 2012;13(1):171. doi:10.1186/1471-2474-13-171
15. Stoffel K, Zderic I, Gras F, et al. Biomechanical Evaluation of the Femoral Neck System in Unstable Pauwels III Femoral Neck Fractures: A Comparison with the Dynamic Hip Screw and Cannulated Screws. *J Orthop Trauma*. 2017;31(3):131-137. doi:10.1097/BOT.0000000000000739
16. DePuy Synthes Static Cut Through Rotation Test in bone Foam with Cannulated Screws . 11/8/2018. Windchill Document # 0000277853.
17. DePuy Synthes Performance Memo: FNS Design and Procedure Comparison. 5/3/2018. Windchill # 0000274963
18. He C, Lu Y, Wang Q, et al. Comparison of the clinical efficacy of a femoral neck system versus cannulated screws in the treatment of femoral neck fracture in young adults. *BMC Musculoskeletal Disord*. 2021;22(1):994. doi:10.1186/s12891-021-04888-0
19. Cintean R, Pankratz C, Hofmann M, Gebhard F, Schütze K. Early Results in Non-Displaced Femoral Neck Fractures Using the Femoral Neck System. *Geriatr Orthop Surg Rehabil*. 2021;12:21514593211050150. doi:10.1177/21514593211050153
20. Zhou XQ, Li ZQ, Xu RJ, et al. Comparison of Early Clinical Results for Femoral Neck System and Cannulated Screws in the Treatment of Unstable Femoral Neck Fractures. *Orthop Surg*. 2021;13(6):1802-1809. doi:10.1111/os.13098
21. Stoffel K, Michelitsch C, Arora R, et al. Clinical performance of the Femoral Neck System within 1 year in 125 patients with acute femoral neck fractures, a prospective observational case series. *Arch Orthop Trauma Surg*. Published online December 2, 2022. doi:10.1007/s00402-022-04686-w
22. Kazley JM, Banerjee S, Abousayed MM, Rosenbaum AJ. Classifications in Brief: Garden Classification of Femoral Neck Fractures. *Clin Orthop*. 2018;476(2):441-445. doi:10.1007/s11999.0000000000000066
23. DePuy Synthes Test Report Static Cut Through Rotation Test in Bone Foam. 3/6/2019. Windchill Document #0000165855
24. Vazquez O, Gamulin A, Hannouche D, Belaieff W. Osteosynthesis of non-displaced femoral neck fractures in the elderly population using the femoral neck system (FNS): short-term clinical and radiological outcomes. *J Orthop Surg*. 2021;16(1):477. doi:10.1186/s13018-021-02622-z
25. Nibe Y, Matsumura T, Takahashi T, Kubo T, Matsumoto Y, Takeshita K. A comparison between the femoral neck system and other implants for elderly patients with femoral neck fracture: A preliminary report of a newly developed implant. *J Orthop Sci Off J Jpn Orthop Assoc*. 2022;27(4):876-880. doi:10.1016/j.jjos.2021.04.016
26. Rogmark C, Carlsson A, Johnell O, Sembo I. Costs of internal fixation and arthroplasty for displaced femoral neck fractures: a randomized study of 68 patients. *Acta Orthop Scand*. 2003;74(3):293-298. doi:10.1080/00016470310014210
27. Palmer SJ, Parker MJ, Hollingworth W. The cost and implications of reoperation after surgery for fracture of the hip. *J Bone Joint Surg Br*. 2000;82(6):864-866. doi:10.1302/0301-620x.82b6.9974
28. Vanderkarr M, Ruppenkamp J, Lerner J, Holy CE, Sparks C. Reoperations in intramedullary fixation of pertrochanteric hip fractures. *J Med Econ*. 2019;22(7):706-712. doi:10.1080/13696998.2019.1600526
29. Patel S, Kumar V, Baburaj V, Dhillon MS. The use of the femoral neck system (FNS) leads to better outcomes in the surgical management of femoral neck fractures in adults compared to fixation with cannulated screws: A systematic review and meta-analysis. *Eur J Orthop Surg Traumatol Orthop Traumatol*. Published online October 6, 2022. doi:10.1007/s00590-022-03407-8



Manufactured or distributed by:

Synthes USA Products, LLC
1302 Wrights Lane Est
West Chester, PA 19380

Synthes USA, LLC
1101 Synthes Avenue
Monument, CO 80132

Synthes GmbH
Luzernstrasse 21
4528 Zuchwil, Switzerland

jnmedtech.com

Please refer to the Instructions For Use for a complete list of indications, contraindications, warnings and precautions.

© DePuy Synthes 2023. All rights reserved.
248676-230519 DSUS