**STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional**  
 Knotless Tissue Control Device  
**Antibacterial**  
 Dyed and Undyed Monofilament Synthetic Absorbable Device

**DESCRIPTION**

The STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Knotless Tissue Control Device is an antibacterial monofilament, synthetic absorbable device prepared from a copolymer of glycolide and ε-caprolactone. The device contains IRGACARE®\* MP (triclosan), a broad spectrum antibacterial agent, at no more than 2360 µg/m. The device is available in a dyed and undyed version. The pigment for the dyed version is D&C Violet No. 2. Poliglecaprone 25 copolymer has been found to be nonpyrogenic and elicits only a slight tissue reaction during absorption.

The STRATAFIX™ Spiral MONOCRYL™ Plus Knotless Tissue Control Bidirectional Device Design consists of barbed suture material, armed with a surgical needle on each end. The device also contains an unbarbed center transition zone that facilitates the initiation of the device use. The STRATAFIX™ Spiral MONOCRYL™ Plus Device barbs allow for tissue approximation without the need to tie surgical knots.

While the formation of barbs in the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device reduces the tensile strength relative to non-barbed suture material of the same size, tying of knots in non-barbed suture materials also reduces their effective strength. For this reason, the strength of the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device can be compared to USP knot strength of non-barbed sutures. The actual diameter of the non-barbed section fiber is one size greater than the designated size with a maximum overage of 0.1 mm.

The USP and EU Pharmacopoeia sizes of the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device are further defined in Table 1.

**Table 1. Diameter Comparison**

USP DEVICE SIZE DESIGNATION	EU PHARMACOPOEIA DEVICE SIZE (Metric / Ph. Eur.) DESIGNATION	STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional	
		USP	Metric / Ph. Eur.
0	3.5	0	3.5
2-0	3	2-0	3
3-0	2	3-0	2
4-0	1.5	4-0	1.5

**TENSILE STRENGTH**

STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device straight tensile strength meets the knot tensile strength for a USP and EU Pharmacopoeia poliglecaprone 25 device of the equivalent size as shown in Table 2.

**Table 2. Tensile Strength Comparison**

USP DEVICE SIZE DESIGNATION	EU PHARMACOPOEIA DEVICE SIZE (Metric / Ph. Eur.) DESIGNATION	Device Minimum Knot Tensile Strength	
		USP (kgf)	Metric / Ph. Eur. (N)
0	3.5	3.90	39.0
2-0	3	2.68	26.8
3-0	2	1.77	17.5
4-0	1.5	0.95	9.5

**ACTIONS**

The results of implantation studies in animals using the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device indicate that approximately 62% of its original strength remains after 7 days implantation and approximately 27% of its original tensile strength at 14 days post implantation. All of the original tensile strength is lost by 21 days post implantation. The absolute strength remaining 14 days post implantation meets or exceeds that historically observed with plain or chromic surgical gut devices (Table 3).

**Table 3. In Vivo Tensile Strength Profile**

DAYS IMPLANTATION	APPROXIMATE % ORIGINAL STRENGTH REMAINING
7 Days	62% of USP knot tensile strength
14 Days	27% of USP knot tensile strength

Data obtained from implantation studies show that absorption of STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device is essentially complete by 91 days.

Using *in vitro* studies, STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device with IRGACARE®\* MP (triclosan) has been shown to inhibit colonization of the device by *Staphylococcus aureus*, *Staphylococcus epidermidis*, methicillin-resistant *S. aureus*, methicillin-resistant *S. epidermidis*, *Escherichia coli*, and *Klebsiella pneumoniae*, which are microorganisms known to be associated with surgical site infections. Animal studies have demonstrated that STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device inhibits bacterial colonization of the device after direct *in vivo* challenge with bacteria.

**INDICATIONS**

STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device is indicated for use in soft tissue approximation where the use of absorbable sutures is appropriate. STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device is not indicated in cardiovascular or neurological tissues, microsurgery or ophthalmic surgery.

**CONTRAINDICATIONS**

The STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device is not to be used where extended approximation of tissues under stress is required and is not to be used in conjunction with or for fixation of prosthetic devices (e.g., heart valves or synthetic grafts) that are nonabsorbable in nature. This device, being absorbable, should not be used where extended approximation of tissue under stress is required, such as in fascia. The STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device should not be used in patients with known allergic reactions to IRGACARE®\* MP (triclosan).

**WARNINGS**

Users should be familiar with surgical procedures and techniques before employing the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device for wound closure, as the risk of wound separation may vary with the site of application and the suture material used. Physicians should consider the *in vivo* performance (under ACTIONS section) when selecting a device for use in patients. The use of this device may be inappropriate in elderly, malnourished or debilitated patients, or in patients suffering from conditions which may delay wound healing.

As with any foreign body, prolonged contact of any device with salt solutions, such as those found in the urinary or biliary tracts may result in calculus formation. As an absorbable device, STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device may act transiently as a foreign body. Acceptable surgical practice should be followed for the management of contaminated or infected wounds.

As this is an absorbable device, the use of supplemental nonabsorbable sutures should be considered by the surgeon in the closure of the sites which may undergo expansion, stretching or distention, or which may require additional support.

The safety and effectiveness of the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device has not been established for use in fascial closures (including abdominal wall, thoracic and extremity fascial closures), gastrointestinal anastomoses, cardiovascular tissue, neural tissue, osseous tissue, tendinous tissue, ophthalmic surgery or for use in microsurgery.

Care should be undertaken to avoid leaving barbed suture ends adjacent to the peritoneum in extra-peritoneal tissue closure.

Small bowel obstruction (SBO): including volvulus, bowel infarction, and significant morbidity, have been reported due to barbs or barbed suture ends hooking onto adjacent small bowel and/or mesentery, such as in peritoneal closure.

Avoid contacting STRATAFIX™ Device and associated needles with other materials (e.g., surgical gauze, drapes, etc.) in the surgical field to prevent ensnaring on the barbs. If the barbs catch, carefully pull the material in the opposite direction of the needle to disengage it from the barbs.

Do not resterilize/reuse. Reuse of this device (or portions of this device) may create a risk of product degradation, which may result in device failure and/or contamination, which may lead to infection or transmission of blood-borne pathogens to patients and users. Discard open packages and unused devices. Discard packaging that is damaged or opened. Do not use after expiration date.

## PRECAUTIONS

The STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device contains bidirectionally oriented barbs to anchor tissues and does not require knots to approximate opposing edges of an incision. Tying of knots on the barbed section of the device will damage the barbs and potentially reduce the device tensile strength and barb effectiveness. For the bidirectional forces to be created and for the device to function properly, both sides of the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device must be engaged in tissue. Additionally, when completing placement, additional backstitch or passes through tissue to the end of the incision or wound are required to lock the device in place.

Devices that remain in skin longer than 7 days may cause localized irritation and should be snipped off or removed, if indicated. Subcuticular devices should be placed as deeply as possible to minimize the erythema and induration normally associated with absorption.

Care should be taken to avoid damage when handling. Avoid crushing or crimping the device material with surgical instruments, such as needle holders and forceps. When removing the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device out of the package take care not to have the barbs catch on one another. Do not attempt to remove memory in the polymer by running fingers down the suture material as this can damage the barbs.

Infections, erythema, foreign body reactions, transient inflammatory reactions and in rare instances wound separation are typical or foreseeable risks associated with any suture and hence are also potential complications associated with the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device.

Acceptable surgical practice should be followed with respect to drainage and closure of infected wounds. To avoid damaging needle points and swage areas, grasp the needle in an area one-third (1/3) to one-half (1/2) of the distance from the swaged end to the point. Reshaping needles may cause them to lose strength and be less resistant to bending and breaking. Users should exercise caution when handling surgical needles to avoid inadvertent needle sticks. Discard used needles in designated biohazard "sharps" containers.

## STORAGE

Store at or below 30°C (86°F). Do not use after expiration date.

## ADVERSE REACTIONS

Adverse effects associated with the use of synthetic absorbable devices include tissue separation, failure to provide adequate wound support in closure of the sites where expansion, stretching, or distension occurs, infection, minimal acute inflammatory tissue reaction, localized irritation when skin sutures are left in place for greater than 7 days, device extrusion and delayed absorption in tissue with poor blood supply, and transitory local irritation at the wound site, as well as allergic reaction to IRGACARE™ MP (triclosan). Broken needles may result in extended or additional surgeries or residual foreign bodies. Inadvertent needle sticks with contaminated surgical needles may result in the transmission of blood-borne pathogens.

## MRI SAFETY INFORMATION

STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device is MR Safe.

## INSTRUCTIONS FOR USE

STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device is designed to be used in continuous suture patterns. STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device is intended to be used without anchoring knots to begin or terminate the device line.

Use as required per surgical procedure.

1. To deploy the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device, the first end should be pulled through the tissue until the center transition zone has reached the tissue. Estimation of the center transition zone can be aided by taking a single bite or pass of tissue, then aligning the two needles together. Take two loose passes through the tissue with one arm of the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device.
2. Then take two loose passes with the other arm of the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device. After two passes through the tissue have been completed per each arm, each of the strands can then be grasped and the tissue can be approximated to the desired tension. Continue approximating the incision utilizing continuous stitching technique, while adjusting to the desired tension with each pass. Care should be taken to utilize the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device on the barbed segments only. Do not attempt to approximate incisions using the non-barbed segments near the needles, as the barbs of the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device are required for successful tissue approximation.
3. For a STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device with straight needles, the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device can be deployed at right angles to the incision. Starting at one edge of the incision, the straight needle is deployed at right angles to the incision along the desired tissue plane and the needle is advanced through the tissue until the desired amount of tissue has been engaged, then the needle can exit at the desired location remote from the incision edge. The same procedure is repeated with the second arm of the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device.

4. After successful deployment of both arms of the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device and needle exits, the center transition zone of the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device should be roughly in the center of the incision. The ends of the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device where they exit remote to the incision can now be grasped and the tissue can be pushed together to facilitate the soft tissue approximation of the desired tension.
5. To complete a subcuticular closure, take at least one pass in the reverse direction. Then pass the needle across the incision and take a split thickness bite perpendicular to the incision and exit the skin. To complete the closure in all other tissue layers, take at least two passes in the reverse direction.
6. Then gently pull on the free ends of the device and cut flush with the surface of the tissue.

To remove the entire STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device: Cut the STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device at the midline or center transition zone between opposing barbed segments and then pull laterally to remove the device segments. Multiple cuts may be necessary to remove longer segments from the tissue.


## STERILITY

STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device is sterilized by ethylene oxide gas. Do not resterilize. Do not use if package is opened or damaged. Discard opened, unused sutures. Do not use after expiration date.

## HOW SUPPLIED

STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device is available sterile, dyed and undyed in various lengths and with a USP diameter size 0 through 4-0 (metric sizes 3.5 through 1.5). Each STRATAFIX™ Spiral MONOCRYL™ Plus Bidirectional Device is double armed with various needle sizes.

## SYMBOLS USED ON LABELING

	Caution
	Do not re-use
	Do not resterilize
	Do not use if package is damaged
	Caution: Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
	Manufacturer
	Sterilized using ethylene oxide
	Batch code
	Use-by date
	Upper limit of temperature
	Packaging unit
	Catalogue number

\*Registered Trademark of BASF Group.

**ETHICON**  
PART OF THE  FAMILY OF COMPANIES

  
Ethicon, LLC  
475 C Street  
Los Frailes Industrial Park, Suite 401  
Guaynabo, Puerto Rico 00969  
USA  
1-877-ETHICON  
+1-513-337-6928

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© Ethicon, Inc. 2021