

DEVICE SPECS

Ethicon Barbed Sutures



STRATAFIX™ Symmetric PDS™ Plus
Knotless Tissue Control Device



STRATAFIX™ Spiral PDS™ Plus
Knotless Tissue Control Device



STRATAFIX™ Spiral MONOCRYL™ Plus
Knotless Tissue Control Device

Codes	See catalog	See catalog	See catalog
Color	Violet	Violet	Violet or Undyed
Antibacterial Coating	Yes	Yes	Yes
Material/Composition	Polydioxanone, Triclosan	Polydioxanone, Triclosan	Poliglecaprone 25, Triclosan
Suture Structure	Monofilament, Barbed	Monofilament, Barbed	Monofilament, Barbed
Absorbable	Yes	Yes	Yes
Absorption Rate	210 days	210 days	91-119 days
Directionality	Unidirectional	Unidirectional	Bidirectional or Unidirectional
Available Needle Options	Taperpoint, Tapercut Surgical Needle, Precision Point Reverse Cutting, Reverse Cutting, Blunt Point Ethiguard	Taperpoint, Tapercut Surgical Needle, Precision Point Reverse Cutting, Reverse Cutting, Blunt Point Ethiguard	Taperpoint, Reverse Cutting, Precision Point Reverse Cutting
Special Options	MultiPass™	MultiPass™	MultiPass™
Pledget Options	No	No	No
Strands per Pack	1	1	1
Quantity per Box	12	12	12