

Let's put it
all together.
With ease
and efficiency.

Their lives. Your skills. Our products.

DESIGNED
FOR OR
EFFICIENCY,
AND EASE
OF USE

All screws can
be used with all
plates within each
Matrix system.




A wide selection
of implants to
address each
individual
surgical needs.



All implants also in
sterile version.



Self-retaining
screws and
blades are
optimized for
ease of use and
handling.



Plates, screws
and instruments
are color-coded to
aid identification
and improve OR
communication.



One screwdriver blade for all screws within each Matrix system for reduced instrument change time.

INNOVATIVE SOLUTIONS FOCUSED ON PATIENT OUTCOMES

Configurable compact modules consolidate and simplify implants and instrument selection.*



A range of preformed implants to support anatomical restoration for complex indications.



Nano case provides configurable, small and cost effective set solution.



Maxillofacial Matrix system integrates fully with TRUMATCH® CMF Personalized Solutions.



REDUCED INVENTORY TO RATIONALISE RESOURCE ALLOCATION


MatrixWAVE MMF system combines strength and simplicity to achieve maxillomandibular fixation.



*Extensive configuration sets are also available.



MatrixMANDIBLE™
Plating System



MatrixMIDFACE™
Plating System



MatrixORTHOGNATHIC™
Plating System



MatrixWAVE™ MMF

MAXILLOFACIAL Matrix PLATFORM:

Successful maxillofacial surgery means restoring physical function, and rebuilding the geography and aesthetic beauty of the human face. Whether correcting a structural defect, rebuilding facial features following extensive resection or repairing massive facial trauma, skilled surgeons working with the right tools, plates and screws, achieve results that patients deserve.

The Maxillofacial Matrix platform has been developed under medical guidance of surgeons to deliver on the highest clinical and hospital needs and standards for life-changing outcomes. As a comprehensive plates and screws system for all major maxillofacial procedures, the Matrix portfolio is approved by the AO Foundation, and has become one of the most widely used maxillofacial surgery systems worldwide.

Johnson & Johnson Medical Pty Ltd t/a DePuy Synthes
1-5 Khartoum Road, North Ryde, NSW 2113, Australia

Johnson & Johnson (NZ) Ltd
570 Mt Wellington Highway, Mt Wellington,
Auckland 1060, New Zealand

106610-190130 June 2019



Synthes GmbH
Einmattstrasse
4436 Oberdorf
Switzerland

Tel: +41 61 965 61 11

Fax: +41 61 965 66 00

Not all products are currently available in all markets. This publication is not intended for distribution in the USA.



0123