

## DEVICE SPECS

# ECHELON ENDOPATH™ Staple Line Reinforcement



**ECHELON ENDOPATH™ Staple Line Reinforcement**

<b>Codes</b>	ECH60R
<b>Compatible with</b>	GST60W, GST60B, GST60D, GST60G, GST60T
<b>Description</b>	Staple line buttress
<b>Size</b>	60mm
<b>Attachment Method</b>	Attach with temporary adhesive
<b>Material Composition</b>	Polyglactin 910 (same material as VICRYL® suture) and 2x polydioxanone (PDO) film
<b>Material Architecture</b>	Laminated knit: designed at an optimized angle to increase strength
<b>Total Thickness (2 layers)</b>	0.3mm
<b>Absorption Time</b>	17 weeks
<b>Surface Area of Implanted Buttress Material</b>	270cm <sup>2</sup>
<b>Qty./Box</b>	6