

DEVICE SPECS

# Ethicon Laparoscopic Stapler Portfolio



**ECHELON™ 3000 Stapling System**



**ECHELON™+ Stapling System**



**ECHELON FLEX™ ENDOPATH™ Staplers**

	<b>ECHELON™ 3000 Stapling System</b>	<b>ECHELON™+ Stapling System</b>	<b>ECHELON FLEX™ ENDOPATH™ Staplers</b>
<b>Codes</b>	ECH45C, ECH45S, ECH45L, ECH60C, ECH60S, ECH60L	PCEE45A, PSEE45A, PLEE45A, PCEE60A, PSEE60A, PLEE60A	SC45A, EC45A, EC45AL, SC60A, EC60A, LONG60A
<b>Associated Reloads</b>	GST45W, GST60W, GST45B, GST60B, GST45D, GST60D, GST45G, GST60G, GST45T, GST60T	GST45W, GST60W, GST45B, GST60B, GST45D, GST60D, GST45G, GST60G, GST45T, GST60T	GST45W, GST60W, GST45B, GST60B, GST45D, GST60D, GST45G, GST60G, GST45T, GST60T
<b>Description</b>	45mm and 60mm laparoscopic staplers	45mm and 60mm laparoscopic staplers	45mm and 60mm laparoscopic staplers
<b>Firing</b>	Powered	Powered	Manual
<b>Articulation</b>	Powered	Manual	Manual
<b>Shaft Length</b>	Compact 28cm, Standard 34cm, Long 44cm	Compact 28cm, Standard 34cm, Long 44cm	Compact 28cm, Standard 34cm, Long 44cm
<b>Degrees of Articulation</b>	>110°	88°	90°
<b>Staple Rows</b>	6	6	6
<b>Max Firings</b>	12	12	12
<b>Staple Form</b>	B-Form	B-Form	B-Form
<b>Qty./Box</b>	3	3	3