

Stapling

Laparoscopic Staplers



ETHICON™ 4000 Stapler System

ECHELON™ 3000 Stapling System

ECHELON™+ Stapling System

ECHELON FLEX™ ENDOPATH™ Staplers

	ETHICON™ 4000 Stapler System	ECHELON™ 3000 Stapling System	ECHELON™+ Stapling System	ECHELON FLEX™ ENDOPATH™ Staplers
Codes	EC3D60C, EC3D60S, EC3D60L	ECH45C, ECH45S, ECH45L, ECH60C, ECH60S, ECH60L	PCEE45A, PSEE45A, PLEE45A, PCEE60A, PSEE60A, PLEE60A	SC45A, EC45A, EC45AL, SC60A, EC60A, LONG60A
Associated Reloads	ER60W, ER60B, ER60G, ER60T	GST45W, GST60W, GST45B, GST60B, GST45D, GST60D, GST45G, GST60G, GST45T, GST60T	GST45W, GST60W, GST45B, GST60B, GST45D, GST60D, GST45G, GST60G, GST45T, GST60T	GST45W, GST60W, GST45B, GST60B, GST45D, GST60D, GST45G, GST60G, GST45T, GST60T, ECR45M, ECR60M
Description	60mm laparoscopic staplers	45mm and 60mm laparoscopic staplers	45mm and 60mm laparoscopic staplers	45mm and 60mm laparoscopic staplers
Firing	Powered	Powered	Powered	Manual
Articulation	Powered	Powered	Manual	Manual
Shaft Length	Compact 28cm, Standard 34cm, Long 44cm	Compact 28cm, Standard 34cm, Long 44cm	Compact 28cm, Standard 34cm, Long 44cm	Compact 28cm, Standard 34cm, Long 44cm
Degrees of Articulation	>110°	>110°	88°	--
Staple Rows	6	6	6	6
Max Firings	12	12	12	12
Staple Form	3D and 2D B-Form	B-Form	B-Form	B-Form
Qty./Box	3	3	3	3