

REIMAGINE NAVIGATION




TruDi[®] Curette Navigation Instrument

Introducing TruDi[®] Curette, Accuracy for Every Dissection

- Distal tip sensors enable plug-and-play solution
- Customizable to help reach targeted anatomy
- Gripping features assist with tactile feel during dissection
- Strong enough to meet your clinical needs

120mm



65°

90°


Acclarent[®]
PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

TruDi® Curette

Navigation Instrument

Precise Dissection Around Critical Anatomy

- Consistent accuracy for ENT and ENT Skull Base procedures
- Sharp enough to accurately dissect bone, tissue, and cartilage
- Low profile design improves maneuverability and dissection in tight anatomy



TruDi® Curette enables the advanced functionality of TruDi® software features like **TruView**, **Planning Points**, **Beacon**, and **FAM**

Catalog code	Product name	Size	Packaging unit
TDC0005	TruDi® Curette - 0 Degree Tip Angle	5"	5PK

Important Safety Information: Acclarent Navigation Technology is intended for use by or under the direction of a physician. Devices guided by the Acclarent Navigation Technology have described risks including the potential for unintended injury to the orbits or CNS as a consequence of inaccurate navigation. Prior to use, it is important to read the Instructions for Use and to understand the registration process and the contraindications, warnings and precautions described for the system.

Caution: Federal (US) law restricts the sale, distribution or use of these devices to, by or on the order of a physician. Third party trademarks used herein are trademarks of their respective owners. This content is intended for audience within the US only.

Acclarent, Inc.
Irvine, CA 92618

© Acclarent, Inc. 2020. All Rights Reserved. 151257-200827

Acclarent®
PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES