

# DERMABOND ADVANCED®

## Topical Skin Adhesive

### A protective barrier that adds strength and inhibits bacteria

- When used in addition to sutures, DERMABOND ADVANCED Adhesive was shown ex vivo to add 75% more strength to the wound closure than sutures alone<sup>1</sup>
- DERMABOND ADVANCED Adhesive demonstrated in vitro inhibition of gram-positive bacteria (MRSA and MRSE) and gram-negative bacteria (E Coli)<sup>2\*</sup>
- Provides a microbial barrier with 99% protection in vitro for at least 72 hours against organisms commonly responsible for SSIs<sup>2\*†</sup>

### A patented, thicker formula with a unique blend of plasticizers and an initiator for precise control<sup>§</sup>

- Plasticizers increase flexibility and 3-dimensional strength for secure wound closure and patient comfort<sup>3,4</sup>
- Initiator enables consistent setting time across various conditions and skin types<sup>5,6</sup>



Only  
**1**  
Layer Required<sup>7</sup>

**2X**  
Faster Setting  
Time (95 seconds)<sup>6†</sup>

**4X**  
More Viscous<sup>8†</sup>

**8**  
US Patients<sup>9</sup>

\*Clinical significance is unknown.

†As compared to High Viscosity DERMABOND® Topical Skin Adhesive.

‡*Staphylococcus epidermidis*, *Staphylococcus aureus*,  
*Escherichia coli*, *Enterococcus faecium* and *Pseudomonas aeruginosa*.

§Compared to DERMABOND® Mini Topical Skin Adhesive

**ETHICON**  
PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

For more information, contact your  
Ethicon representative or call 1-800-255-2500.

**For complete indications, contraindications, precautions, warnings,  
and adverse reactions, please reference full package insert.**

# DERMABOND ADVANCED Adhesive

## Ordering Information

### An applicator ergonomically designed for improved handling and easy coverage

- Angled tip for both fine and wide line application - up to 15cm<sup>10</sup>
- Innovative applicator designed to easily activate and avoid clogging
- Innovative design provides greater control and ability to see product level



Product Code	Description	Size	Packaging
DNX12	DERMABOND ADVANCED Adhesive <ul style="list-style-type: none"> <li>• Unique formulation with an easy-to-use applicator for precise control for wide and fine line applications</li> </ul>	0.7 mL	12 Applicators per box
DNX6	DERMABOND ADVANCED Adhesive <ul style="list-style-type: none"> <li>• Unique formulation with an easy-to-use applicator for precise control for wide and fine line applications</li> </ul>	0.7 mL	6 Applicators per box

For more information, contact your representative or call 1-800-255-2500.

**For complete indications, contraindications, precautions, warnings, and adverse reactions, please reference full package insert.**

**References:** 1. Protocol to determine the strength benefits when adding DERMABOND ADVANCED™ Topical Skin Adhesive to a wound or incision using suture as the primary closure, July 2011. 2. Bhende S, Rothenburger S, Spangler DJ, Dito M. In vitro assessment of microbial barrier properties of DERMABOND® Topical Skin Adhesive. Surg Infect (Larchmt). 2002;3(3):251-257. 3. Material Specification 53MS-0000681. 4. Singer A et al. Evaluation and management of traumatic lacerations. New England Journal of Medicine. 1997;337:1142-1148. 5. ADAPTIV Part PM62028, MWI-0000480. 6. ADAPTIV Report SUBOTH0004215. 7. DERMABOND ADVANCED Adhesive Instructions for Use. 8. FG-0000104 of 49cP, FG-0000101 of 214cP. 9. LAB0012187(v3). 10. 09PDO44 Final Protocol Report: Design Verification-Application Length and Width Evaluation of DERMABOND(R) NX Topical Skin Adhesive, 2010.

©2021 Ethicon, Inc. All rights reserved. 028909-210902

**Please add DERMABOND ADVANCED® Topical Skin Adhesive to my preference card:**

Doctor Name: \_\_\_\_\_

Signature: \_\_\_\_\_

Procedures: \_\_\_\_\_

All Procedures