



ANSPACH® Power Tools

G1 Dissection Tools

Reference Chart




Minimal Access Attachments

Combine the prefix and suffix from the chart below to determine the dissection tool Catalog Number.
Example: MA10-1SB

Catalog Number Suffix

Catalog Number Prefix



ATTACHMENT TYPE

FLUTED BALL

Catalog Number Suffix	Description	MA10-	MA15-	MA19-
1SB	1 mm Fluted Ball	●		
15SB	1.5 mm Fluted Ball	●		
2SB	2 mm Fluted Ball	●	●	●
25SB	2.5 mm Fluted Ball	●		
3SB	3 mm Fluted Ball	●	●	●

DIAMOND BALL

Catalog Number Suffix	Description	MA10-	MA15-	MA19-
05SD	0.5 mm Diamond Ball	●		
08SD	0.8 mm Diamond Ball	●		
15SD	1.5 mm Diamond Ball		●	●
2SD	2 mm Diamond Ball	●	●	●
3SD	3 mm Diamond Ball	●	●	●
4D	4 mm Diamond Ball	●	●	●
4DC	4 mm Coarse Diamond Ball		●	

MATCHSTICK

Catalog Number Suffix	Description	MA10-	MA15-	MA19-
8N-18	1.8 mm Fluted Matchstick	●	●	●
8N	2.3 mm Fluted Matchstick	●	●	●
8NS	3 mm Fluted Matchstick	●	●	●
8NSDC	3 mm Coarse Diamond Matchstick	●	●	●

OCM-G1 OTOLOGIC Curved Micro Attachment, G1

DIAMOND BALL CURVED BURRS

Burr Diameter	5.6 cm Sleeve Length	7.3 cm Sleeve Length
0.5 mm	OCM5-05SD	OCM7-05SD
0.6 mm	OCM5-06SD	OCM7-06SD
0.7 mm	OCM5-07SD	OCM7-07SD
0.8 mm	OCM5-08SD	OCM7-08SD
1.0 mm	OCM5-1SD	OCM7-1SD
1.5 mm	OCM5-15SD	OCM7-15SD
2.0 mm	OCM5-2SD	OCM7-2SD

FLUTED BALL CURVED BURRS

Burr Diameter	5.6 cm Sleeve Length	7.3 cm Sleeve Length
1.0 mm	OCM5-1SB	OCM7-1SB
1.5 mm	OCM5-15SB	OCM7-15SB
2.0 mm	OCM5-2SB	OCM7-2SB

Please also refer to the package insert(s) or other labeling associated with the devices identified in this brochure.

CAUTION: Federal Law restricts these devices to sale by or on the order of a physician.

Some devices listed in this brochure may not have been licensed in accordance with Canadian law and may not be for sale in Canada. Please contact your sales consultant for items approved for sale in Canada.

Not all products may currently be available in all markets.



Manufacturer
The Anspach Effort, Inc.
4500 Riverside Drive
Palm Beach Gardens, FL 33410
USA
To Order (USA): (800) 327-6887
To Order (Canada): (855) 946-8999

EC REP Synthes GmbH

Eimattstrasse 3
4436 Oberdorf
Switzerland
Tel: +41 61 965 61 11
Fax: +41 61 965 66 00

www.jnjmedicaldevices.com