



EG1[®] High Speed Electric System

Precision and Power



ANSPACH® EG1 High Speed System

Precise

- Up to 80,000 rpm
- High cutting precision
- New dissection tools and coupling mechanism
- Minimal start-up kick

Easy to Use


- Place and lock attachment
- Push to lock dissection tool
- User friendly, stackable console
- Compatible with automated washing

Powerful

- 30% more power than the ANSPACH® XMAX® and EMAX® 2 Plus Systems
- Higher cutting speed for comfort and precision
- Air cooled motor for demanding applications



EG1® High Speed Electric System is intended for cutting and shaping bone including the spine and cranium.



Versatile System

- Neuro applications
- Spine applications
- Otology/Neurotology applications



Dedicated Console

- Stackable
- Optional irrigation



Dedicated Cleaning Basket

- Automated washing and sterilization
- Simple to load and unload

Ordering Information

ANSPACH® Power Tools Electric System

| | |
|--------------|---|
| EG1A | G1 Electric System Handpiece |
| SC3001 | G1 System Console without Irrigation |
| SC3002 | G1 System Console with Irrigation |
| E-FP | Electric Systems Foot Control |
| E-FP-DIR | Electric Systems Foot Control with Direction |
| E-FP-DIR/IRR | Electric Systems Foot Control with Direction and Irrigation |

ANSPACH Power Tools Electric System Accessories

| | |
|-------------|--|
| G1-BSKT-BTM | G1 Cleaning and Sterilization Basket, bottom |
| G1-BSKT-INS | G1 Cleaning and Sterilization Basket, insert |
| G1-BSKT-LID | G1 Cleaning and Sterilization Basket, lid |
| G1-BSKT-SLV | G1 Cleaning and Sterilization Basket, sleeve |

ANSPACH Power Tools Irrigation Accessories

| | |
|----------------|---|
| 05.001.178.01S | Sterile Tubing for Irrigation System |
| HOSECLIPS-G1 | Irrigation Tube Hoseclips for EG1A, Quantity 6 |
| IRR-CLIP-10 | Irrigation Clip for SHORT-G1 and QD8-G1 |
| IRR-CLIP-30 | Irrigation Clip for SHORT-HD-G1 |
| IRR-CLIP-40 | Irrigation Clip for MEDIUM-HD-G1 and LONG-HD-G1 |
| IRR-CLIP-50 | Irrigation Clip for MEDIUM-G1, LONG-G1, MIA16-G1, and QD11-G1 |

ANSPACH Power Tools Attachments

| | |
|--------------|---|
| SHORT-G1 | 5 cm Short Attachment |
| MEDIUM-G1 | 8 cm Medium Attachment |
| LONG-G1 | 11 cm Long Attachment |
| MIA16-G1 | 16 cm Minimally Invasive Attachment |
| SHORT-HD-G1 | 5 cm Heavy Duty Short Attachment |
| MEDIUM-HD-G1 | 8 cm Heavy Duty Medium Attachment |
| LONG-HD-G1 | 12.4 cm Heavy Duty Long Attachment |
| CRANI-A-G1 | 6.5 cm Adult Craniotome |
| CRANI-P-G1 | 6.5 cm Pediatric Craniotome |
| CRANI-L-G1 | 7.5 cm Large Craniotome |
| QD8-G1 | 8 cm Quick Disconnect Angle Attachment |
| QD8-S-G1 | 7.5 cm Quick Disconnect Angle Attachment |
| QD11-G1 | 11 cm Quick Disconnect Angle Attachment |
| QD11-S-G1 | 10.5 cm Quick Disconnect Angle Attachment |
| QD14-G1 | 14 cm Quick Disconnect Angle Attachment |
| QD14-S-G1 | 13.5 cm Quick Disconnect angle Attachment |
| OCM-G1 | Otologic Curved Micro (OCM) Attachment |
| CSR60-G1 | Perforator Driver with Hudson End |

ANSPACH Power Tools Minimal Access Attachment

| | |
|-----------|------------------|
| MA-D20-G1 | 20° Angle Driver |
|-----------|------------------|

ANSPACH Power Tools Bearing Sleeves

| | |
|---------|--------------------------------------|
| MA-15S | Bearing Sleeve, 15 cm Straight |
| MA-15C | Bearing Sleeve, 15 cm Curved |
| MA-15ST | Bearing Sleeve, 15 cm Straight Taper |
| MA-19ST | Bearing Sleeve, 19 cm Straight Taper |
| MA-10S | Bearing Sleeve, 10 cm Straight |
| MA-10C | Bearing Sleeve, 10 cm Curved |

Limited Warranty and Disclaimer: *DePuy Synthes* products are sold with a limited warranty to the original purchaser against defects in workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.

Please see the package insert(s), surgical technique guide(s), or other product labeling associated with the devices contained in this brochure for information on indications, warnings, precautions, and instructions for use.

CAUTION: Federal Law restricts these devices to sale by or on the order of a physician.

Some devices listed in this brochure may not have been licensed in accordance with Canadian law and may not be for sale in Canada. Please contact your sales consultant for items approved for sale in Canada.

Not all products are currently available in all markets.



PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES



Manufacturer

The Anspach Effort, Inc.

4500 Riverside Drive
Palm Beach Gardens, FL 33410
USA

To Order (USA): (800) 327-6887

To Order (Canada): (855) 946-8999

www.depuyssynthes.com

EC REP Synthes GmbH

Eimattstrasse 3
4436 Oberdorf
Switzerland

Tel. +41 61 965 61 11

Fax +41 61 965 66 00