


Sterile Tube Ordering Information

INTRODUCING STERILE TUBE
FROM DEPUY SYNTHES

DESIGNED TO SAVE
YOU SPACE* AND TIME



Implants Overview for Trauma



2.0 mm Screws

Cortex Screw Stardrive® Ø 2.0 mm	page 3
Locking Screw Stardrive® Ø 2.0 mm	page 3
VA Locking Screw Stardrive® Ø 2.0 mm	page 4



2.4 mm Screws

Cortex Screw Stardrive® Ø 2.4 mm	pages 4, 5
Locking Screw Stardrive® Ø 2.4 mm	page 5
VA Locking Screw Stardrive® Ø 2.4 mm	page 6



2.7 mm Screws

Cortex Screw Stardrive® Ø 2.7 mm	page 7
Low Profile Metaphyseal Screw Stardrive® Ø 2.7 mm	page 7
Locking Screw Stardrive® Ø 2.7 mm	page 8
VA Locking Screw Stardrive® Ø 2.7 mm	page 8



3.5 mm Screws

Cortex Screw Stardrive® Ø 3.5 mm	page 9
Locking Screw Stardrive® Ø 3.5 mm	page 9
Locking Screw Stardrive® perforated Ø 3.5 mm	page 10
VA Locking Screw Stardrive® Ø 3.5 mm	page 10



4.0 mm Screws

Cancellous Bone Screw Ø 4.0 mm	page 11
--------------------------------------	---------

Note: This Usage Guide does not include information necessary for selection and use of the implants. Please refer to the corresponding Surgical Technique for all necessary implant specific information.

Ordering Information

Cortex Screw Stardrive® Ø 2.0 mm, self-tapping, sterile tube

Stainless Steel	Titanium (TiCP)	Length (mm)
201.356TS	401.356TS	6
201.357TS	401.357TS	7
201.358TS	401.358TS	8
201.359TS	401.359TS	9
201.360TS	401.360TS	10
201.361TS	401.361TS	11
201.362TS	401.362TS	12
201.363TS	401.363TS	13
201.364TS	401.364TS	14
201.365TS	401.365TS	16
201.366TS	401.366TS	18
201.368TS	401.368TS	20
201.370TS	401.370TS	22
201.372TS	401.372TS	24
201.374TS	401.374TS	26
201.375TS	401.375TS	28
201.376TS	401.376TS	30
201.377TS	401.377TS	30
201.378TS	401.378TS	32
201.379TS	401.379TS	34
201.380TS	401.380TS	36
201.381TS	401.381TS	38
201.382TS	401.382TS	40

Locking Screw Stardrive® Ø 2.0 mm, self-tapping, sterile tube

Stainless Steel	Titanium (TiCP)	Length (mm)
201.876TS	401.876TS	6
201.877TS	401.877TS	7
201.878TS	401.878TS	8
201.879TS	401.879TS	9
201.880TS	401.880TS	10
201.881TS	401.881TS	11
201.882TS	401.882TS	12
201.883TS	401.883TS	13
201.884TS	401.884TS	14
201.886TS	401.886TS	16
201.888TS	401.888TS	18
201.890TS	401.890TS	20
201.892TS	401.892TS	22
201.894TS	401.894TS	24
201.896TS	401.896TS	26
201.898TS	401.898TS	28
201.900TS	401.900TS	30

VA Locking Screw Stardrive®
Ø 2.0 mm, self-tapping, sterile tube

Stainless Steel	Titanium (TiCP)	Length (mm)
02.130.306TS	04.130.306TS	6
02.130.307TS	04.130.307TS	7
02.130.308TS	04.130.308TS	8
02.130.309TS	04.130.309TS	9
02.130.310TS	04.130.310TS	10
02.130.311TS	04.130.311TS	11
02.130.312TS	04.130.312TS	12
02.130.313TS	04.130.313TS	13
02.130.314TS	04.130.314TS	14
02.130.315TS	04.130.315TS	15
02.130.316TS	04.130.316TS	16
02.130.318TS	04.130.318TS	18
02.130.320TS	04.130.320TS	20
02.130.322TS	04.130.322TS	22
02.130.324TS	04.130.324TS	24
02.130.326TS	04.130.326TS	26
02.130.328TS	04.130.328TS	28
02.130.330TS	04.130.330TS	30

Cortex Screw Stardrive®
Ø 2.4 mm, self-tapping, sterile tube

Stainless Steel	Titanium (TiCP)	Length (mm)
201.756TS	401.756TS	6
201.757TS	401.757TS	7
201.758TS	401.758TS	8
201.759TS	401.759TS	9
201.760TS	401.760TS	10
201.761TS	401.761TS	11
201.762TS	401.762TS	12
201.763TS	401.763TS	13
201.764TS	401.764TS	14
201.766TS	401.766TS	16
201.768TS	401.768TS	18
201.770TS	401.770TS	20
201.772TS	401.772TS	22
201.774TS	401.774TS	24
201.776TS	401.776TS	26
201.778TS	401.778TS	28
201.780TS	401.780TS	30
201.782TS	401.782TS	32
201.784TS	401.784TS	34
201.786TS	401.786TS	36
201.788TS	401.788TS	38
201.790TS	401.790TS	40

Cortex Screw Stardrive®
Ø 2.4 mm, self-tapping, sterile tube

Stainless Steel	Titanium (TiCP)	Length (mm)
02.210.942TS	04.210.942TS	42
02.210.944TS	04.210.944TS	44
02.210.946TS	04.210.946TS	46
02.210.948TS	04.210.948TS	48
02.210.950TS	04.210.950TS	50
02.210.952TS	04.210.952TS	52
02.210.954TS	04.210.954TS	54
02.210.956TS	04.210.956TS	56

Locking Screw Stardrive®
Ø 2.4 mm, self-tapping, sterile tube

Stainless Steel	Titanium (TiCP)	Length (mm)
212.806TS	412.806TS	6
212.807TS	412.807TS	7
212.808TS	412.808TS	8
212.809TS	412.809TS	9
212.810TS	412.810TS	10
212.811TS	412.811TS	11
212.812TS	412.812TS	12
212.813TS	412.813TS	13
212.814TS	412.814TS	14
212.816TS	412.816TS	16
212.818TS	412.818TS	18
212.820TS	412.820TS	20
212.822TS	412.822TS	22
212.824TS	412.824TS	24
212.826TS	412.826TS	26
212.828TS	412.828TS	28
212.830TS	412.830TS	30

VA Locking Screw Stardrive®
Ø 2.4 mm, self-tapping, sterile tube

Stainless Steel	Titanium (TiCP)	Length (mm)
02.210.108TS	04.210.108TS	8
02.210.110TS	04.210.110TS	10
02.210.112TS	04.210.112TS	12
02.210.114TS	04.210.114TS	14
02.210.116TS	04.210.116TS	16
02.210.118TS	04.210.118TS	18
02.210.120TS	04.210.120TS	20
02.210.122TS	04.210.122TS	22
02.210.124TS	04.210.124TS	24
02.210.126TS	04.210.126TS	26
02.210.128TS	04.210.128TS	28
02.210.130TS	04.210.130TS	30
02.210.106TS	04.210.106TS	6
02.210.107TS	04.210.107TS	7
02.210.109TS	04.210.109TS	9
02.210.111TS	04.210.111TS	11
02.210.113TS	04.210.113TS	13

Stainless Steel	Titanium (TiCP)	Length (mm)
02.210.132TS	04.210.132TS	32
02.210.134TS	04.210.134TS	34
02.210.136TS	04.210.136TS	36
02.210.138TS	04.210.138TS	38
02.210.140TS	04.210.140TS	40
02.210.142TS	04.210.142TS	42
02.210.144TS	04.210.144TS	44
02.210.146TS	04.210.146TS	46
02.210.148TS	04.210.148TS	48
02.210.150TS	04.210.150TS	50
02.210.152TS	04.210.152TS	52
02.210.154TS	04.210.154TS	54
02.210.156TS	04.210.156TS	56

Cortex Screw Stardrive® Ø 2.7 mm, self-tapping, sterile tube

Stainless Steel	Titanium (TiCP)	Length (mm)
202.866TS	402.866TS	6
202.868TS	402.868TS	8
202.870TS	402.870TS	10
202.872TS	402.872TS	12
202.874TS	402.874TS	14
202.876TS	402.876TS	16
202.878TS	402.878TS	18
202.880TS	402.880TS	20
202.882TS	402.882TS	22
202.884TS	402.884TS	24
202.886TS	402.886TS	26
202.888TS	402.888TS	28
202.890TS	402.890TS	30
202.892TS	402.892TS	32
202.894TS	402.894TS	34
202.896TS	402.896TS	36
202.898TS	402.898TS	38
202.900TS	402.900TS	40
202.962TS	402.962TS	42
202.963TS	402.963TS	44
202.964TS	402.964TS	45
202.965TS	402.965TS	46
202.966TS	402.966TS	48
202.967TS	402.967TS	50
202.968TS	402.968TS	55

Low Profile Metaphyseal Screw Stardrive® Ø 2.7 mm, self-tapping, sterile tube

Stainless Steel	Titanium (TiCP)	Length (mm)
02.118.510TS	04.118.510TS	10
02.118.512TS	04.118.512TS	12
02.118.514TS	04.118.514TS	14
02.118.516TS	04.118.516TS	16
02.118.518TS	04.118.518TS	18
02.118.520TS	04.118.520TS	20
02.118.522TS	04.118.522TS	22
02.118.524TS	04.118.524TS	24
02.118.526TS	04.118.526TS	26
02.118.528TS	04.118.528TS	28
02.118.530TS	04.118.530TS	30
02.118.532TS	04.118.532TS	32
02.118.534TS	04.118.534TS	34
02.118.536TS	04.118.536TS	36
02.118.538TS	04.118.538TS	38
02.118.540TS	04.118.540TS	40
02.118.542TS	04.118.542TS	42
02.118.544TS	04.118.544TS	44
02.118.546TS	04.118.546TS	46
02.118.548TS	04.118.548TS	48
02.118.550TS	04.118.550TS	50
02.118.552TS	04.118.552TS	52
02.118.554TS	04.118.554TS	54
02.118.556TS	04.118.556TS	56

**Locking Screw Stardrive® Ø 2.7 mm
(head LCP 2.4), self-tapping,
sterile tube**

Stainless Steel	Titanium (TiCP)	Length (mm)
202.206TS	402.206TS	6
202.207TS	402.207TS	7
202.208TS	402.208TS	8
202.209TS	402.209TS	9
202.210TS	402.210TS	10
202.211TS	402.211TS	11
202.212TS	402.212TS	12
202.213TS	402.213TS	13
202.214TS	402.214TS	14
202.216TS	402.216TS	16
202.218TS	402.218TS	18
202.220TS	402.220TS	20
202.222TS	402.222TS	22
202.224TS	402.224TS	24
202.226TS	402.226TS	26
202.228TS	402.228TS	28
202.230TS	402.230TS	30
202.232TS	402.232TS	32
202.234TS	402.234TS	34
202.236TS	402.236TS	36
202.238TS	402.238TS	38
202.240TS	402.240TS	40
202.242TS	402.242TS	42
202.244TS	402.244TS	44
202.246TS	402.246TS	46
202.248TS	402.248TS	48
202.250TS	402.250TS	50
202.255TS	402.255TS	55

**VA Locking Screw Stardrive® Ø 2.7 mm
(head 2.4), self-tapping, sterile tube**

Stainless Steel	Titanium (TiCP)	Length (mm)
02.211.010TS	04.211.010TS	10
02.211.012TS	04.211.012TS	12
02.211.014TS	04.211.014TS	14
02.211.016TS	04.211.016TS	16
02.211.018TS	04.211.018TS	18
02.211.020TS	04.211.020TS	20
02.211.022TS	04.211.022TS	22
02.211.024TS	04.211.024TS	24
02.211.026TS	04.211.026TS	26
02.211.028TS	04.211.028TS	28
02.211.030TS	04.211.030TS	30
02.211.032TS	04.211.032TS	32
02.211.034TS	04.211.034TS	34
02.211.036TS	04.211.036TS	36
02.211.038TS	04.211.038TS	38
02.211.040TS	04.211.040TS	40
02.211.042TS	04.211.042TS	42
02.211.044TS	04.211.044TS	44
02.211.046TS	04.211.046TS	46
02.211.048TS	04.211.048TS	48
02.211.050TS	04.211.050TS	50
02.211.052TS	04.211.052TS	52
02.211.054TS	04.211.054TS	54
02.211.056TS	04.211.056TS	56

Cortex Screw Stardrive® Ø 3.5 mm, self-tapping, sterile tube

Stainless Steel	Titanium (TiCP)	Length (mm)
02.200.010TS	04.200.010TS	10
02.200.012TS	04.200.012TS	12
02.200.014TS	04.200.014TS	14
02.200.016TS	04.200.016TS	16
02.200.018TS	04.200.018TS	18
02.200.020TS	04.200.020TS	20
02.200.022TS	04.200.022TS	22
02.200.024TS	04.200.024TS	24
02.200.026TS	04.200.026TS	26
02.200.028TS	04.200.028TS	28
02.200.030TS	04.200.030TS	30
02.200.032TS	04.200.032TS	32
02.200.034TS	04.200.034TS	34
02.200.036TS	04.200.036TS	36
02.200.038TS	04.200.038TS	38
02.200.040TS	04.200.040TS	40
02.200.042TS	04.200.042TS	42
02.200.044TS	04.200.044TS	44
02.200.045TS	04.200.045TS	45
02.200.046TS	04.200.046TS	46
02.200.048TS	04.200.048TS	48
02.200.050TS	04.200.050TS	50
02.200.055TS	04.200.055TS	55

Locking Screw Stardrive® Ø 3.5 mm, self-tapping, sterile tube

Stainless Steel	Titanium (TiCP)	Length (mm)
212.101TS	412.101TS	10
212.102TS	412.102TS	12
212.103TS	412.103TS	14
212.104TS	412.104TS	16
212.105TS	412.105TS	18
212.106TS	412.106TS	20
212.107TS	412.107TS	22
212.108TS	412.108TS	24
212.109TS	412.109TS	26
212.110TS	412.110TS	28
212.111TS	412.111TS	30
212.112TS	412.112TS	32
212.113TS	412.113TS	34
212.114TS	412.114TS	35
212.115TS	412.115TS	36
212.116TS	412.116TS	38
212.117TS	412.117TS	40
212.118TS	412.118TS	42
212.119TS	412.119TS	45
212.120TS	412.120TS	48
212.121TS	412.121TS	50
212.122TS	412.122TS	52
212.123TS	412.123TS	55
212.134TS	412.134TS	44
212.136TS	412.136TS	46
02.212.054TS	04.212.054TS	54
02.212.056TS	04.212.056TS	56

Locking Screw Stardrive Ø 3.5 mm, perforated, sterile tube

Stainless Steel	Titanium (TiCP)	Length (mm)
02.125.124TS	04.125.124TS	24
02.125.126TS	04.125.126TS	26
02.125.128TS	04.125.128TS	28
02.125.130TS	04.125.130TS	30
02.125.132TS	04.125.132TS	32
02.125.134TS	04.125.134TS	34
02.125.136TS	04.125.136TS	36
02.125.138TS	04.125.138TS	38
02.125.140TS	04.125.140TS	40
02.125.142TS	04.125.142TS	42
02.125.144TS	04.125.144TS	44
02.125.146TS	04.125.146TS	46
02.125.148TS	04.125.148TS	48
02.125.150TS	04.125.150TS	50
02.125.152TS	04.125.152TS	52
02.125.154TS	04.125.154TS	54

VA Locking Screw Stardrive® Ø 3.5 mm, sterile tube

Stainless Steel	Titanium (TiCP)	Length (mm)
02.127.110TS	N/A	10
02.127.112TS	N/A	12
02.127.114TS	N/A	14
02.127.116TS	N/A	16
02.127.118TS	N/A	18
02.127.120TS	N/A	20
02.127.122TS	N/A	22
02.127.124TS	N/A	24
02.127.126TS	N/A	26
02.127.128TS	N/A	28
02.127.130TS	N/A	30
02.127.132TS	N/A	32
02.127.134TS	N/A	34
02.127.136TS	N/A	36
02.127.138TS	N/A	38
02.127.140TS	N/A	40
02.127.142TS	N/A	42
02.127.144TS	N/A	44
02.127.146TS	N/A	46
02.127.148TS	N/A	48
02.127.150TS	N/A	50
02.127.152TS	N/A	52
02.127.154TS	N/A	54
02.127.156TS	N/A	56

**Cancellous Bone Screw Ø 4.0 mm,
fully threaded, hex, sterile tube**

Stainless Steel	Titanium (TiCP)	Length (mm)
206.010TS	406.010TS	10
206.012TS	406.012TS	12
206.014TS	406.014TS	14
206.016TS	406.016TS	16
206.018TS	406.018TS	18
206.020TS	406.020TS	20
206.022TS	406.022TS	22
206.024TS	406.024TS	24
206.026TS	406.026TS	26
206.028TS	406.028TS	28
206.030TS	406.030TS	30
206.032TS	406.032TS	32
206.035TS	406.035TS	35
206.036TS	N/A	36
206.040TS	406.040TS	40
206.045TS	406.045TS	45
206.050TS	406.050TS	50
206.055TS	406.055TS	55

**Cancellous Bone Screw Ø 4.0 mm,
short thread, hex, sterile tube**

Stainless Steel	Titanium (TiCP)	Length (mm)
207.010TS	407.010TS	10/5
207.012TS	407.012TS	12/5
207.014TS	407.014TS	14/5
207.016TS	407.016TS	16/6
207.018TS	407.018TS	18/7
207.020TS	407.020TS	20/8
207.022TS	407.022TS	22/9
207.024TS	407.024TS	24/10
207.026TS	407.026TS	26/12
207.028TS	407.028TS	28/14
207.030TS	407.030TS	30/14
207.035TS	407.035TS	35/14
207.040TS	407.040TS	40/14
207.045TS	407.045TS	45/15
207.050TS	407.050TS	50/15
207.055TS	407.055TS	55/16



DePuy Synthes

PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

People inspired™



Synthes GmbH
Eimattstrasse 3
4436 Oberdorf
Switzerland
Tel: +41 61 965 61 11
Fax: +41 61 965 66 00
www.depuySynthes.com

This publication is not intended for distribution
in the USA.

All surgical techniques are available as PDF files at
www.depuySynthes.com/ifu

CE
0123