

Thin. Nimble. And **incredibly precise***.



Introducing the
**ECHELON FLEX™ Powered Vascular
Stapler with Advanced
Placement Tip**

Enables the **most precise***
placement on fragile vessels

- **26% narrower curved, blunt tip anvil,**¹ providing superior maneuverability around vascular structures
- **26% thinner shaft,**¹ offering the greatest angle of reach in the tight intercostal space
- **11% greater manual articulation in each direction,**² allowing more flexibility during final placement

A narrower anvil, thinner shaft and more articulation than Endo GIA™ Curved Tip Reload with Tri-Staple™ Technology.

Enhanced stability on fragile vessels

- **83% reduction in tip movement**³ during firing for less movement during transection, potentially resulting in less trauma


Examples of applicable procedures

Designed specifically for transection of fragile and/or difficult-to-reach vessels

- **Thoracic**—lobectomy, segmentectomy, pneumonectomy
- **Solid organ**—nephrectomy, donor nephrectomy, hepatectomy

Product components

- ECHELON FLEX™ Powered Vascular Stapler with Advanced Placement Tip
- ENDOPATH ECHELON™ Vascular White Reload for Advanced Placement Tip



ECHELON FLEX™ Powered Vascular Stapler

CODE	DESCRIPTION	QUANTITY PER SALES UNIT
PVE35A	35mm Standard Articulating Endoscopic Linear Cutter,* 320mm	3

*Instrument does not contain a reload.



ENDOPATH ECHELON™ Vascular White Reload for Advanced Placement Tip (35mm, 4row)

CODE	COLOR	DESCRIPTION	OPEN STAPLE HEIGHT	STAPLE ROWS	CLOSED STAPLE HEIGHT	QUANTITY PER SALES UNIT
VASECR35	White	White reload intended for use in vascular/thin tissue	2.5mm	4	1.0mm	12

To learn more, or for more information on our products please call customer service in Australia on 02 9815 4278 or 1800 252 194 and in New Zealand on 0800 803 988

AUSTRALIA: Johnson & Johnson Medical Pty Ltd, 1-5 Khartoum Road, North Ryde NSW 2113

NEW ZEALAND: Johnson & Johnson (NZ) Ltd, 507 Mt Wellington Hwy, Mt Wellington. Auckland 1060 New Zealand

©2016 Ethicon, Inc. All rights reserved. 088177-180227 - April 2020

For complete product details, see Instructions for Use.

The third-party trademarks used herein are trademarks of their respective owners.