


Megadyne™

MEGADYNE™ LLETZ Loops provide
**a full line of electrodes for
histological tissue excisions**



Easy-to-use, versatile LLETZ Loop Electrodes

- Color-coded for easy size and shape recognition
- Disposable
- Cut, coag and blend capabilities
- Fine tungsten wire loop
- Stable “Y” construction
- Full insulation over all but the distal wire loop
- Fits most standard electrocautery pencils






MEGADYNE™ LLETZ Loop Electrodes are available in nine color-coded shapes and sizes

Ordering information

Loop Electrodes

DESCRIPTION	CATALOG NUMBER	QUANTITY PER CASE
 10mm wide x 7mm deep with 11.8cm shaft orange	0410	12
 10mm wide x 5mm deep with 11.8cm shaft purple	0420	12
 15mm wide x 6mm deep with 11.8cm shaft grey	0430	12
 20mm wide x 10mm deep with 11.8cm shaft red	0440	12
 10mm wide x 10mm deep with 11.8cm shaft yellow	0450	12
 15mm wide x 12mm deep with 11.8cm shaft green	0460	12
 20mm wide x 12mm deep with 11.8cm shaft white	0470	12
 20mm wide x 15mm deep with 11.8cm shaft blue	0480	12
 20mm wide x 8mm deep with 11.8cm shaft tan	0490	12

Accessories

DESCRIPTION	CATALOG NUMBER	QUANTITY PER CASE
 5-inch (12.7cm) Ball Electrode	0009	12
 Rocker Switch Pencil with E-Z Clean Electrosurgical Tip	0030H	50
 Button Switch Pencil with E-Z Clean Electrosurgical Tip	0035H	50