

Megadyne™

MEGADYNE™

Resposable Laparoscopic Electrode System



MEGADYNE™ Resposable Laparoscopic Electrode System

Disposable MegaTIPs and Reusable Indicator Shafts


The MEGADYNE™ Resposable Laparoscopic Electrode System provides an economic solution for monopolar energy during laparoscopy.

Disposable MegaTIPs

- 9 different tip options to meet your surgical needs
- Coated with PTFE to aid in the easy removal of eschar with a damp gauze or sponge

Reusable Indicator Shaft










- Provides a clear indication when replacement is needed
 - 2 layers of insulation: Black outer layer and yellow inner layer
 - Insulation damage (cuts, nicks and scrapes) is evident when yellow is visible on the shaft. This indicates the shaft should be discarded and replaced.
- Compatible with steam sterilization







For complete product details, see Instructions for Use.

Ordering information


MEGADYNE™ Disposable MegaTIPs

DESCRIPTION		CATALOG NUMBER	QUANTITY PER CASE
 MEGADYNE	L-Wire Disposable MegaTIP	0600	6
 MEGADYNE	L-Wire Modified Disposable MegaTIP	0600M	6
 MEGADYNE	J-Wire Disposable MegaTIP	0605	6
 MEGADYNE	Curved Spatula Disposable MegaTIP	0618	6
 MEGADYNE	Curved Blade Disposable MegaTIP	0619	6
 MEGADYNE	L-Hook Disposable MegaTIP	0620	6
 MEGADYNE	L-Hook Modified Disposable MegaTIP	0620M	6
 MEGADYNE	J-Hook Disposable MegaTIP	0621	6
 MEGADYNE	J-Hook Modified Disposable MegaTIP	0621M	6

MEGADYNE™ Reusable Indicator Shafts

DESCRIPTION		CATALOG NUMBER	QUANTITY PER CASE
 32cm Reusable Indicator Shaft	32cm Reusable Indicator Shaft	0690	2
 32cm Reusable Indicator Shaft with split stem	32cm Reusable Indicator Shaft with split stem	0690S	2
 38cm Reusable Indicator Shaft	38cm Reusable Indicator Shaft	0695	2
 38cm Reusable Indicator Shaft with split stem	38cm Reusable Indicator Shaft with split stem	0695S	2

Accessory

DESCRIPTION		CATALOG NUMBER	QUANTITY PER CASE
 Reusable Foot-Control Cord (10-foot)	Reusable Foot-Control Cord (10-foot)	0075	1