


Megadyne™

Designed to **enhance** **patient protection**

MEGADYNE™ MEGA SOFT™
Reusable Patient Return Electrodes



ETHICON
PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

Shaping
the future
of surgery

MEGADYNE™ MEGA SOFT™ Reusable Patient Return Electrodes eliminate the need for small disposable sticky return pads

Versatile and easy-to-use patient return electrodes designed for a wide variety of electrosurgical procedures

Enhanced patient protection

- No adhesives—eliminates adhesive-related injuries
- Radiolucent and latex-free
- Can be used with jewelry, tattoos, piercings and implants—they will not concentrate energy and increase temperature

75+ million procedures




0 Pad-site burns¹

Ease of use

- No skin prep required and no shaving
- MEGADYNE™ MEGA SOFT™ Universal and Universal Plus pads can be used for all patient sizes²
- No skin contact required: May be used with a sheet and draw sheet between patient and pad
- Can be used with patient warming devices³

Compatible with



Most isolated electrosurgical generators⁴

Reduced cost and waste with increased efficiencies

- Reduces costs vs. sticky pads⁵
- Meets three of five facility green initiatives⁶
- Reduces SKUs and ordering hassles⁵
- Reduces waste and waste-disposal costs⁵
- Saves time: Less than half the steps in the operating room vs. sticky pads⁷




Reusable for 24 months⁵



Environmentally friendly

¹ April 2018, per internal research conducted by Ethicon. Pad site burn is defined as a burn that occurs between the pad and patient and does not include alternate site or chemical burns. For a complete list of warning and precautions please refer to the Instructions for Use. (C2681) ² As per the IFU the device is intended to be used as a 24-month reusable patient return electrode for patients from .35 kg (0.8lb) and up (no maximum restriction). (C2630) ³ Refer to MEGA SOFT™ and respective patient warming device IFUs. (C2411) ⁴ Per internal research based on IEC 5th edition. Consult the IFU of the electrosurgical generator being used for further information, warnings and contraindications. (C2500) ⁵ Compared to traditional "sticky" patient return electrode. (C2381) (C2341) ⁶ Kwakye G, Brat GA, Makary MA. Green surgical practices for health care. Arch Surg 2011;146:131-6. ⁷ Summary of AORN guidelines for adhesive return electrode preparation, placement and removal and capacitive return electrode preparation, placement and removal.

A full line of MEGADYNE™ MEGA SOFT™ Reusable Patient Return Electrodes for your OR

	MEGA SOFT™	MEGA SOFT™ Universal Plus	MEGA SOFT™ Universal
Product			
Pad size	117 cm (46") long x 51 cm (20") wide	91 cm (36") long x 51 cm (20") wide	91 cm (36") long x 51 cm (20") wide
Pad thickness	1.25 cm (1/2")	0.635 cm (1/4")	0.3 cm (1/8")
Key features	Pressure reduction	Patient comfort, light-weight, pressure reduction	Versatile and lightweight
Patient size	11.5 kg (25 lbs.) and up	350 grams (0.8 lbs.) and up Pediatrics to adults.	350 grams (0.8 lbs.) and up Pediatrics to adults.
Pad weight	8.4 kg (18.5 lbs.)	3.6 kg (8 lbs.)	1.8 kg (4 lbs.)
Pad use	One-sided	One-sided	Two-sided
Lightweight		✓	✓
Flexible		✓	✓
Pressure reduction	✓	✓	
Comfort	✓	✓	
Max life	24 months	24 months	24 months
Product codes & cord options	0830 Single cord option. 0835 Dual cord option for use with two generators.	0847 Single cord option. 0848 Dual cord option for use with two generators.	0845 Single cord option. 0846 Dual cord option for use with two generators.

Designed for a wide variety of electrosurgical procedures



Flexible for easy handling in the OR.

MEGA SOFT™ Universal Plus is available in a dual-corded pad option for use with two generators. A compatibility cable is also available for use with Contact Quality Monitoring isolated generators and robots.

Ordering information

Please contact Johnson & Johnson Customer Service Department.

- Australia on 02 9815 4278 or 1800 252 194
- New Zealand on 0800 803 988

Product Codes

- 0848 Mega Soft™ Universal Plus Single Cord
- 0847 Mega Soft™ Universal Plus Dual Cord
- 0846 Mega Soft™ Universal Single Cord
- 0845 Mega Soft™ Universal Dual Cord
- 0840 Mega Soft™ Paediatric Single Cord
- 0835 Mega Soft™ Dual Cord
- 0830 Mega Soft™ Single Cord
- M2K01 Mega Soft™ Reusable Cable, Standard Connector 2.4m
- M2K02 Mega Soft™ Reusable Cable, Standard Connector 4.4m
- M2K06 Mega Soft™ Reusable Cable, Argon Beam Connector 2.4m
- M2K07 Mega Soft™ Replacement Cable
- M2K08 Mega Soft™ Reusable Cable, Compatibility Connector 2.4m
- M2K09 Mega Soft™ Reusable Cable, Compatibility Connector 4.4m
- PKIT001 Mega Soft™ Repair Kit

For customer support, contact your local Ethicon Product Specialist.