


# RIGIDLOOP® Cortical Fixation System



## A Complete Family of Cortical Button Devices


- Comprehensive product offering to address multiple sports medicine procedures
- Strong, dependable cortical fixation<sup>1</sup>
- Accommodates various tunnel sizes



# The RIGIDLOOP® Cortical System Family of Devices Delivers Full-Scale Reliability and Versatility



This versatile titanium button comes pre-loaded with a white leading suture and green flipping suture. It can be used for fixation of soft tissue to bone in shoulder, elbow and knee procedures.



The fixed loop button offers strong, dependable fixation strength while minimizing graft laxity.<sup>1</sup> It comes in multiple loop sizes ranging from 15mm – 60mm. An oversized bailout button is also available.



The adjustable button delivers simplicity and performance with infinite loop tensioning. It has the lowest displacement vs. other tested adjustable cortical implants.<sup>2</sup>

Limited Warranty and Disclaimer: DePuy Mitek, Inc. products are sold with a limited warranty to the original purchaser against defects in workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.

WARNING: In the USA, this product has labeling limitations. See package insert for complete information.

CAUTION: USA Law restricts these devices to sale by or on the order of a physician.

## Ordering Information

SKU	Implants (Sterile)
232447	RIGIDLOOP Adjustable Cortical System – STND
232448	RIGIDLOOP Adjustable Cortical System – LONG
232449	RIGIDLOOP Adjustable Cortical System – XL

SKU	Instruments (Sterile)
232453	Combo Beath/Drill Pin – 4.5mm
232460	Arrow Tip Combo Beath/Drill Pin - 4.5mm
218034	RIGIDLOOP Adjustable Disposables Kit*

\*Includes: Malleable graft retractor, ruler, marker, tibial cannula, 4.5mm combo beath pin/reamer, 14" drill pin, 15" Nitinol guide wire.

SKU	Implants (Sterile)
232008	RIGIDLOOP Titanium Button
232009	RIGIDLOOP 15mm Fixed Loop
232011	RIGIDLOOP 20mm Fixed Loop
232014	RIGIDLOOP 25mm Fixed Loop
232016	RIGIDLOOP 30mm Fixed Loop
232017	RIGIDLOOP 35mm Fixed Loop
232018	RIGIDLOOP 40mm Fixed Loop
232019	RIGIDLOOP 45mm Fixed Loop
232021	RIGIDLOOP 50mm Fixed Loop
232022	RIGIDLOOP 55mm Fixed Loop
232023	RIGIDLOOP 60mm Fixed Loop
232034	RIGIDLOOP XL Button

SKU	Instruments
232037	RIGIDLOOP ACL Disposables Kit (sterile)*
232025	RIGIDLOOP Beath Pin
232028	RIGIDLOOP 4.5 mm Cannulated Reamer
232030	RIGIDLOOP Standard Depth Gauge
232031	RIGIDLOOP Variable Depth Gauge
232032	RIGIDLOOP Instrument Caddy
232033	RIGIDLOOP Silicon Mat
215451	Instrument tray lid
232036	Tray and Lid (BOM)

\*Includes RIGIDLOOP custom beath pin, 4.5mm Reamer, Standard Depth Gauge, Guidewire, Malleable graft retractor, Cannula, Ruler, Marker

## References:

1. DHF-102035-MEMO14 &16
2. Adjustable loop ACL suspension devices demonstrate less reliability in terms of reproducibility and irreversible displacement," Ahmad, S.S., Hirschmann, M.T., Voumard, B. et al. Knee Surg Sports Traumatology Arthroscopy (2017)



**Medos International SÀRL**  
Chemin-Blanc 38  
CH-2400 Le Locle  
Switzerland

US Representative:  
**DePuy Mitek, Inc.**  
325 Paramount Drive  
Raynham, MA 02767  
T. +1 (800) 382-4682

**Johnson & Johnson Medical Pty Ltd**  
t/a DePuy Synthes  
1-5 Khartoum Road  
North Ryde NSW 2113  
Australia

