

NO METAL LEFT BEHIND

# HEALIX ADVANCE™ Self-Punching Anchor



Single-step anchor insertion with the  
HEALIX Advance Anchor™ Technology you can trust.  
All without implanted metal.

The HEALIX ADVANCE™ SP Anchor delivers a self-punching solution that demonstrates improved usability and versatility over previous generations of anchors.<sup>5</sup>

## Usability

- Minimal steps for anchor insertion
- Locating pre-made drill holes is no longer necessary

## Versatility

Full compatibility with Mitek Sports Medicine Sutures

### DYNACORD® Suture

- Self-tensioning when needed, to assist in the control of post-operative suture laxity<sup>2</sup>
- 6 x Capacity (5.5mm anchor)

### PERMATAPE® Suture

- 100% Ultra-high molecular weight polyethylene (UHMWPE) flat braid for strength and handling<sup>3,4</sup>
- 3 x Capacity (5.5mm anchor)<sup>3,4</sup>

## HEALIX ADVANCE™ Anchor Technology<sup>1</sup>

- Multithread configuration engages both cortical and cancellous bone more effectively than single-thread designs
- Cannulation to channel blood flow to surface



## Product Codes

### HEALIX ADVANCE™ SELF-PUNCHING Anchors

228055	4.9mm HEALIX ADVANCE SP PEEK Anchor
228056	5.5mm HEALIX ADVANCE SP PEEK Anchor

### HEALIX ADVANCE™ Anchor Instruments

222314	4.5mm HEALIX ADVANCE Anchor Awl
222315	5.5/6.5mm HEALIX ADVANCE Anchor Awl
222316	4.5mm HEALIX ADVANCE Anchor Cortical Awl/Tap
222317	5.5/6.5mm HEALIX ADVANCE Anchor Cortical Awl/Tap
222329	4.5mm HEALIX ADVANCE Anchor Full Tap
222328	5.5/6.5mm HEALIX ADVANCE Anchor Full Tap

## Reference

1. Internal test report RR-101380-R Rev A ACS
2. Internal Test Report 103394861
3. Internal Test Report 103275220
4. Internal Test Report 103294799
5. Internal Test Report 103522100

Developed in collaboration with Thomas Knapp, M.D., Santa Monica, CA

Limited Warranty and Disclaimer: DePuy Synthes Mitek Sports Medicine products are sold with a limited warranty to the original purchaser against defects in workmanship and materials. Any other express or implied warranties, including warranties of merchantability or fitness, are hereby disclaimed.

WARNING: Please also refer to the package insert(s) or other labeling associated with the devices identified in this brochure for additional information.

CAUTION: The third-party trademarks used herein are the trademarks of their respective owners.



Johnson & Johnson Medical Pty Ltd  
t/a DePuy Synthes  
1-5 Khartoum Road  
North Ryde NSW 2113  
Australia

Medos International SÀRL  
Chemin-Blanc 38  
CH-2400 Le Locle  
Switzerland

US Representative:  
DePuy Mitek, Inc.  
325 Paramount Drive  
Raynham, MA 02767  
T. +1 (800) 382-4682