


LIGACLIP®

Outstanding **control** and **performance**

Endoscopic Clip Appliers



Features

- 360° rotation knob/shaft – designed for easy rotation and one-handed use
- Rigid and low glare shaft – made of a material that minimizes reflection and minimizes instrument movement during firing
- Fully automatic clip loading system for a smooth, controlled one-step firing sequence
- Anti-clip drop ratchet mechanism
- Last clip lockout with indicator bar

Product Overview

- The LIGACLIP® Endoscopic Clip Applier is a sterile single-patient-use instrument
- Options include 8.8mm closed clip length for ER320 and 11.0mm for ER420
- Devices are preloaded with 20 titanium clips with lateral and transverse grooves to resist clip dislodgement

Ordering Information

Product Code	Description	Clip Size	Shaft Length	Trocar Size	Color Code	Qty/box
ER320	LIGACLIP® Endoscopic Rotating Multiple-Clip Applier	Medium/Large	29cm	10mm/11mm	Gray	6
ER420	LIGACLIP® Endoscopic Rotating Multiple-Clip Applier	Large	34cm	12mm	Green	6

How to Order

All purchase orders are made to Johnson & Johnson Health Care Systems, Inc. (JJHCS) or your authorized distributor.

Electronic Ordering Options

(Note: Placing orders electronically avoids minimum order fees for your hospital.)

Johnson & Johnson Gateway:
www.jnjgateway.com/commerce

For questions about your order, visit the website or call 1-866-JNJ-GATE (565-4283).

Global Healthcare Exchange (GHX):
www.ghx.com

For questions about your order, visit the website or call 1-800-YOUR-GHX (968-7449).

Electronic Data Interchange (EDI):
JJHCS EDI Help Line: 1-800-262-2888

Nonelectronic/Manual Ordering Options

Call JJHCS at 1-800-255-2500 between 8:30 a.m. and 9:00 p.m. Eastern Standard Time or fax your order to 1-732-562-2212.

Customer support

For product use assistance, clinical guidelines, service and repair, emergency assistance, copy of a 510(k) clearance letter or complaints, please contact our Customer Support Center at **customersupport@eesus.jnj.com** or by calling **1-877-ETHICON** (384-4266). Our Support Center is staffed 24 hours a day, 7 days a week by qualified nurses to answer your product-related questions.