

KINCISE™
SURGICAL AUTOMATED SYSTEM

The power of precision in your hands



Attune®
Knee System
Revision Solutions

Another Innovation for the
KINCISE™ Surgical Automated System

VELYS™
Digital Surgery

The KINCISE™ System ATTUNE® Revision Broach Adapter, along with the ATTUNE® Revision Knee System, is designed to automate the sleeve-preparation process for revision knee surgery.

DePuy Synthes
PART OF THE JOHNSON & JOHNSON FAMILY OF COMPANIES

The KINCISE™ Surgical Automated System— Transforming Orthopaedic Care

Constant and Consistent Application of Energy. Augments Surgical Skills. Aids in Reducing Surgeon Fatigue and Work-Related Injuries.

The KINCISE™ System connects directly to the ATTUNE® Revision Knee System Tibial and Femoral Broaches to impact the broach into the bone until the desired preparation is made. For the surgical technique, please refer to the ATTUNE® Revision Knee System Rotating Platform and Fixed Bearing (Femoral Broaching) Surgical Technique guide, appendix A.



The KINCISE™ Surgical Automated System in conjunction with the KINCISE™ System ATTUNE® Revision Broach Adapter can now be used for ATTUNE® Revision RP tibial and femoral sleeve bone preparation.

KINCISE™ ATTUNE® Revision Broach Adapter

| Part No. | Product Description |
|-------------|---|
| 2544-01-108 | KINCISE™ System ATTUNE® Revision Broach Adapter |
| 2544-01-106 | KINCISE™ System ATTUNE® Adapters Tray |
| 2004-20-933 | TruMatch™ TRIAL SIZE LID |

KINCISE™ Surgical Automated System

| Part No. | Product Description |
|-------------|---|
| 1000-00-101 | KINCISE™ Automated Surgical Impactor |
| 1002-00-102 | KINCISE™ Rechargeable Lithium-ion Battery |
| 1003-00-101 | KINCISE™ 4-Port Battery Charger |
| ME03 | KINCISE™ Sterilization Container/Tray - Surgical Impactor |
| ME04 | KINCISE™ Sterilization Container/Tray - Battery |
| ME03KIT | KINCISE™ Surgical Impactor Tray Disposable Kit |
| ME04KIT | KINCISE™ Battery Tray Disposable Kit |

Contact your DePuy Synthes representative or visit jnmedicaldevices.com

For more information on the KINCISE™ System, visit KINCISEsystem.com

Please refer to the IFU (Instructions for Use) for a complete list of indications, contraindications, precautions and warnings. For further information on DePuy Synthes products, please contact your local DePuy Synthes representative.

Please refer to the package insert(s) or other labeling associated with the devices identified in this sales document for additional information.

CAUTION: Federal Law restricts these devices to sale by or on the order of a physician. Some devices listed in this sales aid may not have been licensed in accordance with Canadian law and may not be for sale in Canada. Please contact your sales consultant for items approved for sale in Canada.

Not all products may currently be available in all markets.

www.depuyssynthes.com

The third party trademarks used herein are the trademarks of their respective owners.