

Inventory Control Form

2.0 mm Quick Insertion Screw (QIS) and 2.0, 2.5, 3.0, 3.5, 4.0 mm Cannulated Compression Headless Screws (CCHS)

Case Information (eg, Case ID):

DO NOT PUT ANY
PATIENT IDENTIFIABLE
INFORMATION IN
THIS FIELD

To order (USA): 800-523-0322
To order (Canada): 844-243-4321

Date: _____
Hospital: _____
Surgeon: _____
Procedure: _____

Implants The screws are made of titanium alloy (TAV or TAN)

2.0 mm Quick Insertion Screw (TAN)



Length	Part Number
<input type="checkbox"/>	11 mm 04.228.511
<input type="checkbox"/>	12 mm 04.228.512
<input type="checkbox"/>	13 mm 04.228.513
<input type="checkbox"/>	14 mm 04.228.514
<input type="checkbox"/>	15 mm 04.228.515
<input type="checkbox"/>	16 mm 04.228.516
<input type="checkbox"/>	17 mm 04.228.517
<input type="checkbox"/>	18 mm 04.228.518

2.5 mm Cannulated Compression Headless Screw (TAV)



Length	Short Thread Part Number	Length	Long Thread Part Number
<input type="checkbox"/>	10 mm 04.333.110	<input type="checkbox"/>	N/A
<input type="checkbox"/>	12 mm 04.333.112	<input type="checkbox"/>	N/A
<input type="checkbox"/>	14 mm 04.333.114	<input type="checkbox"/>	N/A
<input type="checkbox"/>	16 mm 04.333.116	<input type="checkbox"/>	N/A
<input type="checkbox"/>	18 mm 04.333.118	<input type="checkbox"/>	N/A
<input type="checkbox"/>	20 mm 04.333.120	<input type="checkbox"/>	20 mm 04.334.120
<input type="checkbox"/>	22 mm 04.333.122	<input type="checkbox"/>	22 mm 04.334.122
<input type="checkbox"/>	24 mm 04.333.124	<input type="checkbox"/>	24 mm 04.334.124
<input type="checkbox"/>	26 mm 04.333.126	<input type="checkbox"/>	26 mm 04.334.126
<input type="checkbox"/>	28 mm 04.333.128	<input type="checkbox"/>	28 mm 04.334.128
<input type="checkbox"/>	30 mm 04.333.130	<input type="checkbox"/>	30 mm 04.334.130
<input type="checkbox"/>	32 mm 04.333.132	<input type="checkbox"/>	32 mm 04.334.132
<input type="checkbox"/>	34 mm 04.333.134	<input type="checkbox"/>	34 mm 04.334.134
<input type="checkbox"/>	36 mm 04.333.136	<input type="checkbox"/>	36 mm 04.334.136
<input type="checkbox"/>	38 mm 04.333.138	<input type="checkbox"/>	38 mm 04.334.138
<input type="checkbox"/>	40 mm 04.333.140	<input type="checkbox"/>	40 mm 04.334.140

2.0 mm Cannulated Compression Headless Screw (TAV)[◇]



Length	Short Thread Part Number	Length	Long Thread Part Number
<input type="checkbox"/>	10 mm 04.333.010	<input type="checkbox"/>	N/A
<input type="checkbox"/>	12 mm 04.333.012	<input type="checkbox"/>	N/A
<input type="checkbox"/>	14 mm 04.333.014	<input type="checkbox"/>	N/A
<input type="checkbox"/>	16 mm 04.333.016	<input type="checkbox"/>	N/A
<input type="checkbox"/>	18 mm 04.333.018	<input type="checkbox"/>	N/A
<input type="checkbox"/>	20 mm 04.333.020	<input type="checkbox"/>	20 mm 04.334.020
<input type="checkbox"/>	22 mm 04.333.022	<input type="checkbox"/>	22 mm 04.334.022
<input type="checkbox"/>	24 mm 04.333.024	<input type="checkbox"/>	24 mm 04.334.024
<input type="checkbox"/>	26 mm 04.333.026	<input type="checkbox"/>	26 mm 04.334.026
<input type="checkbox"/>	28 mm 04.333.028	<input type="checkbox"/>	28 mm 04.334.028
<input type="checkbox"/>	30 mm 04.333.030	<input type="checkbox"/>	30 mm 04.334.030

Thread Length
Short Thread: approx. 25% Total Length or 4 mm
Long Thread: approx. 40% Total Length

Thread Length
Short Thread: approx. 25% Total Length or 4 mm
Long Thread: approx. 40% Total Length

[◇] All implants are also available packed in sterile tubes. Add suffix "TS" to part number
Sterile part availability might be different depending on the country or region.

Inventory Control Form

2.0 mm Quick Insertion Screw (QIS) and 2.0, 2.5, 3.0, 3.5, 4.0 mm Cannulated Compression Headless Screws (CCHS)

Implants The screws are made of titanium alloy (TAV)

3.0 mm Cannulated Compression Headless Screw[◇]



Length	Short Thread Part Number		Length	Long Thread Part Number
<input type="checkbox"/>	10 mm	04.333.210	<input type="checkbox"/>	N/A
<input type="checkbox"/>	12 mm	04.333.212	<input type="checkbox"/>	N/A
<input type="checkbox"/>	14 mm	04.333.214	<input type="checkbox"/>	N/A
<input type="checkbox"/>	16 mm	04.333.216	<input type="checkbox"/>	N/A
<input type="checkbox"/>	18 mm	04.333.218	<input type="checkbox"/>	N/A
<input type="checkbox"/>	20 mm	04.333.220	<input type="checkbox"/>	20 mm 04.334.220
<input type="checkbox"/>	22 mm	04.333.222	<input type="checkbox"/>	22 mm 04.334.222
<input type="checkbox"/>	24 mm	04.333.224	<input type="checkbox"/>	24 mm 04.334.224
<input type="checkbox"/>	26 mm	04.333.226	<input type="checkbox"/>	26 mm 04.334.226
<input type="checkbox"/>	28 mm	04.333.228	<input type="checkbox"/>	28 mm 04.334.228
<input type="checkbox"/>	30 mm	04.333.230	<input type="checkbox"/>	30 mm 04.334.230
<input type="checkbox"/>	32 mm	04.333.232	<input type="checkbox"/>	32 mm 04.334.232
<input type="checkbox"/>	34 mm	04.333.234	<input type="checkbox"/>	34 mm 04.334.234
<input type="checkbox"/>	36 mm	04.333.236	<input type="checkbox"/>	36 mm 04.334.236
<input type="checkbox"/>	38 mm	04.333.238	<input type="checkbox"/>	38 mm 04.334.238
<input type="checkbox"/>	40 mm	04.333.240	<input type="checkbox"/>	40 mm 04.334.240

Thread Length
 Short Thread: approx. 25% Total Length or 4 mm
 Long Thread: approx. 40% Total Length

3.5 mm Cannulated Compression Headless Screw (TAV)[◇]



Length	Short Thread Part Number		Length	Long Thread Part Number
<input type="checkbox"/>	14 mm	04.333.314	<input type="checkbox"/>	N/A
<input type="checkbox"/>	16 mm	04.333.316	<input type="checkbox"/>	N/A
<input type="checkbox"/>	18 mm	04.333.318	<input type="checkbox"/>	N/A
<input type="checkbox"/>	20 mm	04.333.320	<input type="checkbox"/>	N/A
<input type="checkbox"/>	22 mm	04.333.322	<input type="checkbox"/>	N/A
<input type="checkbox"/>	24 mm	04.333.324	<input type="checkbox"/>	24 mm 04.334.324
<input type="checkbox"/>	26 mm	04.333.326	<input type="checkbox"/>	26 mm 04.334.326
<input type="checkbox"/>	28 mm	04.333.328	<input type="checkbox"/>	28 mm 04.334.328
<input type="checkbox"/>	30 mm	04.333.330	<input type="checkbox"/>	30 mm 04.334.330
<input type="checkbox"/>	32 mm	04.333.332	<input type="checkbox"/>	32 mm 04.334.332
<input type="checkbox"/>	34 mm	04.333.334	<input type="checkbox"/>	34 mm 04.334.334
<input type="checkbox"/>	36 mm	04.333.336	<input type="checkbox"/>	36 mm 04.334.336
<input type="checkbox"/>	38 mm	04.333.338	<input type="checkbox"/>	38 mm 04.334.338
<input type="checkbox"/>	40 mm	04.333.340	<input type="checkbox"/>	40 mm 04.334.340
<input type="checkbox"/>	42 mm	04.333.342	<input type="checkbox"/>	42 mm 04.334.342
<input type="checkbox"/>	44 mm	04.333.344	<input type="checkbox"/>	44 mm 04.334.344
<input type="checkbox"/>	46 mm	04.333.346	<input type="checkbox"/>	46 mm 04.334.346
<input type="checkbox"/>	48 mm	04.333.348	<input type="checkbox"/>	48 mm 04.334.348
<input type="checkbox"/>	50 mm	04.333.350	<input type="checkbox"/>	50 mm 04.334.350

Thread Length
 Short Thread: approx. 25% Total Length or 4 mm
 Long Thread: approx. 40% Total Length

[◇] All implants are also available packed in sterile tubes. Add suffix "TS" to part number
 Sterile part availability might be different depending on the country or region.

Inventory Control Form

2.0 mm Quick Insertion Screw (QIS) and 2.0, 2.5, 3.0, 3.5, 4.0 mm Cannulated Compression Headless Screws (CCHS)

Implants The screws are made of titanium alloy (TAV)

4.0 mm Cannulated Compression Headless Screw (TAV)[◇]



Length	Short Thread Part Number	Thread Length	Long Thread Part Number
<input type="checkbox"/>	14 mm 04.333.414	N/A	
<input type="checkbox"/>	16 mm 04.333.416	N/A	
<input type="checkbox"/>	18 mm 04.333.418	N/A	
<input type="checkbox"/>	20 mm 04.333.420	N/A	
<input type="checkbox"/>	22 mm 04.333.422	N/A	
<input type="checkbox"/>	24 mm 04.333.424	24 mm	04.334.424
<input type="checkbox"/>	26 mm 04.333.426	26 mm	04.334.426
<input type="checkbox"/>	28 mm 04.333.428	28 mm	04.334.428
<input type="checkbox"/>	30 mm 04.333.430	30 mm	04.334.430
<input type="checkbox"/>	32 mm 04.333.432	32 mm	04.334.432
<input type="checkbox"/>	34 mm 04.333.434	34 mm	04.334.434
<input type="checkbox"/>	36 mm 04.333.436	36 mm	04.334.436
<input type="checkbox"/>	38 mm 04.333.438	38 mm	04.334.438
<input type="checkbox"/>	40 mm 04.333.440	40 mm	04.334.440
<input type="checkbox"/>	42 mm 04.333.442	42 mm	04.334.442
<input type="checkbox"/>	44 mm 04.333.444	44 mm	04.334.444
<input type="checkbox"/>	46 mm 04.333.446	46 mm	04.334.446
<input type="checkbox"/>	48 mm 04.333.448	48 mm	04.334.448
<input type="checkbox"/>	50 mm 04.333.450	50 mm	04.334.450
<input type="checkbox"/>	52 mm 04.333.452	52 mm	04.334.452
<input type="checkbox"/>	54 mm 04.333.454	54 mm	04.334.454
<input type="checkbox"/>	56 mm 04.333.456	56 mm	04.334.456
<input type="checkbox"/>	58 mm 04.333.458*	58 mm	04.334.458*
<input type="checkbox"/>	60 mm 04.333.460*	60 mm	04.334.460*

Thread Length

Short Thread: approx. 25% Total Length or 5 mm

Long Thread: approx. 40% Total Length

[◇] Implants are also available sterile packaged. Add suffix "S" to part number for packaging in peel pouch (see lengths with * above) or "TS" for packaging in sterile tube. Sterile part availability might be different depending on the country or region.

Instruments

Instruments for 2.0 mm Screws ●

03.333.000 0.8 mm Guide Wire, L 100 mm Trocar Tip
 03.333.000S* (6 pk)

03.333.100 1.6 mm Cannulated Drill Bit, QC
 03.333.100S*

03.333.200 2.0 mm Cannulated Countersink, QC
 03.333.200S*

Instruments for 2.5/3.0 mm Screws ● ●

03.333.001 1.1 mm Guide Wire, L 150 mm, Trocar Tip
 03.333.001S* (6 pk)

03.333.101 2.0 mm Cannulated Drill Bit, QC
 03.333.101S*

03.333.201 2.5/3.0 mm Cannulated Countersink, QC
 03.333.201S*

Instruments for 3.5/4.0 mm Screws ● ●

03.333.002 1.4 mm Guide Wire, L 150 mm, Trocar Tip
 03.333.002S* (6 pk)

03.333.102 2.7 mm Cannulated Drill Bit, QC
 03.333.102S*

03.333.202 3.5/4.0 mm Cannulated Countersink, QC
 03.333.202S*

*Part numbers with a trailing 'S' designate the sterile part number. Sterile part availability might be different depending on the country or region.

