

DEVICE SPECS

Ethicon Coated VICRYL™ Sutures



**Coated VICRYL™ Plus
Antibacterial (polyglactin 910) Suture**



**Coated VICRYL™
(polyglactin 910) Suture**

Codes	See catalog	See catalog
Color	Violet or Undyed	Violet or Undyed
Antibacterial Coating	Yes	No
Material/Composition	Polyglactin 910, Triclosan	Polyglactin 910
Suture Structure	Braided	Braided
Absorbable	Yes	Yes
Absorption Rate	56-70 days	56-70 days
Directionality	Unidirectional	Unidirectional
Available Needle Options	Taperpoint, Tapercut, Cutting Edge Reverse, Cutting Edge Prime Reverse/Conventional, Blunt Point Ethiguard	Taperpoint, Tapercut, Cutting Edge Reverse/Conventional, Cutting Edge Spatula Reverse/Conventional, Cutting Edge Prime Reverse/Conventional, Cutting Edge Spatulated Center Point, Blunt Point Ethiguard, Blunt Point
Special Options	MutliPass™, CONTROL RELEASE™	MutliPass™, CONTROL RELEASE™
Pledget Options	No	No
Strands per Pack	1, 3, 4, 6, 8, or 12	1, 3, 4, 5, 6, 8, 10, or 12
Quantity per Box	12, 24, or 36	12, 24, or 36