

DERMABOND ADVANCED®

Topical Skin Adhesive

A protective barrier that adds strength and inhibits bacteria

- When used in addition to sutures, DERMABOND ADVANCED Adhesive was shown ex vivo to add 75% more strength to the wound closure than sutures alone¹
- DERMABOND ADVANCED Adhesive demonstrated in vitro inhibition of gram-positive bacteria (MRSA and MRSE) and gram-negative bacteria (E Coli)^{2*}
- Provides a microbial barrier with 99% protection in vitro for at least 72 hours against organisms commonly responsible for SSIs^{2*†}

A patented, thicker formula with a unique blend of plasticizers and an initiator for precise control[§]

- Plasticizers increase flexibility and 3-dimensional strength for secure wound closure and patient comfort^{3,4}
- Initiator enables consistent setting time across various conditions and skin types^{5,6}



Only
1
Layer Required⁷

2X
Faster Setting
Time (95 seconds)^{6†}

4X
More Viscous^{8‡}

8
US Patents⁹

ETHICON
PART OF THE *Johnson & Johnson* FAMILY OF COMPANIES

*Clinical significance is unknown.
†As compared to High Viscosity DERMABOND® Topical Skin Adhesive.
‡*Staphylococcus epidermidis, Staphylococcus aureus, Escherichia coli, Enterococcus faecium and Pseudomonas aeruginosa.*
§Compared to DERMABOND® Mini Topical Skin Adhesive

For more information, contact your Ethicon representative or call 1-800-255-2500.

For complete indications, contraindications, precautions, warnings, and adverse reactions, please reference full package insert.

DERMABOND ADVANCED Adhesive

Ordering Information

An applicator ergonomically designed for improved handling and easy coverage

- Angled tip for both fine and wide line application - up to 15cm¹⁰
- Innovative applicator designed to easily activate and avoid clogging
- Innovative design provides greater control and ability to see product level



Product Code	Description	Size	Packaging
DNX12	DERMABOND ADVANCED Adhesive <ul style="list-style-type: none"> • Unique formulation with an easy-to-use applicator for precise control for wide and fine line applications 	0.7 mL	12 Applicators per box
DNX6	DERMABOND ADVANCED Adhesive <ul style="list-style-type: none"> • Unique formulation with an easy-to-use applicator for precise control for wide and fine line applications 	0.7 mL	6 Applicators per box

For more information, contact your representative or call 1-800-255-2500.

For complete indications, contraindications, precautions, warnings, and adverse reactions, please reference full package insert.

References: 1. Protocol to determine the strength benefits when adding DERMABOND ADVANCED™ Topical Skin Adhesive to a wound or incision using suture as the primary closure, July 20112. 2. Bhende S, Rothenburger S, Spangler DJ, Dito M. In vitro assessment of microbial barrier properties of DERMABOND® Topical Skin Adhesive. Surg Infect (Larchmt). 2002;3(3):251-257. 3. Material Specification 53MS-0000681. 4. Singer A et al. Evaluation and management of traumatic lacerations. New England Journal of Medicine. 1997;337:1142-1148. 5. ADAPTIV Part PM62028, MWI-0000480. 6. ADAPTIV Report SUBOTH0004215. 7. DERMABOND ADVANCED Adhesive Instructions for Use. 8. FG-0000104 of 49cP, FG-0000101 of 214cP. 9. LAB0012187(v3). 10. 09PDO44 Final Protocol Report: Design Verification-Application Length and Width Evaluation of DERMABOND(R) NX Topical Skin Adhesive, 2010.

©2021 Ethicon, Inc. All rights reserved. 028909-210902

Please add DERMABOND ADVANCED® Topical Skin Adhesive to my preference card:

Doctor Name: _____

Signature: _____

Procedures: _____

All Procedures