

DEVICE SPECS

Ethicon Topical Skin Adhesives



DERMABOND™ PRINEO™
Skin Closure System



DERMABOND ADVANCED™
Topical Skin Adhesive



DERMABOND™ Mini
Topical Skin Adhesive

Codes	CLR222US, CLR422US, CLR602US	DNX6, DNX12	DHVM12
Sizes Available	22cm (CLR222US) 42cm (CLR422US) 60cm (CLR602US)	--	--
Adhesive Volume	2.9 mL (CLR222US) 3.8 mL (CLR422US, CLR602US)	0.7 mL	0.36 mL
Formulation	Mesh, Highly purified 2-OCA monomer	Highly purified 2-OCA monomer	Highly purified 2-OCA monomer
Application	Single layer over mesh	Single layer	≥ 2 layers
Drying Time	60 seconds	1.7 minutes	2.4 minutes
Microbial Barrier	Yes	Yes	Yes
Quantity per Box	1 adhesive applicator, and: 1 × 22cm mesh patch (CLR222US) 2 × 22cm mesh patches (CLR422US) 1 × 60cm mesh dispenser (CLR602US)	6 (DNX6) 12 (DNX12)	12