

SPORTS MEDICINE

PRODUCT CATALOG 2023




TABLE OF CONTENTS

Arthroscopic Equipment & Visualization	3
<i>Visualization</i>	4
<i>Fluid Management & Tissue Debridement</i>	9
<i>Radiofrequency</i>	12
Shoulder	15
<i>Rotator Cuff</i>	
<i>Suture Anchors</i>	16
<i>All-Suture Anchors</i>	26
<i>Knotless Anchors</i>	28
<i>Proximal Biceps</i>	
<i>Anchors & Screws</i>	30
<i>Labrum</i>	
<i>Suture Anchors</i>	32
<i>All-Suture Anchors</i>	35
<i>Glenoid</i>	
<i>Latarjet Instability Systems</i>	36
<i>Sutures & Suture Passers</i>	39
<i>Graft Preparation</i>	46
<i>Arthroscopic Instruments & Accessories</i>	47
Knee	51
<i>Ligament</i>	
<i>Cortical Fixation</i>	52
<i>Interference Screws</i>	54
<i>Screw and Sheath Fixation</i>	59
<i>Cross Pin Fixation</i>	61
<i>Graft Preparation</i>	62
<i>Meniscus</i>	64
<i>Cartilage</i>	66
<i>Tibial Realignment</i>	67
<i>Arthroscopic Instruments & Accessories</i>	68

Foot & Ankle	77
<i>Ankle & Tarsals</i>	
<i>Suture Anchors</i>	78
<i>All-Suture Anchors</i>	87
<i>Knotless Anchors</i>	88
<i>Metatarsals & Phalanges</i>	
<i>Suture Anchors</i>	89
<i>Sutures & Suture Passers</i>	93
<i>Graft Preparation</i>	98
<i>Arthroscopic Instruments & Accessories</i>	99
Hand, Wrist & Elbow	100
<i>Carpals & Phalanges</i>	
<i>Suture Anchors</i>	101
<i>Distal Biceps</i>	
<i>Suture Anchors</i>	104
<i>Sutures</i>	105
<i>Arthroscopic Instruments & Accessories</i>	109
Hip	110
<i>Labrum</i>	
<i>Anchors</i>	111
<i>Arthroscopic Instruments & Accessories</i>	112
Suture Portfolio	113
<i>Sutures & Suture Passers</i>	114
Allograft Tissue*	121
<i>Tendons</i>	122
<i>Cartilage</i>	125
<i>Acellular Dermis / Scaffolds</i>	126
<i>Structural Bone Grafts</i>	130
Procedure Systems	136
Non-Surgical Solutions	158

*Allograft tissue distributed by Parametrics Medical



Arthroscopic Equipment & Visualization



Arthroscopic Equipment

↳ Visualization

Visualization Equipment



Video Equipment

CODE	DESCRIPTION
242302	PUREVUE™ 4K Camera Control System with Integrated Light Source, Photo & Video Capture and Wireless Adaptor
242401	Camera Head with C-mount, Autoclavable
242288*	Video Cart with 110V Integrated Power Source and CO2 Tank Holder
242451	Video Cart Power Cord (US)
242284	Primary Multi-directional Display Mount
242518	Secondary Multi-directional Display Mount
242519**	Rotatable Tablet Arm
242285	Secondary Display Cart
218170	USB Keyboard (US)
242370	Sony UP-DR80MD Printer
242371	Sony 8.5 x 11 Ink & Paper RF UPC-R81MD

*Video cart does not include primary display mount; please order separately. [CODE: 242284].

**Tablet arm does not include mobile device-specific mounting bracket.

Software Accessories

CODE	DESCRIPTION
242327	DICOM Activation Key
242329	PUREVUE™ 4K Software Expansion Kit

***Current customers with EVO 4K will need to have upgrade kit (24232) ordered to upgrade to AIO software.



Arthroscopic Equipment

↳ Visualization

Visualization Equipment

Displays and Accessories

CODE	DESCRIPTION
242671	Medical Display, NDS Radiance Ultra 32" HD
242672	Medical Display, NDS Radiance Ultra 32" 4K/UHD
242394	Zerewire G2 Transmitter & Receiver Kit
242395	Zerewire G2 Access Kit - Pair: Radiance Ultra 32"
242396	Zerewire G2 Receiver Only
242397	Zerewire G2 Transmitter Only
242398	Zerewire G2 Access Kit - Single: Radiance Ultra 32"



Arthroscopic Equipment

↳ Visualization

Visualization Equipment

Couplers, Light Guides & Adaptors

CODE	DESCRIPTION
242437	Coupler, Fixed Focal, 19 mm, C-Mount, Autoclavable
242435	Coupler, Fixed Focal, 22 mm, C-Mount, Autoclavable
242436	Coupler, Zoom, C-Mount, Autoclavable
242386	Light Guide, OLYM-ACMI, 3.5 mm X 3.0 m
242389	Light Guide, UNIV, 3.5 mm X 3.0 m
242387	Light Guide, OLYM-ACMI, 5.0 mm X 3.0 m
242390	Light Guide, UNIV, 5.0 mm X 3.0 m
242358	Adaptor ACMI Light Source Universal Tip
242359	Adaptor Olympus Light Source Universal Tip
242360	Adaptor Storz Light Source Universal Tip
242361	Adaptor Wolf Light Source Universal Tip
242362	Adaptor Storz/Olympus Instrument Universal Tip
242363	Adaptor Wolf Instrument Universal Tip
242364	Adaptor ACMI Instrument Universal Tip

All couplers listed above are compatible with the DePuy Synthes PUREVUE™ Visualization System.

Sterilization Trays

CODE	DESCRIPTION
242335	Visualization Tray (Square: 11.3"L x 11.3"W x 2.2"T)
242336	Visualization Tray Lid (Square)
242337	Visualization Scope Tray (Rectangular: 19.8"L x 9.4"W x 2.2"T)
242338*	Visualization Tray Lid (Rectangular)
242339	Visualization Tray (Rectangular: 19.8"L x 9.4"W x 2.2"T)
212331	Laparoscope Tray
242332	Laparoscope Tray Lid
242307**	Camera Head System Tray, Stackable
242309**	Endoscope Tray, Stackable
242313***	Camera and Endoscope Tray Lid, Stackable



*Tray Lid 242338 is compatible with CODES 242337 and 242339

**Tray Lid 242313 is compatible with CODES 242307 and 242309

***Trays 242307 & 242309 are interlocking with each other, only requiring one lid if interlocking feature is used.



Arthroscopic Equipment

↳ Visualization

Arthroscopes

Scope Type	Arthroscope		Sheath		Obturator	Trocar	
	CODE	DESCRIPTION	Single Stopcock	Dual Stopcock			
Mitek Lock Push-Button Style	Shoulder, Knee & Hip	242017	HD Arthroscope, 4.0 mm, 0 deg, 167 mm	242033	242024	242258	242259
		242018	HD Arthroscope, 4.0 mm, 30 deg, 167 mm	242034	242025	242258	242259
		242019	HD Arthroscope, 4.0 mm, 45 deg, 167 mm	242035	242026	242258	242259
		242021	HD Arthroscope, 4.0 mm, 70 deg, 167 mm	242036	242027	242258	242259
		242022	HD Arthroscope, 4.0 mm, 30 deg, 270 mm	242037	242028	242254	242255
		242023	HD Arthroscope, 4.0 mm, 70 deg, 270 mm	242038	242029	242254	242255
		242136	4K Arthroscope, 0°, 4 mm x 167 mm	242031	242032	242143	242139
		242137	4K Arthroscope, 30°, 4 mm x 167 mm	242031	242032	242143	242139
		242138	4K Arthroscope, 70°, 4 mm x 167 mm	242031	242032	242143	242139
		242041	Hip Arthroscope, 30°, 3.5 mm x 205 mm	242061	242062	242063	242064
		242042	Hip Arthroscope, 70°, 3.5 mm x 205 mm	242061	242062	242063	242064
		242094	Arthroscope, 0°, 1.9 mm x 60 mm	242071	242072	242073	242074
		242095	Arthroscope, 30°, 1.9 mm x 60 mm	242071	242072	242073	242074
		242105	Arthroscope, 0°, 2.7 mm x 75 mm	242083	242084	242085	242086
		242107	Arthroscope, 30°, 2.7 mm x 75 mm	242083	242084	242085	242086
		242109	Arthroscope, 70°, 2.7 mm x 75 mm	242083	242084	242085	242086
C-Mount Arthroscopes with Mitek Lock	Shoulder, Knee, Hip & Small Joint	242267	4K C-Mount Arthroscope, 0°, 4 mm x 167 mm	242031	242032	242143	242139
		242268	4K C-Mount Arthroscope, 30°, 4 mm x 167 mm	242031	242032	242143	242139
		242269	4K C-Mount Arthroscope, 70°, 4 mm x 167 mm	242031	242032	242143	242139
		242271	C-Mount Hip Arthroscope, 30°, 3.5 mm x 205 mm	242061	242062	242063	242064
		242272	C-Mount Hip Arthroscope, 70°, 3.5 mm x 205 mm	242061	242062	242063	242064
		242273	C-Mount Arthroscope, 0°, 1.9 mm x 60 mm	242071	242072	242073	242074
		242274	C-Mount Arthroscope, 30°, 1.9 mm x 60 mm	242071	242072	242073	242074
		242275	C-Mount Arthroscope, 0°, 2.7 mm x 75 mm	242083	242084	242085	242086
		242276	C-Mount Arthroscope, 30°, 2.7 mm x 75 mm	242083	242084	242085	242086
		242277	C-Mount Arthroscope, 70°, 2.7 mm x 75 mm	242083	242084	242085	242086
242266	C-Mount Wrench	NA	NA	NA	NA		
Rotating Lock Compatible with Storz Style	Small Joint	242047	HD Arthroscope, 1.9 mm, 0 deg, 60 mm	242561	242523	242556	242555
		242048	HD Arthroscope, 1.9 mm, 30 deg, 60 mm	242562	242524	242556	242555
		242055	HD Arthroscope, 2.7 mm, 0 deg, 75 mm	242565	242528	242560	242559
		242056	HD Arthroscope, 2.7 mm, 30 deg, 75 mm	242566	242529	242560	242559
		242058	HD Arthroscope, 2.7 mm, 70 deg, 75 mm	242567	242540	242560	242559
	Shoulder, Knee & Hip	242005	HD Arthroscope, 4.0 mm, 0 deg, 175 mm	N/A	242547	242585	242584
		242006	HD Arthroscope, 4.0 mm, 30 deg, 175 mm	N/A	242548	242585	242584
		242008	HD Arthroscope, 4.0 mm, 45 deg, 175 mm	N/A	242550	242585	242584
		242009	HD Arthroscope, 4.0 mm, 70 deg, 175 mm	N/A	242551	242585	242584



Arthroscopic Equipment

↳ Visualization

Visualization Equipment

Insufflation and Accessories

CODE	DESCRIPTION
222701	Insufflator 50L w/Power Cord
222702	Tube Set for Insufflation, Single-Use (20 Pieces per each)
222703	Heated Tube Set for Insufflation, Single-Use (10 Pieces per each)
222705	DISS Adapter for Universal Gas Connector
222780	Power Cord US, 2.5 m
222755	Hose for CO2 Gas Bottle Pin, 1.5 m
222779	Hose for CO2 Central Gas DISS:DISS, 3 m (requires DISS Adaptor 222705)



Laparoscopes

CODE	DESCRIPTION
242632	HD Laparoscope, 5.5 mm, Eye Piece, 0 deg, 300 mm
242633	HD Laparoscope, 5.5 mm, Eye Piece, 30 deg, 300 mm
242634	HD Laparoscope, 5.5 mm, Eye Piece, 45 deg, 300 mm
242635	HD Laparoscope, 10 mm, Eye Piece, 0 deg, 330 mm
242636	HD Laparoscope, 10 mm, Eye Piece, 30 deg, 330 mm
242637	HD Laparoscope, 10 mm, Eye Piece, 45 deg, 330 mm
242638	HD Laparoscope, 5.5 mm, Eye Piece, 0 deg, 430 mm
242639	HD Laparoscope, 5.5 mm, Eye Piece, 30 deg, 430 mm
242640	HD Laparoscope, 5.5 mm, Eye Piece, 45 deg, 430 mm
242652	HD Laparoscope, 10 mm, Eye Piece, 0 deg, 450 mm
242653	HD Laparoscope, 10 mm, Eye Piece, 30 deg, 450 mm
242654	HD Laparoscope, 10 mm, Eye Piece, 45 deg, 450 mm



Insufflation Tubing

CODE	DESCRIPTION	QTY
242702	Standard Tube Set for Insufflation, Single-Use	20/PK
242703	Heated Tube Set for Insufflation, Single-Use	10/PK



Arthroscopic Equipment

↳ Fluid Management & Tissue Debridement

Pumps and Instruments

CODE	DESCRIPTION	QTY
284004	FMS VUE Fluid Management & Tissue Debridement System	1/each
283573	Foot Pedal (FMS VUE System)	1/each
283551	Remote Control (FMS VUE System)	1/each



Tubing

CODE	DESCRIPTION	QTY
284508	Inflow Tubing (FMS VUE System)	1/each
284610	Outflow Tubing without One-Way Valve	1/each
284649	Outflow Tubing with One-Way Valve	1/each
281103	Intermediary Tubing without One-Way Valve	1/each
281142	Intermediary Tubing with One-Way Valve	1/each



CODE	DESCRIPTION	QTY
242671	Medical Display, NDS Radiance Ultra 32" HD	1/each
242672	Medical Display, NDS Radiance Ultra 32" 4K/UHD	1/each
242394	Zerewire G2 Transmitter & Receiver kit	1/each
242395	Zerewire G2 Access Kit - Pair: Radiance Ultra 32"	1/each
242396	Zerewire G2 Receiver Only	1/each
242397	Zerewire G2 Transmitter Only	1/each
242398	ZeroWire G2 Access Kit - Single: Radiance Ultra 32"	1/each

Arthroscopic Equipment

↳ Fluid Management & Tissue Debridement














Shaver Handpiece

CODE	DESCRIPTION
283812	MICRO Tornado II Hand Controlled Shaver Handpiece
283712	MICRO Tornado II No-Button Shaver Handpiece



Sterile Single-Use Blades & Burrs

All blades and burrs are 13cm long and packaged 5 per box

BLADE/BURR COLOR RING	PRODUCT	CODE	DESCRIPTION	
COMBO	 ○	OMNICUT RESECTION BLADE A combination blade excellent for soft-tissue and bone debridement.	283345 4.2 mm (for use with MICRO Tornado Handpiece) 283346 5.2 mm (for use with MICRO Tornado Handpiece)	
	 ●	ULTRA-AGGRESSIVE BLADE PLUS Ultra-aggressive all-purpose blade ideal for thick-tissue removal.	283329 3.5 mm (for use with MICRO Tornado Handpiece) 283429 4.0 mm (for use with MICRO Tornado Handpiece) 283729 5.0 mm (for use with MICRO Tornado Handpiece)	
BLADES	 ●	AGGRESSIVE BLADE PLUS Aggressive meniscal trimming, plica, and synovia removal.	283319 3.5 mm (for use with MICRO Tornado Handpiece) 283419 4.0 mm (for use with MICRO Tornado Handpiece) 283719 5.0 mm (for use with MICRO Tornado Handpiece)	
	 ●	MULTI BLADE PLUS Great for meniscal trimming, joint debridement, and general tissue removal.	283439 4.0 mm (for use with MICRO Tornado Handpiece) 283739 5.0 mm (for use with MICRO Tornado Handpiece)	
	 ●	FULL-RADIUS BLADE PLUS Excellent for general debridement.	283309 3.5 mm (for use with MICRO Tornado Handpiece) 283409 4.0 mm (for use with MICRO Tornado Handpiece) 283709 5.0 mm (for use with MICRO Tornado Handpiece)	
	 ●	CURVED ULTRA-AGGRESSIVE BLADE PLUS Ultra-aggressive all-purpose blade ideal for thick-tissue removal.	283422 4.2 mm (for use with MICRO Tornado Handpiece)	
	 ●	CURVED AGGRESSIVE BLADE PLUS Aggressive meniscal trimming, plica, and synovia removal.	283412 4.2 mm (for use with MICRO Tornado Handpiece)	
	 ●	CURVED FULL-RADIUS BLADE PLUS Excellent for general debridement.	283403 4.2 mm (for use with MICRO Tornado Handpiece)	
	 ●	CURVED MULTI BLADE PLUS Great for meniscal trimming, joint debridement, and general tissue removal	283432 4.2 mm (for use with MICRO Tornado Handpiece)	
	BURRS	 ●	BARREL TORNADO BURR PLUS Ultra-aggressive burr for bone resection during acromioplasty and notchplasty.	283489 4.0 mm (for use with MICRO Tornado Handpiece) 283889 5.0 mm (for use with MICRO Tornado Handpiece)
		 ●	BARREL BURR PLUS Ideal for resecting bone and cartilage during acromioplasty and notchplasty.	283469 4.0 mm (for use with MICRO Tornado Handpiece) 283869 5.0 mm (for use with MICRO Tornado Handpiece)
		 ●	ROUND TORNADO BURR PLUS Ultra-aggressive bone, cartilage, osteochondral debridement, and osteophyte resection.	283479 4.0 mm (for use with MICRO Tornado Handpiece) 283879 5.0 mm (for use with MICRO Tornado Handpiece)
 ●		ROUND BURR PLUS Bone cartilage, osteochondral debridement, and osteophyte resection.	283459 4.0 mm (for use with MICRO Tornado Handpiece) 283859 5.0 mm (for use with MICRO Tornado Handpiece)	







Arthroscopic Equipment

↳ Fluid Management & Tissue Debridement











Sterile Single-Use Mini Blades & Mini Burrs

All mini blades and mini burrs are 8 cm long and packaged 5 per box

BLADE/BURR COLOR RING		PRODUCT	CODE	DESCRIPTION
BLADES		ULTRA-AGGRESSIVE BLADE PLUS Ultra-aggressive all-purpose blade ideal for thick-tissue removal.	283229 283629	2.7 mm Mini (for use with MICRO Tornado Handpiece) 3.5 mm Mini (for use with MICRO Tornado Handpiece)
		AGGRESSIVE BLADE PLUS Aggressive meniscal trimming, plica, and synovia removal.	283619	3.5 mm Mini (for use with MICRO Tornado Handpiece)
BURRS		FULL-RADIUS BLADE PLUS Excellent for general debridement.	283209 283609	2.7 mm Mini (for use with MICRO Tornado Handpiece) 3.5 mm Mini (for use with MICRO Tornado Handpiece)
		MINI ROUND TORNADO BURR PLUS Ultra-aggressive bone, cartilage, osteochondral debridement, and osteophyte resection.	283055 283255	2.0 mm Mini (for use with MICRO Tornado Handpiece) 3.0 mm Mini (for use with MICRO Tornado Handpiece)

Sterile Single-Use XL Blades & XL Burrs

All XL blades and XL burrs are 19 cm long and packaged 5 per box

BLADE/BURR COLOR RING		PRODUCT	CODE	DESCRIPTION
BLADES		ULTRA-AGGRESSIVE BLADE PLUS Ultra-aggressive all-purpose blade ideal for thick-tissue removal.	283428	4.2 mm XL (for use with MICRO Tornado Handpiece)
		AGGRESSIVE BLADE PLUS Aggressive meniscal trimming, plica, and synovia removal.	283418	4.2 mm XL (for use with MICRO Tornado Handpiece)
		MULTI BLADE PLUS Great for meniscal trimming, joint debridement, and general tissue removal.	283431	4.2 mm XL (for use with MICRO Tornado Handpiece)
		FULL-RADIUS BLADE PLUS Excellent for general debridement.	283408	4.2 mm XL (for use with MICRO Tornado Handpiece)
		CURVED ULTRA-AGGRESSIVE BLADE PLUS Ultra-aggressive all-purpose blade ideal for thick-tissue removal.	283421	4.2 mm XL (for use with MICRO Tornado Handpiece)
		CURVED AGGRESSIVE BLADE PLUS Aggressive meniscal trimming, plica, and synovia removal	283411	4.2 mm XL (for use with MICRO Tornado Handpiece)
		CURVED FULL-RADIUS BLADE PLUS Excellent for general debridement.	283404	4.2 mm XL (for use with MICRO Tornado Handpiece)
		CURVED MULTI BLADE PLUS Great for meniscal trimming, joint debridement, and general tissue removal.	283433	4.2 mm XL (for use with MICRO Tornado Handpiece)
BURRS		XL ROUND TORNADO BURR PLUS Ultra-aggressive bone, cartilage, osteochondral debridement, and osteophyte resection.	283478 283878	4.0 mm XL (for use with MICRO Tornado Handpiece) 5.0 mm XL (for use with MICRO Tornado Handpiece)
		XL ROUND BURR PLUS Bone cartilage, osteochondral debridement, and osteophyte resection.	283458 283858	4.0 mm XL (for use with MICRO Tornado Handpiece) 5.0 mm XL (for use with MICRO Tornado Handpiece)



Arthroscopic Equipment

↳ Radiofrequency

RF System and Accessories



CODE	DESCRIPTION
225024	VAPRVUE System Generator
227214	VAPRVUE System Wireless Footswitch
225023	VAPRVUE System Wired Footswitch
225002	VAPR Handpiece

Arthroscopic Equipment

↳ Radiofrequency

Sterile Single-Use Mini Blades & Mini Burrs

All mini blades and mini burrs are 8 cm long and packaged 5 per box

Sterile Single-Use Electrodes

VAPR ELECTRODES	PRODUCT	CODE (HC): Hand Controls (2P): Two-piece Electrode	MODE	DEFAULT				MAXIMUM			
				ABLATE/CUT		DES/COAG		ABLATE/CUT		DES/COAG	
				Basic	Advanced (watts)	Basic	Advanced (watts)	Basic	Advanced (watts)	Basic	Advance (watts)
	TRIPOLAR 90	225028 (HC)	CP	4	220	4	100	10	380	6	160
	COOLPULSE Curve	225031	CP	4	200	3	80	5	260	4	100
	COOLPULSE Curve - XL	225032	CP	4	200	3	80	5	260	4	100
	ARCTIC Dome (Articulating)	225029	CP	4	200	3	80	5	260	4	100
	ARCTIC Chisel (Articulating)	225030									
	COOLPULSE 90	228146 228147 (HC)	CP	4	220	3	80	10	380	6	160
	Premiere 90	227204	V3	7	240	5	120	8	260	6	160
	S90	225370 228370 (HC)									
	Premiere 50	227504 228504 (HC)	V2	6	160	2	45	7	200	3	60
	S50	227355	V2	3	60	2	45	4	90	3	60



SUCTION ELECTRODES




Arthroscopic Equipment

↳ Radiofrequency

Sterile Single-Use Electrodes

VAPR ELECTRODES	PRODUCT	CODE (HC): Hand Controls (2P): Two-piece Electrode	MODE	DEFAULT				MAXIMUM			
				ABLATE/CUT		DES/COAG		ABLATE/CUT		DES/COAG	
				Basic	Advanced (watts)	Basic	Advanced (watts)	Basic	Advanced (watts)	Basic	Advanced (watts)
	Straight End Effect	225303	V2								
	Flexible End Effect	225314		4	90	3	60	5	130	4	90
	Angled End Effect	225304									
	2.3 mm Side Effect	225201	V2	3	60	2	45	4	90	3	60

Shoulder



Shoulder

↳ Rotator Cuff

↳ Suture Anchors

HEALIX ADVANCE™ Anchors with DYNATAPE™ Sutures



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
280723	4.5 mm HEALIX ADVANCE BR Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280735	4.5 mm HEALIX ADVANCE BR Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each
280724	5.5 mm HEALIX ADVANCE BR Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280736	5.5 mm HEALIX ADVANCE BR Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each
280725	6.5 mm HEALIX ADVANCE BR Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280737	6.5 mm HEALIX ADVANCE BR Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each

PEEK

CODE	DESCRIPTION	QTY
280717	4.5 mm HEALIX ADVANCE PEEK Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280729	4.5 mm HEALIX ADVANCE PEEK Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each
280718	5.5 mm HEALIX ADVANCE PEEK Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280730	5.5 mm HEALIX ADVANCE PEEK Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each
280719	6.5 mm HEALIX ADVANCE PEEK Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280731	6.5 mm HEALIX ADVANCE PEEK Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each



Shoulder

↳ Rotator Cuff

↳ Suture Anchors

HEALIX ADVANCE™ Anchors with DYNACORD™ Sutures



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
222001	4.5 mm HEALIX ADVANCE BR Anchor with #2 DYNACORD Suture	1/each
222004	5.5 mm HEALIX ADVANCE BR Anchor with #2 DYNACORD Suture	1/each
222011	6.5 mm HEALIX ADVANCE BR Anchor with #2 DYNACORD Suture	1/each
222000	4.5 mm HEALIX ADVANCE BR Anchor with #2 DYNACORD Suture (with needles)	1/each
222006	5.5 mm HEALIX ADVANCE BR Anchor with #2 DYNACORD Suture (with needles)	1/each
222015	6.5 mm HEALIX ADVANCE BR Anchor with #2 DYNACORD Suture (with needles)	1/each
222002	4.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 DYNACORD Suture	1/each
222005	5.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 DYNACORD Suture	1/each
222013	6.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 DYNACORD Suture	1/each

PEEK

CODE	DESCRIPTION	QTY
222028	4.5 mm HEALIX ADVANCE PEEK Anchor with #2 DYNACORD Suture	1/each
222032	5.5 mm HEALIX ADVANCE PEEK Anchor with #2 DYNACORD Suture	1/each
222036	6.5 mm HEALIX ADVANCE PEEK Anchor with #2 DYNACORD Suture	1/each
222027	4.5 mm HEALIX ADVANCE PEEK Anchor with #2 DYNACORD Suture (with needles)	1/each
222034	5.5 mm HEALIX ADVANCE PEEK Anchor with #2 DYNACORD Suture (with needles)	1/each
222038	6.5 mm HEALIX ADVANCE PEEK Anchor with #2 DYNACORD Suture (with needles)	1/each
222029	4.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 DYNACORD Suture	1/each
222033	5.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 DYNACORD Suture	1/each
222037	6.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 DYNACORD Suture	1/each



Shoulder

↳ Rotator Cuff

↳ Suture Anchors

HEALIX ADVANCE™ Anchors with ORTHOCORD™ Sutures

BIOCRYL RAPIDE™ Biocomposite Material



CODE	DESCRIPTION	QTY
222295	4.5 mm HEALIX ADVANCE BR Anchor with #2 ORTHOCORD Suture	1/each
222297	5.5 mm HEALIX ADVANCE BR Anchor with #2 ORTHOCORD Suture	1/each
222300	6.5 mm HEALIX ADVANCE BR Anchor with #2 ORTHOCORD Suture	1/each
222319	4.5 mm HEALIX ADVANCE BR Anchor Pack with #2 ORTHOCORD Suture	5/pack
222323	5.5 mm HEALIX ADVANCE BR Anchor Pack with #2 ORTHOCORD Suture	5/pack
222299	5.5 mm HEALIX ADVANCE BR Anchor with #2 ORTHOCORD Suture (with needles)	1/each
222302	6.5 mm HEALIX ADVANCE BR Anchor with #2 ORTHOCORD Suture (with needles)	1/each
222296	4.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 ORTHOCORD Suture	1/each
222298	5.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 ORTHOCORD Suture	1/each
222301	6.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 ORTHOCORD Suture	1/each
222320	4.5 mm HEALIX ADVANCE BR Anchor Pack with Triple Strand #2 ORTHOCORD Suture	5/pack
222324	5.5 mm HEALIX ADVANCE BR Anchor Pack with Triple Strand #2 ORTHOCORD Suture	5/pack

PEEK

CODE	DESCRIPTION	QTY
222303	4.5 mm HEALIX ADVANCE PEEK Anchor with #2 ORTHOCORD Suture	1/each
222305	5.5 mm HEALIX ADVANCE PEEK Anchor with #2 ORTHOCORD Suture	1/each
222309	6.5 mm HEALIX ADVANCE PEEK Anchor with #2 ORTHOCORD Suture	1/each
222321	4.5 mm HEALIX ADVANCE PEEK Anchor Pack with #2 ORTHOCORD Suture	5/pack
222325	5.5 mm HEALIX ADVANCE PEEK Anchor Pack with #2 ORTHOCORD Suture	5/pack
222308	5.5 mm HEALIX ADVANCE PEEK Anchor with #2 ORTHOCORD Suture (with needles)	1/each
222313	6.5 mm HEALIX ADVANCE PEEK Anchor with #2 ORTHOCORD Suture (with needles)	1/each
222304	4.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 ORTHOCORD Suture	1/each
222307	5.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 ORTHOCORD Suture	1/each
222311	6.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 ORTHOCORD Suture	1/each
222322	4.5 mm HEALIX ADVANCE PEEK Anchor Pack with Triple Strand #2 ORTHOCORD Suture	5/pack
222318	5.5 mm HEALIX ADVANCE PEEK Anchor Pack with Triple Strand #2 ORTHOCORD Suture	5/pack



Shoulder

↳ Rotator Cuff

↳ Suture Anchors

HEALIX ADVANCE™ Anchors with PERMATAPE™ Sutures



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
222359	4.5 mm HEALIX ADVANCE BR Anchor with PERMATAPE Suture (Blue)	1/each
222360	4.5 mm HEALIX ADVANCE BR Anchor with PERMATAPE Suture (Blue/White)	1/each
222361	5.5 mm HEALIX ADVANCE BR Anchor with PERMATAPE Suture (Blue)	1/each
222362	5.5 mm HEALIX ADVANCE BR Anchor with PERMATAPE Suture (Blue/White)	1/each
222376	6.5 mm HEALIX ADVANCE BR Anchor with PERMATAPE Suture (Blue)	1/each
222377	6.5 mm HEALIX ADVANCE BR Anchor with PERMATAPE Suture (Blue/White)	1/each

PEEK

CODE	DESCRIPTION	QTY
222363	4.5 mm HEALIX ADVANCE PEEK Anchor with PERMATAPE Suture (Blue)	1/each
222364	4.5 mm HEALIX ADVANCE PEEK Anchor with PERMATAPE Suture (Blue/White)	1/each
222365	5.5 mm HEALIX ADVANCE PEEK Anchor with PERMATAPE Suture (Blue)	1/each
222366	5.5 mm HEALIX ADVANCE PEEK Anchor with PERMATAPE Suture (Blue/White)	1/each
222378	6.5 mm HEALIX ADVANCE PEEK Anchor with PERMATAPE Suture (Blue)	1/each
222379	6.5 mm HEALIX ADVANCE PEEK Anchor with PERMATAPE Suture (Blue/White)	1/each

Shoulder

↳ Rotator Cuff

↳ Suture Anchors

HEALIX ADVANCE™ Anchors with PERMACORD™ Sutures



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
223128	4.5 mm HEALIX ADVANCE BR Anchor with #2 PERMACORD Suture	1/each
223130	5.5 mm HEALIX ADVANCE BR Anchor with #2 PERMACORD Suture	1/each
223132	6.5 mm HEALIX ADVANCE BR Anchor with #2 PERMACORD Suture	1/each
223129	4.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 PERMACORD Suture	1/each
223131	5.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 PERMACORD Suture	1/each
223133	6.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 PERMACORD Suture	1/each

PEEK

CODE	DESCRIPTION	QTY
223134	4.5 mm HEALIX ADVANCE PEEK Anchor with #2 PERMACORD Suture	1/each
223136	5.5 mm HEALIX ADVANCE PEEK Anchor with #2 PERMACORD Suture	1/each
223138	6.5 mm HEALIX ADVANCE PEEK Anchor with #2 PERMACORD Suture	1/each
223135	4.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 PERMACORD Suture	1/each
223137	5.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 PERMACORD Suture	1/each
223139	6.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 PERMACORD Suture	1/each



Shoulder

↳ Rotator Cuff

↳ Suture Anchors

HEALIX ADVANCE™ Anchor Instruments



HEALIX ADVANCE™ AWL



HEALIX ADVANCE™ AWL-TAP

CODE	DESCRIPTION	QTY
222314	HEALIX ADVANCE Anchor Awl, 4.5 mm	1/each
222315	HEALIX ADVANCE Anchor Awl, 5.5 mm/6.5 mm	1/each
222316	HEALIX ADVANCE Anchor Cortical Awl/Tap, 4.5 mm	1/each
222317	HEALIX Anchor Cortical Awl/Tap, 5.5 mm/6.5 mm	1/each
222329	HEALIX ADVANCE Anchor Full Tap, 4.5 mm	1/each
222328	HEALIX ADVANCE Anchor Full Tap, 5.5 mm/6.5 mm	1/each
222327	HEALIX ADVANCE Anchor Drill Bit, 3.5 mm (DO NOT USE with HEALIX ADVANCE Drill Guide - 222375)	1/each
222326	HEALIX ADVANCE Anchor Drill Bit, 4 mm (DO NOT USE with HEALIX ADVANCE Drill Guide - 222375)	1/each
222358	HEALIX ADVANCE Hard Bone Kit (3.5 mm Drill Bit, Free Needle Loop, 4.5 mm Anchor Tap)	1/each
222375	HEALIX ADVANCE Open Drill Guide*	1/each

*The HEALIX ADVANCE™ Drill Guide can only be used with the 3.5 mm HEALIX ADVANCE™ drill bit of the Hard Bone Kit.

Shoulder

↳ Rotator Cuff

↳ Suture Anchors

HEALIX TI™ Anchors with DYNATAPE™ Suture



Titanium

CODE	DESCRIPTION	QTY
280807	4.5 mm HEALIX TI Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280813	4.5 mm HEALIX TI Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each
280808	5.5 mm HEALIX TI Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280814	5.5 mm HEALIX TI Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each
280809	6.5 mm HEALIX TI Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280815	6.5 mm HEALIX TI Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each

HEALIX TI™ Anchors with DYNACORD™ Suture



Titanium

CODE	DESCRIPTION	QTY
222051	4.5 mm HEALIX TI Anchor with #2 DYNACORD Suture	1/each
222054	5.5 mm HEALIX TI Anchor with #2 DYNACORD Suture	1/each
222058	6.5 mm HEALIX TI Anchor with #2 DYNACORD Suture	1/each
222055	5.5 mm HEALIX TI Anchor with #2 DYNACORD Suture (with needles)	1/each
222059	6.5 mm HEALIX TI Anchor with #2 DYNACORD Suture (with needles)	1/each
222052	4.5 mm HEALIX TI Anchor with Triple Strand #2 DYNACORD Suture	1/each
222056	5.5 mm HEALIX TI Anchor with Triple Strand #2 DYNACORD Suture	1/each
222061	6.5 mm HEALIX TI Anchor with Triple Strand #2 DYNACORD Suture	1/each



Shoulder

↳ Rotator Cuff

↳ Suture Anchors

HEALIX TI™ Anchors with ORTHOCORD™ Suture



Titanium

CODE	DESCRIPTION	QTY
222253	4.5 mm HEALIX TI Anchor with #2 ORTHOCORD Suture	1/each
222245	5.5 mm HEALIX TI Anchor with #2 ORTHOCORD Suture	1/each
222249	6.5 mm HEALIX TI Anchor with #2 ORTHOCORD Suture	1/each
222242	5.5 mm HEALIX TI Anchor with #2 ORTHOCORD Suture (with needles)	1/each
222248	6.5 mm HEALIX TI Anchor with #2 ORTHOCORD Suture (with needles)	1/each
222271	4.5 mm HEALIX TI Anchor with Triple Strand #2 ORTHOCORD Suture	1/each
222243	5.5 mm HEALIX TI Anchor with Triple Strand #2 ORTHOCORD Suture	1/each
222241	6.5 mm HEALIX TI Anchor with Triple Strand #2 ORTHOCORD Suture	1/each

HEALIX TI™ Anchors with PERMACORD™ Suture



Titanium

CODE	DESCRIPTION	QTY
223140	4.5 mm HEALIX TI Anchor with #2 PERMACORD Suture	1/each
223150	5.5 mm HEALIX TI Anchor with #2 PERMACORD Suture	1/each
223144	6.5 mm HEALIX TI Anchor with #2 PERMACORD Suture	1/each
223141	4.5 mm HEALIX TI Anchor with Triple Strand #2 PERMACORD Suture	1/each
223151	5.5 mm HEALIX TI Anchor with Triple Strand #2 PERMACORD Suture	1/each
223145	6.5 mm HEALIX TI Anchor with Triple Strand #2 PERMACORD Suture	1/each



Shoulder

↳ Rotator Cuff

↳ Suture Anchors

HEALIX TRANSTEND™ Implant System



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
222265	HEALIX TRANSTEND BR Anchor with ORTHOCORD Suture 3.4 mm	1/each
222261	HEALIX TRANSTEND BR Implant System (1 BR Anchor, PERCANNULA System)	1/each
222263	HEALIX TRANSTEND BR Implant System (2 BR Anchors, PERCANNULA System)	1/each
222338	HEALIX TRANSTEND BR Implant System (1 BR, Violet, 1 BR, Blue, PERCANNULA System)	1/each
222337	HEALIX TRANSTEND BR Dual Strand Implant System (1 BR DS, PERCANNULA System)	1/each

PEEK

CODE	DESCRIPTION	QTY
222276	HEALIX TRANSTEND PEEK Anchor with ORTHOCORD Suture 3.4 mm	1/each
222278	HEALIX TRANSTEND PEEK Implant System (1 PEEK Anchor, PERCANNULA System)	1/each
222280	HEALIX TRANSTEND PEEK Implant System (2 PEEK Anchors, PERCANNULA System)	1/each
222340	HEALIX TRANSTEND PEEK Implant System (1 PEEK, Violet, 1 PEEK, Blue, PERCANNULA)	1/each
222339	HEALIX TRANSTEND PEEK Dual Strand Implant System (1 PEEK DS, PERCANNULA System)	1/each

Titanium

CODE	DESCRIPTION	QTY
222260	HEALIX TRANSTEND TI Anchor with ORTHOCORD Suture 2.9 mm	1/each
222266	HEALIX TRANSTEND TI Implant System (1 TI Anchor, PERCANNULA System)	1/each
222268	HEALIX TRANSTEND TI Implant System (2 TI Anchors, PERCANNULA System)	1/each

Instruments

CODE	DESCRIPTION	QTY
222274	HEALIX PERCANNULA System 4.0 mm	1/each
222273	HEALIX TRANSTEND Anchor Awl/Tap, 3.4 mm (for BR and PEEK Anchors)	1/each
222281	HEALIX TRANSTEND Anchor Drill Bit, 2.2 mm (optional for TI Anchor)	1/each



Shoulder ↳ Rotator Cuff ↳ Suture Anchors

GII QUICKANCHOR™ Suture Anchors



GII QUICKANCHOR Suture Anchor with #2 ORTHOCORD Suture (with CP-2 needles)

Anchors

CODE	DESCRIPTION	QTY
222983	GII QUICKANCHOR Suture Anchor with #2 ORTHOCORD Suture (with CP-2 needles)	1/each



Instruments

CODE	DESCRIPTION	QTY
211093	GII Drill Bit, 2.4 mm	1/each
211502	Arthroscopic GII Drill Bit	1/each
213207	GII/SUPER/RC/GLS Slotted Drill Guide	1/each
214561	Arthroscopic Sharkmouth Drill Guide	1/each
214562	Arthroscopic Sawtooth Drill Guide	1/each
212094	GII Inserter	1/each
214542	Arthroscopic T-Handle Reamer	1/each
211093	GII Drill Bit, 2.4 mm	1/each



Super QUICKANCHOR™ Suture Anchors

Anchors

CODE	DESCRIPTION	QTY
222984	Super QUICKANCHOR™ Plus with ORTHOCORD™ Suture (CP-2 Needles)	1/each



Instruments

CODE	DESCRIPTION	QTY
211113	SUPER/RC Drill Bit (Open), 2.9 mm	1/each
213207	GII/SUPER/RC Slotted Drill Guide (Open)	1/each
211513	Arthroscopic SUPER/RC Drill Guide, Arthroscopic Drill, 2.9 mm	1/each



Shoulder

↳ Rotator Cuff

↳ All-Suture Anchors

VERSALOOPTM All-Suture Anchors



Anchors

CODE	DESCRIPTION	QTY
210165	1.5 mm VERSALOOPTM Anchor with a Size #2 Suture	1/each
210166	1.5 mm VERSALOOPTM Anchor with a 1.5 mm Tape	1/each
210167	1.8 mm VERSALOOPTM Anchor with 2 Size #2 Sutures	1/each
210168	1.8 mm VERSALOOPTM Anchor with 2 1.5 mm Tapes	1/each
210169	2.5 mm VERSALOOPTM Anchor with 2 Size #2 Sutures	1/each
210170	2.5 mm VERSALOOPTM Anchor with 2 1.5 mm Tapes	1/each
210171	2.5 mm VERSALOOPTM Anchor with 3 Size #2 Sutures	1/each

Shoulder

↳ Rotator Cuff

↳ All-Suture Anchors

VERSALOOPTM Instruments



1.8 mm Family of Instruments



1.5 mm 25° guide, 1.8 mm 12.5° guide & 2.5 mm Straight Guide

Obturator

CODE	DESCRIPTION	QTY
210205	VERSALOOPTM Obturator for 1.5 mm Anchor	1/each
210213	VERSALOOPTM Obturator for 1.8 mm Anchor	1/each
210214	VERSALOOPTM Obturator for 2.5 mm Anchor	1/each

Curved Instruments

CODE	DESCRIPTION	QTY
210195	VERSALOOPTM 12.5° Guide for 1.5 mm Anchor, Sawtooth	1/each
210196	VERSALOOPTM 12.5° Guide for 1.8 mm Anchor, Sawtooth	1/each
210197	VERSALOOPTM 25° Guide for 1.5 mm Anchor, Sawtooth	1/each
210198	VERSALOOPTM 25° Guide for 1.8 mm Anchor, Sawtooth	1/each
210199	VERSALOOPTM 1.5 mm Disposable Drill	1/each
210203	VERSALOOPTM 1.8 mm Disposable Drill	1/each

Straight Instruments

CODE	DESCRIPTION	QTY
210175	VERSALOOPTM Straight Guide for 1.5 mm Anchor, Fishmouth	1/each
210176	VERSALOOPTM Straight Guide for 1.8 mm Anchor, Fishmouth	1/each
210177	VERSALOOPTM Straight Guide for 2.5 mm Anchor, Fishmouth	1/each
210178	VERSALOOPTM Straight Guide for 1.5 mm Anchor, Sawtooth	1/each
210179	VERSALOOPTM Straight Guide for 1.8 mm Anchor, Sawtooth	1/each
210180	VERSALOOPTM Straight Guide for 2.5 mm Anchor, Sawtooth	1/each
210181	VERSALOOPTM 1.5 mm Reusable Straight Drill	1/each
210182	VERSALOOPTM 1.8 mm Reusable Straight Drill	1/each
210184	VERSALOOPTM 2.5 mm Reusable Straight Drill	1/each
210207	VERSALOOPTM 2.5 mm Reusable Awl	1/each



Shoulder ↳ Rotator Cuff ↳ Knotless Anchors

HEALIX ADVANCE™ Self-Punching (SP) Anchors



BIOCRYL™ Biocomposite Material

CODE	DESCRIPTION	QTY
222422	4.9 mm HEALIX ADVANCE Self-Punching Biocomposite Anchor	1/each
222424	5.5 mm HEALIX ADVANCE Self-Punching Biocomposite Anchor	1/each
222426	6.5 mm HEALIX ADVANCE Self-Punching Biocomposite Anchor	1/each

PEEK

CODE	DESCRIPTION	QTY
222423	4.9 mm HEALIX ADVANCE Self-Punching PEEK Anchor	1/each
222425	5.5 mm HEALIX ADVANCE Self-Punching PEEK Anchor	1/each
222427	6.5 mm HEALIX ADVANCE Self-Punching PEEK Anchor	1/each

HEALIX ADVANCE™ Knotless Anchors



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
222581	4.75 mm HEALIX ADVANCE Knotless BR Anchor (Double Kite)	1/each
222582	5.5 mm HEALIX ADVANCE Knotless BR Anchor (Double Kite)	1/each
222583	6.5 mm HEALIX ADVANCE Knotless BR Anchor (Double Kite)	1/each

PEEK

CODE	DESCRIPTION	QTY
222584	4.75 mm HEALIX ADVANCE Knotless PEEK Anchor (Double Kite)	1/each
222585	5.5 mm HEALIX ADVANCE Knotless PEEK Anchor (Double Kite)	1/each
222586	6.5 mm HEALIX ADVANCE Knotless PEEK Anchor (Double Kite)	1/each



Shoulder

↳ Rotator Cuff

↳ Knotless Anchors

VERSALOK™ Anchors



PEEK

CODE	DESCRIPTION	QTY
210818	VERSALOK PEEK Anchor (PEEK Sleeve) with #2 ORTHOCORD Suture	1/each

Titanium

CODE	DESCRIPTION	QTY
210808	VERSALOK Titanium Anchor (PEEK Sleeve) with #2 ORTHOCORD Suture	1/each



Instruments

CODE	DESCRIPTION	QTY
214710	VERSALOK Anchor Deployment Gun	1/each
214711	VERSALOK Anchor Awl	1/each
215440	VERSALOK Caddy	1/each

Shoulder

↳ Proximal Biceps

↳ Anchors & Screws

TIGHT-N™ Tendon Docking Anchors



Anchors

CODE	DESCRIPTION	QTY
208881	TIGHT-N Tendon Anchor, Small, 5.5 mm	1/each
208882	TIGHT-N Tendon Anchor, Medium, 7.0 mm	1/each
208883	TIGHT-N Tendon Anchor, Large, 8.5 mm	1/each

Instruments

CODE	DESCRIPTION	QTY
208886	TIGHT-N Tendon Sizer and Biceps Guide	1/each
208887	TIGHT-N Small (5.5 mm) Cannulated Reamer	1/each
208888	TIGHT-N Med (7.0 mm) Cannulated Reamer	1/each
208889	TIGHT-N Large (8.5 mm) Cannulated Reamer	1/each
208891	TIGHT-N Small (5.5 mm) Self-Centering Reamer	1/each
208892	TIGHT-N Med (7.0 mm) Self-Centering Reamer	1/each
208893	TIGHT-N Large (8.5 mm) Self-Centering Reamer	1/each

MILAGRO™ Interference Screws

BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
231814	MILAGRO BR Screw, 5 mm x 12 mm	1/each
231803	MILAGRO BR Screw, 5 mm x 23 mm	1/each
231804	MILAGRO BR Screw, 5 mm x 30 mm	1/each
231813	MILAGRO BR Screw, 6 mm x 12 mm	1/each
231807	MILAGRO BR Screw, 6 mm x 23 mm	1/each
231809	MILAGRO BR Screw, 6 mm x 30 mm	1/each
231811	MILAGRO BR Screw, 7 mm x 15 mm	1/each
231812	MILAGRO BR Screw, 8 mm x 15 mm	1/each




Shoulder

↳ Proximal Biceps

↳ Anchors & Screws

MILAGRO™ Instruments

CODE	DESCRIPTION	QTY
219974	Tendon Fork, Large (5.5 mm wide)	1/each
219621	Tendon Fork, Small (4 mm wide)	1/each
219448	Modular Driver, 12 mm	1/each
219447	Modular Driver, 15 mm	1/each
219973	Modular Driver, 23 mm	1/each
219446	Modular Driver, 5 mm x 23 mm and 6 mm x 23 mm Screws	1/each
219438	Modular Driver, 5 mm x 30 mm and 6 mm x 30 mm Screws	1/each
219445	Tendon Sizer, 4 mm to 9 mm	1/each
219971	Guidewire Cap	1/each
254729	Calibrated Drill Pin, 2.4 mm x 14" (35.6 cm)	1/each
232426	Sterile Tenodesis Kit	1/each
219440	Acorn Reamer, 5 mm	1/each
219441	Acorn Reamer, 5.5 mm	1/each
219442	Fully Fluted Reamer, 5 mm	1/each
219443	Fully Fluted Reamer, 5.5 mm	1/each
214015	Guidewire, Sterile, Stainless Steel, 1.1 mm x 14" (35.6 cm)	1/each
214152	Guidewire, Sterile, Stainless Steel, 1.6 mm x 14" (35.6 cm)	1/each
211321	Trocar Tip Guidewire, Non-Sterile, Stainless Steel, 1.6 mm x 14" (35.6 cm)	1/each



1.6 mm GUIDEWIRE

Note: If doing a Biceps Tenodesis with 1.6 mm guidewire, surgeon must use 23 mm tenodesis driver (219973).
The 1.1 mm guidewire must be used when using a MILAGRO™ Advance screw.

RIGIDLOOP™ Titanium Buttons

The RIGIDLOOP™ Titanium button is preloaded with a white UHMWPE (ultra high molecular weight polyethylene) leading utility suture and a green #5 ETHIBOND trailing suture.

CODE	DESCRIPTION	QTY
232008	RIGIDLOOP Titanium Button	1/each
232034	RIGIDLOOP XL Button	1/each



Shoulder ↳ Labrum ↳ Suture Anchors

GRYPHON™ PROKNOT™ Suture Anchors

BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
210820	GRYPHON PROKNOT BR Anchor with PERMACORD Suture	1/each



PEEK

CODE	DESCRIPTION	QTY
210821	GRYPHON PROKNOT PEEK Anchor with PERMACORD Suture	1/each

Instruments

CODE	DESCRIPTION	QTY
210824	GRYPHON PROKNOT Anchor Slotted Sawtooth Guide	1/each
210826	GRYPHON PROKNOT Anchor Slotted Fishmouth Guide	1/each



GRYPHON™ Anchors with DYNACORD™ Suture



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
210120	GRYPHON-P BR Anchor with #2 DYNACORD Suture	1/each
210121	GRYPHON-P BR Anchor with Dual Strand #2 DYNACORD Suture	1/each

PEEK

CODE	DESCRIPTION	QTY
210122	GRYPHON-P PEEK Anchor with #2 DYNACORD Suture	1/each
210123	GRYPHON-P PEEK Anchor with Dual Strand #2 DYNACORD Suture	1/each



Shoulder ↳ Labrum ↳ Suture Anchors

GRYPHON™ Anchors with ORTHOCORD™ Suture



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
210814	GRYPHON-P BR Anchor with #2 ORTHOCORD Suture	1/each
210813	GRYPHON-P BR Anchor with Dual Strand #2 ORTHOCORD Suture	1/each
210812	GRYPHON-T BR Anchor with #2 ORTHOCORD Suture	1/each
210811	GRYPHON-T BR Anchor with Dual Strand #2 ORTHOCORD Suture	1/each

PEEK

CODE	DESCRIPTION	QTY
210031	GRYPHON-P PEEK Anchor with #2 ORTHOCORD Suture	1/each
210030	GRYPHON-P PEEK Anchor with Dual Strand #2 ORTHOCORD Suture	1/each

GRYPHON™ Anchors with PERMACORD™ Suture



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
223146	GRYPHON-P BR Anchor with #2 PERMACORD Suture	1/each
223147	GRYPHON-P BR Anchor with Dual Strand #2 PERMACORD Suture	1/each

PEEK

CODE	DESCRIPTION	QTY
223148	GRYPHON-P PEEK Anchor with #2 PERMACORD Suture	1/each
223149	GRYPHON-P PEEK Anchor with Dual Strand #2 PERMACORD Suture	1/each



Shoulder

↳ Labrum

↳ Suture Anchors

GRYPHON™ Anchor Instruments

CODE	DESCRIPTION	QTY
211035	GRYPHON Anchor Tap (for use with Threaded Anchor only)	1/each
211036	GRYPHON Anchor Drill Bit, 2.4 mm	1/each
210824	GRYPHON PROKNOT Anchor Slotted Sawtooth Guide	1/each
211039	GRYPHON Anchor Sawtooth Guide	1/each
210826	GRYPHON PROKNOT Anchor Slotted Fishmouth Guide	1/each
211037	GRYPHON Anchor Fishmouth Guide	1/each
211038	GRYPHON Anchor Sharp Obturator	1/each
211041	GRYPHON Anchor Blunt Obturator	1/each
210835	Knot Manipulator	1/each
211601	GRYPHON Anchor Hard Bone Drill Bit, 2.5 mm	1/each



LUPINE™ Anchors

BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
210711	LUPINE BR Anchor with ORTHOCORD Suture	1/each
210712	LUPINE BR Anchor with Dual Strand ORTHOCORD Suture	1/each



Instruments

CODE	DESCRIPTION	QTY
211033	LUPINE Anchor Drill Bit	1/each
211042	LUPINE Anchor Spiral Fluted Drill Bit	1/each
214639	LUPINE Anchor Guide Obturator	1/each
213816	LUPINE Anchor Dovetail Guide	1/each
213817	LUPINE Anchor Fishmouth Guide	1/each
213818	LUPINE Anchor Hybrid Fishmouth Guide	1/each
213819	LUPINE Anchor Sawtooth Drill Guide	1/each



Dovetail



Fishmouth



Hybrid Fishmouth



Sawtooth



Shoulder ↳ Labrum ↳ All-Suture Anchors

VERSALOOPTM All-Suture Anchors

Anchors

CODE	DESCRIPTION	QTY
210165	1.5 mm VERSALOOPTM Anchor with a Size #2 Suture	1/each
210166	1.5 mm VERSALOOPTM Anchor with a 1.5 mm Tape	1/each
210167	1.8 mm VERSALOOPTM Anchor with 2 Size #2 Sutures	1/each
210168	1.8 mm VERSALOOPTM Anchor with 2 1.5 mm Tapes	1/each



Obturator

CODE	DESCRIPTION	QTY
210205	VERSALOOPTM Obturator for 1.5 mm Anchor	1/each
210213	VERSALOOPTM Obturator for 1.8 mm Anchor	1/each

Curved Instruments

CODE	DESCRIPTION	QTY
210195	VERSALOOPTM 12.5° Guide for 1.5 mm Anchor, Sawtooth	1/each
210196	VERSALOOPTM 12.5° Guide for 1.8 mm Anchor, Sawtooth	1/each
210197	VERSALOOPTM 25° Guide for 1.5 mm Anchor, Sawtooth	1/each
210198	VERSALOOPTM 25° Guide for 1.8 mm Anchor, Sawtooth	1/each
210199	VERSALOOPTM 1.5 mm Disposable Drill	1/each
210203	VERSALOOPTM 1.8 mm Disposable Drill	1/each



1.8 mm Family of Instruments

Straight Instruments

CODE	DESCRIPTION	QTY
210175	VERSALOOPTM Straight Guide for 1.5 mm Anchor, Fishmouth	1/each
210176	VERSALOOPTM Straight Guide for 1.8 mm Anchor, Fishmouth	1/each
210178	VERSALOOPTM Straight Guide for 1.5 mm Anchor, Sawtooth	1/each
210179	VERSALOOPTM Straight Guide for 1.8 mm Anchor, Sawtooth	1/each
210181	VERSALOOPTM 1.5 mm Reusable Straight Drill	1/each
210182	VERSALOOPTM 1.8 mm Reusable Straight Drill	1/each



1.5 mm 25° guide, 1.8 mm 12.5° guide & 2.5 mm Straight Guide



Shoulder ↳ Glenoid ↳ Latarjet Instability Systems

LATARJET EXPERIENCE Instability Shoulder System



LATARJET EXPERIENCE Screw and Top Hat

Implants

CODE	DESCRIPTION	QTY
288222	Sterile Latarjet Screw, 28 mm	1/each
288223	Sterile Latarjet Screw, 30 mm	1/each
288224	Sterile Latarjet Screw, 32 mm	1/each
288225	Sterile Latarjet Screw, 34 mm	1/each
288226	Sterile Latarjet Screw, 36 mm	1/each
288227	Sterile Latarjet Screw, 38 mm	1/each
288228	Sterile Latarjet Screw, 40 mm	1/each
288229	Sterile Latarjet Screw, 42 mm	1/each
288230	Sterile Latarjet Screw, 44 mm	1/each
288231	Sterile Top Hats (contains 2 each)	2/each
288200	Inline Coracoid Drill Guide	1/each

Note: screw size is 4.5 mm



LATARJET Procedure Kit



BRISTOW Procedure Kit

Disposable Kits

CODE	DESCRIPTION	QTY
288238	LATARJET Disposable Kit	1/each
288239	BRISTOW Disposable Kit	1/each



Shoulder ↳ Glenoid ↳ Latarjet Instability Systems

LATARJET EXPERIENCE Instability Shoulder System

Instruments

CODE	DESCRIPTION	QTY
		
288101	LATARJET Small Switch Stick, 4 mm (with ruler)	1/each
288240	LATARJET Large Switch Stick, 4.75 mm	1/each
		
288200	Inline Coracoid Drill Guide	1/each
		
288201	Coracoid Drill, 3.2 mm	1/each
		
288202	Coracoid Top Hat Drill, 3.2 mm–4.8 mm	1/each
		
288203	Top Hat Tap	1/each
		
288204	Sharp Curved Osteotome	1/each
		
288205	Subscapularis Channeler	1/each
		
288206	Coracoid Offset Pin	1/each
288207	Glenoid K-Wire, 1.5 mm x 380 mm	1/each
		
288208	Glenoid Offset Pin	1/each
		
288209	Coracoid Screw, 3.5 mm	1/each

CODE	DESCRIPTION	QTY
		
288210	Reference Tube	1/each
		
288211	LATARJET Combo Screwdriver	1/each
		
288212	LATARJET Humeral Head Retractor	1/each
		
288213	LATARJET Hohmann Retractor	1/each
		
288214	Glenoid Double Prong Lever	1/each
		
288215	Subscapularis Spreader	1/each
		
288221	Glenoid Guide	1/each
		
285300	LATARJET Glenoid Long Drill Bit, 3.2 mm	1/each
		
288242	LATARJET Small Probe	1/each
288217	LATARJET Case Lid	1/each
288218	LATARJET Tray	1/each
288219	LATARJET Case	1/each



Shoulder ↳ Glenoid ↳ Latarjet Instability Systems

LATARJET EXPERIENCE Instability Shoulder System

Instruments






CODE	DESCRIPTION	QTY
285365	LATARJET Coracoid Drill Guide	1/each
285260	LATARJET Coracoid Short K-Wire, 150 mm	1/each
285270	LATARJET Curved Osteotome	1/each
285300	LATARJET Glenoid Drill Bit, 3.2 mm	1/each
211607	U-Shaped Cannula, Medium	1/each
211610	U-Shaped Cannula Obturator, Medium	1/each
211608	U-Shaped Cannula, Large	1/each
211611	U-Shaped Cannula Obturator, Large	1/each








Shoulder


↳ Sutures & Suture Passers

Needle Guide


27 mm	13 mm	26 mm	27 mm	13 mm
				
CP-2 Needle (Reverse Cutting)	P-3 Needle (Reverse Cutting)	SH Needle (Taper Point)	CT-2 Needle (Taper Point)	C-1 Needle (Taper Point)

3/8 circle 17 mm	1/2 circle 17 mm	1/2 circle 17 mm	1/2 circle 22 mm	1/2 circle 17 mm
				
V-4 Needle (TAPERCUT™)	V-5 Needle (TAPERCUT™)	OS-6 36 mm (Reverse Cutting)	MO-6 27 mm (Taper Point)	RB-1 Needle (Taper Point)


Taper Point
Designed for soft, easy penetrated tissue.



Reverse Cutting
Cutting edge on outer curve. Designed for tough, difficult-to-penetrate tissue.



TAPERCUT™ Technology
Cutting tip, taper body designed for tough tissue, like two needles in one.



Free Needles



CODE	DESCRIPTION	QTY
223173	Free Needle Loop (OS-6)	12/pack



Shoulder


↳ Sutures & Suture Passers

DYNATAPE™ Sutures



CODE	DESCRIPTION	QTY
280825	2.5 mm DYNATAPE Suture Pack, White, without needles	6/pack
280828	2.5 mm DYNATAPE Suture Pack, Black, without needles	6/pack
223173	Free Needle Loop OS-6	12/pack

DYNACORD™ Sutures



CODE	DESCRIPTION	QTY
222067	Size #2 DYNACORD Suture Pack, Blue, without Needle	12/pack
222065	Size #2 DYNACORD Suture Pack, Blue, with OS-6 Needle	12/pack
222066	Size #2 DYNACORD Suture Pack, Blue, with MO-7 Needle	12/pack
222068	Size #2 DYNACORD Suture Pack, Striped, without Needle	12/pack
222069	Size #2 DYNACORD Suture Pack, Striped/Blue, without Needle	12/pack
222071	Size #2 DYNACORD Suture Pack, Striped/Blue, with MO-7 Needle	12/pack
222073	Size #2 DYNACORD Suture Pack, Striped/Blue, with OS-6 Needle	12/pack

Shoulder ↳ Sutures & Suture Passers

ORTHOCORD™ Sutures



Size #2 Sutures

CODE	DESCRIPTION	QTY
223102	ORTHOCORD Suture #2, Violet, with MO-6 1/2 Circle, Taper Point Needle, 26 mm	12/pack
223103	ORTHOCORD Suture #2, Violet, with OS-6 1/2 Circle, Reverse Cutting Needle, 36 mm	12/pack
223104	ORTHOCORD Suture #2, Violet, with MO-7 1/2 Circle, Taper Point Needle, 22 mm	12/pack
223105	ORTHOCORD Suture #2, Violet, without Needles	12/pack
223111	ORTHOCORD Suture #2, Blue, without Needles	12/pack
223113	ORTHOCORD Suture #2, Violet/Blue, without Needles	12/pack
223114	ORTHOCORD Suture #2, Violet/Blue, with MO-7 1/2 Circle, Taper Point Needle, 22 mm	12/pack
223115	ORTHOCORD Suture #2, Violet/Blue, with MO-6 1/2 Circle, Taper Point Needle, 26 mm	12/pack
223116	ORTHOCORD Suture #2, Violet/Blue, with OS-6 1/2 Circle, Reverse Cutting Needle, 36 mm	12/pack

Small Size Sutures

CODE	DESCRIPTION	QTY
223121	#4-0 ORTHOCORD Suture Loop 51cm Taper Point, 3/8 Circle, 17.9 mm	12/pack
223122	#2-0 ORTHOCORD Suture 46cm Taper Point, 1/2 Circle, 26.5 mm	12/pack
223123	#4-0 ORTHOCORD Suture 46cm Diamond Point, 3/8 Circle, 18.7 mm	12/pack
223124	#0 ORTHOCORD Suture 96.5cm Taper Point, 1/2 Circle, 22.2 mm	12/pack
223125	#3-0 ORTHOCORD Suture 46cm Taper Point, 3/8 Circle, 15 mm	12/pack
223126	#2-0 ORTHOCORD Suture 46cm Taper Point, 3/8 Circle, 17.9 mm	12/pack
223127	#4-0 ORTHOCORD Suture 46cm Taper Point, 3/8 Circle, 2.3 mm	12/pack



Shoulder

↳ Sutures & Suture Passers


2.5 mm PERMATAPE™ Sutures



CODE	DESCRIPTION	QTY
223165	2.5 mm PERMATAPE Suture, Blue Pack	6/pack
223166	2.5 mm PERMATAPE Suture, White Pack	6/pack
223169	2.5 mm PERMATAPE Suture, Blue	1/each
223170	2.5 mm PERMATAPE Suture, White	1/each

1.3 mm PERMATAPE™ Sutures

CODE	DESCRIPTION	QTY
280772	1.3 mm PERMATAPE Suture, Blue, without Needles	12/pack
280773	1.3 mm PERMATAPE Suture, Blue, OS6 Needles	12/pack
280774	1.3 mm PERMATAPE Suture, Blue, MO-7 Needles	12/pack
280775	1.3 mm PERMATAPE Suture, White-Blue, without Needles	12/pack
280776	1.3 mm PERMATAPE Suture, White-Blue, OS6 Needles	12/pack
280777	1.3 mm PERMATAPE Suture, White-Blue, MO-7 Needles	12/pack



PERMALOOP™ Sutures

CODE	DESCRIPTION	QTY
280936	1.3 mm PERMALOOP Suture Straight Needle	12/pack
208938	1.3 mm PERMALOOP Suture MO-5 Needle	12/pack



Shoulder

↳ Sutures & Suture Passers

EXPRESSEW™ III Suture Passers



CODE	DESCRIPTION	QTY
214124	EXPRESSEW III with Hook	1/each
214140	EXPRESSEW III without Hook	1/each
288233	EXPRESSEW III with AutoCapture+	1/each
288235	Retention plate with Loader	1/each
288234	Retention plate with Loader	5/pack
214141	EXPRESSEW III Needles	5/pack

Shoulder

↳ Sutures & Suture Passers

IDEAL™ Suture Graspers & Shuttles



IDEAL Suture Graspers

CODE	DESCRIPTION	QTY
251720	IDEAL Suture Grasper, 15°	1/each
251721	IDEAL Suture Grasper, 30°	1/each
251722	IDEAL Suture Grasper, 45°	1/each
251723	IDEAL Suture Grasper, 60°	1/each

IDEAL Suture Shuttles

CODE	DESCRIPTION	QTY
251001	IDEAL Suture Shuttle, 25° Left	1/each
251002	IDEAL Suture Shuttle, 25° Right	1/each
251003	IDEAL Suture Shuttle, 45° Left	1/each
251004	IDEAL Suture Shuttle, 45° Right	1/each
251005	IDEAL Suture Shuttle, 90° Up	1/each
251006	IDEAL Suture Shuttle, Straight Hook	1/each
251007	IDEAL Suture Shuttle, Straight Crescent	1/each



CHIA PERCPASSER™ Suture Passer

CODE	DESCRIPTION	QTY
214101	CHIA PERCPASSER Suture Passer	5/pack



Shoulder

↳ Sutures & Suture Passers

Suture Leaders

CODE	DESCRIPTION	QTY
214608	Suture Leader, 20° Up	1/each
214609	Suture Leader, 30° Right	1/each
214610	Suture Leader, 30° Left	1/each
214611	Suture Leader, 45° Right	1/each
214612	Suture Leader, 45° Left	1/each
214613	Suture Leader, 70° Right	1/each
214614	Suture Leader, 70° Left	1/each
214628	Suture Leader, 180° Right	1/each
214629	Suture Leader, 180° Left	1/each



Shoulder ↳ Graft Preparation

SPEEDTRAP™ Graft Preparation System

CODE	DESCRIPTION	QTY
223746	SPEEDTRAP System, White, 20 mm	1/each
223751	SPEEDTRAP System, Green/White, 20 mm	1/each
223747	SPEEDTRAP System, White, 30 mm	1/each
223748	SPEEDTRAP System, Green/White, 30 mm	1/each
223749	SPEEDTRAP System Pack, White, 30 mm	4/pack
223750	SPEEDTRAP System Pack, Green/White, 30 mm	4/pack



PERMALOOP™ Sutures

CODE	DESCRIPTION	QTY
280936	1.3 mm PERMALOOP Suture Straight Needle	12/pack
208938	1.3 mm PERMALOOP Suture MO-5 needle	12/pack



Graft Preparation System



CODE	DESCRIPTION	QTY
215810	Graft Preparation System Tray	1/each
219961	Graft Preparation Board	1/each
219962	Slider Unit (requires 2 each for system)	1/each
219963	Tensioning Clamp	1/each
219964	Graft Clamp (requires 2 each for system)	1/each
219965	Suture Clamp	1/each
219444	Graft Sizer, 5-12 mm (white color)	1/each
219966	Graft Sizer, 6-12 mm (amber color)	1/each
219967	Drill Guide	1/each
219968	Cutting Board	1/each



Shoulder

↳ Arthroscopic Instruments & Accessories

Arthroscopic Instruments

Clever Hook / Sixter

CODE	DESCRIPTION	QTY
214640	Clever Hook, Right	1/each
214641	Clever Hook, Left	1/each
214642	Sixter, Right	1/each
214643	Sixter, Left	1/each

General


CODE	DESCRIPTION	QTY
214600	Penetrating Grasper, Straight	1/each
214601	Penetrating Grasper, 35° Up	1/each
214602	Penetrating Grasper, 45° Left	1/each
214603	Penetrating Grasper, 45° Right	1/each
214604	Combo Grasper	1/each
214575	Arthroscopic Suture Grasper	1/each
214627	Sliding Suture Cutter	1/each
270120	Grasper / Grabber	1/each
214626	Suture Manipulator	1/each
214646	CordCutter (2 parts)	1/each
210263	Flush Suture Cutter	1/each
214606	Tissue Grasper w/ Ratchet	1/each
214607	Forked Tissue Grasper	1/each
214615	Knot Manipulator, Full Loop (open)	1/each
210835	PROKNOT Knot Manipulator (closed)	1/each
214616	Crochet Hook	1/each
214617	Calibrated Probe	1/each
214618	Switching Stick, 4.5 mm	1/each
214502	Skid Rod with Handle, #2-0 or #0 Suture	1/each
214605	Suture Scissors	1/each
214681	Switching Stick Semi Blunt, 4 mm	1/each
214619	Rasp, 20° Up, Bend	1/each
214620	Rasp, 20° Down, Bend	1/each



Penetrating Grasper



Combo Grasper




CordCutter



Flush Suture Cutter



Knot Manipulator



Close

Open



Shoulder

↳ Arthroscopic Instruments & Accessories

Arthroscopic Instruments

General (continued)

CODE	DESCRIPTION	QTY
214621	Tissue Liberator, Atraumatic Up	1/each
214622	Tissue Liberator, Atraumatic Down	1/each
214623	Tissue Liberator, Blade Up	1/each
214624	Tissue Liberator, Blade Down	1/each
214008	Atraumatic Soft Tissue Liberator	1/each
214529	Arthroscopic Periosteal Elevator	1/each

Flex Soft Cannulas

CODE	DESCRIPTION	QTY
210252	Flex Soft Cannula 8x50 mm	5/pack
210253	Flex Soft Cannula 10x50 mm	5/pack
210254	Flex Soft Cannula 12x50 mm	5/pack
210242	Flex Soft Cannula 10x30 mm	5/pack



Shoulder

↳ Arthroscopic Instruments & Accessories

Clear Cannulas

CODE	DESCRIPTION	QTY
214104	Clear Cannula, Threaded, 5.5 mm x 55 mm	5/pack
214106	Clear Cannula, Smooth, 5.5 mm x 75 mm	5/pack
214107	Clear Cannula, Smooth, Distal Rib, 5.5 mm x 75 mm	5/pack
214108	Clear Cannula, Threaded, 5.5 mm x 75 mm	5/pack
214110	Clear Cannula, Threaded, 5.5 mm x 90 mm	5/pack
214114	Clear Cannula, Smooth, 7 mm x 75 mm	5/pack
214115	Clear Cannula, Smooth, Distal Rib, 7 mm x 75 mm	5/pack
214116	Clear Cannula, Threaded, 7 mm x 75 mm	5/pack
214112	Clear Cannula, Partially Threaded, 7 mm x 110 mm	5/pack
214118	Clear Cannula, Threaded, 8.5 mm x 55 mm	5/pack
214120	Clear Cannula, Threaded, 8.5 mm x 75 mm	5/pack
214122	Clear Cannula, Threaded, 8.5 mm x 90 mm	5/pack
214103	5.5 mm x 55 mm Reusable Obturator	1/each
214105	5.5 mm x 75 mm Reusable Obturator	1/each
214109	5.5 mm x 90 mm Reusable Obturator	1/each
214113	7 mm x 75 mm Reusable Obturator	1/each
214111	7 mm x 110 mm Reusable Obturator	1/each
214117	8.5 mm x 55 mm Reusable Obturator	1/each
214119	8.5 mm x 75 mm Reusable Obturator	1/each
214121	8.5 mm x 90 mm Reusable Obturator	1/each
214123	4 mm Switching Stick, Blunt/Sharp	1/each
214681	4 mm Switching Stick, Blunt/Blunt	1/each



Shoulder


↳ Arthroscopic Instruments & Accessories

Sterilization Trays

CODE	DESCRIPTION	QTY
215401	Shoulder Procedural Tray (Outer Case)	1/each
215402	Suture Management Sterilization Case	1/each
215403	RC/Instability Fixation Tray (Top Shelf In Shoulder Tray)	1/each
215406	Arthroscopic Instability Sterilization Tray (10 Slot)	1/each
215407	Knot Pusher Sterilization Tray	1/each
215408	Knot Pusher Sterilization Tray Pin Mat	1/each
215409	Suture Leader Caddy	1/each
215410	Suture Leader Caddy Pin Mat	1/each




Knee



Knee ↳ Ligament ↳ Cortical Fixation

RIGIDLOOP™ Adjustable Cortical Fixation System



Implants

CODE	DESCRIPTION	QTY
232447	RIGIDLOOP Adjustable System Standard Loop Implant	1/each
232448	RIGIDLOOP Adjustable System Long Loop Implant	1/each
232449	RIGIDLOOP Adjustable System XL Button	1/each
232480	RIGIDLOOP T Adjustable System Standard Button	1/each
232481	RIGIDLOOP T Adjustable System XL Button	1/each
232485	RIGIDLOOP Adjustable BTB Button	1/each
232486	RIGIDLOOP Adjustable BTB XL Button	1/each
232484	RIGIDLOOP Suture Loop	1/each
232482	RIGIDLOOP Clip Cortical Button (Standard)	1/each
232483	RIGIDLOOP Clip Cortical Button (XL)	1/each



Instruments

CODE	DESCRIPTION	QTY
232453	RIGIDLOOP Adjustable System Spade Tip Combo Reamer/Passing Pin, 4.5 mm	1/each
232460	RIGIDLOOP Adjustable System Spade Tip Combo Reamer/Passing Pin, 4.5 mm (Longer Tip)	1/each
218034	RIGIDLOOP Adjustable Disposable Kit (2.4 mm x 14" Drill Tip Passing Pin, 4.5 mm Spade Tip Combo Reamer/Passing Pin, 1.1 mm x 15" Nitinol Guidewire, Threaded Cannula, Malleable Graft Retractor, Marking Pen, Ruler)	1/each



Knee ↳ Ligament ↳ Cortical Fixation

RIGIDLOOP™ Cortical Fixation System

Implants

CODE	DESCRIPTION	QTY
232008	RIGIDLOOP Titanium Button	1/each
232009	RIGIDLOOP Fixed System Loop Implant, 15 mm	1/each
232011	RIGIDLOOP Fixed System Loop Implant, 20 mm	1/each
232014	RIGIDLOOP Fixed System Loop Implant, 25 mm	1/each
232016	RIGIDLOOP Fixed System Loop Implant, 30 mm	1/each
232017	RIGIDLOOP Fixed System Loop Implant, 35 mm	1/each
232018	RIGIDLOOP Fixed System Loop Implant, 40 mm	1/each
232019	RIGIDLOOP Fixed System Loop Implant, 45 mm	1/each
232021	RIGIDLOOP Fixed System Loop Implant, 50 mm	1/each
232022	RIGIDLOOP Fixed System Loop Implant, 55 mm	1/each
232023	RIGIDLOOP Fixed System Loop Implant, 60 mm	1/each
232034	RIGIDLOOP XL Button	1/each



RIGIDLOOP Titanium Button



RIGIDLOOP Fixed System Loop Implant



RIGIDLOOP XL Button

Instruments

CODE	DESCRIPTION	QTY
232025	RIGIDLOOP System Trocar Tip Passing Pin, 2.4 mm x 18"	1/each
232028	RIGIDLOOP System Cannulated Reamer, 4.5 mm	1/each
232030	RIGIDLOOP System Standard Depth Gauge	1/each
232031	RIGIDLOOP System Variable Depth Gauge	1/each
232037	RIGIDLOOP Disposable Kit (2.4 mm x 18" Trocar Tip Passing Pin, 4.5 mm Cannulated Reamer, Depth Gauge, 1.1 mm x 15" Nitinol Guidewire, Threaded Cannula, Malleable Graft Retractor, Marking Pen, Ruler)	1/each
232032	RIGIDLOOP System Instrument Caddy	1/each
232033	RIGIDLOOP System Silicone Mat	1/each
215451	Large Sterilization Tray Lid	1/each



Standard Depth Gauge



Variable Depth Gauge



RIGIDLOOP Disposable Kit



Instrument Caddy and Lid



Knee ↳ Ligament ↳ Interference Screws



MILAGRO™ ADVANCE Interference Screws

BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
231816	MILAGRO ADVANCE BR Interference Screw, 7 mm x 23 mm	1/each
231831	MILAGRO ADVANCE BR Interference Screw, 7 mm x 30 mm	1/each
231817	MILAGRO ADVANCE BR Interference Screw, 8 mm x 23 mm	1/each
231821	MILAGRO ADVANCE BR Interference Screw, 8 mm x 30 mm	1/each
231818	MILAGRO ADVANCE BR Interference Screw, 9 mm x 23 mm	1/each
231822	MILAGRO ADVANCE BR Interference Screw, 9 mm x 30 mm	1/each
231826	MILAGRO ADVANCE BR Interference Screw, 9 mm x 35 mm	1/each
231819	MILAGRO ADVANCE BR Interference Screw, 10 mm x 23 mm	1/each
231823	MILAGRO ADVANCE BR Interference Screw, 10 mm x 30 mm	1/each
231827	MILAGRO ADVANCE BR Interference Screw, 10 mm x 35 mm	1/each
231824	MILAGRO ADVANCE BR Interference Screw, 11 mm x 30 mm	1/each
231828	MILAGRO ADVANCE BR Interference Screw, 11 mm x 35 mm	1/each
231832	MILAGRO ADVANCE BR Interference Screw, 12 mm x 30 mm	1/each
231829	MILAGRO ADVANCE BR Interference Screw, 12 mm x 35 mm	1/each

PEEK

CODE	DESCRIPTION	QTY
231856	MILAGRO ADVANCE PEEK Interference Screw, 7 mm x 23 mm	1/each
231866	MILAGRO ADVANCE PEEK Interference Screw, 7 mm x 30 mm	1/each
231857	MILAGRO ADVANCE PEEK Interference Screw, 8 mm x 23 mm	1/each
231867	MILAGRO ADVANCE PEEK Interference Screw, 8 mm x 30 mm	1/each
231858	MILAGRO ADVANCE PEEK Interference Screw, 9 mm x 23 mm	1/each
231868	MILAGRO ADVANCE PEEK Interference Screw, 9 mm x 30 mm	1/each
231859	MILAGRO ADVANCE PEEK Interference Screw, 10 mm x 23 mm	1/each
231869	MILAGRO ADVANCE PEEK Interference Screw, 10 mm x 30 mm	1/each
231870	MILAGRO ADVANCE PEEK Interference Screw, 11 mm x 30 mm	1/each
231871	MILAGRO ADVANCE PEEK Interference Screw, 12 mm x 30 mm	1/each



Knee ↳ Ligament ↳ Interference Screws

MILAGRO™ ADVANCE Instruments



Universal Modular Driver



Quick-Connect Ratcheting Handle

Instruments

CODE	DESCRIPTION	QTY
215452	Medium Sterilization Tray	1/each
215453	Medium Sterilization Tray Lid	1/each
219622	5 mm MILAGRO Screw Tap (for use with Quick-Connect)	1/each
219623	6 mm MILAGRO Screw Tap (for use with Quick-Connect)	1/each
219496	MILAGRO ADVANCE Screw Tap, 7-8 mm	1/each
219497	MILAGRO ADVANCE Screw Tap, 9-10 mm	1/each
219498	MILAGRO ADVANCE Screw Tap, 11-12 mm	1/each
219500	MILAGRO ADVANCE Screw Universal Modular Driver	1/each
254847	Quick-Connect Ratcheting Handle	1/each
254715	Grey Ratcheting Handle with Jacob's Chuck	1/each
219960	Modular Driver, 23 mm	1/each
219958	Modular Driver, 30 mm	1/each
219959	Modular Driver, 35 mm	1/each
270105	30 mm Hex Driver with Teardrop Handle	1/each
219952	Pathmaker Tunnel Notcher	1/each

Knee ↳ Ligament ↳ Interference Screws

Small Size MILAGRO™ Interference Screws



Implants

CODE	DESCRIPTION	QTY
231814	MILAGRO Bioabsorbable Screw, 5 mm x 12 mm	1/each
231803	MILAGRO Bioabsorbable Screw, 5 mm x 23 mm	1/each
231804	MILAGRO Bioabsorbable Screw, 5 mm x 30 mm	1/each
231813	MILAGRO Bioabsorbable Screw, 6 mm x 12 mm	1/each
231807	MILAGRO Bioabsorbable Screw, 6 mm x 23 mm	1/each
231809	MILAGRO Bioabsorbable Screw, 6 mm x 30 mm	1/each
231811	MILAGRO Bioabsorbable Screw, 7 mm x 15 mm	1/each
231812	MILAGRO Bioabsorbable Screw, 8 mm x 15 mm	1/each



Universal Modular Driver



Quick-Connect Ratcheting Handle

Instruments

CODE	DESCRIPTION	QTY
219448	Modular Driver – Normal Cannulation, 12 mm	1/each
219447	Modular Driver – Large Cannulation, 15 mm	1/each
219446	Modular Driver (5/6 mm Screws), 23 mm	1/each
219438	Modular Driver (5/6 mm Screws), 30 mm	1/each
219960	Modular Driver (Standard Length), 23 mm	1/each
219950	Starter Tap	1/each
219959	Modular Driver, 35 mm	1/each
219957	Modular Driver (24cm Length), 23 mm	1/each



Knee ↳ Ligament ↳ Interference Screws

PROFILE™ Titanium Interference Screws



Titanium

CODE	DESCRIPTION	QTY
230320	PROFILE Titanium Interference Screw, 7 mm x 20 mm	1/each
230325	PROFILE Titanium Interference Screw, 7 mm x 25 mm	1/each
230330	PROFILE Titanium Interference Screw, 7 mm x 30 mm	1/each
230420	PROFILE Titanium Interference Screw, 8 mm x 20 mm	1/each
230425	PROFILE Titanium Interference Screw, 8 mm x 25 mm	1/each
230430	PROFILE Titanium Interference Screw, 8 mm x 30 mm	1/each
230520	PROFILE Titanium Interference Screw, 9 mm x 20 mm	1/each
230525	PROFILE Titanium Interference Screw, 9 mm x 25 mm	1/each
230530	PROFILE Titanium Interference Screw, 9 mm x 30 mm	1/each



Instruments

CODE	DESCRIPTION	QTY
219203	Cannulated Hex Driver, 3.5 mm (for 7 mm and 8 mm Titanium Screws)	1/each
219204	Cannulated Hex Driver, 4 mm (for 9 mm, 11 mm and 13 mm Titanium Screws)	1/each

Knee ↳ Ligament ↳ Interference Screws

Screw and Washer System



Implants

CODE	DESCRIPTION	QTY	CODE	DESCRIPTION	QTY
220225	Cancellous Screw, 6.5 mm x 25 mm	1/each	220450	Cortical Screw, 4.5 mm x 50 mm	1/each
220230	Cancellous Screw, 6.5 mm x 30 mm	1/each	220452	Cortical Screw, 4.5 mm x 52 mm	1/each
220235	Cancellous Screw, 6.5 mm x 35 mm	1/each	220454	Cortical Screw, 4.5 mm x 54 mm	1/each
220240	Cancellous Screw, 6.5 mm x 40 mm	1/each	220456	Cortical Screw, 4.5 mm x 56 mm	1/each
220245	Cancellous Screw, 6.5 mm x 45 mm	1/each	220458	Cortical Screw, 4.5 mm x 58 mm	1/each
220250	Cancellous Screw, 6.5 mm x 50 mm	1/each	220460	Cortical Screw, 4.5 mm x 60 mm	1/each
220255	Cancellous Screw, 6.5 mm x 55 mm	1/each	220462	Cortical Screw, 4.5 mm x 62 mm	1/each
220426	Cortical Screw, 4.5 mm x 26 mm	1/each	220464	Cortical Screw, 4.5 mm x 64 mm	1/each
220428	Cortical Screw, 4.5 mm x 28 mm	1/each	220466	Cortical Screw, 4.5 mm x 66 mm	1/each
220430	Cortical Screw, 4.5 mm x 30 mm	1/each	220468	Cortical Screw, 4.5 mm x 68 mm	1/each
220432	Cortical Screw, 4.5 mm x 32 mm	1/each	220470	Cortical Screw, 4.5 mm x 70 mm	1/each
220434	Cortical Screw, 4.5 mm x 34 mm	1/each	220472	Cortical Screw, 4.5 mm x 72 mm	1/each
220436	Cortical Screw, 4.5 mm x 36 mm	1/each	220474	Cortical Screw, 4.5 mm x 74 mm	1/each
220438	Cortical Screw, 4.5 mm x 38 mm	1/each	220476	Cortical Screw, 4.5 mm x 76 mm	1/each
220440	Cortical Screw, 4.5 mm x 40 mm	1/each	221014	Spiked Washer, 14 mm (OD)	1/each
220442	Cortical Screw, 4.5 mm x 42 mm	1/each	221017	Spiked Washer, 17 mm (OD)	1/each
220444	Cortical Screw, 4.5 mm x 44 mm	1/each	221020	Spiked Washer, 20 mm (OD)	1/each
220446	Cortical Screw, 4.5 mm x 46 mm	1/each	221114	Flat Washer, 14 mm (OD)	1/each
220448	Cortical Screw, 4.5 mm x 48 mm	1/each	221117	Flat Washer, 17 mm (OD)	1/each

Instruments

CODE	DESCRIPTION	QTY
211024	Drill Bit, 3.2 mm	1/each
219114	New Depth Gauge (Screw & Washer)	1/each
219124	Hex Driver, 2.5 mm	1/each
219200	Cancellous Tap, 6.5 mm	1/each
219400	Cortical Tap, 4.5 mm	1/each
215024	Screw and Washer Implant Tray	1/each
215124	Screw and Washer System Base & Tray	1/each



Knee ↳ Ligament ↳ Screw & Sheath Fixation

INTRAFIX™ ADVANCE Tibial Fastener System

BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
254800	INTRAFIX ADVANCE BR Screw and Small Sheath, 7 mm x 23 mm	1/each
254806	INTRAFIX ADVANCE BR Screw and Small Sheath, 7 mm x 30 mm	1/each
254801	INTRAFIX ADVANCE BR Screw and Small Sheath, 8 mm x 23 mm	1/each
254807	INTRAFIX ADVANCE BR Screw and Small Sheath, 8 mm x 30 mm	1/each
254802	INTRAFIX ADVANCE BR Screw and Large Sheath, 9 mm x 23 mm	1/each
254808	INTRAFIX ADVANCE BR Screw and Large Sheath, 9 mm x 30 mm	1/each
254803	INTRAFIX ADVANCE BR Screw and Large Sheath, 10 mm x 23 mm	1/each
254809	INTRAFIX ADVANCE BR Screw and Large Sheath, 10 mm x 30 mm	1/each



PEEK

CODE	DESCRIPTION	QTY
254818	INTRAFIX ADVANCE Polypropylene Sheath, Small (7-8 mm), 23 mm	1/each
254825	INTRAFIX ADVANCE Polypropylene Sheath, Small (7-8 mm), 30 mm	1/each
254819	INTRAFIX ADVANCE Polypropylene Sheath, Large (9-10 mm), 23 mm	1/each
254826	INTRAFIX ADVANCE Polypropylene Sheath, Large (9-10 mm), 30 mm	1/each
254827	INTRAFIX ADVANCE Polypropylene Sheath, XL (11-12 mm), 30 mm	1/each
254814	INTRAFIX ADVANCE PEEK Screw, 7 mm x 23 mm	1/each
254820	INTRAFIX ADVANCE PEEK Screw, 7 mm x 30 mm	1/each
254815	INTRAFIX ADVANCE PEEK Screw, 8 mm x 23 mm	1/each
254821	INTRAFIX ADVANCE PEEK Screw, 8 mm x 30 mm	1/each
254816	INTRAFIX ADVANCE PEEK Screw, 9 mm x 23 mm	1/each
254822	INTRAFIX ADVANCE PEEK Screw, 9 mm x 30 mm	1/each
254817	INTRAFIX ADVANCE PEEK Screw, 10 mm x 23 mm	1/each
254823	INTRAFIX ADVANCE PEEK Screw, 10 mm x 30 mm	1/each
254824	INTRAFIX ADVANCE PEEK Screw, 11 mm x 30 mm	1/each
254841	INTRAFIX ADVANCE PEEK Screw, 12 mm x 30 mm	1/each



Knee ↳ Ligament ↳ Screw & Sheath Fixation

INTRAFIX™ ADVANCE Tibial Fastener System Instruments

Instruments

CODE	DESCRIPTION	QTY
254828	INTRAFIX ADVANCE Small Trial/Dilator	1/each
254829	INTRAFIX ADVANCE Large Trial/Dilator	1/each
254830	INTRAFIX ADVANCE XL Trial/Dilator	1/each
254831	INTRAFIX ADVANCE Sheath Inserter Small, 23 mm	1/each
254832	INTRAFIX ADVANCE Sheath Inserter Large, 23 mm	1/each
254833	INTRAFIX ADVANCE Sheath Inserter Small, 30 mm	1/each
254834	INTRAFIX ADVANCE Sheath Inserter Large, 30 mm	1/each
254835	INTRAFIX ADVANCE Sheath Inserter XL, 30 mm	1/each
254836	INTRAFIX ADVANCE Hex Head Screwdriver	1/each
254838	INTRAFIX ADVANCE Tensioner	1/each
254605	Tie Tensioner	1/each
254844	Guidewire 12" (30.5 cm)	1/each
254842	INTRAFIX ADVANCE Slap Hammer	1/each
254847	INTRAFIX ADVANCE Ratchet Handle	1/each
254839	INTRAFIX ADVANCE Sterilization Case (Outer)	1/each
254845	INTRAFIX ADVANCE Sterilization Tray 1	1/each
254846	INTRAFIX ADVANCE Sterilization Tray 2	1/each
254848	INTRAFIX ADVANCE Sterilization Lid	1/each



INTRAFIX ADVANCE Sheath Inserters





Tie Tensioner

Knee ↳ Ligament ↳ Cross Pin Fixation

RIGIDFIX™ Cross Pin System


CODE	DESCRIPTION	QTY
210133	RIGIDFIX PLA ST Femoral Cross Pin Kit, 3.3 mm	1/each
210816	RIGIDFIX BIOCRYL Material ST Femoral Cross Pin Kit, 3.3 mm	1/each
210127	RIGIDFIX PLA BTB Femoral Cross Pin Kit, 2.7 mm	1/each
210815	RIGIDFIX BIOCRYL Material BTB Femoral Cross Pin Kit, 2.7 mm	1/each
213701	RIGIDFIX System ST Guide Frame	1/each
213702	RIGIDFIX System BTB Guide Frame	1/each
213713	RIGIDFIX System BTB Stepped Trocar, Short	1/each
213703	RIGIDFIX System BTB Stepped Trocar, Long	1/each
213704	RIGIDFIX System Stepped Pin Inserter Rod	1/each
270118	RIGIDFIX System Femoral Rod, 6 mm (ST Only)	1/each
270119	RIGIDFIX System Femoral Rod, 6.5 mm (ST Only)	1/each
213707	RIGIDFIX System Femoral Rod, 7 mm (ST Only)	1/each
213717	RIGIDFIX System Femoral Rod, 7.5 mm (ST Only)	1/each
213708	RIGIDFIX System Femoral Rod, 8 mm (ST Only)	1/each
213718	RIGIDFIX System Femoral Rod, 8.5 mm (ST Only)	1/each
213709	RIGIDFIX System Femoral Rod, 9 mm	1/each
213719	RIGIDFIX System Femoral Rod, 9.5 mm	1/each
213710	RIGIDFIX System Femoral Rod, 10 mm	1/each
213720	RIGIDFIX System Femoral Rod, 10.5 mm	1/each
213711	RIGIDFIX System Femoral Rod, 11 mm	1/each
213721	RIGIDFIX System Femoral Rod, 11.5 mm	1/each
213712	RIGIDFIX System Femoral Rod, 12 mm	1/each
213705	RIGIDFIX System Mallet	1/each
270127	Hex Driver with Teardrop Handle, 23 mm	1/each
270128	Hex Driver with Teardrop Handle, 35 mm	1/each
213706	RIGIDFIX System Sleeve Removal Tool	1/each
215422	RIGIDFIX System Sterilization Tray	1/each
215423	ST RIGIDFIX System Caddy	1/each
215424	BTB RIGIDFIX System Caddy	1/each
215425	Combo RIGIDFIX System Caddy	1/each
213600	RIGIDFIX System/Tibial Rod Replacement Thumb Screw	1/each



Knee ↳ Graft Preparation

SPEEDTRAP™ Graft Preparation System

CODE	DESCRIPTION	QTY
223746	SPEEDTRAP System, White, 20 mm	1/each
223751	SPEEDTRAP System, Green/White, 20 mm	1/each
223747	SPEEDTRAP System, White, 30 mm	1/each
223748	SPEEDTRAP System, Green/White, 30 mm	1/each
223749	SPEEDTRAP System Pack, White, 30 mm	4/pack
223750	SPEEDTRAP System Pack, Green/White, 30 mm	4/pack



PERMALOOP™ Sutures

CODE	DESCRIPTION	QTY
280936	1.3 mm PERMALOOP Suture Straight Needle	12/pack
208938	1.3 mm PERMALOOP Suture MO-5 needle	12/pack



Graft Preparation System



CODE	DESCRIPTION	QTY
215810	Graft Preparation System Tray	1/each
219961	Graft Preparation Board	1/each
219962	Slider Unit (requires 2 each for system)	1/each
219963	Tensioning Clamp	1/each
219964	Graft Clamp (requires 2 each for system)	1/each
219965	Suture Clamp	1/each
219444	Graft Sizer, 5-12 mm (white color)	1/each
219966	Graft Sizer, 6-12 mm (amber color)	1/each
219967	Drill Guide	1/each
219968	Cutting Board	1/each



Knee

↳ Graft Preparation

ACL Graft Knives

CODE	DESCRIPTION	QTY
232108	ACL Graft Knife, 8 mm (Sterile, Single-Use)	10/pack
232109	ACL Graft Knife, 9 mm (Sterile, Single-Use)	10/pack
232110	ACL Graft Knife, 10 mm (Sterile, Single-Use)	10/pack
232111	ACL Graft Knife, 11 mm (Sterile, Single-Use)	10/pack
232112	ACL Graft Knife, 12 mm (Sterile, Single-Use)	10/pack



Knee ↳ Meniscus

TRUESPAN™ Meniscal Repair System



PLGA


CODE	DESCRIPTION	QTY
228160	TRUESPAN PLGA 0° Needle	1/each
228161	TRUESPAN PLGA 12° Needle	1/each
228162	TRUESPAN PLGA 24° Needle	1/each

PEEK

CODE	DESCRIPTION	QTY
228150	TRUESPAN PEEK 0° Needle	1/each
228151	TRUESPAN PEEK 12° Needle	1/each
228152	TRUESPAN PEEK 24° Needle	1/each

Meniscal Repair System

#2-0 ORTHOCORD Suture with swedged-on double-armed 13-inch stainless steel needles (228144).



#2-0 ORTHOCORD (3 free strands) with a single 10-inch nitinol needle with wire eyelet (228145).

Meniscal Needles (Sterile)

CODE	DESCRIPTION	QTY
228144	#2-0 ORTHOCORD Sutures with double-armed meniscal needles	10/pack
228145	#2-0 ORTHOCORD Sutures (3 free strands) with single meniscal needle with wire loop	10/pack



Knee

↳ Meniscus

Meniscal Repair System



Instruments

CODE	DESCRIPTION	QTY
228302	Arthroscopic Pusher/Cutter	1/each
228301	Arthroscopic Pusher	1/each
232024	Malleable Graft Retractor (non-sterile)	5/pack
253110	Meniscal Repair System Rasp, Straight Up with 45° Face	1/each
253111	Meniscal Repair System Rasp, Straight Up with 90° Face	1/each
253112	Meniscal Repair System Rasp, Straight with 90° Face	1/each
253113	Meniscal Calibration Probe Straight	1/each
215426	Meniscal Repair System Caddy	1/each
215427	System Caddy Pin Mat	1/each

Knee ↳ Cartilage

COR™ Precision Targeting Cartilage Repair System



CODE	DESCRIPTION	QTY
252107	COR System Disposable Kit, 4 mm	1/each
252108	COR System Disposable Kit, 6 mm	1/each
252109	COR System Disposable Kit, 8 mm	1/each
252110	COR System Disposable Kit, 10 mm	1/each
252111	COR System Disposable Kit, 6 mm with Perpendicularity	1/each
252112	COR System Disposable Kit, 8 mm with Perpendicularity	1/each
252113	COR System Disposable Kit, 10 mm with Perpendicularity	1/each
252307	COR System Plunger/Sizer 4 mm	1/each
252313	COR System Plunger/Sizer 6 mm	1/each
252319	COR System Plunger/Sizer 8 mm	1/each

Knee ↳ Tibial Realignment

TRACKER™ AMZ Guide System



CODE	DESCRIPTION	QTY
219465	TRACKER System Anchor Pin, 2 mm x 45 mm	1/each
219466	TRACKER System Anchor Pin, 2 mm x 64 mm	1/each
219467	TRACKER System Anchor Pin, 2 mm x 77 mm	1/each
215187	TRACKER System Guide Tray	1/each
219452	TRACKER System Cutting Block	1/each
219453	TRACKER System Cutting Block with Pivot	1/each
219454	TRACKER System Retractor	1/each
219455	TRACKER System Slope Selector Arm	1/each
219456	TRACKER System Anchor Pin Puller	1/each
219457	TRACKER System Bone Clamp	1/each
219458	TRACKER System Positioning Card	1/each
219460	TRACKER System Drill Bit, 2 mm x 6.5"	1/each
219461	TRACKER System Drill Bit, 3.2 mm x 7.0"	1/each
219462	TRACKER System Drill Bit, 4.5 mm x 6.5"	1/each

Knee


↳ Arthroscopic Instruments & Accessories

TWISTR™ Retrograde Reamer and Cruciate+ Instruments



TWISTR™ Retrograde Reamer and Disposables

CODE	DESCRIPTION	QTY
232000	TWISTR Retrograde Reamer, 6-12 mm	1/each
232301	TWISTR ACL Disposables Kit (2.4 mm x 13" Outside-In Forked Tip Passing Pin, 2.4 mm x 14" Drill Tip Passing Pin, 1.1 mm x 15" Nitinol Guidewire, Malleable Graft Retractor, Marking Pen, Ruler)	1/each



Cruciate+ Instruments

CODE	DESCRIPTION	QTY
219001	Side-Loading Guide Carriage	1/each
219002	Retrograde Reamer Bullet Obturator	1/each
219003	Beath Pin Bullet	1/each
219004	Retrograde Reamer Bullet	1/each
219005	Insert for RR Bullet, 2.4 mm	1/each
219006	ACL Tibial Point Aimer	1/each
219007	ACL Tibial Capture Aimer	1/each
219008	ACL/PCL Femoral Point Aimer	1/each
219009	ACL/PCL Femoral Capture Aimer	1/each
219011	PCL Tibial Point Aimer	1/each
219012	PCL Tibial Disc Aimer	1/each
219013	PCL Rasp	1/each
219015	Cannulated PCL Elevator	1/each
219017	Pin Depth Gage	1/each
219018	Outside-In Suture Passing Pin	1/each
215203	Sterilization Tray	1/each
215205	Sterilization Tray Lid	1/each



Knee

↳ Arthroscopic Instruments & Accessories

ACL Reconstruction Instruments

General Instruments

CODE	DESCRIPTION	QTY
219352	Pin Puller	1/each
215105	Screw Insertion Tray	1/each
215106	Silicon Mat-Screw Insertion Tray	1/each
215158	Large Container System	1/each

Guidewires and Guidepins

CODE	DESCRIPTION	QTY
254729	Drill Tip Calibrated Passing Pin, 2.4 mm X 14" (Non-Sterile)	1/each
254514	Round Tip Guidewire, 1.1 mm X 15" (Sterile) (Nitinol)	1/each
214152	Trocar Tip Guidewire, 1.6 mm x14" (Sterile) (Stainless Steel)	1/each
211321	Trocar Tip Guidewire, 1.6mm x 14" (Non-Sterile) (Stainless Steel)	5/pack
214015	Trocar Tip Guidewire, 1.1 mm X 14" (Sterile) (Stainless Steel)	1/each

Tibial Aimer

CODE	DESCRIPTION	QTY
219301	Two Prong Tibial Guide/Ratchet (3 parts, inclusive system)	1/each
219390	Tibial Guide Bullet	1/each
219391	Tibial Guide Rail	1/each
219392	Tibial Guide - Elbow Aiming Arm	1/each
219393	Tibial Guide - Tip Aiming Arm	1/each



Tibial Aimer

Knee

↳ Arthroscopic Instruments & Accessories

ACL Reconstruction Instruments

Tendon Harvester

CODE	DESCRIPTION	QTY
232002	Pigtail Tendon Peeler	1/each
254730	Tendon Stripper, Closed	1/each
232004	Tendon Harvester (used with 232005)	1/each
232005	Sterile Tendon Harvester Disposable Cutting Tip (Used with 232004)	1/each



Additional Guidewires

CODE	DESCRIPTION	QTY
254417	Round Tip Guidewire, 1.1 mm X 9" (Non-Sterile) (Nitinol)	6/pack
211340	Round Tip Guidewire 1.1 mm X 14" (Non-Sterile) (Stainless Steel)	1/each
211332	Round Tip Guidewire, 2.4 mm X 14" (Non-Sterile) (Nitinol)	5/pack
211322	Round Tip Guidewire, 1.6 mm X 16" (Non-Sterile) (Nitinol)	5/pack
211320	Round Tip Guidewire, 1.6 mm X 14" (Non-Sterile) (Stainless Steel)	5/pack
211330	Round Tip Guidewire, 2.4 mm X 14" (Non-Sterile) (Stainless Steel)	5/pack
254843	Round Tip Guidewire, 1.1 mm X 12" (Sterile) (Nitinol)	6/pack
254844	Round Tip Guidewire, 1.1 mm X 12" (Sterile) (Nitinol)	1/each

ACL Disposable Kits

CODE	DESCRIPTION	QTY
232300	ACL Disposable Kit (2.4 mm x 18" Bayonet Tip Passing Pin, 2.4 mm x 14" Drill Tip Passing Pin, 1.1 mm x 15" Nitinol Guidewire, Threaded Cannula, Malleable Graft Retractor, Marking Pen, Ruler)	1/each



Knee

↳ Arthroscopic Instruments & Accessories

Reamers



Fully Fluted Reamers (Reusable)

CODE	DESCRIPTION	QTY
219442	Fully Fluted Reamer, 5.0 mm	1/each
219443	Fully Fluted Reamer, 5.5 mm	1/each
219346	Fully Fluted Reamer, 6.0 mm	1/each
219553	Fully Fluted Reamer, 6.5 mm	1/each
219347	Fully Fluted Reamer, 7.0 mm	1/each
219554	Fully Fluted Reamer, 7.5 mm	1/each
219348	Fully Fluted Reamer, 8.0 mm	1/each
219555	Fully Fluted Reamer, 8.5 mm	1/each
219349	Fully Fluted Reamer, 9.0 mm	1/each
219556	Fully Fluted Reamer, 9.5 mm	1/each
219340	Fully Fluted Reamer, 10.0 mm	1/each
219557	Fully Fluted Reamer, 10.5 mm	1/each
219341	Fully Fluted Reamer, 11.0 mm	1/each
219558	Fully Fluted Reamer, 11.5 mm	1/each
219342	Fully Fluted Reamer, 12.0 mm	1/each
219343	Fully Fluted Reamer, 13.0 mm	1/each
215417	Fluted Reamer Caddy	1/each
215418	Reamer Pin Mat	1/each

Acorn Reamers (Reusable)

CODE	DESCRIPTION	QTY
219440	Acorn Reamer, 5.0 mm	1/each
219441	Acorn Reamer, 5.5 mm	1/each
219336	Acorn Reamer, 6.0 mm	1/each
219547	Acorn Reamer, 6.5 mm	1/each
219337	Acorn Reamer, 7 mm	1/each
219548	Acorn Reamer, 7.5 mm	1/each
219338	Acorn Reamer, 8.0 mm	1/each
219549	Acorn Reamer, 8.5 mm	1/each
219339	Acorn Reamer, 9.0 mm	1/each
219550	Acorn Reamer, 9.5 mm	1/each
219335	Acorn Reamer, 10.0 mm	1/each
219551	Acorn Reamer, 10.5 mm	1/each
219331	Acorn Reamer, 11.0 mm	1/each
219552	Acorn Reamer, 11.5 mm	1/each
219332	Acorn Reamer, 12.0 mm	1/each
219333	Acorn Reamer, 13.0 mm	1/each
215416	Acorn Reamer Caddy	1/each
215418	Reamer Pin Mat	1/each

Knee

↳ Arthroscopic Instruments & Accessories

Reamers



Fully Fluted Reamers (Sterile, Single-Use)

CODE	DESCRIPTION	QTY
232413	Sterile Fully Fluted Reamer, 6 mm	1/each
232414	Sterile Fully Fluted Reamer, 6.5 mm	1/each
232415	Sterile Fully Fluted Reamer, 7 mm	1/each
232416	Sterile Fully Fluted Reamer, 7.5 mm	1/each
232417	Sterile Fully Fluted Reamer, 8 mm	1/each
232418	Sterile Fully Fluted Reamer, 8.5 mm	1/each
232419	Sterile Fully Fluted Reamer, 9 mm	1/each
232420	Sterile Fully Fluted Reamer, 9.5 mm	1/each
232421	Sterile Fully Fluted Reamer, 10 mm	1/each
232422	Sterile Fully Fluted Reamer, 10.5 mm	1/each
232423	Sterile Fully Fluted Reamer, 11 mm	1/each
232424	Sterile Fully Fluted Reamer, 11.5 mm	1/each
232425	Sterile Fully Fluted Reamer, 12 mm	1/each

Acorn Reamers (Sterile, Single-Use)

CODE	DESCRIPTION	QTY
232400	Sterile Acorn Reamer, 6 mm	1/each
232401	Sterile Acorn Reamer, 6.5 mm	1/each
232402	Sterile Acorn Reamer, 7 mm	1/each
232403	Sterile Acorn Reamer, 7.5 mm	1/each
232404	Sterile Acorn Reamer, 8 mm	1/each
232405	Sterile Acorn Reamer, 8.5 mm	1/each
232406	Sterile Acorn Reamer, 9 mm	1/each
232407	Sterile Acorn Reamer, 9.5 mm	1/each
232408	Sterile Acorn Reamer, 10 mm	1/each
232409	Sterile Acorn Reamer, 10.5 mm	1/each
232410	Sterile Acorn Reamer, 11 mm	1/each
232411	Sterile Acorn Reamer, 11.5 mm	1/each
232412	Sterile Acorn Reamer, 12 mm	1/each

Knee

↳ Arthroscopic Instruments & Accessories

Low Profile Reamers (Non-Sterile)












Reusable

CODE	DESCRIPTION	QTY
219624	Low Profile Reamer, 6 mm	1/each
219625	Low Profile Reamer, 6.5 mm	1/each
219626	Low Profile Reamer, 7 mm	1/each
219627	Low Profile Reamer, 7.5 mm	1/each
219628	Low Profile Reamer, 8 mm	1/each
219629	Low Profile Reamer, 8.5 mm	1/each
219630	Low Profile Reamer, 9 mm	1/each
219631	Low Profile Reamer, 9.5 mm	1/each
219632	Low Profile Reamer, 10 mm	1/each
219633	Low Profile Reamer, 10.5 mm	1/each
219634	Low Profile Reamer, 11 mm	1/each
219635	Low Profile Reamer, 11.5 mm	1/each
219636	Low Profile Reamer, 12 mm	1/each
219637	Low Profile Reamer, 13 mm	1/each
215416	Acorn Reamer Caddy	1/each
215418	Reamer Pin Mat	1/each

Knee

↳ Arthroscopic Instruments & Accessories

Notchplasty Instruments

CODE	DESCRIPTION	QTY
		
219303	Osteotome w/ Scale	1/each
219304	Mallet	1/each
		
219305	Arthroscopic Gauge	1/each
		
219306	Ring Curette, 4 mm	1/each
219307	Ring Curette, 7 mm	1/each
		
219308	Cup Curette	1/each
		
219302	Flat Rasp	1/each
		
219309	Curved Rasp	1/each
		
219311	Graft Harvest Gouge	1/each
		
219313	Femoral Awl, 45° Round	1/each
		
219314	Tibial Awl, 90° Round	1/each

Knee

↳ Arthroscopic Instruments & Accessories

Femoral Aimers



Transtibial Femoral Aimer



AM Femoral Aimer

Femoral Aimers

CODE	DESCRIPTION	QTY
219355	Offset Femoral Aimer, 5 mm	1/each
219356	Offset Femoral Aimer, 6 mm	1/each
219357	Offset Femoral Aimer, 7 mm	1/each
219358	Offset Femoral Aimer, 8 mm	1/each

Anteromedial Femoral Aimers

CODE	DESCRIPTION	QTY
219435	AM Femoral Aimer, 4.5 mm	1/each
219434	AM Femoral Aimer, 5 mm	1/each
254672	AM Femoral Aimer, 5.5 mm	1/each
254686	AM Femoral Aimer, 6.5 mm	1/each
254687	AM Femoral Aimer, 7.5 mm	1/each

ACL Accessories

CODE	DESCRIPTION	QTY
219403	Tunnel Dilator, 6 mm	1/each
219404	Tunnel Dilator, 6.5 mm	1/each
219405	Tunnel Dilator, 7 mm	1/each
219406	Tunnel Dilator, 7.5 mm	1/each
219407	Tunnel Dilator, 8 mm	1/each
219408	Tunnel Dilator, 8.5 mm	1/each
219409	Tunnel Dilator, 9 mm	1/each
219410	Tunnel Dilator, 9.5 mm	1/each
219411	Tunnel Dilator, 10 mm	1/each
219412	Tunnel Dilator, 10.5 mm	1/each
219413	Tunnel Dilator, 11 mm	1/each
219414	Tunnel Dilator, 11.5 mm	1/each
219415	Tunnel Dilator, 12 mm	1/each
215419	Tunnel Dilator Tray	1/each
215451	Large Sterilization Tray Lid	1/each
215700	Tunnel Dilator Sterilization Tray	1/each



Tunnel Dilator Tray


Knee

↳ Revision Instruments

Revision Instruments

CODE	DESCRIPTION	QTY
219303	5 mm Osteotome with Millimeter Scale	1/each
219306	Ring Curette, 4 mm	1/each
219307	Ring Curette, 7 mm	1/each
219308	Cup Curette	1/each
219302	Flat Rasp	1/each
219309	Curved Rasp	1/each
219381	Screw-Out Tip, Small Hex	1/each
219304	Large Mallet	1/each
219305	Arthroscopic Gouge	1/each
219311	Graft Harvest Gouge	1/each
219313	Femoral Awl, 45° Round	1/each
219314	Tibial Awl, 90° Round	1/each

Foot & Ankle



Foot & Ankle

↳ Ankle & Tarsals

↳ Suture Anchors

HEALIX ADVANCE™ Anchors with DYNATAPE™ Sutures



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
280723	4.5 mm HEALIX ADVANCE BR Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280735	4.5 mm HEALIX ADVANCE BR Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each
280724	5.5 mm HEALIX ADVANCE BR Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280736	5.5 mm HEALIX ADVANCE BR Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each
280725	6.5 mm HEALIX ADVANCE BR Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280737	6.5 mm HEALIX ADVANCE BR Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each

PEEK

CODE	DESCRIPTION	QTY
280717	4.5 mm HEALIX ADVANCE PEEK Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280729	4.5 mm HEALIX ADVANCE PEEK Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each
280718	5.5 mm HEALIX ADVANCE PEEK Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280730	5.5 mm HEALIX ADVANCE PEEK Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each
280719	6.5 mm HEALIX ADVANCE PEEK Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280731	6.5 mm HEALIX ADVANCE PEEK Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each



Foot & Ankle

↳ Ankle & Tarsals

↳ Suture Anchors

HEALIX ADVANCE™ Anchors with DYNACORD™ Sutures



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
222001	4.5 mm HEALIX ADVANCE BR Anchor with #2 DYNACORD Suture	1/each
222004	5.5 mm HEALIX ADVANCE BR Anchor with #2 DYNACORD Suture	1/each
222011	6.5 mm HEALIX ADVANCE BR Anchor with #2 DYNACORD Suture	1/each
222000	4.5 mm HEALIX ADVANCE BR Anchor with #2 DYNACORD Suture (with needles)	1/each
222006	5.5 mm HEALIX ADVANCE BR Anchor with #2 DYNACORD Suture (with needles)	1/each
222015	6.5 mm HEALIX ADVANCE BR Anchor with #2 DYNACORD Suture (with needles)	1/each
222002	4.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 DYNACORD Suture	1/each
222005	5.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 DYNACORD Suture	1/each
222013	6.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 DYNACORD Suture	1/each

PEEK

CODE	DESCRIPTION	QTY
222028	4.5 mm HEALIX ADVANCE PEEK Anchor with #2 DYNACORD Suture	1/each
222032	5.5 mm HEALIX ADVANCE PEEK Anchor with #2 DYNACORD Suture	1/each
222036	6.5 mm HEALIX ADVANCE PEEK Anchor with #2 DYNACORD Suture	1/each
222027	4.5 mm HEALIX ADVANCE PEEK Anchor with #2 DYNACORD Suture (with needles)	1/each
222034	5.5 mm HEALIX ADVANCE PEEK Anchor with #2 DYNACORD Suture (with needles)	1/each
222038	6.5 mm HEALIX ADVANCE PEEK Anchor with #2 DYNACORD Suture (with needles)	1/each
222029	4.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 DYNACORD Suture	1/each
222033	5.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 DYNACORD Suture	1/each
222037	6.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 DYNACORD Suture	1/each



Foot & Ankle

↳ Ankle & Tarsals

↳ Suture Anchors

HEALIX ADVANCE™ Anchors with ORTHOCORD™ Sutures



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
222295	4.5 mm HEALIX ADVANCE BR Anchor with #2 ORTHOCORD Suture	1/each
222297	5.5 mm HEALIX ADVANCE BR Anchor with #2 ORTHOCORD Suture	1/each
222300	6.5 mm HEALIX ADVANCE BR Anchor with #2 ORTHOCORD Suture	1/each
222319	4.5 mm HEALIX ADVANCE BR Anchor Pack with #2 ORTHOCORD Suture	5/pack
222323	5.5 mm HEALIX ADVANCE BR Anchor Pack with #2 ORTHOCORD Suture	5/pack
222299	5.5 mm HEALIX ADVANCE BR Anchor with #2 ORTHOCORD Suture (with needles)	1/each
222302	6.5 mm HEALIX ADVANCE BR Anchor with #2 ORTHOCORD Suture (with needles)	1/each
222296	4.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 ORTHOCORD Suture	1/each
222298	5.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 ORTHOCORD Suture	1/each
222301	6.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 ORTHOCORD Suture	1/each
222320	4.5 mm HEALIX ADVANCE BR Anchor Pack with Triple Strand #2 ORTHOCORD Suture	5/pack
222324	5.5 mm HEALIX ADVANCE BR Anchor Pack with Triple Strand #2 ORTHOCORD Suture	5/pack

PEEK

CODE	DESCRIPTION	QTY
222303	4.5 mm HEALIX ADVANCE PEEK Anchor with #2 ORTHOCORD Suture	1/each
222305	5.5 mm HEALIX ADVANCE PEEK Anchor with #2 ORTHOCORD Suture	1/each
222309	6.5 mm HEALIX ADVANCE PEEK Anchor with #2 ORTHOCORD Suture	1/each
222321	4.5 mm HEALIX ADVANCE PEEK Anchor Pack with #2 ORTHOCORD Suture	5/pack
222325	5.5 mm HEALIX ADVANCE PEEK Anchor Pack with #2 ORTHOCORD Suture	5/pack
222308	5.5 mm HEALIX ADVANCE PEEK Anchor with #2 ORTHOCORD Suture (with needles)	1/each
222313	6.5 mm HEALIX ADVANCE PEEK Anchor with #2 ORTHOCORD Suture (with needles)	1/each
222304	4.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 ORTHOCORD Suture	1/each
222307	5.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 ORTHOCORD Suture	1/each
222311	6.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 ORTHOCORD Suture	1/each
222322	4.5 mm HEALIX ADVANCE PEEK Anchor Pack with Triple Strand #2 ORTHOCORD Suture	5/pack
222318	5.5 mm HEALIX ADVANCE PEEK Anchor Pack with Triple Strand #2 ORTHOCORD Suture	5/pack



Foot & Ankle

↳ Ankle & Tarsals

↳ Suture Anchors

HEALIX ADVANCE™ Anchors with PERMATAPE™ Sutures



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
222359	4.5 mm HEALIX ADVANCE BR Anchor with PERMATAPE Suture (Blue)	1/each
222360	4.5 mm HEALIX ADVANCE BR Anchor with PERMATAPE Suture (Blue/White)	1/each
222361	5.5 mm HEALIX ADVANCE BR Anchor with PERMATAPE Suture (Blue)	1/each
222362	5.5 mm HEALIX ADVANCE BR Anchor with PERMATAPE Suture (Blue/White)	1/each
222376	6.5 mm HEALIX ADVANCE BR Anchor with PERMATAPE Suture (Blue)	1/each
222377	6.5 mm HEALIX ADVANCE BR Anchor with PERMATAPE Suture (Blue/White)	1/each

PEEK

CODE	DESCRIPTION	QTY
222363	4.5 mm HEALIX ADVANCE PEEK Anchor with PERMATAPE Suture (Blue)	1/each
222364	4.5 mm HEALIX ADVANCE PEEK Anchor with PERMATAPE Suture (Blue/White)	1/each
222365	5.5 mm HEALIX ADVANCE PEEK Anchor with PERMATAPE Suture (Blue)	1/each
222366	5.5 mm HEALIX ADVANCE PEEK Anchor with PERMATAPE Suture (Blue/White)	1/each
222378	6.5 mm HEALIX ADVANCE PEEK Anchor with PERMATAPE Suture (Blue)	1/each
222379	6.5 mm HEALIX ADVANCE PEEK Anchor with PERMATAPE Suture (Blue/White)	1/each

Foot & Ankle

↳ Ankle & Tarsals

↳ Suture Anchors

HEALIX ADVANCE™ Anchors with PERMACORD™ Sutures



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
223128	4.5 mm HEALIX ADVANCE BR Anchor with #2 PERMACORD Suture	1/each
223130	5.5 mm HEALIX ADVANCE BR Anchor with #2 PERMACORD Suture	1/each
223132	6.5 mm HEALIX ADVANCE BR Anchor with #2 PERMACORD Suture	1/each
223129	4.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 PERMACORD Suture	1/each
223131	5.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 PERMACORD Suture	1/each
223133	6.5 mm HEALIX ADVANCE BR Anchor with Triple Strand #2 PERMACORD Suture	1/each

PEEK


CODE	DESCRIPTION	QTY
223134	4.5 mm HEALIX ADVANCE PEEK Anchor with #2 PERMACORD Suture	1/each
223136	5.5 mm HEALIX ADVANCE PEEK Anchor with #2 PERMACORD Suture	1/each
223138	6.5 mm HEALIX ADVANCE PEEK Anchor with #2 PERMACORD Suture	1/each
223135	4.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 PERMACORD Suture	1/each
223137	5.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 PERMACORD Suture	1/each
223139	6.5 mm HEALIX ADVANCE PEEK Anchor with Triple Strand #2 PERMACORD Suture	1/each



Foot & Ankle

- ↳ Ankle & Tarsals
- ↳ Suture Anchors

HEALIX ADVANCE™ Anchor Instruments



HEALIX ADVANCE™ AWL



HEALIX ADVANCE™ AWL-TAP

CODE	DESCRIPTION	QTY
222314	HEALIX ADVANCE Anchor Awl, 4.5 mm	1/each
222315	HEALIX ADVANCE Anchor Awl, 5.5 mm/6.5 mm	1/each
222316	HEALIX ADVANCE Anchor Cortical Awl/Tap, 4.5 mm	1/each
222317	HEALIX Anchor Cortical Awl/Tap, 5.5 mm/6.5 mm	1/each
222329	HEALIX ADVANCE Anchor Full Tap, 4.5 mm	1/each
222328	HEALIX ADVANCE Anchor Full Tap, 5.5 mm/6.5 mm	1/each
222327	HEALIX ADVANCE Anchor Drill Bit, 3.5 mm (DO NOT USE with HEALIX ADVANCE Drill Guide - 222375)	1/each
222326	HEALIX ADVANCE Anchor Drill Bit, 4 mm (DO NOT USE with HEALIX ADVANCE Drill Guide - 222375)	1/each
222358	HEALIX ADVANCE Hard Bone Kit (3.5 mm Drill Bit, Free Needle Loop, 4.5 mm Anchor Tap)	1/each
222375	HEALIX ADVANCE Open Drill Guide*	1/each

*The HEALIX ADVANCE™ Drill Guide can only be used with the 3.5 mm HEALIX ADVANCE™ drill bit of the Hard Bone Kit.

Foot & Ankle

- ↳ Ankle & Tarsals
- ↳ Suture Anchors

HEALIX TI™ Anchors with DYNATAPE™ Suture



Titanium

CODE	DESCRIPTION	QTY
280807	4.5 mm HEALIX TI Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280813	4.5 mm HEALIX TI Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each
280808	5.5 mm HEALIX TI Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280814	5.5 mm HEALIX TI Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each
280809	6.5 mm HEALIX TI Anchor with 2.5 mm DYNATAPE Suture (White)	1/each
280815	6.5 mm HEALIX TI Anchor with 2.5 mm DYNATAPE Suture (Black)	1/each

HEALIX TI™ Anchors with DYNACORD™ Suture



Titanium

CODE	DESCRIPTION	QTY
222051	4.5 mm HEALIX TI Anchor with #2 DYNACORD Suture	1/each
222054	5.5 mm HEALIX TI Anchor with #2 DYNACORD Suture	1/each
222058	6.5 mm HEALIX TI Anchor with #2 DYNACORD Suture	1/each
222055	5.5 mm HEALIX TI Anchor with #2 DYNACORD Suture (with needles)	1/each
222059	6.5 mm HEALIX TI Anchor with #2 DYNACORD Suture (with needles)	1/each
222052	4.5 mm HEALIX TI Anchor with Triple Strand #2 DYNACORD Suture	1/each
222056	5.5 mm HEALIX TI Anchor with Triple Strand #2 DYNACORD Suture	1/each
222061	6.5 mm HEALIX TI Anchor with Triple Strand #2 DYNACORD Suture	1/each



Foot & Ankle

↳ Ankle & Tarsals

↳ Suture Anchors

HEALIX TI™ Anchors with ORTHOCORD™ Suture



Titanium

CODE	DESCRIPTION	QTY
222253	4.5 mm HEALIX TI Anchor with #2 ORTHOCORD Suture	1/each
222245	5.5 mm HEALIX TI Anchor with #2 ORTHOCORD Suture	1/each
222249	6.5 mm HEALIX TI Anchor with #2 ORTHOCORD Suture	1/each
222242	5.5 mm HEALIX TI Anchor with #2 ORTHOCORD Suture (with needles)	1/each
222248	6.5 mm HEALIX TI Anchor with #2 ORTHOCORD Suture (with needles)	1/each
222271	4.5 mm HEALIX TI Anchor with Triple Strand #2 ORTHOCORD Suture	1/each
222243	5.5 mm HEALIX TI Anchor with Triple Strand #2 ORTHOCORD Suture	1/each
222241	6.5 mm HEALIX TI Anchor with Triple Strand #2 ORTHOCORD Suture	1/each

HEALIX TI™ Anchors with PERMACORD™ Suture



Titanium

CODE	DESCRIPTION	QTY
223140	4.5 mm HEALIX TI Anchor with #2 PERMACORD Suture	1/each
223150	5.5 mm HEALIX TI Anchor with #2 PERMACORD Suture	1/each
223144	6.5 mm HEALIX TI Anchor with #2 PERMACORD Suture	1/each
223141	4.5 mm HEALIX TI Anchor with Triple Strand #2 PERMACORD Suture	1/each
223151	5.5 mm HEALIX TI Anchor with Triple Strand #2 PERMACORD Suture	1/each
223145	6.5 mm HEALIX TI Anchor with Triple Strand #2 PERMACORD Suture	1/each



Foot & Ankle

↳ Ankle & Tarsals

↳ Suture Anchors

HEALIX TRANSTEND™ Implant System



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
222265	HEALIX TRANSTEND BR Anchor with ORTHOCORD Suture 3.4 mm	1/each
222261	HEALIX TRANSTEND BR Implant System (1 BR Anchor, PERCANNULA System)	1/each
222263	HEALIX TRANSTEND BR Implant System (2 BR Anchors, PERCANNULA System)	1/each
222338	HEALIX TRANSTEND BR Implant System (1 BR, Violet, 1 BR, Blue, PERCANNULA System)	1/each
222337	HEALIX TRANSTEND BR Dual Strand Implant System (1 BR DS, PERCANNULA System)	1/each

PEEK

CODE	DESCRIPTION	QTY
222276	HEALIX TRANSTEND PEEK Anchor with ORTHOCORD Suture 3.4 mm	1/each
222278	HEALIX TRANSTEND PEEK Implant System (1 PEEK Anchor, PERCANNULA System)	1/each
222280	HEALIX TRANSTEND PEEK Implant System (2 PEEK Anchors, PERCANNULA System)	1/each
222340	HEALIX TRANSTEND PEEK Implant System (1 PEEK, Violet, 1 PEEK, Blue, PERCANNULA)	1/each
222339	HEALIX TRANSTEND PEEK Dual Strand Implant System (1 PEEK DS, PERCANNULA System)	1/each

Titanium

CODE	DESCRIPTION	QTY
222260	HEALIX TRANSTEND TI Anchor with ORTHOCORD Suture 2.9 mm	1/each
222266	HEALIX TRANSTEND TI Implant System (1 TI Anchor, PERCANNULA System)	1/each
222268	HEALIX TRANSTEND TI Implant System (2 TI Anchors, PERCANNULA System)	1/each

Instruments

CODE	DESCRIPTION	QTY
222274	HEALIX PERCANNULA System 4.0 mm	1/each
222273	HEALIX TRANSTEND Anchor Awl/Tap, 3.4 mm (for BR and PEEK Anchors)	1/each
222281	HEALIX TRANSTEND Anchor Drill Bit, 2.2 mm (optional for TI Anchor)	1/each



Foot & Ankle

↳ Ankle & Tarsals

↳ All-Suture Anchors

VERSALOOPTM All-Suture Anchors

Anchors

CODE	DESCRIPTION	QTY
210165	1.5 mm VERSALOOPTM Anchor with a Size #2 Suture	1/each
210166	1.5 mm VERSALOOPTM Anchor with a 1.5 mm Tape	1/each
210167	1.8 mm VERSALOOPTM Anchor with 2 Size #2 Sutures	1/each
210168	1.8 mm VERSALOOPTM Anchor with 2 1.5 mm Tapes	1/each



Obturator

CODE	DESCRIPTION	QTY
210205	VERSALOOPTM Obturator for 1.5 mm Anchor	1/each
210213	VERSALOOPTM Obturator for 1.8 mm Anchor	1/each

Curved Instruments

CODE	DESCRIPTION	QTY
210195	VERSALOOPTM 12.5° Guide for 1.5 mm Anchor, Sawtooth	1/each
210196	VERSALOOPTM 12.5° Guide for 1.8 mm Anchor, Sawtooth	1/each
210197	VERSALOOPTM 25° Guide for 1.5 mm Anchor, Sawtooth	1/each
210198	VERSALOOPTM 25° Guide for 1.8 mm Anchor, Sawtooth	1/each
210199	VERSALOOPTM 1.5 mm Disposable Drill	1/each
210203	VERSALOOPTM 1.8 mm Disposable Drill	1/each



1.8 mm Family of Instruments

Straight Instruments

CODE	DESCRIPTION	QTY
210175	VERSALOOPTM Straight Guide for 1.5 mm Anchor, Fishmouth	1/each
210176	VERSALOOPTM Straight Guide for 1.8 mm Anchor, Fishmouth	1/each
210178	VERSALOOPTM Straight Guide for 1.5 mm Anchor, Sawtooth	1/each
210179	VERSALOOPTM Straight Guide for 1.8 mm Anchor, Sawtooth	1/each
210181	VERSALOOPTM 1.5 mm Reusable Straight Drill	1/each
210182	VERSALOOPTM 1.8 mm Reusable Straight Drill	1/each



1.5 mm 25° guide, 1.8 mm 12.5° guide & 2.5 mm Straight Guide



Foot & Ankle

↳ Ankle & Tarsals

↳ Knotless Anchors

HEALIX ADVANCE™ Self-Punching (SP) Anchors



BIOCRYL™ Biocomposite Material

CODE	DESCRIPTION	QTY
222422	4.9 mm HEALIX ADVANCE Self-Punching Biocomposite Anchor	1/each
222424	5.5 mm HEALIX ADVANCE Self-Punching Biocomposite Anchor	1/each
222426	6.5 mm HEALIX ADVANCE Self-Punching Biocomposite Anchor	1/each

PEEK

CODE	DESCRIPTION	QTY
222423	4.9 mm HEALIX ADVANCE Self-Punching PEEK Anchor	1/each
222425	5.5 mm HEALIX ADVANCE Self-Punching PEEK Anchor	1/each
222427	6.5 mm HEALIX ADVANCE Self-Punching PEEK Anchor	1/each

HEALIX ADVANCE™ Knotless Anchors



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
222581	4.75 mm HEALIX ADVANCE Knotless BR Anchor (Double Kite)	1/each
222582	5.5 mm HEALIX ADVANCE Knotless BR Anchor (Double Kite)	1/each
222583	6.5 mm HEALIX ADVANCE Knotless BR Anchor (Double Kite)	1/each

PEEK

CODE	DESCRIPTION	QTY
222584	4.75 mm HEALIX ADVANCE Knotless PEEK Anchor (Double Kite)	1/each
222585	5.5 mm HEALIX ADVANCE Knotless PEEK Anchor (Double Kite)	1/each
222586	6.5 mm HEALIX ADVANCE Knotless PEEK Anchor (Double Kite)	1/each



Foot & Ankle

↳ Metatarsals & Phalanges

↳ Suture Anchors

GRYPHON™ Anchors with DYNACORD™ Suture



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
210120	GRYPHON-P BR Anchor with #2 DYNACORD Suture	1/each
210121	GRYPHON-P BR Anchor with Dual Strand #2 DYNACORD Suture	1/each

PEEK

CODE	DESCRIPTION	QTY
210122	GRYPHON-P PEEK Anchor with #2 DYNACORD Suture	1/each
210123	GRYPHON-P PEEK Anchor with Dual Strand #2 DYNACORD Suture	1/each

GRYPHON™ Anchors with ORTHOCORD™ Suture

BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
210814	GRYPHON-P BR Anchor with #2 ORTHOCORD Suture	1/each
210813	GRYPHON-P BR Anchor with Dual Strand #2 ORTHOCORD Suture	1/each
210812	GRYPHON-T BR Anchor with #2 ORTHOCORD Suture	1/each
210811	GRYPHON-T BR Anchor with Dual Strand #2 ORTHOCORD Suture	1/each

PEEK

CODE	DESCRIPTION	QTY
210031	GRYPHON-P PEEK Anchor with #2 ORTHOCORD Suture	1/each
210030	GRYPHON-P PEEK Anchor with Dual Strand #2 ORTHOCORD Suture	1/each



Foot & Ankle

↳ Metatarsals & Phalanges

↳ Suture Anchors

GRYPHON™ Anchors with PERMACORD™ Suture



BIOCRYL RAPIDE™ Biocomposite Material

CODE	DESCRIPTION	QTY
223146	GRYPHON-P BR Anchor with #2 PERMACORD Suture	1/each
223147	GRYPHON-P BR Anchor with Dual Strand #2 PERMACORD Suture	1/each

PEEK

CODE	DESCRIPTION	QTY
223148	GRYPHON-P PEEK Anchor with #2 PERMACORD Suture	1/each
223149	GRYPHON-P PEEK Anchor with Dual Strand #2 PERMACORD Suture	1/each

GRYPHON™ Instruments



Instruments

CODE	DESCRIPTION	QTY
211035	GRYPHON Anchor Tap (for use with Threaded Anchor only)	1/each
211036	GRYPHON Anchor Drill Bit, 2.4 mm	1/each
210824	GRYPHON PROKNOT Anchor Slotted Sawtooth Guide	1/each
211039	GRYPHON Anchor Sawtooth Guide	1/each
210826	GRYPHON PROKNOT Anchor Slotted Fishmouth Guide	1/each
211037	GRYPHON Anchor Fishmouth Guide	1/each
211038	GRYPHON Anchor Sharp Obturator	1/each
211041	GRYPHON Anchor Blunt Obturator	1/each
210835	Knot Manipulator	1/each
211601	GRYPHON Anchor Hard Bone Drill Bit, 2.5 mm	1/each



Foot & Ankle

↳ Metatarsals & Phalanges

↳ Suture Anchors

Super QUICKANCHOR™ Suture Anchor



Anchor

CODE	DESCRIPTION	QTY
222984	Super QUICKANCHOR Plus with ORTHOCORD Suture (CP-2 Needles)	1/each

Instruments

CODE	DESCRIPTION	QTY
211113	SUPER/RC Drill Bit (Open), 2.9 mm	1/each
213207	GII/SUPER/RC Slotted Drill Guide (Open)	1/each
211513	Arthroscopic SUPER/RC Drill Guide, Arthroscopic Drill, 2.9 mm	1/each

GII QUICKANCHOR™ Suture Anchor



Anchor

CODE	DESCRIPTION	QTY
222983	GII QUICKANCHOR Suture Anchor with #2 ORTHOCORD Suture (with CP-2 needles)	1/each

Instruments



CODE	DESCRIPTION	QTY
211093	GII Drill Bit, 2.4 mm	1/each
211502	Arthroscopic GII Drill Bit	1/each
213207	GII/SUPER/RC/GLS Slotted Drill Guide	1/each
214561	Arthroscopic Sharkmouth Drill Guide	1/each
214562	Arthroscopic Sawtooth Drill Guide	1/each
212094	GII Inserter	1/each
214542	Arthroscopic T-Handle Reamer	1/each



Foot & Ankle

↳ Metatarsals & Phalanges

↳ Suture Anchors

MINI QUICKANCHOR™ Suture Anchors



Anchors

CODE	DESCRIPTION	QTY
212035	MINI QUICKANCHOR PLUS with #2/0 ORTHOCORD Suture (V-5) Drill Bit, 2 mm	1/each
212038	MINI QUICKANCHOR PLUS with #0 ORTHOCORD Suture (OS-2) Drill Bit, 2 mm	1/each



Instruments

CODE	DESCRIPTION	QTY
211083	MINI Drill Bit, 2.1 mm	1/each
213072	MINI Drill Guide	1/each

MINILOK™ QUICKANCHOR™ Suture Anchors








CODE	DESCRIPTION	QTY
212849	MINILOK QUICKANCHOR with #0 ORTHOCORD Suture (OS-2) with Drill Bit	1/each
212853	MINILOK QUICKANCHOR with #2/0 ORTHOCORD Suture (V-5) with Drill Bit	1/each
212854	MINILOK QUICKANCHOR with #2/0 ORTHOCORD Suture (RB-1) with Drill Bit	1/each








Foot & Ankle


↳ Sutures & Suture Passers

Needle Guide


27 mm	13 mm	26 mm	27 mm	13 mm
				
CP-2 Needle (Reverse Cutting)	P-3 Needle (Reverse Cutting)	SH Needle (Taper Point)	CT-2 Needle (Taper Point)	C-1 Needle (Taper Point)

3/8 circle 17 mm	1/2 circle 17 mm	1/2 circle 17 mm	1/2 circle 22 mm	1/2 circle 17 mm
				
V-4 Needle (TAPERCUT™)	V-5 Needle (TAPERCUT™)	OS-6 36 mm (Reverse Cutting)	MO-6 27 mm (Taper Point)	RB-1 Needle (Taper Point)


Taper Point
Designed for soft, easy penetrated tissue.



Reverse Cutting
Cutting edge on outer curve. Designed for tough, difficult-to-penetrate tissue.



TAPERCUT™ Technology
Cutting tip, taper body designed for tough tissue, like two needles in one.



Free Needles



CODE	DESCRIPTION	QTY
223173	Free Needle Loop (OS-6)	12/each



Foot & Ankle

↳ Sutures & Suture Passers

DYNATAPE™ Sutures



CODE	DESCRIPTION	QTY
280825	2.5 mm DYNATAPE Suture Pack, White, without needles	6/pack
280828	2.5 mm DYNATAPE Suture Pack, Black, without needles	6/pack
223173	Free Needle Loop OS-6	12/pack

DYNACORD™ Sutures



CODE	DESCRIPTION	QTY
222067	Size #2 DYNACORD Suture Pack, Blue, without Needle	12/pack
222065	Size #2 DYNACORD Suture Pack, Blue, with OS-6 Needle	12/pack
222066	Size #2 DYNACORD Suture Pack, Blue, with MO-7 Needle	12/pack
222068	Size #2 DYNACORD Suture Pack, Striped, without Needle	12/pack
222069	Size #2 DYNACORD Suture Pack, Striped/Blue, without Needle	12/pack
222071	Size #2 DYNACORD Suture Pack, Striped/Blue, with MO-7 Needle	12/pack
222073	Size #2 DYNACORD Suture Pack, Striped/Blue, with OS-6 Needle	12/pack

Foot & Ankle

↳ Sutures & Suture Passers

ORTHOCORD™ Sutures



Size #2 Sutures

CODE	DESCRIPTION	QTY
223102	ORTHOCORD Suture #2, Violet, with MO-6 1/2 Circle, Taper Point Needle, 26 mm	12/pack
223103	ORTHOCORD Suture #2, Violet, with OS-6 1/2 Circle, Reverse Cutting Needle, 36 mm	12/pack
223104	ORTHOCORD Suture #2, Violet, with MO-7 1/2 Circle, Taper Point Needle, 22 mm	12/pack
223105	ORTHOCORD Suture #2, Violet, without Needles	12/pack
223111	ORTHOCORD Suture #2, Blue, without Needles	12/pack
223113	ORTHOCORD Suture #2, Violet/Blue, without Needles	12/pack
223114	ORTHOCORD Suture #2, Violet/Blue, with MO-7 1/2 Circle, Taper Point Needle, 22 mm	12/pack
223115	ORTHOCORD Suture #2, Violet/Blue, with MO-6 1/2 Circle, Taper Point Needle, 26 mm	12/pack
223116	ORTHOCORD Suture #2, Violet/Blue, with OS-6 1/2 Circle, Reverse Cutting Needle, 36 mm	12/pack

Small Size Sutures

CODE	DESCRIPTION	QTY
223121	#4-0 ORTHOCORD Suture Loop 51cm Taper Point, 3/8 Circle, 17.9 mm	12/pack
223122	#2-0 ORTHOCORD Suture 46cm Taper Point, 1/2 Circle, 26.5 mm	12/pack
223123	#4-0 ORTHOCORD Suture 46cm Diamond Point, 3/8 Circle, 18.7 mm	12/pack
223124	#0 ORTHOCORD Suture 96.5cm Taper Point, 1/2 Circle, 22.2 mm	12/pack
223125	#3-0 ORTHOCORD Suture 46cm Taper Point, 3/8 Circle, 15 mm	12/pack
223126	#2-0 ORTHOCORD Suture 46cm Taper Point, 3/8 Circle, 17.9 mm	12/pack
223127	#4-0 ORTHOCORD Suture 46cm Taper Point, 3/8 Circle, 2.3 mm	12/pack



Foot & Ankle

↳ Sutures & Suture Passers


2.5 mm PERMATAPE™ Sutures



CODE	DESCRIPTION	QTY
223165	2.5 mm PERMATAPE Suture, Blue Pack	6/pack
223166	2.5 mm PERMATAPE Suture, White Pack	6/pack
223169	2.5 mm PERMATAPE Suture, Blue	1/each
223170	2.5 mm PERMATAPE Suture, White	1/each

1.3 mm PERMATAPE™ Sutures

CODE	DESCRIPTION	QTY
280772	1.3 mm PERMATAPE Suture, Blue, without Needles	12/pack
280773	1.3 mm PERMATAPE Suture, Blue, OS6 Needles	12/pack
280774	1.3 mm PERMATAPE Suture, Blue, MO-7 Needles	12/pack
280775	1.3 mm PERMATAPE Suture, White-Blue, without Needles	12/pack
280776	1.3 mm PERMATAPE Suture, White-Blue, OS6 Needles	12/pack
280777	1.3 mm PERMATAPE Suture, White-Blue, MO-7 Needles	12/pack



PERMALOOP™ Sutures

CODE	DESCRIPTION	QTY
280936	1.3 mm PERMALOOP Suture Straight Needle	12/pack
208938	1.3 mm PERMALOOP Suture MO-5 Needle	12/pack

Foot & Ankle

↳ Sutures & Suture Passers

IDEAL™ Suture Graspers & Shuttles



IDEAL Suture Graspers

CODE	DESCRIPTION	QTY
251720	IDEAL Suture Grasper, 15°	1/each
251721	IDEAL Suture Grasper, 30°	1/each
251722	IDEAL Suture Grasper, 45°	1/each
251723	IDEAL Suture Grasper, 60°	1/each

IDEAL Suture Shuttles

CODE	DESCRIPTION	QTY
251001	IDEAL Suture Shuttle, 25° Left	1/each
251002	IDEAL Suture Shuttle, 25° Right	1/each
251003	IDEAL Suture Shuttle, 45° Left	1/each
251004	IDEAL Suture Shuttle, 45° Right	1/each
251005	IDEAL Suture Shuttle, 90° Up	1/each
251006	IDEAL Suture Shuttle, Straight Hook	1/each
251007	IDEAL Suture Shuttle, Straight Crescent	1/each



CHIA PERCPASSER™ Suture Passer

CODE	DESCRIPTION	QTY
214101	CHIA PERCPASSER Suture Passer	5/each




Foot & Ankle

↳ Graft Preparation

SPEEDTRAP™ Graft Preparation System

CODE	DESCRIPTION	QTY
223746	SPEEDTRAP System, White, 20 mm	1/each
223751	SPEEDTRAP System, Green/White, 20 mm	1/each
223747	SPEEDTRAP System, White, 30 mm	1/each
223748	SPEEDTRAP System, Green/White, 30 mm	1/each
223749	SPEEDTRAP System Pack, White, 30 mm	4/each
223750	SPEEDTRAP System Pack, Green/White, 30 mm	4/each




Foot & Ankle


↳ Arthroscopic Instruments & Accessories

STRADA™ Small Joint Instruments

CODE	DESCRIPTION	QTY
224200	Grasper, 2.5 mm x 90 mm	1/each
224201	Straight Punch Blunt, 2.1 mm x 90 mm	1/each
224202	Right Punch Blunt, 2.1 mm x 90 mm	1/each
224203	Left Punch Blunt, 2.1 mm x 90 mm	1/each
224204	Down Punch Blunt, 2.1 mm x 90 mm	1/each
224205	Large Punch Blunt, 2.7 mm x 90 mm	1/each
224206	Straight Probe Marked, 1 mm x 65 mm	1/each
224207	30° Probe Marked, 1 mm x 65 mm	1/each
224208	90° Probe Marked, 1 mm x 65 mm	1/each
224209	Scissors, 2.5 mm x 90 mm	1/each
224210	Small Joint Shear Pin	1/each
215201	Small Joint Arthroscopy Tray	1/each



Hand, Wrist & Elbow



Hand, Wrist & Elbow

- ↳ Carpals & Phalanges
- ↳ Suture Anchors

MINI QUICKANCHOR™ Suture Anchors



Anchors


CODE	DESCRIPTION	QTY
212035	MINI QUICKANCHOR PLUS with #2/0 ORTHOCORD Suture (V-5) Drill Bit, 2 mm	1/each
212038	MINI QUICKANCHOR PLUS with #0 ORTHOCORD Suture (OS-2) Drill Bit, 2 mm	1/each



Instruments

CODE	DESCRIPTION	QTY
211083	MINI Drill Bit, 2.1 mm	1/each
213072	MINI Drill Guide	1/each

MICRO QUICKANCHOR™ Suture Anchors




CODE	DESCRIPTION	QTY
212843	MICRO QUICKANCHOR PLUS with #3/0 ORTHOCORD Suture (V-4) Drill Bit, 1.3 mm	1/each
212041	MICRO QUICKANCHOR PLUS with #4/0 ORTHOCORD Suture (C-1) Drill Bit, 1.3 mm	1/each
212042	MICRO QUICKANCHOR PLUS with #4/0 ORTHOCORD Suture (P-3) Drill Bit, 1.3 mm	1/each

Hand, Wrist & Elbow


- ↳ Carpals & Phalanges
- ↳ Suture Anchors

MICROFIX™ QUICKANCHOR™ Suture Anchors



CODE	DESCRIPTION	QTY
212859	MICROFIX Anchor with #3/0 ORTHOCORD Suture (V-4) with Drill Bit	1/each
212045	MICROFIX Anchor with #4/0 ORTHOCORD Suture (J-1) with Drill Bit	1/each
212044	MICROFIX Anchor with #4/0 ORTHOCORD Suture (C-1) with Drill Bit	1/each

MINILOK™ QUICKANCHOR™ Suture Anchors



CODE	DESCRIPTION	QTY
212849	MINILOK Anchor with #0 ORTHOCORD Suture (OS-2) with Drill Bit	1/each
212853	MINILOK Anchor with #2/0 ORTHOCORD Suture (V-5) with Drill Bit	1/each
212854	MINILOK Anchor with #2/0 ORTHOCORD Suture (RB-1) with Drill Bit	1/each

GII QUICKANCHOR™ Suture Anchor



Anchors

CODE	DESCRIPTION	QTY
222983	GII QUICKANCHOR Suture Anchor with #2 ORTHOCORD Suture (with CP-2 needles)	1/each



Hand, Wrist & Elbow

- ↳ Carpals & Phalanges
- ↳ Suture Anchors

GII QUICKANCHOR™ Suture Anchor

Instruments

CODE	DESCRIPTION	QTY
211093	GII Drill Bit, 2.4 mm	1/each
211502	Arthroscopic GII Drill Bit	1/each
213207	GII/SUPER/RC/GLS Slotted Drill Guide	1/each
214561	Arthroscopic Sharkmouth Drill Guide	1/each
214562	Arthroscopic Sawtooth Drill Guide	1/each
212094	GII Inserter	1/each
214542	Arthroscopic T-Handle Reamer	1/each



Super QUICKANCHOR™ Suture Anchor

Anchor

CODE	DESCRIPTION	QTY
222984	Super QUICKANCHOR Plus with ORTHOCORD Suture (CP-2 Needles)	1/each



Instruments

CODE	DESCRIPTION	QTY
211113	SUPER/RC Drill Bit (Open), 2.9 mm	1/each
213207	GII/SUPER/RC Slotted Drill Guide (Open)	1/each
211513	Arthroscopic SUPER/RC Drill Guide, Arthroscopic Drill, 2.9 mm	1/each



Hand, Wrist & Elbow

↳ Distal Biceps

↳ Suture Anchors

GII QUICKANCHOR™ Suture Anchor

Anchor

CODE	DESCRIPTION	QTY
222983	GII QUICKANCHOR Suture Anchor with #2 ORTHOCORD Suture (with CP-2 needles)	1/each



Instruments

CODE	DESCRIPTION	QTY
211093	GII Drill Bit, 2.4 mm	1/each
211502	Arthroscopic GII Drill Bit	1/each
213207	GII/SUPER/RC/GLS Slotted Drill Guide	1/each
214561	Arthroscopic Sharkmouth Drill Guide	1/each
214562	Arthroscopic Sawtooth Drill Guide	1/each
212094	GII Inserter	1/each
214542	Arthroscopic T-Handle Reamer	1/each



Super QUICKANCHOR™ Suture Anchor

Anchor

CODE	DESCRIPTION	QTY
222984	Super QUICKANCHOR Plus with ORTHOCORD Suture (CP-2 Needles)	1/each








Instruments






CODE	DESCRIPTION	QTY
211113	SUPER/RC Drill Bit (Open), 2.9 mm	1/each
213207	GII/SUPER/RC Slotted Drill Guide (Open)	1/each
211513	Arthroscopic SUPER/RC Drill Guide, Arthroscopic Drill, 2.9 mm	1/each




Hand, Wrist & Elbow ↳ Sutures

Needle Guide


27 mm	13 mm	26 mm	27 mm	13 mm
				
CP-2 Needle (Reverse Cutting)	P-3 Needle (Reverse Cutting)	SH Needle (Taper Point)	CT-2 Needle (Taper Point)	C-1 Needle (Taper Point)

3/8 circle 17 mm	1/2 circle 17 mm	1/2 circle 17 mm	1/2 circle 22 mm	1/2 circle 17 mm
				
V-4 Needle (TAPERCUT™)	V-5 Needle (TAPERCUT™)	OS-6 36 mm (Reverse Cutting)	MO-6 27 mm (Taper Point)	RB-1 Needle (Taper Point)


Taper Point
Designed for soft, easy
penetrated tissue.



Reverse Cutting
Cutting edge on outer curve.
Designed for tough,
difficult-to-penetrate tissue.



TAPERCUT™ Technology
Cutting tip, taper body designed for
tough tissue, like two needles in one.



Free Needles



CODE	DESCRIPTION	QTY
223173	Free Needle Loop (OS-6)	12/pack



Hand, Wrist & Elbow ↳ Sutures

DYNATAPE™ Sutures



CODE	DESCRIPTION	QTY
280825	2.5 mm DYNATAPE Suture Pack, White, without needles	6/pack
280828	2.5 mm DYNATAPE Suture Pack, Black, without needles	6/pack
223173	Free Needle Loop OS-6	12/pack

DYNACORD™ Sutures



CODE	DESCRIPTION	QTY
222067	Size #2 DYNACORD Suture Pack, Blue, without Needle	12/pack
222065	Size #2 DYNACORD Suture Pack, Blue, with OS-6 Needle	12/pack
222066	Size #2 DYNACORD Suture Pack, Blue, with MO-7 Needle	12/pack
222068	Size #2 DYNACORD Suture Pack, Striped, without Needle	12/pack
222069	Size #2 DYNACORD Suture Pack, Striped/Blue, without Needle	12/pack
222071	Size #2 DYNACORD Suture Pack, Striped/Blue, with MO-7 Needle	12/pack
222073	Size #2 DYNACORD Suture Pack, Striped/Blue, with OS-6 Needle	12/pack

Hand, Wrist & Elbow ↳ Sutures

ORTHOCORD™ Sutures



Size #2 Sutures

CODE	DESCRIPTION	QTY
223102	ORTHOCORD Suture #2, Violet, with MO-6 1/2 Circle, Taper Point Needle, 26 mm	12/pack
223103	ORTHOCORD Suture #2, Violet, with OS-6 1/2 Circle, Reverse Cutting Needle, 36 mm	12/pack
223104	ORTHOCORD Suture #2, Violet, with MO-7 1/2 Circle, Taper Point Needle, 22 mm	12/pack
223105	ORTHOCORD Suture #2, Violet, without Needles	12/pack
223111	ORTHOCORD Suture #2, Blue, without Needles	12/pack
223113	ORTHOCORD Suture #2, Violet/Blue, without Needles	12/pack
223114	ORTHOCORD Suture #2, Violet/Blue, with MO-7 1/2 Circle, Taper Point Needle, 22 mm	12/pack
223115	ORTHOCORD Suture #2, Violet/Blue, with MO-6 1/2 Circle, Taper Point Needle, 26 mm	12/pack
223116	ORTHOCORD Suture #2, Violet/Blue, with OS-6 1/2 Circle, Reverse Cutting Needle, 36 mm	12/pack

Small Size Sutures

CODE	DESCRIPTION	QTY
223121	#4-0 ORTHOCORD Suture Loop 51cm Taper Point, 3/8 Circle, 17.9 mm	12/pack
223122	#2-0 ORTHOCORD Suture 46cm Taper Point, 1/2 Circle, 26.5 mm	12/pack
223123	#4-0 ORTHOCORD Suture 46cm Diamond Point, 3/8 Circle, 18.7 mm	12/pack
223124	#0 ORTHOCORD Suture 96.5cm Taper Point, 1/2 Circle, 22.2 mm	12/pack
223125	#3-0 ORTHOCORD Suture 46cm Taper Point, 3/8 Circle, 15 mm	12/pack
223126	#2-0 ORTHOCORD Suture 46cm Taper Point, 3/8 Circle, 17.9 mm	12/pack
223127	#4-0 ORTHOCORD Suture 46cm Taper Point, 3/8 Circle, 2.3 mm	12/pack



Hand, Wrist & Elbow ↳ Sutures


2.5 mm PERMATAPE™ Sutures



CODE	DESCRIPTION	QTY
223165	2.5 mm PERMATAPE Suture, Blue Pack	6/pack
223166	2.5 mm PERMATAPE Suture, White Pack	6/pack
223169	2.5 mm PERMATAPE Suture, Blue	1/each
223170	2.5 mm PERMATAPE Suture, White	1/each

1.3 mm PERMATAPE™ Sutures

CODE	DESCRIPTION	QTY
280772	1.3 mm PERMATAPE Suture, Blue, without Needles	12/pack
280773	1.3 mm PERMATAPE Suture, Blue, OS6 Needles	12/pack
280774	1.3 mm PERMATAPE Suture, Blue, MO-7 Needles	12/pack
280775	1.3 mm PERMATAPE Suture, White-Blue, without Needles	12/pack
280776	1.3 mm PERMATAPE Suture, White-Blue, OS6 Needles	12/pack
280777	1.3 mm PERMATAPE Suture, White-Blue, MO-7 Needles	12/pack



PERMALOOP™ Sutures

CODE	DESCRIPTION	QTY
280936	1.3 mm PERMALOOP Suture Straight Needle	12/pack
208938	1.3 mm PERMALOOP Suture MO-5 Needle	12/pack




Hand, Wrist & Elbow


↳ Arthroscopic Instruments & Accessories

STRADA™ Small Joint Instruments

CODE	DESCRIPTION	QTY
224200	Grasper, 2.5 mm x 90 mm	1/each
224201	Straight Punch Blunt, 2.1 mm x 90 mm	1/each
224202	Right Punch Blunt, 2.1 mm x 90 mm	1/each
224203	Left Punch Blunt, 2.1 mm x 90 mm	1/each
224204	Down Punch Blunt, 2.1 mm x 90 mm	1/each
224205	Large Punch Blunt, 2.7 mm x 90 mm	1/each
224206	Straight Probe Marked, 1 mm x 65 mm	1/each
224207	30° Probe Marked, 1 mm x 65 mm	1/each
224208	90° Probe Marked, 1 mm x 65 mm	1/each
224209	Scissors, 2.5 mm x 90 mm	1/each
224210	Small Joint Shear Pin	1/each
215201	Small Joint Arthroscopy Tray	1/each



Hip




Hip ↳ Labrum ↳ Anchors

VERSALOOPTM All-Suture Anchors

Anchors

CODE	DESCRIPTION	QTY
210215	1.5 mm VERSALOOPTM Anchor XL with a Size #2 Suture	1/each
210216	1.5 mm VERSALOOPTM Anchor XL with a 1.5 mm Tape	1/each



Instruments

CODE	DESCRIPTION	QTY
210217	VERSALOOPTM XL Straight Guide for 1.5 mm Anchor XL, Sawtooth	1/each
210218	VERSALOOPTM XL 12.5° Guide for 1.5 mm Anchor XL, Sawtooth	1/each
210219	VERSALOOPTM XL 1.5 mm Reusable Straight Drill	1/each
210221	VERSALOOPTM XL 1.5 mm Disposable Drill	1/each
210222	VERSALOOPTM XL Cannulated Obturator for 1.5 mm Anchor XL	1/each



GRYPHONTM Suture Anchors

BIOCRYL RAPIDETM Biocomposite Material

CODE	DESCRIPTION	QTY
210811	GRYPHON Threaded BR with Dual Strand #2 ORTHOCORD Sutures	1/each
210812	GRYPHON Threaded BR with #2 ORTHOCORD Suture	1/each
210813	GRYPHON Push-In BR with Dual Strand #2 ORTHOCORD Sutures	1/each
210814	GRYPHON Push-In BR with ORTHOCORD Suture	1/each



Instruments

CODE	DESCRIPTION	QTY
213042	GRYPHON Hip Guide Sawtooth Drill Guide	1/each
213044	GRYPHON Hip Guide Fishmouth Drill Guide	1/each
219102	GRYPHON Anchor Hip Obturator	1/each
219103	GRYPHON Anchor Slotted Hip Guide	1/each
211036	GRYPHON Anchor Drill Bit, 2.4 mm	1/each



Hip ↳ Arthroscopic Instruments & Accessories

Hip Procedural Instrument System

CODE	DESCRIPTION	QTY
214009	Hip Suture Passer Up	1/each
214013	Hip Suture Passer Straight	1/each
214111	7 x 110 mm Reusable Obturator	1/each
214131	Hip Probe 90°	1/each
214132	5.5 Hip Cannulated SwitchStick	1/each
214133	6.0 Hip Cannulated SwitchStick	1/each
214134	4.0 Hip Cannulated SwitchStick	1/each
214135	5.5 Hip Slotted Cannula	1/each
214136	6.0 Hip Slotted Cannula	1/each
214137	Hip Biter 15° Up	1/each
214138	Hip CordCutter	1/each
214139	Hip Biter Straight	1/each
214142	Hip Curette 20°	1/each
214143	Hip Grasper 20° Left	1/each
214144	Hip Grasper 20° Right	1/each
214145	Hip Knot Manipulator	1/each
214146	Hip Suture Manipulator	1/each
214147	Hip Grasper Straight	1/each
214148	Hip TissLib Blade 20° down	1/each




CODE	DESCRIPTION	QTY
214151	Handle, Hip SwitchStick	1/each
214153	5.0 Hip Cannulated SwitchStick	1/each
214154	4.5 Hip Cannulated SwitchStick	1/each
215504	Hip Slotted Cannula Tray	1/each
215505	Hip Tissue & Suture Management Tray	1/each
214109	5.5 mm x 90 mm Reusable Obturator	1/each
214121	8.5 mm x 90 mm Reusable Obturator	1/each
213705	RIGIDFIX™ System Mallet	1/each
215506	Hip Entry Tray Pin Mat Accessory	1/each
215507	Hip Entry & Tissue Removal Tray	1/each
215508	Hip Labral Repair Tray	1/each
215509	Hip Sterilization Case	1/each
215510	Lid, Hip Sterilization Case	1/each

Hip Access Kits

CODE	DESCRIPTION
260012	Hip Access Kit No Blade (17 g x 6" Spinal Needle, 17 g x 7" Spinal Needle, 0.049" x 17.5" Round Tip Guidewire, Marking Pen, Ruler, 30 mL Syringe)
260009	Hip Access Kit Curved Blade (Curved Blade, 16 g x 7" Spinal Needle, 16 g x 6" Spinal Needle, Two 1.25 mm x 17.5" Round Tip Guidewires, Marking Pen, Ruler, 30 mL Syringe)
260011	Hip Access Kit Straight Blade (Straight Blade, 16g x 7" Spinal Needle, 16 g x 6" Spinal Needle, Two 1.25 mm x 17.5" Round Tip Guidewires, Marking Pen, Ruler, 30 mL Syringe)
260003	Hip Blade Curved
260010	Hip Blade Straight













Suture Portfolio




Sutures & Suture Passers

Needle Guide


27 mm	13 mm	26 mm	27 mm	13 mm
				
CP-2 Needle (Reverse Cutting)	P-3 Needle (Reverse Cutting)	SH Needle (Taper Point)	CT-2 Needle (Taper Point)	C-1 Needle (Taper Point)

3/8 circle 17 mm	1/2 circle 17 mm	1/2 circle 17 mm	1/2 circle 22 mm	1/2 circle 17 mm
				
V-4 Needle (TAPERCUT™)	V-5 Needle (TAPERCUT™)	OS-6 36 mm (Reverse Cutting)	MO-6 27 mm (Taper Point)	RB-1 Needle (Taper Point)


Taper Point
Designed for soft, easy
penetrated tissue.



Reverse Cutting
Cutting edge on outer curve.
Designed for tough,
difficult-to-penetrate tissue.



TAPERCUT™ Technology
Cutting tip, taper body designed for
tough tissue, like two needles in one.



Free Needles



CODE	DESCRIPTION	QTY
223173	Free Needle Loop (OS-6)	12/pack



Sutures & Suture Passers

DYNATAPE™ Sutures



CODE	DESCRIPTION	QTY
280825	2.5 mm DYNATAPE Suture Pack, White, without needles	6/pack
280828	2.5 mm DYNATAPE Suture Pack, Black, without needles	6/pack
223173	Free Needle Loop OS-6	12/pack

DYNACORD™ Sutures



CODE	DESCRIPTION	QTY
222067	Size #2 DYNACORD Suture Pack, Blue, without Needle	12/pack
222065	Size #2 DYNACORD Suture Pack, Blue, with OS-6 Needle	12/pack
222066	Size #2 DYNACORD Suture Pack, Blue, with MO-7 Needle	12/pack
222068	Size #2 DYNACORD Suture Pack, Striped, without Needle	12/pack
222069	Size #2 DYNACORD Suture Pack, Striped/Blue, without Needle	12/pack
222071	Size #2 DYNACORD Suture Pack, Striped/Blue, with MO-7 Needle	12/pack
222073	Size #2 DYNACORD Suture Pack, Striped/Blue, with OS-6 Needle	12/pack

Sutures & Suture Passers

ORTHOCORD™ Sutures



Size #2 Sutures

CODE	DESCRIPTION	QTY
223102	ORTHOCORD Suture #2, Violet, with MO-6 1/2 Circle, Taper Point Needle, 26 mm	12/pack
223103	ORTHOCORD Suture #2, Violet, with OS-6 1/2 Circle, Reverse Cutting Needle, 36 mm	12/pack
223104	ORTHOCORD Suture #2, Violet, with MO-7 1/2 Circle, Taper Point Needle, 22 mm	12/pack
223105	ORTHOCORD Suture #2, Violet, without Needles	12/pack
223111	ORTHOCORD Suture #2, Blue, without Needles	12/pack
223113	ORTHOCORD Suture #2, Violet/Blue, without Needles	12/pack
223114	ORTHOCORD Suture #2, Violet/Blue, with MO-7 1/2 Circle, Taper Point Needle, 22 mm	12/pack
223115	ORTHOCORD Suture #2, Violet/Blue, with MO-6 1/2 Circle, Taper Point Needle, 26 mm	12/pack
223116	ORTHOCORD Suture #2, Violet/Blue, with OS-6 1/2 Circle, Reverse Cutting Needle, 36 mm	12/pack

Small Size Sutures

CODE	DESCRIPTION	QTY
223121	#4-0 ORTHOCORD Suture Loop 51cm Taper Point, 3/8 Circle, 17.9 mm	12/pack
223122	#2-0 ORTHOCORD Suture 46cm Taper Point, 1/2 Circle, 26.5 mm	12/pack
223123	#4-0 ORTHOCORD Suture 46cm Diamond Point, 3/8 Circle, 18.7 mm	12/pack
223124	#0 ORTHOCORD Suture 96.5cm Taper Point, 1/2 Circle, 22.2 mm	12/pack
223125	#3-0 ORTHOCORD Suture 46cm Taper Point, 3/8 Circle, 15 mm	12/pack
223126	#2-0 ORTHOCORD Suture 46cm Taper Point, 3/8 Circle, 17.9 mm	12/pack
223127	#4-0 ORTHOCORD Suture 46cm Taper Point, 3/8 Circle, 2.3 mm	12/pack



Sutures & Suture Passers


2.5 mm PERMATAPE™ Sutures



CODE	DESCRIPTION	QTY
223165	2.5 mm PERMATAPE Suture, Blue Pack	6/pack
223166	2.5 mm PERMATAPE Suture, White Pack	6/pack
223169	2.5 mm PERMATAPE Suture, Blue	1/each
223170	2.5 mm PERMATAPE Suture, White	1/each

1.3 mm PERMATAPE™ Sutures

CODE	DESCRIPTION	QTY
280772	1.3 mm PERMATAPE Suture, Blue, without Needles	12/pack
280773	1.3 mm PERMATAPE Suture, Blue, OS6 Needles	12/pack
280774	1.3 mm PERMATAPE Suture, Blue, MO-7 Needles	12/pack
280775	1.3 mm PERMATAPE Suture, White-Blue, without Needles	12/pack
280776	1.3 mm PERMATAPE Suture, White-Blue, OS6 Needles	12/pack
280777	1.3 mm PERMATAPE Suture, White-Blue, MO-7 Needles	12/pack



PERMALOOP™ Sutures

CODE	DESCRIPTION	QTY
280936	1.3 mm PERMALOOP Suture Straight Needle	12/pack
208938	1.3 mm PERMALOOP Suture MO-5 Needle	12/pack



Sutures & Suture Passers

EXPRESSEW™ III Suture Passers



CODE	DESCRIPTION	QTY
214124	EXPRESSEW III with Hook	1/each
214140	EXPRESSEW III without Hook	1/each
288233	EXPRESSEW III with AutoCapture+	1/each
288235	Retention plate with Loader	1/each
288234	Retention plate with Loader	5/each
214141	EXPRESSEW III Needles	5/each

Sutures & Suture Passers

IDEAL™ Suture Graspers & Shuttles



IDEAL Suture Graspers

CODE	DESCRIPTION	QTY
251720	IDEAL Suture Grasper, 15°	1/each
251721	IDEAL Suture Grasper, 30°	1/each
251722	IDEAL Suture Grasper, 45°	1/each
251723	IDEAL Suture Grasper, 60°	1/each

IDEAL Suture Shuttles

CODE	DESCRIPTION	QTY
251001	IDEAL Suture Shuttle, 25° Left	1/each
251002	IDEAL Suture Shuttle, 25° Right	1/each
251003	IDEAL Suture Shuttle, 45° Left	1/each
251004	IDEAL Suture Shuttle, 45° Right	1/each
251005	IDEAL Suture Shuttle, 90° Up	1/each
251006	IDEAL Suture Shuttle, Straight Hook	1/each
251007	IDEAL Suture Shuttle, Straight Crescent	1/each



CHIA PERCPASSER™ Suture Passer

CODE	DESCRIPTION	QTY
214101	CHIA PERCPASSER Suture Passer	5/each




Sutures & Suture Passers

Suture Leaders

CODE	DESCRIPTION	QTY
214608	Suture Leader, 20° Up	1/each
214609	Suture Leader, 30° Right	1/each
214610	Suture Leader, 30° Left	1/each
214611	Suture Leader, 45° Right	1/each
214612	Suture Leader, 45° Left	1/each
214613	Suture Leader, 70° Right	1/each
214614	Suture Leader, 70° Left	1/each
214628	Suture Leader, 180° Right	1/each
214629	Suture Leader, 180° Left	1/each



Allograft Tissue



Allograft ↳ Tendons

We offer an extensive portfolio of allograft tendons used in numerous clinical applications, including but not limited to: ACL reconstruction, MPFL reconstruction, and lateral ankle stabilization. Parametrics Medical offers three unique processing options for allograft tendons (Coll-e-Strong™, Lo-RAD™, and No-RAD™).

Coll-e-Strong®

Coll-e-Strong allografts are Parametrics Medical's trademarked electron beam sterilized tendons. This innovative processing methodology achieves a sterility assurance level (SAL) of 10^{-6} without the use of gamma irradiation. Coll-e-Strong grafts provide Parametrics Medical's best option of sterility, safety, and strength for surgeons and patients.

Lo-RAD™

Lo-RAD allografts are Parametrics Medical's trademarked low-dose gamma-irradiated tendons. These allografts are processed with an average exposure rate of 1.25 megarads which is among the lowest level of irradiation on the market. The Lo-RAD processing methodology achieves an SAL of 10^{-6} .

No-RAD™

No-Rad allografts are Parametrics Medical's trademarked aseptically processed tendons. No irradiation is used in processing or after packaging. The No-Rad processing methodology achieves an SAL of 10^{-3} .



**ACL Reconstruction
with BTB**

Allograft ↳ Tendons

Achilles

CODE	DESCRIPTION
PACH1001E	Preshaped Achilles – Coll-e-Strong 10 mm
PREAT2000	Preshaped Achilles – Lo-RAD 10 mm
NPREAT2000	Preshaped Achilles – No-RAD 10 mm
ACH1001E	Achilles Tendon with bone – Coll-e-Strong
AT2000	Achilles Tendon with bone – Lo-RAD
NAT2000	Achilles Tendon with bone – No-RAD
ACH1001EN	Achilles Tendon without bone – Coll-e-Strong
AT2000N	Achilles Tendon without bone – Lo-RAD
NAT2000N	Achilles Tendon without bone – No-RAD



Bone-Tendon-Bone

CODE	DESCRIPTION
PS10501E	Preshaped Hemi BTB – Coll-e-Strong 10 mm
PRE2039	Preshaped Hemi BTB – Lo-RAD 10 mm
NPRE2039	Preshaped Hemi BTB – No-RAD 10 mm
BTB3001E	Hemi BTB – Coll-e-Strong
PTH2039	Hemi BTB – Lo-RAD
NPTH2039	Hemi BTB – No-RAD
PTW2040	Whole BTB – Lo-RAD
NPTW2040	Whole BTB – No-RAD
PTWQ2040	Whole BTB with Extensors – Lo-RAD
NPTWQ2040	Whole BTB with Extensors – No-RAD



Quadricep Tendons




CODE	DESCRIPTION
QT1903E	Quad Tendon with Bone Block – Coll-e-Strong
QT1903	Quad Tendon with Bone Block – Lo-RAD
NQT1903	Quad Tendon with Bone Block – No-RAD



Allograft ↳ Tendons

Non-Bone Tendons

CODE	DESCRIPTION
ATIB2001E	Anterior Tibialis – Coll-e-Strong
AT2052	Anterior Tibialis – Lo-RAD
NAT2052	Anterior Tibialis – No-RAD
PTIB4001E	Posterior Tibialis – Coll-e-Strong
PT2053	Posterior Tibialis – Lo-RAD
NPT2053	Posterior Tibialis – No-RAD
PL8001E	Peroneus Tendon – Coll-e-Strong
PL2043	Peroneus Tendon – Lo-RAD
NPL2043	Peroneus Tendon – No-RAD
ST6001E	Semitendinosus – Coll-e-Strong
ST2044	Semitendinosus – Lo-RAD
NST2044	Semitendinosus – No-RAD
GR9001E	Gracilis – Coll-e-Strong
GR2033	Gracilis – Lo-RAD
NGR2033	Gracilis – No-RAD
STG7001E	Semitendinosus/Gracilis – Coll-e-Strong
STG2055	Semitendinosus/Gracilis – Lo-RAD
NSTG2055	Semitendinosus/Gracilis – No-RAD
FT2050	Flexor Tendon
PL2315E	Plantaris (Single Strand) – Coll-e-Strong
PL2315	Plantaris (Single Strand) – Lo-RAD
NPL2315	Plantaris (Single Strand) – No-RAD



Allograft

↳ Cartilage

Costal Cartilage

Costal Cartilage is hyaline cartilage from the rib that may be used in numerous clinical applications, including but not limited to: wrist, nasal, and ear procedures.

CODE	DESCRIPTION
OSCC20905	Costal Cartilage, Hydrated ≤5 cm
OSCC20958	Costal Cartilage, Hydrated 5 cm-8 cm
OSCC20908	Costal Cartilage, Hydrated >8 cm



Distal Tibia with Cartilage

Distal tibia allograft may be used to treat shoulder instability due to glenoid bone loss. The distal tibia's lateral portion matches the patient's glenoid because it has similar curvature, bone density, and cartilage.


CODE	DESCRIPTION
OSTD2070RFC	Distal Tibia with Cartilage, Right, Frozen, Irradiated
OSTD2070LFC	Distal Tibia with Cartilage, Left, Frozen, Irradiated



Meniscus

The biomechanical role of the meniscus is load distribution and stability in the knee joint. Parametrics Medical's allograft meniscus with hemi plateau bone block is used in meniscal allograft transplant procedures and one of the few viable treatment options for the young meniscectomized knee. Our allograft meniscus are aseptically processed. No irradiation is used in processing or after packaging.

CODE	DESCRIPTION
NMN3051	Meniscus, Left Medial, Frozen, Non-Irradiated
NMN3053	Meniscus, Left Lateral, Frozen, Non-Irradiated
NMN3052	Meniscus, Right Medial, Frozen, Non-Irradiated
NMN3054	Meniscus, Right Lateral, Frozen, Non-Irradiated




Allograft

↳ Acellular Dermis / Scaffolds

DERMIS ON DEMAND™ Allograft

CODE	DESCRIPTION
DOD15453	Dermis on Demand, Freeze-dried, 15 mm x 4.5 mm x 3 mm

Rotator Cuff Repair
with Augmentation of
DERMIS ON DEMAND™
Allograft




Allograft

↳ Acellular Dermis / Scaffolds

Coll-e-Derm™ RT

Coll-e-Derm RT is a prehydrated human acellular dermal matrix (ADM) that retains angiogenin and collagen type IV. Angiogenin and collagen type IV may play a key role in supporting revascularization.¹ Coll-e-Derm RT may be used in numerous clinical applications, including but not limited to: superior capsule reconstruction and Achilles repair.

**Superior Capsule
Reconstruction
with Coll-e-Derm RT**



¹ Testing performed by an independent laboratory. Data on file, Aziyo Biologics. Animal and bench testing results may not necessarily be indicative of clinical performance.



Allograft

↳ Acellular Dermis / Scaffolds

Coll-e-Derm™ **RT**, Meshed (0.4 – 0.8 mm thick)

CODE	DESCRIPTION
ACDMHY0404RT	Coll-e-Derm RT Patch, Meshed, 4x4 cm
ACDMHY0408RT	Coll-e-Derm RT Patch, Meshed, 4x8 cm

Coll-e-Derm™ **RT**, Thin (0.63 – 1.48 mm thick)

CODE	DESCRIPTION
ADTHY0101RT	Coll-e-Derm RT Patch, Thin, 1x1 cm
ADTHY0102RT	Coll-e-Derm RT Patch, Thin, 1x2 cm
ADTHY0104RT	Coll-e-Derm RT Patch, Thin, 1x4 cm
ADTHY0202RT	Coll-e-Derm RT Patch, Thin, 2x2 cm
ADTHY0204RT	Coll-e-Derm RT Patch, Thin, 2x4 cm
ADTHY54RT	Coll-e-Derm RT Patch, Thin, 5x4 cm



Coll-e-Derm™ **RT**, Medium (0.90 – 1.99 mm thick)

CODE	DESCRIPTION
ADMHY0307RT	Coll-e-Derm RT Patch, Medium, 3x7 cm
ADMHY0407RT	Coll-e-Derm RT Patch, Medium, 4x7 cm
ADMHY0412RT	Coll-e-Derm RT Patch, Medium, 4x12 cm
ADMHY0416RT	Coll-e-Derm RT Patch, Medium, 4x16 cm
ADMHY54RT	Coll-e-Derm RT Patch, Medium, 5x4 cm
ADMHY105RT	Coll-e-Derm RT Patch, Medium, 5x10 cm

Coll-e-Derm™ **RT**, SCR (3.00 – 3.50 mm thick)

CODE	DESCRIPTION
SCRHY473RT	Coll-e-Derm RT Patch, SCR, 4x7 cm
SCRHY543RT	Coll-e-Derm RT Patch, SCR, 5x4 cm



Allograft

↳ Acellular Dermis / Scaffolds

Fascia Lata

Fascia Lata is a fibrous covering of muscle recovered from the lateral aspect of the thigh used in various clinical applications such as periodontal surgery, dural defect repair, hip labral reconstruction, and superior capsule reconstruction. These allografts are gamma-irradiated in their final packaging to achieve an SAL of 10^{-6} .

CODE	DESCRIPTION
FAS2061	Fascia, Freeze Dried, 1-24 cm ²
FAS2062	Fascia, Freeze Dried, 25-49 cm ²
FAS2063	Fascia, Freeze Dried, 50-100 cm ²
FAS2064	Fascia, Freeze Dried, 101-150 cm ²
FAS2065	Fascia, Freeze Dried, 151-200 cm ²
FAS2066	Fascia, Freeze Dried, 201-250 cm ²



Allograft

↳ Structural Bone Grafts

Parametrics Medical offers a wide selection of structural bone grafts. These allografts are gamma-irradiated in their final packaging to achieve an SAL of 10^{-6} .

Backfill Plug, 20-25 mm L (for use with COR™ Precision Targeting System)

Cancellous backfill plugs, designed for use with the COR™ Precision Targeting System, are used to fill autograft harvest sites created during the repair of osteochondral defects. The osteoconductive backfill plugs act as a scaffold for the growth of natural bone and may reduce the risk of donor site morbidity. The backfill plugs are gamma-irradiated in their final packaging to achieve an SAL of 10^{-6} .


CODE	DESCRIPTION
OSCOR2004	Backfill Plug, Ambient, Hydrated, 4 mm
OSCOR2006	Backfill Plug, Ambient, Hydrated, 6 mm
OSCOR2008	Backfill Plug, Ambient, Hydrated, 8 mm
OSCOR2010	Backfill Plug, Ambient, Hydrated, 10 mm



Bone Dowel

Allograft bone dowels are commonly used for ACL revision tunnel backfill. The cylindrical graft with cortical face, cancellous body, and a pilot hole is designed to provide an immediate scaffold that may help restore the open void and promotes natural bone growth. The bone dowels are gamma-irradiated in their final packaging to achieve an SAL of 10^{-6} .

CODE	DESCRIPTION
OSCLD2007-9	Bone Dowel, Freeze Dried, 9 mm
OSCLD2007-10	Bone Dowel, Freeze Dried, 10 mm
OSCLD2007-12	Bone Dowel, Freeze Dried, 12 mm
OSCLD2007-14	Bone Dowel, Freeze Dried, 14 mm
OSCLD2007-16	Bone Dowel, Freeze Dried, 16 mm




Allograft

↳ Structural Bone Grafts

Parametrics Medical offers a wide selection of structural bone grafts. These allografts are gamma-irradiated in their final packaging to achieve an SAL of 10^{-6} .


Tri-Cortical Block

CODE	DESCRIPTION
OSITB2035	Ilium, Tri-Cortical Block, Freeze Dried, 1.0 cm
OSITB2036	Ilium, Tri-Cortical Block, Freeze Dried, 1.1-1.4 cm
OSITB2037	Ilium, Tri-Cortical Block, Freeze Dried, 1.5-1.9 cm
OSITB2038	Ilium, Tri-Cortical Block, Freeze Dried, 2.0-2.9 cm
OSITB2039	Ilium, Tri-Cortical Block, Freeze Dried, 3.0-3.9 cm
OSITB2040	Ilium, Tri-Cortical Block, Freeze Dried, 4.0-4.9 cm
OSITB2041	Ilium, Tri-Cortical Block, Freeze Dried, 5.0+ cm



Tri-Cortical Strip

CODE	DESCRIPTION
OSITS2040	Ilium, Tri-Cortical Strip, Freeze Dried, Small
OSITS2045	Ilium, Tri-Cortical Strip, Freeze Dried, Medium
OSITS2050	Ilium, Tri-Cortical Strip, Freeze Dried, Large



Ilium

CODE	DESCRIPTION
OSIH2055	Ilium, Hemi, Freeze Dried
OSIH2055F	Ilium, Hemi, Frozen
OSIW2056F	Ilium, Whole, Frozen




Allograft

↳ Structural Bone Grafts

Parametrics Medical offers a wide selection of structural bone grafts. These allografts are gamma-irradiated in their final packaging to achieve an SAL of 10^{-6} .


Cancellous Bone Block

CODE	DESCRIPTION
OSCLB2009-1	Cancellous Bone Block, Freeze Dried, 1.7-2.2x2.5-3.5 cm
OSCLB2009-2	Cancellous Bone Block, Freeze Dried, 1.0-1.5x2.5-3.5 cm



Cortical Strut

CODE	DESCRIPTION
OSCS2002	Cortical Strut, Freeze Dried, All Sizes



Humerus

CODE	DESCRIPTION
OSHHD1148F	Humerus, Distal, Frozen
OSHP2052F	Humerus, Proximal (L,R), Frozen
OSHPT2053F	Humerus, Proximal w/ Soft Tissue, Frozen
OSHS2053	Humerus, Shaft, Freeze Dried
OSHS2053F	Humerus, Shaft, Frozen
OSHH2051F	Humerus, Head, Frozen
OSHW2054F	Humerus, Whole (L,R), Frozen



Radius

CODE	DESCRIPTION
OSRD2064F	Radius, Distal- (L,R), Frozen
OSRP2065F	Radius, Prox- (L,R), Frozen
OSRS2063	Radius, Segment, Freeze Dried



Allograft

↳ Structural Bone Grafts

Parametrics Medical offers a wide selection of structural bone grafts. These allografts are gamma-irradiated in their final packaging to achieve an SAL of 10^{-6} .

Femur

CODE	DESCRIPTION
OSFS2016	Femur Shaft, Freeze Dried, < 11.0 cm
OSFS2016F	Femur Shaft, Frozen, < 11.0 cm
OSFSL2017	Femur Shaft, Freeze Dried, > 11.0 cm
OSFS2025	Femur Shaft, Split One Per Pack, Freeze Dried
OSFS2025F	Femur Shaft, Split One Per Pack, Frozen
OSFD2030F	Femur, Distal (L,R), Frozen
OSFPH2031F	Femur, Prox w/ Head (L,R), Frozen
OSFWH2034F	Femur, Whole w/ Head (L,R), Frozen
OSFP2001	Femur Plate, Freeze Dried, ≥ 15 cm < 20 cm
OSFP2000	Femur Plate, Freeze Dried, ≥ 20 cm
OSFSG2018	Femur Segment, Freeze Dried
OSHFCL2195F	Hemi Femoral Condyle, Lateral, Frozen
OSHFCLM2195F	Hemi Femoral Condyle, Medial, Frozen
OSFH2015	Femoral Head w/o Cartilage, Freeze Dried
OSFH2015F	Femoral Head w/o Cartilage, Frozen



Patella

CODE	DESCRIPTION
OSPW2042	Patellar Wedge, Freeze Dried
OSPW2042F	Patellar Wedge, Frozen



Allograft

↳ Structural Bone Grafts

Parametrics Medical offers a wide selection of structural bone grafts. These allografts are gamma-irradiated in their final packaging to achieve an SAL of 10^{-6} .

Tibia

CODE	DESCRIPTION
OSTD2070F	Tibia, Distal, Frozen
OSTPP1181F	Tibia Proximal w/Patella Ligament (L,R), Frozen
OSTPP1182F	Tibia Proximal w/Tuberosity (L,R), Frozen
OSTPP2073F	Tibia Proximal w/o Patella Ligament, Frozen
OSTD2071	Tibia Shaft, Freeze Dried
OSTD2071F	Tibia Shaft, Frozen
OSTIB2001	Tibia Shaft Split, Freeze Dried
OSTIB2001F	Tibia Shaft Split, Frozen
OSTP2000	Tibia Plate ≥ 20 cm, Freeze Dried
OSTS2074	Tibia Segment, Freeze Dried
OSTW2074F	Tibia, Whole (L,R), Frozen



Fibula

CODE	DESCRIPTION
OSFD2040F	Fibula, Distal (L,R), Frozen
OSFP2041F	Fibula, Proximal (L,R), Frozen
OSFS2029	Fibula, Segment ≤ 5 cm, Freeze Dried
OSFS2030	Fibula, Segment 6-10 cm, Freeze Dried
OSFS2031	Fibula, Segment 11+ cm, Freeze Dried
OSFS2031F	Fibula, Segment 11+ cm, Frozen
OSFW2043F	Fibula Whole (L,R), Frozen





Procedure Systems

200

16

17



Procedure Systems

Table of Contents

Shoulder

Shoulder Procedural Set	138
Rotator Cuff	
<i>HEALIX ADVANCE™ Anchor Instrument Set</i>	139
<i>VERSALoop™ 2.5mm Instrument Kit</i>	139
<i>GII™ Anchor Arthroscopic Instrument Set</i>	140
<i>VERSALOK™ Suture Anchor Instrument Set</i>	140
<i>EXPRESSEW™ III / EXPRESSEW™ III with Autocapture+ Flexible Suture Passer Set</i>	140
Biceps	
<i>Biceps Repair Set</i>	141
Labrum	
<i>VERSALoop™ 1.5mm & 1.8mm Instrument Kit</i>	142
<i>GRYPHON™ PROKNOT™ / GRYPHON™ Suture Anchor Instrument Set</i>	142
<i>Clever Hook & Sixter Instrument Set</i>	143
Glenoid	
<i>LATARJET EXPERIENCE Shoulder Instability System Instrument Set</i>	144



Knee

Ligament	
<i>MILAGRO™ ADVANCE / MILAGRO™ BR Interference Screw Instrument Sets</i>	145
<i>PROFILE™ Titanium Interference Screw Set</i>	145
<i>INTRAFIX™ ADVANCE Instrument Set</i>	146
<i>Screw & Washer Instrument Set</i>	147
<i>RIGIDFIX™ BTB Cross Pin System Instrument Set</i>	148
<i>RIGIDFIX™ ST Cross Pin System Instrument Set</i>	148
Graft Prep	
<i>Graft Preparation System Instrument Set</i>	149
Meniscus	
<i>TRUESPAN™ Meniscal Repair System Instrument Set</i>	149
Cartilage	
<i>COR™ Plunger/Sizer Set</i>	149
Tibial Realignment	
<i>TRACKER™ AMZ Guide System Instrument Set</i>	150
Arthroscopic Equipment & Accessories	
<i>Cruciate+ Instrument Set</i>	151
<i>Anteromedial Drilling Set</i>	152
<i>Notchplasty Instrument Set</i>	153
<i>Tunnel Preparation Tray</i>	154
<i>Tunnel Dilator Instrument Set</i>	155
<i>Revision Instrument Set</i>	155

Small Joint

<i>GII™ / Super / Mini Suture Anchor Set</i>	155
<i>HEALIX™ Hard Bone Drill Guide</i>	155
<i>STRADA™ Small Joint Arthroscopy Instrument Set</i>	155

Hip

<i>Hip Procedural Instrument System</i>	156
<i>VERSALoop™ 1.5 mm XL Hip Instrument Kit</i>	156

Procedure Systems

↳ Shoulder

Shoulder Procedural Set

CODE	DESCRIPTION	QTY	CODE	DESCRIPTION	QTY
215401	Shoulder Sterilization Case	1/each	211042	LUPINE™ Anchor Spiral Fluted Drill Bit	1/each
215459	External Sterilization Case Lid	1/each	214617	Calibrated Probe	1/each
215403	RC/Instability Fixation/Sterilization Tray	1/each	215406	Arthroscopic Instability Sterilization Tray (10 slots)	1/each
222314	HEALIX ADVANCE™ Anchor Awl, 4.5 mm	1/each	288233	EXPRESSEW™ III AutoCapture+ Suture Passer	1/each
222315	HEALIX ADVANCE™ Anchor Awl, 5.5 mm/6.5 mm	1/each	214626	Suture Manipulator	1/each
222316	HEALIX ADVANCE™ Anchor Cortical Awl/Tap, 4.5 mm	1/each	270120	Grasper/Grabber	1/each
222317	HEALIX ADVANCE™ Anchor Cortical Awl/Tap, 5.5 mm/6.5 mm	1/each	214600	Penetrating Grasper Straight	1/each
222223	HEALIX™ Anchor Awl	1/each	214601	Penetrating Grasper 35' Up	1/each
219495	HEALIX™ Anchor Awl, 3.4 mm	1/each	214602	Penetrating Grasper 45' Left	1/each
222273	HEALIX TRANSTEND™ Anchor Awl/ Tap, 3.4 mm (BR and PEEK Anchors)	1/each	214603	Penetrating Grasper 45' Right	1/each
222281	HEALIX TRANSTEND™ Anchor Drill	1/each	214604	Combo Grasper	1/each
214681	Bit, 2.2 mm (for TI Anchor) 4 mm Switching Stick, Blunt/Blunt	1/each	214605	Suture Scissors	1/each
214123	Switching Stick 4 mm	1/each	214646	Cordcutter	1/each
214619	Rasp 20 Up Bend	1/each	210263	Flush Suture Cutter	1/each
214620	Rasp 20 Down Bend	1/each	215407	Pathseeker Sterilization Tray	1/each
214623	Tissue Liberator Blade Up	1/each	215408	Pathseeker Sterilization Tray Pin Mat	1/each
214624	Tissue Liberator Blade Down	1/each	214105	5.5 x 75 mm Reusable Obturator	1/each
210824	GRYPHON™ PROKNOT™ Anchor	1/each	214113	7 x 75 mm Reusable Obturator	1/each
210826	GRYPHON™ PROKNOT™ Anchor Slotted Fishmouth Guide	1/each	214119	8.5 x 75 mm Resuable Obturator	1/each
211038	GRYPHON™ Anchor Obturator	1/each	214615	Knot Manipulator Full	1/each
211041	GRYPHON™ Anchor Blunt Obturator	1/each	214616	Crochet Hook	1/each
211036	GRYPHON™ Anchor Drill Bit, 2.4 mm (recommended to order 2)	1/each	210835	Knot Manipulator	1/each
213819	LUPINE™ Anchor Sawtooth Drill Guide	1/each	219304	Large Mallet	1/each
213817	LUPINE™ Anchor Fishmouth Drill Guide	1/each	215410	Additional Pin Mat (at bottom of tray)	1/each
214639	V Anchor Drill Guide Obturator	1/each	214140	EXPRESSEW™ III without Hook	1/each
211033	LUPINE™ Anchor Drill Bit	1/each	214640	Clever Hook Right	1/each
			214641	Clever Hook Left	1/each
			214642	Sixer Right	1/each
			214643	Sixer Left	1/each



Procedure Systems

↳ Shoulder

HEALIX ADVANCE™ Anchor Instrument Set

CODE	DESCRIPTION	QTY
215452	Medium Sterilization Tray	1/each
215453	Medium Sterilization Tray lid	1/each
222314	HEALIX ADVANCE Anchor Awl, 4.5 mm	1/each
222315	HEALIX ADVANCE Anchor Awl, 5.5 mm/6.5 mm	1/each
222316	HEALIX ADVANCE Anchor Cortical Awl/Tap, 4.5 mm	1/each
222329	HEALIX ADVANCE Anchor Full Tap, 4.5 mm	1/each
222328	HEALIX ADVANCE Anchor Full Tap, 5.5 mm/6.5 mm	1/each
222327	HEALIX ADVANCE Anchor Drill Bit, 3.5 mm	1/each
222326	HEALIX ADVANCE Anchor Drill Bit, 4 mm	1/each
222223	HEALIX™ Anchor Awl	1/each
219495	HEALIX Anchor Awl, 3.4 mm	1/each
222273	HEALIX TRANSTEND™ Anchor Awl/Tap, 3.4 mm (for BR and PEEK Anchors)	1/each
222281	HEALIX TRANSTEND Anchor Drill Bit, 2.2 mm (optional for TI Anchor)	1/each

VERSALOOP™ 2.5 mm Instrument Kit

CODE	DESCRIPTION	QTY
215452	Medium Sterilization Tray	1/each
215453	Medium Sterilization Tray Lid	1/each
210177	VERSALOOP Straight Guide for 2.5 mm Anchor, Fishmouth	1/each
210180	VERSALOOP Straight Guide for 2.5 mm Anchor, Sawtooth	1/each
210184	VERSALOOP Reusable Straight Drill for 2.5 mm Anchor	1/each
210207	VERSALOOP Reusable Awl for 2.5 mm Anchor	1/each
210214	VERSALOOP Obturator for 2.5 mm Anchor	1/each



Procedure Systems

↳ Shoulder

GII™ Anchor Arthroscopic Instrument Set

CODE	DESCRIPTION	QTY
215452	Medium Sterilization Tray	1/each
215453	Medium Sterilization Tray Lid	1/each
214528	Arthroscopic Suture Cutter	1/each
211502	Arthroscopic GII Anchor Drill Bit, 2.4 mm	1/each
211513	Arthroscopic Super/RC Drill Bit, 2.9 mm	1/each
212114	Super/RC Inserter	1/each
214541	Arthroscopic GII Anchor Inserter	1/each
214542	Arthroscopic T-Handle Reamer	1/each
214560	Arthroscopic Introducer Drill Guide	1/each
214561	Arthroscopic Sharkmouth Drill Guide	1/each
214562	Arthroscopic Sawtooth Drill Guide	1/each
214529	Arthroscopic Periosteal Elevator	1/each
214575	Arthroscopic Suture Grasper	1/each
214616	Crochet Hook	1/each
214502	Skid Rod with Handle, #2-0 or #0 Suture	1/each
214527	Arthroscopic Knot Pusher	1/each

VERSALOK™ Suture Anchor Instrument Set

CODE	DESCRIPTION	QTY
215452	Medium Sterilization Tray	1/each
215453	Medium Sterilization Tray lid	1/each
214710	VERSALOK Suture Anchor Deployment Gun	1/each
214711	VERSALOK Suture Anchor Awl	1/each

EXPRESSEW™ III / EXPRESSEW™ III w/Autocapture+ Flexible Suture Passer Set

CODE	DESCRIPTION	QTY
215452	Medium Sterilization Tray	1/each
215453	Medium Sterilization Tray lid	1/each
214124	EXPRESSEW III Flexible Suture Passer with Hook	1/each
214140	EXPRESSEW III Flexible Suture Passer without Hook	1/each
288233	EXPRESSEW III AutoCapture+ Flexible Suture Passer	1/each



Procedure Systems

↳ Shoulder

Biceps Repair Set

CODE	DESCRIPTION	QTY
215452	Medium Sterilization Tray	1/each
215453	Medium Sterilization Tray Lid	1/each
254847	Quick-Connect Ratcheting Handle	1/each
219500	MILAGRO™ ADVANCE Screw Universal Modular Driver	1/each
219447	Modular Driver, 15 mm	1/each
219448	Modular Driver, 12 mm	1/each
219446	Modular Driver, 5 mm x 23 mm and 6 mm x 23 mm Screws	1/each
219438	Modular Driver, 5 mm x 30 mm and 6 mm x 30 mm Screws	1/each
219445	Tendon Sizer, 4 mm to 9 mm	1/each
232028	RIGIDLOOP™ System Cannulated Reamer, 4.5 mm	1/each
219440	Acorn Reamer, 5 mm	1/each
219441	Acorn Reamer, 5.5 mm	1/each
219442	Fully Fluted Reamer, 5 mm	1/each
219443	Fully Fluted Reamer, 5.5 mm	1/each
219621	Tendon Fork, Small (4 mm wide)	1/each
219974	Tendon Fork, Large (5.5 mm wide)	1/each
219971	Guidewire Cap	1/each
214617	Calibrated Probe	1/each
214615	Knot Manipulator, Full Loop	1/each
219304	Large Mallet	1/each



Procedure Systems

↳ Shoulder

VERSALOOPTM 1.5mm & 1.8mm Instrument Set

CODE	DESCRIPTION	QTY
215452	Medium Sterilization Tray	1/each
215453	Medium Sterilization Tray Lid	1/each
210175	VERSALOOPTM Straight Guide for 1.5 mm Anchor, Fishmouth	1/each
210176	VERSALOOPTM Straight Guide for 1.8 mm Anchor, Fishmouth	1/each
210178	VERSALOOPTM Straight Guide for 1.5 mm Anchor, Sawtooth	1/each
210179	VERSALOOPTM Straight Guide for 1.8 mm Anchor, Sawtooth	1/each
210181	VERSALOOPTM Reusable Straight Drill for 1.5 mm Anchor	1/each
210182	VERSALOOPTM Reusable Straight Drill for 1.8 mm Anchor	1/each
210195	VERSALOOPTM 12.5° Guide for 1.5 mm Anchor, Sawtooth	1/each
210196	VERSALOOPTM 12.5° Guide for 1.8 mm Anchor, Sawtooth	1/each
210197	VERSALOOPTM 25° Guide for 1.5 mm Anchor, Sawtooth	1/each
210198	VERSALOOPTM 25° Guide for 1.8 mm Anchor, Sawtooth	1/each
210205	VERSALOOPTM Obturator for 1.5 mm Anchor	1/each
210213	VERSALOOPTM Obturator for 1.8 mm Anchor	1/each

GRYPHONTM PROKNOTM / GRYPHONTM Suture Anchor Instrument Set

CODE	DESCRIPTION	QTY
215452	Medium Sterilization Tray	1/each
215453	Medium Sterilization Tray lid	1/each
211035	GRYPHONTM Suture Anchor Tap	1/each
211036	GRYPHONTM Suture Anchor Drill Bit, 2.4 mm	1/each
211037	GRYPHONTM Suture Anchor Fishmouth Guide	1/each
211038	GRYPHONTM Suture Anchor Obturator	1/each
211039	GRYPHONTM Suture Anchor Sawtooth Guide	1/each
211041	GRYPHONTM Suture Anchor Blunt Obturator	1/each
210824	GRYPHONTM PROKNOTM Anchor	1/each
210835	Knot Manipulator	1/each



Procedure Systems

↳ Shoulder

Clever Hook & Sixter Instrument Set

CODE	DESCRIPTION	QTY
215201	Small Arthroscopy Tray	1/each
214640	Clever Hook - Right	1/each
214641	Clever Hook - Left	1/each
214642	Sixter - Right	1/each
214643	Sixter - Left	1/each

Procedure Systems

↳ Shoulder

LATARJET EXPERIENCE Shoulder Instability System Instrument Set

CODE	DESCRIPTION	QTY	CODE	DESCRIPTION	QTY
288219	Latarjet Case	1/each	288210	Reference Tube	1/each
288217	Latarjet Case Lid	1/each	288209	Coracoid Screw, 3.5 mm	1/each
288218	Latarjet Tray	1/each	288214	Glenoid Double-Pronged Lever	1/each
288200	Inline Coracoid Drill Guide	1/each	288212	Humeral Head Retractor	1/each
288206	Coracoid Offset Pin	1/each	288213	Hohmann Retractor	1/each
288201	Coracoid Drill	1/each	288215	Subscapularis Spreader	1/each
288202	Coracoid Top Hat Drill	1/each	214135	5.5 mm Hip Slotted Cannula	1/each
288203	Top Hat Tap	1/each	214136	6 mm Hip Slotted Cannula	1/each
288211	Combo Screwdriver	1/each	285365	Latarjet Coracoid Drill Guide	1/each
288204	Sharp Curved Osteotome	1/each	285260	Latarjet Coracoid Short K-Wire, 150 mm	1/each
288240	Latarjet Switching Stick	1/each	285270	Latarjet Curved Osteotome	1/each
288101	Latarjet Switching Stick, Small	1/each	285300	Latarjet Glenoid Drill Bit, 3.2 mm	1/each
288205	Subscapularis Channeler	1/each	211607	U-Shaped Cannula, Medium	1/each
288221	Glenoid Drill Guide	1/each	211610	U-Shaped Cannula Obturator, Medium	1/each
288208	Glenoid Offset Pin	1/each	211608	U-Shaped Cannula, Large	1/each
288241	Glenoid Drill	1/each	211611	U-Shaped Cannula Obturator, Large	1/each
288242	Small Probe	1/each			

Procedure Systems

↳ Knee

MILAGRO™ ADVANCE / MILAGRO™ BR Interference Screw Instrument Set

CODE	DESCRIPTION	QTY
215452	Medium Sterilization Tray	1/each
215453	Medium Sterilization Tray Lid	1/each
219215	Ratchet Handle	1/each
219950	Starter Tap	1/each
219952	Pathmaker Tunnel Notcher	1/each
219957	Modular Driver - 23 mm	1/each
219958	Modular Driver - 30 mm	1/each
219959	Modular Driver 35 mm	1/each
219496	MILAGRO ADVANCE Screw Tap, 7-8 mm	1/each
219497	MILAGRO ADVANCE Screw Tap, 9-10 mm	1/each

PROFILE™ Titanium Interference Screw Set

CODE	DESCRIPTION	QTY
215450	Large Sterilization Tray	1/each
215451	Large Sterilization Tray Lid	1/each
211332	Round Tip Guidewire, 2.4 mm x 14", Non-Sterile, Nitinol	5/each
211322	Round Tip Guidewire, 1.6 mm x 16", Non-Sterile, Nitinol	5/each
219218	3.5 mm Hex Driver Mod Tip	1/each
219219	4.0 mm Hex Driver Mod Tip	1/each
219215	Rachet Handle	1/each
254511	ACL Tunnel Notcher Shaft,	1/each
254512	Absolute Screw Easy Out Shaft	1/each
254715	ACL Modular Handle, Grey Ratcheting Handle with Jacob's Chuck	1/each



Procedure Systems

↳ Knee

INTRAFIX™ ADVANCE Instrument Set

CODE	DESCRIPTION	QTY
254839	Sterilization Case	1/each
254845	Sterilization Tray 1	1/each
254846	Sterilization Tray 2	1/each
254848	Sterilization Lid	1/each
254828	INTRAFIX ADVANCE Small Trial/Dilator	1/each
254829	INTRAFIX ADVANCE Large Trial/Dilator	1/each
254830	INTRAFIX ADVANCE XL Trial/Dilator	1/each
254831	INTRAFIX ADVANCE Small 23 mm Inserter	1/each
254832	INTRAFIX ADVANCE Large 23 mm Inserter	1/each
254833	INTRAFIX ADVANCE Small 30 mm Inserter	1/each
254834	INTRAFIX ADVANCE Large 30 mm Inserter	1/each
254835	INTRAFIX ADVANCE XL 30 mm Inserter	1/each
254836	INTRAFIX ADVANCE Modular Hex Driver	1/each
254847	Quick-Connect Ratcheting Handle	1/each
254838	Tensioner	1/each
254842	Slap Hammer	1/each
254605	ACL Tie Tensioner	1/each



Procedure Systems

↳ Knee

Screw & Washer Implants and Instrument Set

CODE	DESCRIPTION	QTY	CODE	DESCRIPTION	QTY
215024	Implant Tray	1/each	220442	Cortical Screw, 4.5 mm x 42 mm	1/each
215124	Screw and washer tray - system	1/each	220444	Cortical Screw, 4.5 mm x 44 mm	1/each
211024	Drill Bit, 3.2 mm	1/each	220446	Cortical Screw, 4.5 mm x 46 mm	1/each
254318	3.2 mm GeoFit Drill Guide	1/each	220448	Cortical Screw, 4.5 mm x 48 mm	1/each
219114	New Depth Gauge (Screw & Washer)	1/each	220450	Cortical Screw, 4.5 mm x 50 mm	1/each
219124	Hex Driver, 2.5 mm	1/each	220452	Cortical Screw, 4.5 mm x 52 mm	1/each
219200	Cancellous Tap, 6.5 mm	1/each	220454	Cortical Screw, 4.5 mm x 54 mm	1/each
219400	Cortical Tap, 4.5 mm	1/each	220456	Cortical Screw, 4.5 mm x 56 mm	1/each
220225	Cancellous Screw, 6.5 mm x 25 mm	1/each	220458	Cortical Screw, 4.5 mm x 58 mm	1/each
220230	Cancellous Screw, 6.5 mm x 30 mm	1/each	220460	Cortical Screw, 4.5 mm x 60 mm	1/each
220235	Cancellous Screw, 6.5 mm x 35 mm	1/each	220462	Cortical Screw, 4.5 mm x 62 mm	1/each
220240	Cancellous Screw, 6.5 mm x 40 mm	1/each	220464	Cortical Screw, 4.5 mm x 64 mm	1/each
220245	Cancellous Screw, 6.5 mm x 45 mm	1/each	220466	Cortical Screw, 4.5 mm x 66 mm	1/each
220250	Cancellous Screw, 6.5 mm x 50 mm	1/each	220468	Cortical Screw, 4.5 mm x 68 mm	1/each
220255	Cancellous Screw, 6.5 mm x 55 mm	1/each	220470	Cortical Screw, 4.5 mm x 70 mm	1/each
220426	Cortical Screw, 4.5 mm x 26 mm	1/each	220472	Cortical Screw, 4.5 mm x 72 mm	1/each
220428	Cortical Screw, 4.5 mm x 28 mm	1/each	220474	Cortical Screw, 4.5 mm x 74 mm	1/each
220430	Cortical Screw, 4.5 mm x 30 mm	1/each	220476	Cortical Screw, 4.5 mm x 76 mm	1/each
220432	Cortical Screw, 4.5 mm x 32 mm	1/each	221014	Spiked Washer, 14 mm (OD)	1/each
220434	Cortical Screw, 4.5 mm x 34 mm	1/each	221017	Spiked Washer, 17 mm (OD)	1/each
220436	Cortical Screw, 4.5 mm x 36 mm	1/each	221020	Spiked Washer, 20 mm (OD)	1/each
220438	Cortical Screw, 4.5 mm x 38 mm	1/each	221114	Flat Washer, 14 mm (OD)	1/each
220440	Cortical Screw, 4.5 mm x 40 mm	1/each	221117	Flat Washer, 17 mm (OD)	1/each

Procedure Systems

↳ Knee

RIGIDFIX™ BTB Cross Pin System Instrument Set

CODE	DESCRIPTION	QTY
215422	RIGIDFIX ACL System	1/each
215424	BTB RIGIDFIX Caddy	1/each
215451	Large Sterilization Tray Lid	1/each
213713	RIGIDFIX System BTB Stepped Trocar, Short	1/each
213703	RIGIDFIX System BTB Stepped Trocar, Long	1/each
213702	RIGIDFIX Cross Pin System BTB Guide Frame	1/each
213704	RIGIDFIX Cross Pin System Stepped Pin Inserter Rod	1/each
213705	RIGIDFIX Cross Pin System Mallet	1/each
213706	RIGIDFIX Cross Pin System Sleeve Removal Tool	1/each
213709	RIGIDFIX Cross Pin System Rod, 9 mm	1/each
213710	RIGIDFIX Cross Pin System Rod, 10 mm	1/each
213711	RIGIDFIX Cross Pin System Rod, 11 mm	1/each
213712	RIGIDFIX Cross Pin System Rod, 12 mm	1/each

RIGIDFIX™ ST Cross Pin System Instrument Set

CODE	DESCRIPTION	QTY
215422	RIGIDFIX ACL System	1/each
215423	Rod Caddy	1/each
215451	Large Sterilization Tray Lid	1/each
213701	RIGIDFIX Cross Pin System ST Guide Frame	1/each
213704	RIGIDFIX Cross Pin System Stepped Pin Inserter Rod	1/each
213705	RIGIDFIX Cross Pin System Mallet	1/each
213706	RIGIDFIX System Sleeve Removal Tool	1/each
213707	RIGIDFIX Cross Pin System Rod, 7 mm (ST Only)	1/each
213708	RIGIDFIX Cross Pin System Rod, 8 mm (ST Only)	1/each
213885	RIGIDFIX Cross Pin System Rod, 8.5 mm (ST Only)	1/each
213709	RIGIDFIX Cross Pin System Rod, 9 mm	1/each
213710	RIGIDFIX Cross Pin System Rod, 10 mm	1/each



Procedure Systems

↳ Knee

Graft Preparation System Instrument Set

CODE	DESCRIPTION	QTY
215810	Graft Preparation Sterilization Tray	1/each
219961	Graft Preparation Board	1/each
219962	Slider Unit (Requires 2 each for system)	1/each
219963	Tensioning Clamp	1/each
219964	Graft Clamp (Requires 2 each for system)	1/each
219965	Suture Clamp	1/each
219966	Graft Sizer, 6-12 mm (amber color)	1/each
219444	Graft Sizer, 5-12 mm (white color)	1/each
219967	BTB Drill Guide	1/each
219968	Cutting Board	1/each

TRUESPAN™ Meniscal Repair System Instrument Set

CODE	DESCRIPTION	QTY
215452	Medium Sterilization Tray	1/each
215453	Medium Sterilization Tray Lid	1/each
228301	Arthroscopic Pusher	1/each
228302	Arthroscopic Pusher/Cutter	1/each
253110	Meniscal Repair System Rasp, Straight Up with 45° Face	1/each
253111	Meniscal Repair System Rasp, Straight Up with 90° Face	1/each
253112	Meniscal Repair System Rasp, Straight with 90° Face	1/each
253113	Meniscal Probe	1/each
232024	Malleable Graft Retractor (5PK)	1/each

COR™ Plunger/Sizer Set

CODE	DESCRIPTION	QTY
215454	Small Sterilization Tray	1/each
215455	Small Sterilization Tray Lid	1/each
252307	COR System Plunger, 4 mm	1/each
252313	COR System Plunger, 6 mm	1/each
252319	COR System Plunger, 8 mm	1/each



Procedure Systems

↳ Knee

TRACKER™ AMZ Guide System Instrument Set

CODE	DESCRIPTION	QTY
215187	TRACKER Tray & Lid	1/each
219452	TRACKER AMZ Guide System Cutting Block	1/each
219453	TRACKER AMZ Guide System Cutting Block with Pivot	1/each
219454	TRACKER AMZ Guide System Retractor	1/each
219455	TRACKER AMZ Guide System Slope Selector Arm	1/each
219456	TRACKER AMZ Guide System Anchor Pin Puller	1/each
219457	TRACKER AMZ Guide System Bone Clamp	1/each
219458	TRACKER AMZ Guide System Positioning Card	1/each
219460	TRACKER AMZ Guide System Drill Bit, 2 mm x 6.5"	1/each
219461	TRACKER AMZ Guide System Drill Bit, 3.2 mm x 7.0"	1/each
219462	TRACKER AMZ Guide System Drill Bit, 4.5 mm x 6.5"	1/each
219465	TRACKER AMZ Guide System Anchor Pin, 2 mm x 45 mm	1/each
219466	TRACKER AMZ Guide System Anchor Pin, 2 mm x 64 mm	1/each
219467	TRACKER AMZ Guide System Anchor Pin, 2 mm x 77 mm	1/each

Procedure Systems

↳ Knee

Cruciate+ Instrument Set

CODE	DESCRIPTION	QTY
215203	Cruciate+ Targeting Tray	1/each
215205	Cruciate+ Targeting Tray Lid	1/each
219001	Cruciate+ Guide Carriage	1/each
219006	Cruciate+ ACL Tibial Point	1/each
219007	Cruciate+ ACL Tibial Capture	1/each
219008	Cruciate+ Femoral Point	1/each
219009	Cruciate+ Femoral Capture	1/each
219003	Cruciate+ Beath Pin Bullet	1/each
219004	Cruciate+ Retroreamer Bullet	1/each
219002	Cruciate+ RR Bullet Obturator	1/each
219005	Cruciate+ 2.4 mm Bullet Insert	1/each
215203	Cruciate+ Targeting Tray	1/each
215205	Cruciate+ Targeting Tray Lid	1/each
219011	Cruciate+ PCL Tibial Point	1/each
219012	Cruciate+ PCL Tibial Disc	1/each
219013	Cruciate+ PCL Rasp	1/each
219015	Cruciate+ Cannulated PCL Elevator	1/each
219017	Cruciate+ Pin Depth Gage	1/each
219018	Cruciate+ Suture Passing Pin	1/each

Procedure Systems

↳ Knee

Anteromedial Drilling Set

CODE	DESCRIPTION	QTY
215452	Medium Sterilization Tray	1/each
215453	Medium Sterilization Tray Lid	1/each
215416	Acorn Reamer Caddy	1/each
219624	Low Profile Reamer, 6 mm	1/each
219625	Low Profile Reamer, 6.5 mm	1/each
219626	Low Profile Reamer, 7 mm	1/each
219627	Low Profile Reamer, 7.5 mm	1/each
219628	Low Profile Reamer, 8 mm	1/each
219629	Low Profile Reamer, 8.5 mm	1/each
219630	Low Profile Reamer, 9 mm	1/each
219631	Low Profile Reamer, 9.5 mm	1/each
219632	Low Profile Reamer, 10 mm	1/each
219633	Low Profile Reamer, 10.5 mm	1/each
219634	Low Profile Reamer, 11 mm	1/each
219635	Low Profile Reamer, 11.5 mm	1/each
219636	Low Profile Reamer, 12 mm	1/each
219637	Low Profile Reamer, 13 mm	1/each
219435	Anteromedial Femoral Aimer, 4.5 mm	1/each
219434	Anteromedial Femoral Aimer, 5 mm	1/each
254672	Anteromedial Femoral Aimer, 5.5 mm	1/each
254686	Anteromedial Femoral Aimer, 6.5 mm	1/each



Procedure Systems

↳ Knee

Notchplasty Instrument Set

CODE	DESCRIPTION	QTY
215452	Medium Sterilization Tray	1/each
215453	Medium Sterilization Lid	1/each
215413	Caddy - Notchplasty	1/each
219302	Restore Flat Rasp	1/each
219314	Tibial Awl, 90° Round	1/each
219308	Cup Curette	1/each
219311	Graft Harvest Gouge	1/each
219305	Restore Arthroscopic Gouge	1/each
219304	Large Mallet	1/each
219303	Restore Osteotome with mm Scale	1/each
219307	Ring Curette, 7 mm	1/each
232015	Graft Sizer	1/each

Procedure Systems

↳ Knee

Tunnel Preparation Tray

CODE	DESCRIPTION	QTY
215411	ACL Tunnel Prep Case	1/each
215459	External Case Lid	1/each
215412	Tibial Guide & Graft Harvest Tray	1/each
232002	Pigtail Tendon Peeler	1/each
254730	Tendon Stripper Closed 6050	1/each
219301	Tibial Guide/Ratchet (3 pieces)	1/each
219392	Elbow Aiming Arm	1/each
219393	Tip Aiming Arm	1/each
219312	Back Cutting Burr/Jacob's End	1/each
219316	ACL Patella Retractor	1/each
219444	Graft Sizer, 5-12 mm (white color)	1/each
232024	Malleable Graft Retractor	5/each
219304	Large Mallet	1/each
215414	Notchplasty Pin Mat	1/each
219308	Cup Curette	1/each
219302	Flat Rasp	1/each
219306	Ring Curette, 4 mm	1/each
219313	Femoral Awl, 45° Round	1/each
219303	Osteotome with scale	1/each
219309	Curved Rasp	1/each
215415	Aimer Tray	1/each
215416	Caddy Acorn Reamer Thin	1/each
219336	6 mm Acorn Reamer	1/each
219337	7 mm Acorn Reamer	1/each
219338	8 mm Acorn Reamer Thin Shaft	1/each
219339	9 mm Acorn Reamer Thin Shaft	1/each
219335	10 mm Acorn Reamer Thin Shaft	1/each
219331	11 mm Acorn Reamer Thin Shaft	1/each

CODE	DESCRIPTION	QTY
219332	12 mm Acorn Reamer Thin Shaft	1/each
219333	13 mm Acorn Reamer Thin Shaft	1/each
219547	6.5 mm Acorn Reamer	1/each
219548	7.5 mm Acorn Reamer	1/each
219549	8.5 mm Acorn Reamer	1/each
219550	9.5 mm Acorn Reamer	1/each
219551	10.5 mm Acorn Reamer	1/each
219552	11.5 mm Acorn Reamer	1/each
215417	Caddy for Fully Fluted Reamer	1/each
219346	6 mm Fully Fluted Reamer	1/each
219347	7 mm Fully Fluted Reamer	1/each
219348	8 mm Fully Fluted Reamer	1/each
219349	9 mm Fully Fluted Reamer	1/each
219340	10 mm Fully Fluted Reamer	1/each
219341	11 mm Fully Fluted Reamer	1/each
219342	12 mm Fully Fluted Reamer	1/each
219553	6.5 mm Fully Fluted Reamer	1/each
219554	7.5 mm Fully Fluted Reamer	1/each
219555	8.5 mm Fully Fluted Reamer	1/each
219556	9.5 mm Fully Fluted Reamer	1/each
219557	10.5 mm Fully Fluted Reamer	1/each
219558	11.5 mm Fully Fluted Reamer	1/each
219355	Femoral aimer 5 mm offset	1/each
219356	Femoral aimer 6 mm offset	1/each
219357	Femoral aimer 7 mm offset	1/each
219358	Femoral aimer 8 mm offset	1/each
219352	Pin Puller	1/each



Procedure Systems

↳ Knee

Tunnel Dilator Instrument Set

CODE	DESCRIPTION	QTY
215419	Tunnel Dilator Tray	1/each
215451	Large Sterilization Tray Lid	1/each
219304	Large Mallet	1/each
219403	Tunnel Dilator, 6 mm	1/each
219404	Tunnel Dilator, 6.5 mm	1/each
219405	Tunnel Dilator, 7 mm	1/each
219406	Tunnel Dilator, 7.5 mm	1/each
219407	Tunnel Dilator, 8 mm	1/each
219408	Tunnel Dilator, 8.5 mm	1/each
219409	Tunnel Dilator, 9 mm	1/each
219410	Tunnel Dilator, 9.5 mm	1/each
219411	Tunnel Dilator, 10 mm	1/each
219412	Tunnel Dilator, 10.5 mm	1/each
219413	Tunnel Dilator, 11 mm	1/each
219414	Tunnel Dilator, 11.5 mm	1/each
219415	Tunnel Dilator, 12 mm	1/each
219966	Graft Sizer, 6-12 mm (amber color)	1/each

Revision Instrument Set

CODE	DESCRIPTION	QTY
215452	Medium Sterilization Tray	1/each
215453	Medium Sterilization Tray Lid	1/each
219380	Revision T-Handle	1/each
219381	Screw-Out Tip, Small Hex	1/each
219382	Revision Screw-Out Tip, Large Hex	1/each
219383	Revision Curette, Extra Small	1/each
219384	Revision Curette, Large	1/each
219385	Revision Trephine for 7 mm	1/each
219386	Revision Trephine for 9 mm	1/each



Procedure Systems

↳ Small Joints

GII™/Super/Mini Suture Anchor Set

CODE	DESCRIPTION	QTY
215454	Small Sterilization Tray	1/each
215455	Small Sterilization Tray Lid	1/each
213207	GII Slotted Drill Guide	1/each
211093	GII Drill Bit - 2.4 mm	1/each
211113	Super/RC/GLS Drill Bit - 2.9 mm	1/each
211083	Mini Drill Bit - 2.1 mm	1/each
213072	Mini Drill Guide	1/each

HEALIX™ Hard Bone Drill Guide

CODE	DESCRIPTION	QTY
215454	Small Sterilization Tray	1/each
215455	Small Sterilization Tray Lid	1/each
222375	HEALIX ADVANCE™ Open Drill Guide	1/each

STRADA™ Small Joint Arthroscopy Instrument Set

CODE	DESCRIPTION	QTY
215201	Small Joint Arthroscopy Sterilization Tray	1/each
224200	Small Joint Arthroscopy Grasper, 2.5 mm x 90 mm	1/each
224201	Small Joint Arthroscopy Punch, Straight Blunt, 2.1 mm x 90 mm	1/each
224202	Small Joint Arthroscopy Punch, Right Blunt, 2.1 mm x 90 mm	1/each
224203	Small Joint Arthroscopy Punch, Left Blunt, 2.1 mm x 90 mm	1/each
224204	Small Joint Arthroscopy Punch, Down Blunt, 2.1 mm x 90 mm	1/each
224205	Small Joint Arthroscopy Punch, Large, 2.7 mm x 90 mm	1/each
224206	Small Joint Arthroscopy Punch, Straight, 1 mm x 65 mm	1/each
224207	Small Joint Arthroscopy Probe, 30°, 1 mm x 65 mm	1/each
224208	Small Joint Arthroscopy Probe, 90°, 1 mm x 65 mm	1/each
224209	Small Joint Arthroscopy Scissors, 2.5 mm x 90 mm	1/each



Procedure Systems

↳ Hip

Hip Procedural Instrument System


CODE	DESCRIPTION	QTY	CODE	DESCRIPTION	QTY
214009	Hip Suture Passer Up	1/each	215508	Hip Labral Repair Tray	1/each
214013	Hip Suture Passer Straight	1/each	215509	Hip Sterilization Case	1/each
214111	7 x 110 mm Reusable obturator	1/each	215510	Lid, Hip Sterilization Case	1/each
214131	Hip Probe 90	1/each	214153	5.0 Hip Cannulated SwitchStick	1/each
214132	5.5 Hip Cannulated SwitchStick	1/each	214154	4.5 Hip Cannulated SwitchStick	1/each
214133	6.0 Hip Cannulated SwitchStick	1/each	215504	Hip Slotted Cannula Tray	1/each
214134	4.0 Hip Cannulated SwitchStick	1/each	215505	Hip Tissue & Suture Management Tray	1/each
214135	5.5 Hip Slotted Cannula	1/each	214109	5.5 mm X 90 mm Reusable Obturator	1/each
214136	6.0 Hip Slotted Cannula	1/each	214121	8.5 mm X 90 mm Reusable Obturator	1/each
214137	Hip Biter 15 Up	1/each	219103	GRYPHON™ Anchor Slotted Hip Guide	1/each
214138	Hip CordCutter	1/each	213042	GRYPHON™ Anchor Hip Sawtooth Drill Guide	1/each
214139	Hip Biter Straight	1/each	213044	GRYPHON™ Anchor Hip Fishmouth Drill Guide	1/each
214142	Hip Curette 20	1/each	219102	GRYPHON™ Anchor Hip Obturator	1/each
214143	Hip Grasper 20 Left	1/each	211036	GRYPHON™ Anchor Drill Bit, 2.4 mm	1/each
214144	Hip Grasper 20 Right	1/each	213705	RIGIDFIX™ System Mallet	1/each
214145	Hip Knot Manipulator	1/each	215506	Hip Entry Tray Pin Mat Accessory	1/each
214146	Hip Suture Manipulator	1/each	215507	Hip Entry & Tissue Removal Tray	1/each
214147	Hip Grasper Straight	1/each			
214148	Hip TissLib Blade 20 down	1/each			
214151	Handle, Hip SwitchStick	1/each			

VERSALOOPTM 1.5 mm XL Hip Instrument Kit

CODE	DESCRIPTION	QTY
215452	Medium Sterilization Tray	1/each
215453	Medium Sterilization Tray Lid	1/each
210217	VERSALOOPTM XL Straight Guide for 1.5 mm Anchor XL, Sawtooth	1/each
210218	VERSALOOPTM XL 12.5° Guide for 1.5 mm Anchor XL, Sawtooth	1/each
210219	VERSALOOPTM XL Reusable Straight Drill for 1.5 mm Anchor XL	1/each
210222	VERSALOOPTM XL Cannulated Obturator for 1.5 mm Anchor XL	1/each



Non-Surgical Solutions



Non-Surgical Solutions

MONOVISC™ High Molecular Weight Hyaluronan



CODE	DESCRIPTION	QTY
277515	MONOVISC High Molecular Weight Hyaluronan (Single Injection)	1/each

ORTHOVISC™ High Molecular Weight Hyaluronan




CODE	DESCRIPTION	QTY
277500	ORTHOVISC High Molecular Weight Hyaluronan (Multiple Injection)	1/each

PEAK™ Platelet Rich Plasma (PRP) System



CODE	DESCRIPTION	QTY
278001	PRP 30 ML DISP	1/each
278002	PRP 60 ML DISP	1/each
278010	PRP CENTRIFUGE UNIT	1/each



Contact a Sales Consultant for more information.
Visit [DePuySynthes.com](https://www.depuyorthopaedics.com).

The third-party trademarks used herein are the trademarks of their respective owners.
For further information on the products detailed in the Product Catalog, please refer to the relevant IFUs.



DePuy Mitek, Inc.
325 Paramount Drive
Raynham, MA 02767
T: +1 (800) 382-4682

[depuyorthopaedics.com](https://www.depuyorthopaedics.com)

©DePuy Synthes 2023. All rights reserved.
Version: 1.2 April 2023