

# Inventory Control Form

## QIS AND 2.0 MM

### CANNULATED SCREW SYSTEM PLUS (CSS+)

Case Information (eg, Case ID):

DO NOT PUT ANY  
PATIENT IDENTIFIABLE  
INFORMATION IN  
THIS FIELD

To order (USA): 800-523-0322

To order (Canada): 844-243-4321

Date: \_\_\_\_\_

Hospital: \_\_\_\_\_



Surgeon: \_\_\_\_\_

Procedure: \_\_\_\_\_

## Implants The screws are made of titanium alloy (TAV or TAN)

### 2.0 mm Quick Insertion Screw (TAN)

Length	Part Number
11 mm	04.228.511
12 mm	04.228.512
13 mm	04.228.513
14 mm	04.228.514
15 mm	04.228.515
16 mm	04.228.516
17 mm	04.228.517
18 mm	04.228.518



**Head Height:** 1.65 mm  
**Washer and Head Construct Height:** 1.44 mm

#### Thread Length

**Short Thread:**  
approx. 25% Total Length or 4 mm

**Long Thread:**  
approx. 40% Total Length

### 2.0 mm Cannulated Screw System + (TAV)

Length	Short Thread Part Number	Length	Long Thread Part Number
10 mm	04.353.010	N/A	
11 mm	04.353.011	N/A	
12 mm	04.353.012	N/A	
13 mm	04.353.013	N/A	
14 mm	04.353.014	N/A	
15 mm	04.353.015	N/A	
16 mm	04.353.016	N/A	
17 mm	04.353.017	N/A	
18 mm	04.353.018	N/A	
19 mm	04.353.019	N/A	
20 mm	04.353.020	20 mm	04.354.020
22 mm	04.353.022	22 mm	04.354.022
24 mm	04.353.024	24 mm	04.354.024
26 mm	04.353.026	26 mm	04.354.026
28 mm	04.353.028	28 mm	04.354.028
30 mm	04.353.030	30 mm	04.354.030
32 mm	04.353.032TS	32 mm	04.354.032TS
34 mm	04.353.034TS	34 mm	04.354.034TS
36 mm	04.353.036TS	36 mm	04.354.036TS
38 mm	04.353.038TS	38 mm	04.354.038TS
40 mm	04.353.040TS	40 mm	04.354.040TS

 Sterile Packed Only

#### Washers

Part Number	
04.353.900.05	Washer For 2.0 Cann Screw + 5pk
04.353.900.S	Washer For 2.0 Cann Screw + Sterile

All implants for this diameter are also available sterile packed. Sterile part availability might be different depending on the country or region.




# Inventory Control Form

## QIS AND 2.0 MM

### CANNULATED SCREW SYSTEM PLUS (CSS+)

#### Instruments

##### Instruments for 2.0 mm Screws

Part Number		
<input type="text"/>	03.333.000**	0.8 mm Guide Wire, L 100 mm Trocar Tip
		
<input type="text"/>	03.333.100	1.6 mm Cannulated Drill Bit, QC
<input type="text"/>	03.333.100S*	
		
<input type="text"/>	03.333.210	2.0 mm Cannulated Headed Countersink, QC
<input type="text"/>	03.333.210S*	
		
<b>Screwdriver: T6 StarDrive</b>		
<input type="text"/>	03.333.300	T6 Cannulated Stardrive Screwdriver Shaft/Self-Retaining
		

\*Part numbers with a trailing 'S' or 'TS' designate the sterile part number. Sterile part availability might be different depending on the country or region.

\*\* All guide wires are also available as a four (4) pack. Add ".04" to part number. All guide wires are also available sterile packed (6 pack). Add suffix "S" to part number. Sterile part availability might be different depending on the country or region.

Some devices in this form may not have been licensed in accordance with Canadian law and may not be for sale in Canada. Please contact your sales consultant for items approved for sale in Canada.

CSSPLUSICF001V2 6/22

# Inventory Control Form

## 2.5 and 3.0 MM

### CANNULATED SCREW SYSTEM PLUS (CSS+)

#### Implants The screws are made of titanium alloy (TAV) (cont'd)


##### 2.5 mm Cannulated Screw System +

Length	Short Thread Part Number	Length	Long Thread Part Number
10 mm	04.353.110	N/A	
11 mm	04.353.111	N/A	
12 mm	04.353.112	N/A	
13 mm	04.353.113	N/A	
14 mm	04.353.114	N/A	
15 mm	04.353.115	N/A	
16 mm	04.353.116	N/A	
17 mm	04.353.117	N/A	
18 mm	04.353.118	N/A	
19 mm	04.353.119	N/A	
20 mm	04.353.120	20 mm	04.354.120
22 mm	04.353.122	22 mm	04.354.122
24 mm	04.353.124	24 mm	04.354.124
26 mm	04.353.126	26 mm	04.354.126
28 mm	04.353.128	28 mm	04.354.128
30 mm	04.353.130	30 mm	04.354.130
32 mm	04.353.132	32 mm	04.354.132
34 mm	04.353.134	34 mm	04.354.134
36 mm	04.353.136	36 mm	04.354.136
38 mm	04.353.138	38 mm	04.354.138
40 mm	04.353.140	40 mm	04.354.140
45 mm	04.353.145TS	45 mm	04.354.145TS
50 mm	04.353.150TS	50 mm	04.354.150TS
55 mm	04.353.155TS	55 mm	04.354.155TS

 Sterile Packed Only

##### Washers

Part Number	
04.353.901.05	Washer For 2.5 Cann Screw + 5pk
04.353.901S	Washer For 2.5 Cann Screw + Sterile



All implants for this diameter are also available sterile packed. Sterile part availability might be different depending on the country or region.

##### 3.0 mm Cannulated Screw System +

Length	Short Thread Part Number	Length	Long Thread Part Number
10 mm	04.353.210	N/A	
11 mm	04.353.211	N/A	
12 mm	04.353.212	N/A	
13 mm	04.353.213	N/A	
14 mm	04.353.214	N/A	
15 mm	04.353.215	N/A	
16 mm	04.353.216	N/A	
17 mm	04.353.217	N/A	
18 mm	04.353.218	N/A	
19 mm	04.353.219	N/A	
20 mm	04.353.220	20 mm	04.354.220
22 mm	04.353.222	22 mm	04.354.222
24 mm	04.353.224	24 mm	04.354.224
26 mm	04.353.226	26 mm	04.354.226
28 mm	04.353.228	28 mm	04.354.228
30 mm	04.353.230	30 mm	04.354.230
32 mm	04.353.232	32 mm	04.354.232
34 mm	04.353.234	34 mm	04.354.234
36 mm	04.353.236	36 mm	04.354.236
38 mm	04.353.238	38 mm	04.354.238
40 mm	04.353.240	40 mm	04.354.240
45 mm	04.353.245TS	45 mm	04.354.245TS
50 mm	04.353.250TS	50 mm	04.354.250TS
55 mm	04.353.255TS	55 mm	04.354.255TS

 Sterile Packed Only

##### Washers

Part Number	
04.353.902.05	Washer For 3.0 Cann Screw + 5pk
04.353.902S	Washer For 3.0 Cann Screw + Sterile



All implants for this diameter are also available sterile packed. Sterile part availability might be different depending on the country or region.

# Inventory Control Form

## 2.5 and 3.0 MM


### CANNULATED SCREW SYSTEM PLUS (CSS+)

#### Instruments

Instruments for 2.5/3.0 mm Screws  

Part Number		
<input type="checkbox"/>	03.333.001**	1.1 mm Guide Wire, L 150 mm, Trocar Tip
<input type="checkbox"/>	03.333.011**	1.1 mm Guide Wire, L 150 mm, Threaded Tip
		
<input type="checkbox"/>	03.333.101	2.0 mm Cannulated Drill Bit, QC
<input type="checkbox"/>	03.333.101S*	
		
<input type="checkbox"/>	03.333.211	2.5/3.0 mm Cannulated Headed Countersink, QC
<input type="checkbox"/>	03.333.211S*	
		

Screwdriver: T8 StarDrive

<input type="checkbox"/>	03.333.302	T8 Cannulated Stardrive Screwdriver Shaft/Self-Retaining	
--------------------------	------------	--	--

\*Part numbers with a trailing 'S' or 'TS' designate the sterile part number. Sterile part availability might be different depending on the country or region.

\*\* All guide wires are also available as a four (4) pack. Add ".04" to part number. All guide wires are also available sterile packed (6 pack). Add suffix "S" to part number. Sterile part availability might be different depending on the country or region.

Some devices in this form may not have been licensed in accordance with Canadian law and may not be for sale in Canada. Please contact your sales consultant for items approved for sale in Canada.

CSSPLUSICF001V2 6/22

# Inventory Control Form

## 3.5 MM

### CANNULATED SCREW SYSTEM PLUS (CSS+)

#### Implants The screws are made of titanium alloy (TAV) (cont'd)

##### 3.5 mm Cannulated Screw System +

Length	Short Thread Part Number	Length	Long Thread Part Number	Length	Full Thread Part Number
	10 mm 04.353.310	N/A			10 mm 04.355.310
	12 mm 04.353.312	N/A			12 mm 04.355.312
	14 mm 04.353.314	N/A			14 mm 04.355.314
	16 mm 04.353.316	N/A			16 mm 04.355.316
	18 mm 04.353.318	N/A			18 mm 04.355.318
	20 mm 04.353.320		20 mm 04.354.320		20 mm 04.355.320
	22 mm 04.353.322		22 mm 04.354.322		22 mm 04.355.322
	24 mm 04.353.324		24 mm 04.354.324		24 mm 04.355.324
	26 mm 04.353.326		26 mm 04.354.326		26 mm 04.355.326
	28 mm 04.353.328		28 mm 04.354.328		28 mm 04.355.328
	30 mm 04.353.330		30 mm 04.354.330		30 mm 04.355.330
	32 mm 04.353.332		32 mm 04.354.332		32 mm 04.355.332
	34 mm 04.353.334		34 mm 04.354.334		34 mm 04.355.334
	36 mm 04.353.336		36 mm 04.354.336		36 mm 04.355.336
	38 mm 04.353.338		38 mm 04.354.338		38 mm 04.355.338
	40 mm 04.353.340		40 mm 04.354.340		40 mm 04.355.340
	42 mm 04.353.342		42 mm 04.354.342		42 mm 04.355.342
	44 mm 04.353.344		44 mm 04.354.344		44 mm 04.355.344
	46 mm 04.353.346		46 mm 04.354.346		46 mm 04.355.346
	48 mm 04.353.348		48 mm 04.354.348		48 mm 04.355.348
	50 mm 04.353.350		50 mm 04.354.350		50 mm 04.355.350
	52 mm 04.353.352		52 mm 04.354.352		52 mm 04.355.352
	54 mm 04.353.354		54 mm 04.354.354		54 mm 04.355.354
	56 mm 04.353.356		56 mm 04.354.356		56 mm 04.355.356
	58 mm 04.353.358		58 mm 04.354.358		58 mm 04.355.358
	60 mm 04.353.360		60 mm 04.354.360		60 mm 04.355.360
	65 mm 04.353.365S		65 mm 04.354.365S		65 mm 04.355.365S
	70 mm 04.353.370S		70 mm 04.354.370S		70 mm 04.355.370S
	75 mm 04.353.375S		75 mm 04.354.375S		75 mm 04.355.375S

 Sterile Packed Only

##### Washers

Part Number	
04.353.903.05	Washer For 3.5 Cann Screw + 5pk
04.353.903S	Washer For 3.5 Cann Screw + Sterile



All implants for this diameter are also available sterile packed. Sterile part availability might be different depending on the country or region.

# Inventory Control Form

## 4.0 MM

### CANNULATED SCREW SYSTEM PLUS (CSS+)

#### Implants The screws are made of titanium alloy (TAV) (cont'd)

##### 4.0 mm Cannulated Screw System +

Length		Short Thread Part Number	Length		Long Thread Part Number	Length		Full Thread Part Number
	10 mm	04.353.410	N/A				10 mm	04.355.410
	12 mm	04.353.412	N/A				12 mm	04.355.412
	14 mm	04.353.414	N/A				14 mm	04.355.414
	16 mm	04.353.416	N/A				16 mm	04.355.416
	18 mm	04.353.418	N/A				18 mm	04.355.418
	20 mm	04.353.420		20 mm	04.354.420		20 mm	04.355.420
	22 mm	04.353.422		22 mm	04.354.422		22 mm	04.355.422
	24 mm	04.353.424		24 mm	04.354.424		24 mm	04.355.424
	26 mm	04.353.426		26 mm	04.354.426		26 mm	04.355.426
	28 mm	04.353.428		28 mm	04.354.428		28 mm	04.355.428
	30 mm	04.353.430		30 mm	04.354.430		30 mm	04.355.430
	32 mm	04.353.432		32 mm	04.354.432		32 mm	04.355.432
	34 mm	04.353.434		34 mm	04.354.434		34 mm	04.355.434
	36 mm	04.353.436		36 mm	04.354.436		36 mm	04.355.436
	38 mm	04.353.438		38 mm	04.354.438		38 mm	04.355.438
	40 mm	04.353.440		40 mm	04.354.440		40 mm	04.355.440
	42 mm	04.353.442		42 mm	04.354.442		42 mm	04.355.442
	44 mm	04.353.444		44 mm	04.354.444		44 mm	04.355.444
	46 mm	04.353.446		46 mm	04.354.446		46 mm	04.355.446
	48 mm	04.353.448		48 mm	04.354.448		48 mm	04.355.448
	50 mm	04.353.450		50 mm	04.354.450		50 mm	04.355.450
	52 mm	04.353.452		52 mm	04.354.452		52 mm	04.355.452
	54 mm	04.353.454		54 mm	04.354.454		54 mm	04.355.454
	56 mm	04.353.456		56 mm	04.354.456		56 mm	04.355.456
	58 mm	04.353.458		58 mm	04.354.458		58 mm	04.355.458
	60 mm	04.353.460		60 mm	04.354.460		60 mm	04.355.460
	65 mm	04.353.465S		65 mm	04.354.465S		65 mm	04.355.465S
	70 mm	04.353.470S		70 mm	04.354.470S		70 mm	04.355.470S
	75 mm	04.353.475S		75 mm	04.354.475S		75 mm	04.355.475S

 Sterile Packed Only

##### Washers

Part Number	
04.353.904.05	Washer For 4.0 Cann Screw + 5pk
04.353.904S	Washer For 4.0 Cann Screw + Sterile



**Head Height:** 2.7 mm  
**Washer and Head Construct Height:** 2.29 mm

##### Thread Length

**Short Thread:**  
approx. 25% Total Length or 5 mm

**Long Thread:**  
approx. 40% Total Length

**Full Thread:**  
Total Length minus 5 mm

All implants for this diameter are also available sterile packed. Sterile part availability might be different depending on the country or region.

# Inventory Control Form

## 3.5 and 4.0 MM

### CANNULATED SCREW SYSTEM PLUS (CSS+)

#### Instruments

##### Instruments for 3.5/4.0 mm Screws



##### Part Number

	03.333.002	1.4 mm Guide Wire, L 100 mm Trocar Tip
	03.333.012 **	1.4 mm Guide Wire, L 150 mm, Threaded Tip



	03.333.102	2.7 mm Cannulated Drill Bit, QC
	03.333.102S*	



	03.333.212	3.5/4.0 mm Cannulated Headed Countersink, QC
	03.333.212S*	



##### Screwdriver: T15 StarDrive

	03.333.304	T15 Cannulated StarDrive Screwdriver Shaft/Self-Retaining
--	------------	---



\*Part numbers with a trailing 'S' or 'TS' designate the sterile part number. Sterile part availability might be different depending on the country or region.

\*\* All guide wires are also available as a four (4) pack. Add ".04" to part number. All guide wires are also available sterile packed (6 pack). Add suffix "S" to part number. Sterile part availability might be different depending on the country or region.

Some devices in this form may not have been licensed in accordance with Canadian law and may not be for sale in Canada. Please contact your sales consultant for items approved for sale in Canada.



Manufactured by:  
**Synthes USA, LLC**  
1101 Synthes Avenue  
Monument, CO 80132

Manufactured by:  
**Tyber Medical, LLC**  
83 South Commerce Way, Suite 310  
Bethlehem, PA 18017

To order (USA): 800-523-0322  
To order (Canada): (844)-243-4321

[www.tybermedical.com](http://www.tybermedical.com)

[jnmedicaldevices.com](http://jnmedicaldevices.com)