



DEVICE SPECS

ECHELON CONTOUR™ Curved Cutter

		
	ECHELON CONTOUR™ Curved Cutter	ECHELON CONTOUR™ Curved Cutter
Codes	GCS40G	GCS40B
Associated Reload	GCR40G	GCR40B
Description	40mm Curved Cutter with Green Reload	40mm Curved Cutter with Blue Reload
Device Technology	Curved Head	Curved Head
Size	40mm	40mm
Max Firings	6	6
Staple Rows	4	4
Closed Staple Height	2.0mm	1.5mm
Tissue Thickness	Thick	Regular
Staple Technology	GST + 3D	GST + 3D
Qty./Box	3 (Curved Cutter), 6 (Reloads)	3 (Curved Cutter), 6 (Reloads)