

## DEVICE SPECS

# ENDOPATH® ETS Articulating Linear Cutters



**ENDOPATH™ ETS Articulating Linear Cutters**

<b>Codes</b>	ATS45
<b>Associated Reloads</b>	TR45W, 6R45B, TR45G
<b>Description</b>	45mm ETS Articulating Linear Cutter
<b>Firing</b>	Single Stroke
<b>Articulation</b>	Manual
<b>Shaft Length</b>	34cm
<b>Degrees of Articulation</b>	70°
<b>Staple Rows</b>	4 (when using TR45G) 6 (when using TR45W or 6R45B)
<b>Max Firings</b>	8
<b>Size</b>	45mm
<b>Staple Form</b>	B-Form
<b>Qty./Box</b>	3