

# ENTERPRISE™ Vascular Reconstruction Device and Delivery System

## Patient Information

Your surgeon has implanted the CERENOVUS ENTERPRISE™ Intracranial Stent.

## Device Description

The ENTERPRISE Device is used to treat aneurysms in the head. The device is aimed to help with the placement of another device called microcoils. The microcoils are used to treat aneurysms.

## Device materials

Some people may be allergic to the metal part of the device. Make sure to let your doctor know before your surgery if you have allergies to nickel titanium (Nitinol).

## Information for Safe Use

Be sure to follow your doctors' orders for general safety. Magnetic Resonance Imaging (MRI) is a test used to diagnose certain diseases. Also, it can be used during medical procedures. It is very important to tell your doctors that you have an implanted device before having a test like an MRI. Device implant and MRI testing information can be requested from Cerenovus.

## Patient Information Portal

Any updated information will be provided on our website <https://www.jnjmedicaldevices.com/en-AU/patient-information-leaflets>

Information specific to your implant, including the serial number, unique device identifier, etc., are included on the implant card as well as in patient records kept by your healthcare provider.

## Reporting Adverse Effects

If you wish to report any adverse effects that you believe are a result of your surgery, please speak with your medical team or report the information to Johnson & Johnson Medical Product Safety Department at [productsafetyjmanz@its.jnj.com](mailto:productsafetyjmanz@its.jnj.com) or the Therapeutic Goods Administration at <https://www.tga.gov.au>

For access to this leaflet, please visit: <https://www.jnjmedicaldevices.com/en-AU/patient-information-leaflets>

## Sponsor

Johnson & Johnson Medical Pty Ltd t/a DePuy Synthes  
1-5 Khartoum Rd, North Ryde, NSW 2113  
Tel: 1300 562 711