

DEVICE SPECS

Ethicon Laparoscopic Stapler Portfolio



ECHELON™ 3000 Stapling System



ECHELON™+ Stapling System



ECHELON FLEX™ ENDOPATH™ Staplers

	ECHELON™ 3000 Stapling System	ECHELON™+ Stapling System	ECHELON FLEX™ ENDOPATH™ Staplers
Codes	ECH45C, ECH45S, ECH45L, ECH60C, ECH60S, ECH60L	PCEE45A, PSEE45A, PLEE45A, PCEE60A, PSEE60A, PLEE60A	SC45A, EC45A, EC45AL, SC60A, EC60A, LONG60A
Associated Reloads	GST45W, GST60W, GST45B, GST60B, GST45D, GST60D, GST45G, GST60G, GST45T, GST60T	GST45W, GST60W, GST45B, GST60B, GST45D, GST60D, GST45G, GST60G, GST45T, GST60T	GST45W, GST60W, GST45B, GST60B, GST45D, GST60D, GST45G, GST60G, GST45T, GST60T, ECR45M, ECR60M
Description	45mm and 60mm laparoscopic staplers	45mm and 60mm laparoscopic staplers	45mm and 60mm laparoscopic staplers
Firing	Powered	Powered	Manual
Articulation	Powered	Manual	Manual
Shaft Length	Compact 28cm, Standard 34cm, Long 44cm	Compact 28cm, Standard 34cm, Long 44cm	Compact 28cm, Standard 34cm, Long 44cm
Degrees of Articulation	>110°	88°	--
Staple Rows	6	6	6
Max Firings	12	12	12
Staple Form	B-Form	B-Form	B-Form
Qty./Box	3	3	3