






**DEVICE SPECS**

# Ethicon Laparoscopic Stapler GST Reloads

	 <b>ECHELON ENDOPATH™</b> Reload (White)	 <b>ECHELON ENDOPATH™</b> Reload (Blue)	 <b>ECHELON ENDOPATH™</b> Reload (Gold)	 <b>ECHELON ENDOPATH™</b> Reload (Green)	 <b>ECHELON ENDOPATH™</b> Reload (Black)
<b>Codes</b>	GST45W, GST60W	GST45B, GST60B	GST45D, GST60D	GST45G, GST60G	GST45T, GST60T
<b>Description</b>	Powered and manual stapler reload	Powered and manual stapler reload	Powered and manual stapler reload	Powered and manual stapler reload	Powered and manual stapler reload
<b>Device Technology</b>	Gripping Surface Technology	Gripping Surface Technology	Gripping Surface Technology	Gripping Surface Technology	Gripping Surface Technology
<b>Sizes</b>	45mm, 60mm	45mm, 60mm	45mm, 60mm	45mm, 60mm	45mm, 60mm
<b>Tissue Thickness Range</b>	1.0mm - 2.0mm	1.5mm - 2.4mm	1.8mm - 3.0mm	2.0mm - 3.3mm	2.3mm - 4.0mm
<b>Tissue Thickness Description</b>	Thin	Regular	Regular/Thick	Thick	Very Thick
<b>Closed Staple Height</b>	1.0mm	1.5mm	1.8mm	2.0mm	2.3mm
<b>Open Staple Height</b>	2.6mm	3.6mm	3.8mm	4.1mm	4.2mm
<b>Staple Type</b>	B-Form	B-Form	B-Form	B-Form	B-Form
<b>Qty./Box</b>	12	12	12	12	12