

DEVICE SPECS

# VICRYL RAPIDE™ Suture and Ethicon Gut Sutures



**VICRYL RAPIDE™**  
(polyglactin 910) Suture



**Fast-Absorbing Surgical Gut Suture (Plain)**



**Surgical Gut Suture (Plain)**



**Surgical Gut Suture (Chromic)**

<b>Codes</b>	<a href="#">See catalog</a>	1915G	<a href="#">See catalog</a>	<a href="#">See catalog</a>
<b>Color</b>	Undyed	Yellowish Tan	Yellowish Tan, Brown, or Blue	Blue or Undyed
<b>Antibacterial Coating</b>	No	No	No	No
<b>Material/Composition</b>	Polyglactin 910	Natural, Collagen	Natural, Collagen	Natural, Collagen
<b>Suture Structure</b>	Braided	Monofilament	Monofilament	Monofilament
<b>Absorbable</b>	Yes	Yes	Yes	Yes
<b>Absorption Rate</b>	42 days	21-42 days	70 days	90 days
<b>Directionality</b>	Unidirectional	Unidirectional	Unidirectional	Unidirectional
<b>Available Needle Options</b>	Taperpoint, Tapercut, Cutting Edge Reverse, Cutting Edge Prime Reverse/Conventional	Cutting Edge Conventional	Taperpoint, Tapercut, Cutting Edge Reverse, Cutting Edge Prime Reverse, Cutting Edge Spatula Reverse	Taperpoint, Tapercut, Reverse Cutting, Cutting Edge Reverse, Cutting Edge Ultima Reverse, Cutting Edge Spatula Reverse, Cutting Edge Prime Reverse, Blunt Point
<b>Special Options</b>	MultiPass™	MultiPass™	MultiPass™, CONTROL RELEASE™	MultiPass™, CONTROL RELEASE™
<b>Pledget Options</b>	No	No	No	No
<b>Strands per Pack</b>	1	1	1, 3, 8, or 12	1, 3, 4, 8, or 12
<b>Quantity per Box</b>	12	12	12, 24, or 36	12, 24, or 36