

L107.960

Grosser Fixateur externe

Box01/01

Pos.	Mat.Nr	Description	Qty	Unit	Location	EAN/GTIN-Code
Box01/01		Module 01				

1	390.008	Backe, aufsteckbar, selbsthaltend	10	ST	B01M01	07611819760592
---	---------	-----------------------------------	----	----	--------	----------------



2	68.301.011	Modul für Backen, für Grossen Fixateur externe, zu Nr. 68.301.000	1	ST	B01M01	07611819186880
---	------------	---	---	----	--------	----------------



Pos.	Mat.Nr	Description	Qty	Unit	Location	EAN/GTIN-Code
Box01/01		Module 02				

3	310.370	Spiralbohrer Ø 3.5 mm, Länge 195/170 mm, 2-lippig, für Schnellkupplung	2	ST	B01M02	07611819018730
---	---------	--	---	----	--------	----------------



4	393.910	Verbindungsstab Ø 5.0 mm, Länge 150 mm, Stahl	2	ST	B01M02	07611819025523
---	---------	---	---	----	--------	----------------



5	686.107	Einsatz für Schanzsche Schrauben Ø 5.0 mm, ohne Deckel, ohne Inhalt	1	ST	B01M02	07611819080911
---	---------	---	---	----	--------	----------------



6	686.109	Deckel, zu Nr. 686.107	1	ST	B01M02	07611819080898
---	---------	------------------------	---	----	--------	----------------



Pos.	Mat.Nr	Description	Qty	Unit	Location	EAN/GTIN-Code
Box01/01		Tray 01				

7	390.002	Multi-Pin-Backe, 6 Positionen, gross	4	ST	B01T01	07611819760554
---	---------	--------------------------------------	---	----	--------	----------------



8	390.007	Stab-zu-Stab-Backe	4	ST	B01T01	07611819760585
---	---------	--------------------	---	----	--------	----------------

















9	68.301.000	Vario Case™ für Grossen Fixateur externe, ohne Deckel, ohne Inhalt	1	ST	B01T01	07611819186866
---	------------	--	---	----	--------	----------------



L107.960

Grosser Fixateur externe





Pos.	Mat.Nr	Description	Qty	Unit	Location	EAN/GTIN-Code
10	689.507	Stahldeckel, Grösse 1/1, für Vario Case™	1	ST	B01T01	07611819108097
						
Box01/01 Tray 02						
11	321.160	Ringgabelschlüssel Ø 11.0 mm	1	ST	B01T02	07611819021075
						
12	355.140	Steckschlüssel Ø 11.0 mm, durchbohrt	1	ST	B01T02	07611819023086
						
13	393.100	Universalbohrfutter mit T-Griff	1	ST	B01T02	07611819025196
						
14	393.103	Adapter für Seldrill™ Schanzsche Schraube	2	ST	B01T02	07611819775770
						
15	395.911	Handgriff für Bohrbüchse	1	ST	B01T02	07611819739024
						
16	395.921	Bohrbüchse 6.0/5.0, kurz, mit Gewinde	1	ST	B01T02	07611819739055
						
17	395.923	Bohrbüchse 6.0/5.0, lang, mit Gewinde	1	ST	B01T02	07611819739079
						
18	395.912	Bohrbüchse 5.0/3.5, kurz	1	ST	B01T02	07611819739031
						
19	395.913	Bohrbüchse 5.0/3.5, lang	1	ST	B01T02	07611819739048
						
20	394.055	Stab mit Scharniergelenk für Fixateur externe	1	ST	B01T02	07611819815070
						
21	394.181	Trokar Ø 3.5 mm, kurz	1	ST	B01T02	07611819738997
						
22	394.182	Trokar Ø 3.5 mm, lang	1	ST	B01T02	07611819739000
						
68.301.012		Einsatz 1/1 Nr. 1 zu Vario Case 68.301.000	1	ST	B01T02	
						

L107.960

Grosser Fixateur externe

Box01/01

Tray 03

23	393.500	Stahlrohr Ø 11.0 mm, Länge 100 mm		2	ST	B01T03	07611819025301
24	393.520	Stahlrohr Ø 11.0 mm, Länge 150 mm		2	ST	B01T03	07611819025325
25	393.530	Stahlrohr Ø 11.0 mm, Länge 200 mm		2	ST	B01T03	07611819025332
26	393.540	Stahlrohr Ø 11.0 mm, Länge 250 mm		2	ST	B01T03	07611819025349
27	393.550	Stahlrohr Ø 11.0 mm, Länge 300 mm		2	ST	B01T03	07611819025356
28	393.560	Stahlrohr Ø 11.0 mm, Länge 350 mm		2	ST	B01T03	07611819025363
29	393.570	Stahlrohr Ø 11.0 mm, Länge 400 mm		2	ST	B01T03	07611819025370
30	393.760	Druckspanner, offen		2	ST	B01T03	07611819025455
31	394.800	Kohlefaserstab Ø 11.0 mm, Länge 100 mm		2	ST	B01T03	07611819025905
32	394.820	Kohlefaserstab Ø 11.0 mm, Länge 150 mm		2	ST	B01T03	07611819025929
33	394.830	Kohlefaserstab Ø 11.0 mm, Länge 200 mm		2	ST	B01T03	07611819025936
34	394.840	Kohlefaserstab Ø 11.0 mm, Länge 250 mm		2	ST	B01T03	07611819025943
35	394.850	Kohlefaserstab Ø 11.0 mm, Länge 300 mm		2	ST	B01T03	07611819025950
36	394.860	Kohlefaserstab Ø 11.0 mm, Länge 350 mm		2	ST	B01T03	07611819025967
37	394.870	Kohlefaserstab Ø 11.0 mm, Länge 400 mm		2	ST	B01T03	07611819025974
38	394.970.10	Abschlusskappe f/Stäbe+Rohre (10 ST/1 Pk.)		1	PAK	B01T03	07611819322035
68.301.013	Einsatz 1/1 Nr. 2 zu Vario Case 68.301.000			1	ST	B01T03	07611819186903

L107.960

Grosser Fixateur externe


B01M01



B01M02




B01T01




L107.960

Grosser Fixateur externe

B01T02



B01T03



End of document