






DEVICE SPECS

Ethicon Ligation Portfolio

	 <p>LIGAMAX™5 Endoscopic Multiple Clip Applier</p>	 <p>LIGACLIP™ Endoscopic Rotating Multiple Clip Applier</p>	 <p>LIGACLIP™ Multiple Clip Applier</p>	 <p>LIGACLIP™ Endoscopic Multi-Patient Use Single Clip Appliers</p>	 <p>LIGACLIP™ Multi-Patient Use Single Clip Appliers</p>
Codes	EL5ML	ER320, ER420	MCL20, MCM20, MCS20, MSM20, MCM30	EL214, EL314, EL414	LC307, LC310, LC320, LC407, LC410, LC420, LC430, LX105, LX107, LX110, LX130, LX205, LX207, LX210, LX220, LX230
Usage	Single-patient use	Single-patient use	Single-patient use	Multi-patient use	Multi-patient use
Open or Endoscopic	Endoscopic	Endoscopic	Open	Endoscopic	Open
Length	13in	11.4in (ER320) 13.4in (ER420)	9.375in (MCS20, MSM20) 11in (MCM20, MCM30) 13in (MCL20)	14in.	5.75in (LX105, LX205) 7.5in (LC307, LC407, LX107, LX207) 10.5in (LC310, LC320, LC410, LC420, LC430, LX110, LX130, LX210, LX220, LX230)
Diameter	5mm	10mm (ER320) 12mm (ER420)	--	--	--
Clip Size	Medium/Large	Medium/Large (ER320) Large (ER420)	Small (MCS20) Medium (MCM20, MSM20, MCM30) Large (MCL20)	Medium (EL214) Medium/Large (EL314) Large (EL414)	Small (LX105, LX107, LX110, LX130) Medium (LX205, LX207, LX210, LX220, LX230) Medium/Large (LC307, LC310, LC320) Large (LC407, LC410, LC420, LC430)
Number of Clips	15	20	20 (codes ending in 20) 30 (MCM30)	--	--
Handle Color	--	--	--	--	Blue (LX105, LX107, LX110, LX130) Silver (LX205, LX207, LX210, LX220, LX230) Green (LC307, LC310, LC320) Yellow (LC407, LC410, LC420, LC430)
Qty./Box	3	3	6	1	1