

DEVICE SPECS

Ethicon Nylon Sutures



ETHILON™
Nylon Suture



NUROLON™
Nylon Suture

Codes	See catalog	See catalog
Color	Black, Green, or Undyed	Black or White
Antibacterial Coating	No	No
Material/Composition	Polyamide 6	Polyamide 6.6
Suture Structure	Monofilament	Braided
Absorbable	No	No
Absorption Rate	--	--
Directionality	Unidirectional	Unidirectional
Available Needle Options	Taperpoint, Tapercut, Cutting Edge Reverse/Conventional, Cutting Edge Spatulated Center Point, Cutting Edge Spatula Reverse, Cutting Edge Prime Reverse/Conventional, Cutting Edge Ultima Reverse, Blunt Point	Taperpoint, Cutting Edge Reverse, Cutting Edge Prime Reverse/Conventional
Special Options	MultiPass™, CONTROL RELEASE™	MultiPass™, CONTROL RELEASE™
Pledget Options	No	No
Strands per Pack	1 or 8	1, 2, 3, 5, 6, 8, 10, or 12
Quantity per Box	12, 24, or 36	12, 24, or 36