




DEVICE SPECS

Ethicon PLUS™ Antibacterial Sutures

	 <p>MONOCRYL™ Plus Antibacterial (poliglecaprone 25) Suture</p>	 <p>Coated VICRYL™ Plus Antibacterial (polyglactin 910) Suture</p>	 <p>PDS™ Plus Antibacterial (polydioxanone) Suture</p>
Codes	See catalog	See catalog	See catalog
Color	Violet or Undyed	Violet or Undyed	Violet or Undyed
Antibacterial Coating	Yes	Yes	Yes
Material/Composition	Poliglecaprone 25, Triclosan	Polyglactin 910, Triclosan	Polydioxanone, Triclosan
Suture Structure	Monofilament	Braided	Monofilament
Absorbable	Yes	Yes	Yes
Absorption Rate	90-120 days	56-70 days	182-238 days
Directionality	Unidirectional	Unidirectional or Bidirectional	Unidirectional
Available Needle Options	Taperpoint, Cutting Edge Reverse, Cutting Edge Prime/Conventional, Blunt Point Ethiguard	Taperpoint, Tapercut, Cutting Edge Reverse, Cutting Edge Prime Reverse/Conventional, Blunt Point Ethiguard	Taperpoint, Tapercut, Cutting Edge Reverse, Cutting Edge Prime Reverse/Conventional, Blunt Point Ethiguard, Blunt Point
Special Options	MultiPass™	MutliPass™, CONTROL RELEASE™	MultiPass™
Pledget Options	No	No	No
Strands per Pack	1	1, 3, 4, 6, 8, or 12	1, or 8
Quantity per Box	12 or 36	2, 12, 24, or 36	12, 24, or 36