

PRODUCT SPECS

Ethicon PROLENE™ Sutures



PROLENE™ Polypropylene Suture with HEMO-SEAL™ Technology



PROLENE™ Polypropylene Suture

Codes	See catalog	See catalog
Color	Blue	Blue or Undyed
Antibacterial Coating	No	No
Material/Composition	Polypropylene	Polypropylene
Suture Structure	Monofilament	Monofilament
Absorbable	No	No
Absorption Rate	--	--
Directionality	Unidirectional	Unidirectional
Available Needle Options	Taperpoint	Taperpoint, Taperpoint Cardio, Tapercut, Cutting Edge Reverse, Cutting Edge Prime Reverse/Conventional, Cutting Edge Spatula Reverse/Conventional, Cutting Edge Ultima Reverse, Blunt Point Ethiguard, Blunt Point
Special Options	EVERPOINT™, HEMO-SEAL™	EVERPOINT™, VISI-BLACK™, CONTROL RELEASE™, MultiPass™
Pledget Options	No	Yes
Strands per Pack	1, 2, 4	1, 2, 3, 4, 5, 6, 8, 10
Quantity per Box	12 or 36	6, 12, 24, or 36