





DEVICE SPECS

Ethicon Skin Stapler Portfolio

	 <p>PROXIMATE™ Plus MD Skin Stapler</p>	 <p>PROXIMATE™ Rotating-Head Skin Staplers</p>	 <p>PROXIMATE™ Fixed-Head Skin Staplers</p>	 <p>PROXIMATE™ Skin Staple Extractor</p>
Codes	PMW35	PRR35, PRW35	PXW35	PSX
Description	Sterile, single-patient use skin stapler with 0.58 diameter	Sterile, single-patient use rotating-head skin stapler with 0.53 or 0.58 diameter	Sterile, single-patient use regular fixed-head skin stapler with 0.58 diameter	Sterile, single-patient use squeeze handle staple extractor
Size	Wide	Regular, Wide	Wide	--
Wire Diameter	0.58	0.53, 0.58	0.58	--
Staple Leg Length	3.9	3.9	3.9	--
Qty./Box	6	6	6	12