

DEVICE SPECS

# BIOPATCH™ Protective Disk with CHG

	BIOPATCH™ Protective Disk with CHG	BIOPATCH™ Protective Disk with CHG	BIOPATCH™ Protective Disk with CHG
<b>Code</b>	4150	4151	4152
<b>Disk Diameter</b>	1" (2.5cm)	3/4" (1.9cm)	1" (2.5cm)
<b>Center Hole Diameter</b>	4.0mm	1.5mm	7.0mm
<b>French Size Range</b>	6-12Fr	<6Fr	13-20Fr
<b>Average CHG per disk</b>	92mg	52.5mg	86.8mg
<b>Common Uses</b>	Central Lines, PICC	Peripheral IVs, Huber Needles (port), Arterial Lines, Extended Dwell PIVs, Midlines, PICCs, Pins	Dialysis Catheters, Drains, Sheaths, Cordis Catheters, VAD Drive Lines
<b>Qty. / Box</b>	10	10	10
<b>Qty. / Case</b>	40 (4 boxes)	40 (4 boxes)	40 (4 boxes)