

DEVICE SPECS

Ethicon Circular Stapler Portfolio



ECHELON CIRCULAR™ Powered Stapler



ETHICON™ Circular Stapler, XL Sealed



ETHICON™ Circular Stapler

Codes	CDH23P, CDH25P, CDH29P, CDH31P	ECS21B, ECS25B, ECS29B, ECS33B	CDH21B, CDH25B, CDH29B, CDH33B
Description	Powered Circular Stapler with Gripping Surface Technology and 3D Staples	Circular Stapler, XL Sealed	Circular Stapler
Device Technology	Powered with Adjustable Staple Height	Adjustable Staple Height, XL shaft	Adjustable Staple Height, Standard shaft
Lumen Diameter	23mm, 25mm, 29mm, 31 mm	21mm, 25mm, 29mm, 33mm	21mm, 25mm, 29mm, 33mm
Staple Rows	2	2	2
Open Staple Leg Length	5.2mm	5.2mm	5.2mm
Closed Staple Height	1.5mm-2.2mm	1.5mm-2.2mm	1.5mm-2.2mm
Staple Technology	GST + 3D	B-Form	B-Form
Qty./Box	3	3	3