

DEVICE SPECS

MEGADYNE™ Electrosurgical Pencils Portfolio



	MEGADYNE™ Telescoping Soft Tissue Dissector	MEGADYNE™ Electrocauter Pencils (EZ-CLEAN)	MEGADYNE™ Electrocauter Pencils (Stainless Steel)
Codes	ME7251T, ME725M1T	0030, 0030L, 0030H, 0030HL, 0035, 0035L, 0035H, 0035HL, 0036, 0036H, 0037, 0037H	0038, 0038H, 0039, 0039H, FC39H
Compatible With	MEGEN1	MEGEN1	MEGEN1
Description	Electrosurgical soft tissue dissector with telescoping shaft	Electrosurgical pencils with either: <ul style="list-style-type: none"> • button or rocker switch • no holster or holster • standard or modified E-Z clean tip 	Electrosurgical pencils with either: <ul style="list-style-type: none"> • button switch, rocker switch, or foot control • no holster or holster
Cord Length	10 ft.	10 ft. (all other codes) 15 ft. (0030L, 0030HL, 0035L, 0035HL)	10 ft.
Modality	Monopolar	Monopolar	Monopolar
Tip Configuration	Blade or Modified Blade	Blade or Modified Blade	Blade or Modified Blade
Device Technology	Geometric Electron Modulation (GEM)	PTFE-Coated Blade	Standard Monopolar Electrosurgery
Surgical Approach	Open	Open & Laparoscopic	Open
Qty./Box	6	25 (codes ending in "L") 50 (all other codes)	50